# Vitassay qPCR

## Macrolide Res- 23S rRNA Mutations

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- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>I</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- MIC qPCR Cycler (Bio Molecular Systems)
- LightCycler 480 Instrument II (Roche)
- NEOS-96 Real Time PCR System (Linear)
- QuantStudio™ 5 Real-Time PCR System (Applied Biosystems)
- QuantGene 9600 (BioEr Technology)
- Rotor-Gene® Q (Qiagen) <sup>II</sup>

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## Intended use

Vitassay qPCR Macrolide Res- 23S rRNA Mutations allows the qualitative detection and identification of specific point mutations (conferred by base substitutions in 23S rRNA) implicated in macrolide resistance of *M. genitalium* by Real-Time PCR in genital swabs from patients with confirmed *Mycoplasma genitalium* infection by molecular assays. This product is intended to aid in the potential resistances to macrolide diagnosis, alongside the patient's clinical data, epidemiological risk factors, and other laboratory tests outcomes.

## References

Vitassay qPCR Macrolide Res- 23S rRNA Mutations

4x8-well strip, low profile

7041077

Vitassay qPCR Macrolide Res- 23S rRNA Mutations

4x8-well strip, high profile

7042077

## Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S077/ 7042S077	Macrolide Res- 23S strips low/high profile	-	4 x 8-well strip
7C077	Macrolide Res- 23S Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Additional equipment and material required

- Collection and transport system.

- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

## Summary

Sexually transmitted diseases (STDs) are very common and are transmitted from one person to another through intimate physical contact and sexual activity that includes vaginal, oral, and anal sex. STDs can mostly be prevented by not having sex. STDs do not always cause symptoms, so it is possible to have an infection and not know it.

Over the past few decades, *Mycoplasma genitalium* has become a STD pathogen. The bacterium is a slow growing organism with a genome size of only 580 Kb coding for fewer than 500 genes.

*Mycoplasma* species are found in the mouth, the respiratory and genitourinary tract. *M. genitalium* specifically colonizes the male and female reproductive tract in organs such as the urethra, cervix and fallopian tubes.

Between 40 and 75% of female carriers of STDs lack of specific symptoms. On the other hand, this organism has been linked to a many severe diseases such as cervicitis (vaginal pruritus and discharge, dysuria, irregular bleeding, among others), pelvic inflammatory disease (abdominal and pelvic pain, and abnormal vaginal discharge/bleeding, among others) and urethritis (dysuria, urethral pruritus, and mucopurulent urethral discharge).

Its average prevalence is between 0.5 and 10% in the general population and between 20 and 40% in people with sexually transmitted infections.

Antimicrobial resistance in *M. genitalium* can develop as a result of selective pressure from antibiotic exposure of the organism leading to mutational changes during treatment, or of selection of pre-existing mutations already present. Since no cell wall exists in *Mycoplasmas*, including *M. genitalium*, β-lactam antibiotics and other antibiotics that react with this structure are ineffective. Only a select number of antibiotics such as macrolides are effective against *Mycoplasmas*.

The first-line therapy recommended against this pathogen has been the macrolide azithromycin. However, due to the widespread use of this antibiotic for the treatment of STIs, particularly as a single dose for *Chlamydia trachomatis* and *Neisseria gonorrhoeae*, has probably led to the emergence of macrolide resistance in *M. genitalium* worldwide.

Macrolide resistance is a critical growing problem, with resistance rates currently estimated at 30-100% worldwide. Resistance can develop with a single base mutation, most frequently at position A2058 or A2059 (*Escherichia coli* numbering) of the 23S rRNA, without affecting the lifespan of the organism allowing its transmission.

The prevalence of macrolide resistance has increased rapidly in Europe, particularly in Northern countries (except Sweden) where estimates exceed 50%.

The polymerase chain reaction (PCR) is the method of choice for diagnosis due to the challenges involved in culturing the pathogen and the lack of standardized serological tests for *M. genitalium*. This technique also allows the detection of macrolide resistance due to mutation in the *M. genitalium* gene coding for 23S rRNA. *M. genitalium* and its most common mutations can be simultaneously detected using new gene amplification technologies.

### **Principle of the test**

Vitassay qPCR Macrolide Res- 23S rRNA Mutations is based on the real-time amplification of a conserved region of 23S rRNA gene in which specific point mutations could occur (conferred by base substitutions) implicated in macrolide resistance of *Mycoplasma genitalium*. After DNA extraction, the *M. genitalium* resistant/sensitive to macrolides presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR Macrolide Res- 23S rRNA Mutations is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible inhibition reaction. The amplification of the target sequence is detected through the FAM channel (*M. genitalium* resistant to macrolides) and HEX, VIC, or JOE channel, depending on the equipment used, (*M. genitalium* sensitive to macrolides), whereas the internal control (IC) is detected through the Cy5.

### **Precautions**

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.

- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.
- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## **Procedures**

### **Specimen collection, transport, and conservation**

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR Macrolide Res- 23S rRNA Mutations has been tested in urethral, rectal, endocervical and vaginal swabs collected with a sterile single-use Copan Liquid Amies Elution Swab (ESwab®) with flocked nylon filled with 1 mL of Liquid Amies preservation medium (Copan), as well as in urethral urines collected in sterile tubes

(VACUETTE™, Greiner Bio-One) from patients with confirmed *Mycoplasma genitalium* infection. Other sample types must be validated by the user.

Overall, clinical samples should be collected appropriately in clean containers, correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The urine specimens should be transported at room temperature for up to 1 hour, or at 2-8°C for up to 24 hours. For long term transport (more than 24 hours), we recommend shipping at -20°C or lower. Samples in ESwab® Collection and Preservation System should be transport at 20-25°C for up to 5 days, or at 4°C for up to 7 days. For long term transport (more than 7 days), we recommend shipping at -20°C or lower. The samples can be stored at 2-8°C for up to 24 hours or frozen at -20°C or lower (at -80°C ideally) for long-term conservation, if preservation medium permits. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- IDSA (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94)
- García-Lechuz Moya, J.M., González López, J.J., Orta Mira, N., Sánchez Romero, M.I. (2017). Recogida, transporte y procesamiento general de las muestras en el Laboratorio de Microbiología. 2017. 1b. Sánchez Romero, M.I., (coordinadora). Procedimientos en Microbiología Clínica. Cercenado Mansilla, E., Cantón Moreno, R., (editores). *Sociedad Española de Enfermedades Infecciosas y Microbiología Clínica (SEIMC)*.

## DNA extraction

For nucleic acid isolation from clinical specimens, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

- Invisorb® Spin Universal Kit (Invitex).
- DNA/RNA Extraction Kit (Vazyme).
- MagDEA Dx SV kit, empleando el equipo magLEAD® 12gC instrument (Precision System Science Co.).
- STARMag 96x4 Universal Cartridge Kit (Seegene®), empleando el Microlab® STARlet automatic extraction system (Hamilton).

## Positive control preparation

Reconstitute the lyophilized Macrolide Res- 23S Positive Control (red tube) with 100  $\mu$ L of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

## Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15  $\mu$ L of Resuspension buffer (green tube) and add them into each well.
- Pipette 5  $\mu$ L of DNA sample, negative control (yellow tube) or reconstituted positive control (red tube) and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Polymerase activation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (*M. genitalium* macrolide resistant) and HEX, VIC, or JOE (*M. genitalium* macrolide sensitive) channels, and Cy5 (Internal Control, IC). Depending on the equipment used select the proper detection channel (Attached II).

## **Analysis and interpretation of results**

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, check the internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

### **Positive control**

The positive control used in each run must show an amplification curve ( $Ct \leq 40$ ) in FAM and HEX channel, which validates the reaction.

### **Negative control**

The negative control included in each run must show signal' absence ( $Ct > 40$  or no signal) in FAM and HEX channel, which validates the reaction.

The Internal Control (IC) should show an amplification signal ( $Ct \leq 40$ ) in positive and negative controls wells.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

<i>M. genitalium</i> macrolide resistant (FAM)	<i>M. genitalium</i> macrolide sensitive (HEX)	Internal Control (Cy5)		Interpretation
+	-	+/- <sup>1</sup>	Valid	<i>M. genitalium</i> macrolide resistant DNA detected
-	+	+/- <sup>1</sup>	Valid	<i>M. genitalium</i> macrolide sensitive DNA detected
+	+	+/- <sup>1</sup>	Valid	<i>M. genitalium</i> macrolide resistant and <i>M. genitalium</i> macrolide sensitive DNA detected
-	-	+ <sup>2</sup>	Valid	Targets not detected
-	-	- <sup>2</sup>	Invalid	Test failed <sup>2</sup>

**(+) Positive:** Amplification signal (Ct ≤40)

**(-) Negative:** No amplification signal (Ct >40 or no signal)

<sup>1</sup> Sometimes, the Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

<sup>2</sup> In the case of negative *M. genitalium* macrolide resistant and *M. genitalium* macrolide sensitive target genes detection, IC must show an amplification signal with Ct ≤35. If there is a signal absence or Ct value > 35 of Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

### Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

The Vitassay qPCR Macrolide Res- 23S rRNA Mutations clinical evaluation was performed using a total of 169 urogenital samples (76 urethral, 15 rectal, 60 endocervical, and 6 vaginal swabs and 12 urethral urines), all of which were positive for *Mycoplasma genitalium*. DNA extractions were performed with the STARMag 96x4 Universal Cartridge Kit (Seegene Inc.), using the Microlab® STARlet automatic extraction system (Hamilton) and the thermocycler used was DTprime Real-Time Detection Thermal Cycler (DNA-Technology). The obtained results were compared with the ones from the initial characterization (ResistancePlus® MG Kit (SpeeDx Pty Ltd) used for identification of macrolide resistance in *M. genitalium* positive samples and Allplex™ STI essential assay panel (Seegene Inc.) used for discordant results).

Comparative analysis between the initial characterization and Vitassay qPCR, and after assessing the discrepancies, showed that Vitassay was able to detect 80 true positive and 87 true negative samples for 23S-point mutation (i.e. macrolide resistance).

On the other hand, Vitassay kit reported 87 true positive and 80 true negative samples for wildtype gene (WT) (i.e macrolide sensitivity). In both cases, two discrepancies were found: 2 false negative values for the 23S point mutation and 2 false positive values for WT. Hence, specificity and sensitivity values for 23S point mutation were 1 (0.96-1) and 0.98 (0.91-0.99), respectively, and for WT were 0.98 (0.91-0.99) and 1 (0.96-1), respectively.

In conclusion, Vitassay qPCR Macrolide Res- 23S rRNA Mutations represents an efficient tool for the diagnostic of macrolide resistance-associated mutations in *M. genitalium* positive urogenital samples, as the results show high concordance.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of *M. genitalium* ranging from  $10^7$  to  $10^1$  copies/reaction. This assay has a detection limit of 8 DNA copies/ $\mu$ L for *M. genitalium* resistant to macrolides and 0.2 CFU/ $\mu$ l for *M. genitalium* sensitive to macrolides on vaginal swabs samples (positive rate  $\geq 95\%$ ).

### Analytical specificity

The analytical specificity for the *M. genitalium* detection was tested within the panel of following microorganisms associated to sexually transmitted diseases and with antimicrobial resistance, where no cross-reactivity was observed with any of the species:

Cross-reactivity testing					
<i>Acinetobacter baumannii</i>	-	<i>Enterococcus faecalis</i> resistente a vancomicina A	-	<i>Klebsiella pneumoniae</i> ESBL/ positivo a carbapenemasa (TEM-1(non-ESBL), SHV-1(no-ESBL), CTX-M-2(ESBL), KPC-2)	-
<i>Aspergillus fumigatus</i>	-	<i>Enterococcus faecalis</i> resistente a vancomicina B (equivalente a ATCC® 51299™)*	-	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> serotipo <i>Capsular 2</i>	-
<i>Atopobium vaginae</i>	-	<i>Enterococcus faecium</i> serotipo 11	-	<i>Listeria innocua</i> serotipo 6a	-
<i>Bacteroides fragilis</i>	-	<i>Enterococcus faecium</i> resistente a vancomicina A (LMG16165)	-	<i>Listeria ivanovii</i> subsp. <i>ivanovii</i> serovar 5	-
<i>Candida albicans</i>	-	<i>Enterococcus faecium</i> resistente a vancomicina A (IOWA 1)	-	<i>Listeria monocytogenes</i> serovar 4b	-
<i>Candida dubliniensis</i>	-	<i>Enterococcus faecium</i> resistente a vancomicina B (IOWA 2)	-	<i>Mycoplasma hominis</i>	-
<i>Candida glabrata</i>	-	<i>Enterococcus faecium</i> S (MI12043391) resistente a vancomicina C	-	<i>Neisseria gonorrhoeae</i> cepas <i>St 49226</i> y <i>Lvl Ng PorA</i>	-
<i>Candida krusei</i> ( <i>Issatchenkia orientalis</i> )	-	<i>Enterococcus gallinarum</i> (LMG16289) resistente a vancomicina C	-	<i>Neisseria meningitidis</i> serogrupo A	-
<i>Candida parapsilosis</i>	-	<i>Enterococcus gallinarum</i> (ENT20120142) resistente a vancomicina B y C	-	<i>Proteus mirabilis</i>	-
<i>Candida tropicalis</i>	-	<i>Escherichia coli</i> positivo a carbapenemasa (OXA-244)	-	<i>Pseudomonas aeruginosa</i>	-
<i>Chlamydia trachomatis</i> cepa SW	-	<i>Escherichia coli</i> positivo a carbapenemasa (TEM-1(no-ESBL), IMP-1)	-	<i>Serratia marcescens</i> subsp. <i>marcescens</i>	-
<i>Chlamydia trachomatis</i> genovar F	-	<i>Escherichia coli</i> cepa 0.1285; O18:H7:K1	-	<i>Serratia marcescens</i> positivo a carbapenemasa (OXA-48)	-
<i>Chlamydia trachomatis</i> cepa LGV	-	<i>Gardnerella vaginalis</i>	-	<i>Staphylococcus aureus</i> resistente a metilicina	-

Cross-reactivity testing				
				(variedades N315, ST398 y mecC)
<i>Citrobacter brakii</i> positivo a carbapenemasa (VIM-1)	-	<i>Haemophilus ducreyi</i> cepa Class 1	-	<i>Staphylococcus aureus</i> resistente a meticilina adquirida en la comunidad (oxa <sup>R</sup> , PVL-positivo, spa:t 310)
Complejo <i>Citrobacter freundii</i> positivo a carbapenemasa (KPC-3, VIM-4)	-	<i>Haemophilus influenzae</i> cepa Minn A	-	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
Cytomegalovirus cepa AD-169	-	<i>Helicobacter pylori</i> resistente a claritromicina (A2146G y A2147G)	-	<i>Stenotrophomonas maltophilia</i>
<i>Enterobacter aerogenes</i> serotipo Cloaca B	-	Hepatitis A virus	-	<i>Streptococcus agalactiae</i> cepa Z019
<i>Enterobacter cloacae</i> ESBL/ positivo a carbapenemasa (SHV-12 (ESBL), CTX-M-9 (ESBL), OXA-48)	-	Herpes Simplex Virus 1 cepa MacIntyre	-	<i>Streptococcus pneumoniae</i> cepa Z022
<i>Enterobacter cloacae</i> ESBL/ positivo a carbapenemasa (TEM-1 (non ESBL), SHV-12 (ESBL), CTX-M-15 (ESBL), NDM-1)	-	Herpes Simplex Virus 2 cepa MS	-	<i>Treponema pallidum</i>
Complejo <i>Enterobacter cloacae</i> positivo a carbapenemasa (NDM-7)	-	Human Papillomavirus tipos 16 y 18	-	<i>Trichomonas vaginalis</i>
<i>Enterobacter cloacae</i> serotipo Cloaca A	-	<i>Klebsiella oxytoca</i>	-	<i>Ureaplasma parvum</i>
<i>Enterococcus faecalis</i>	-	<i>Klebsiella pneumoniae</i> positivo a carbapenemasa (SHV-1(no-ESBL), KPC-3, OXA-48)	-	<i>Ureaplasma urealyticum</i>

\*ATCC® 51299™: *Enterococcus faecalis* (Andrewes and Horder) Schleifer and Klipper-Balz, strain NJ-3.

### Analytical reactivity

The Vitassay qPCR Macrolide Res- 23S rRNA Mutations kit's reactivity for *M. genitalium* macrolide-sensitive was evaluated against DNA extracted (as template) from:

*Mycoplasma genitalium* Tully et al., strain M30

Simulated clinical samples positive for *Mycoplasma genitalium*

Simulated clinical samples positive for *Mycoplasma genitalium* strain G37

showing positive results.

The Vitassay qPCR Macrolide Res- 23S rRNA Mutations kit's reactivity for *M. genitalium* macrolide-resistant was evaluated against DNA extracted (as template) from:

*Clinical samples*

Simulated clinical samples positive for *Mycoplasma genitalium* strain M6303

Simulated clinical samples positive for *Mycoplasma genitalium* strain G6593

showing positive results.

### **Compatibles real-time PCR equipment**

Vitassay qPCR Macrolide Res- 23S rRNA Mutations has been validated on the following equipment:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>I</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- MIC qPCR Cyclor (Bio Molecular Systems)
- LightCycler 480 Instrument II (Roche)
- NEOS-96 Real Time PCR System (Linear)
- QuantStudio™ 5 Real-Time PCR System (Applied Biosystems)
- QuantGene 9600 (BioEr Technology)
- Rotor-Gene® Q (Qiagen) <sup>II</sup>

<sup>I</sup>: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

<sup>II</sup>: For Rotor-Gene® Q and MIC qPCR Cyclor the product should be reconstituted following the procedure and transferred into specific tubes of the equipment.

\*If the result in channel FAM is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

## Limitations

- All obtained results must be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from genital swabs (urethral urine, urethral swabs, endocervical swabs, rectal swabs, and vaginal swabs) positive for *M. genitalium* infection by molecular assays. The use of other samples has not been established.
- The correct test performance depends on the sample's quality; proper DNA from clinical specimens must be extracted.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination either by Macrolide Res- 23S Control during its reconstitution, by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including DNA extraction); c) DNA degradation during specimen shipping, storage, and/or preparation; d) pathogen load below the limit of detection for the assay; e) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures; g) mutations or polymorphisms in primer or probe binding regions may affect detection of *M. genitalium* resistant to macrolides and/or *M. genitalium* sensitive to macrolides.
- DNA detection may not indicate the presence of viable and/or infectious bacteria or that *M. genitalium* resistant to macrolides and/or *M. genitalium* sensitive to macrolides is the causative agent for clinical symptoms.
- The negative results do not preclude *M. genitalium* resistant to macrolides and/or *M. genitalium* sensitive to macrolides infection and should not be the sole basis of a patient treatment/management decision. Optimum specimen types for *M. genitalium* resistant to macrolides and/or *M. genitalium* sensitive to macrolides identification and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the bacteria.

- If the patient clinical data, laboratory tests and epidemiological studies suggest possible *M. genitalium* resistant to macrolides and/or *M. genitalium* sensitive to macrolides infection, and other Sexually Transmitted Diseases (STD) have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type, and its previous treatment, among others.

## Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 <sup>(1)</sup>
AriaDx Real-Time PCR System	<b>Applied Biosystems</b>
<b>Applied Biosystems</b>	7300 Real-Time PCR System <sup>(3) (1)</sup>
7500 Fast Real-Time PCR System <sup>(1) (2)</sup>	7500 Real-Time PCR System <sup>(1)</sup>
7500 Fast Dx Real-Time PCR System <sup>(1) (2)</sup>	7900 Real-Time PCR System <sup>(3)</sup>
QuantStudio™ 12K Flex 96-well Fast	ABI PRISM 7000 <sup>(3)</sup>
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7700 <sup>(3)</sup>
QuantStudio™ 7 Flex 96-well Fast	QuantStudio™ 12K Flex 96-well
QuantStudio™ 5 Real-Time PCR System	QuantStudio™ 6 Flex 96-well
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 7 Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System <sup>(3)</sup>	QuantStudio™ 5 Fast Real-Time PCR System
StepOne Plus™ Real-Time PCR System <sup>(3)</sup>	QuantStudio™ 5 Real-Time PCR System
StepOne™ Real-Time PCR System <sup>(4)</sup>	QuantStudio™ 3 Real-Time PCR System <sup>(3)</sup>
ViiA™ 7 Fast Real-Time PCR System	ViiA™ 7 Real-Time PCR System
<b>Azure Biosystems</b>	<b>Agilent Technologies</b>
Azure Cielo 3 <sup>(5)</sup>	Mx3000P™ Real Time PCR System
Azure Cielo 6	Mx3005P™ Real Time PCR System
<b>BIONEER</b>	<b>Analytik Jena</b>
Exicycler™ 96	qTOWER
<b>Bio-Rad</b>	<b>BIONEER</b>
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
Mini Opticon™ Real-Time PCR Detection System <sup>(5)</sup>	<b>BIOER</b>
<b>Roche</b>	QuantGene 9600
LightCycler®480 Real-Time PCR System <sup>(1) (6)</sup>	<b>Bio-Rad</b>
LightCycler®96 Real-Time PCR System <sup>(1)</sup>	CFX96™ Deep Well Real-Time PCR
Cobas z480 Analyzer <sup>(1) (6)</sup>	iCycler iQ™ Real-Time PCR
	iCycler iQ™5 Real-Time PCR
	My iQ™ Real-Time PCR Detection System <sup>(5)</sup>
	My iQ™ 2 Real-Time PCR Detection System <sup>(5)</sup>
<b>Special Formats <sup>(7)</sup></b>	<b>DNA-Technology</b>
<b>Bio Molecular Systems</b>	DTlite Real-Time PCR System <sup>(8)</sup>
Mic Real Time PCR Cyclor	DTprime Real-time Detection Thermal Cyclor <sup>(8)</sup>
<b>Cepheid</b>	<b>Eppendorf</b>
SmartCycler®	Mastercycler™ ep <i>realplex</i>
<b>Precision System Science Co., Ltd.</b>	<b>Qiagen</b>
geneLEAD VIII System	QIAquant 96
<b>Qiagen</b>	
Rotor-Gene® Q	

- (1) A special grid is needed to fit these real-time PCR kits.
- (2) Select Ramp Speed "Standard".
- (3) No ROX caption.
- (4) No Cy5 caption.
- (5) Only FAM and HEX caption.
- (6) Specific compensation color is required.
- (7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, Rotor-Gene® Q or geneLEAD VIII System.
- (8) See Attached III to configure exposure settings.

## Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
<b>Bio-Rad CFX96™</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Applied Biosystems ABI 7500</b>	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Roche Cobas z 480</b>	FAM	465/510	Colour Compensation required
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
<b>Cepheid Smartcycler®</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Agilent Technologies Mx3000P™ Mx 3005P™</b>	FAM	FAM	Passive reference option ROX is none
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Bio Molecular Systems Mic Real Time PCR Cycler</b>	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Agilent Technologies AriaMx</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Qiagen Rotor-Gene® Q</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>BIONEER Exicycler™ 96</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

### Attached III. Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
<b>DTIite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

\*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

## Bibliography/Bibliografía

1. Fernández-Huerta, M., Barberá, M. J., Serra-Pladevall, J., Esperalba, J., Martínez-Gómez, X., Centeno, C., Pich, O. Q., Pumarola, T., & Espasa, M. (2020). *Mycoplasma genitalium* and antimicrobial resistance in Europe: a comprehensive review. *International journal of STD & AIDS*, 31(3), 190–197.
2. Gnanadurai, R., & Fifer, H. (2020). *Mycoplasma genitalium*: A Review. *Microbiology (Reading, England)*, 166(1), 21–29.
3. Heavey E. (2017). *Mycoplasma genitalium*. *Nursing*, 47(7), 61–62.
4. De Carvalho, N. S., Palú, G., & Witkin, S. S. (2020). *Mycoplasma genitalium*, a stealth female reproductive tract. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology*, 39(2), 229–234.

## Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

## Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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