**VA Vitassay**

Real-Time PCR Kits

## Vitassay qPCR MonkeyPox

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Sample	TP	TN	FP	FN	SE	SP	PPV	NPV	LR+	LR-	OA%
<b>Monkey Pox</b>	40	92	0	0	1 (0.95-1)	1 (0.98-1)	1 (0.90-1)	1 (0.95-1)	183.7 (11.57-2916)	0.012 (0.001-0.193)	100 (97.2-100)

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Pruebas de reactividad cruzada					
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Respiratory syncytial virus (RSV) A and B	-	<i>Aspergillus fumigatus</i>	-
Bocavirus	-	<i>Anaplasma marginale</i>	-	<i>Atopobium vaginae</i>	-
<i>Bordetella bronchiseptica</i>	-	<i>Bartonella henselae</i> strain Houston-1	-	<i>Bacteroides fragilis</i>	-
<i>Bordetella holmesii</i>	-	<i>Borrelia hermsii</i>	-	<i>Candida albicans</i>	-
<i>Bordetella parapertussis</i>	-	<i>Borrelia lusitanae</i>	-	<i>Candida dubliniensis</i>	-
<i>Bordetella pertussis</i>	-	<i>Borrelia valaisiana</i>	-	<i>Candida glabrata</i>	-
<i>Chlamydia caviae</i>	-	<i>Borrelia azfeli</i> strain P-Ko/1984	-	<i>Candida krusei</i>	-
<i>Chlamydia psittaci</i> genotypes A and C	-	<i>Borrelia bavariensis</i>	-	<i>Candida parapsilosis</i>	-
<i>Chlamydophila pneumoniae</i> CM-1	-	<i>Borrelia bisetti</i>	-	<i>Candida tropicalis</i>	-
Human coronavirus 229E, OC43, NL63 and HKU1	-	<i>Borrelia burgdorferi</i> strain IRS	-	<i>Chlamydia trachomatis</i> SW	-
MERS-Coronavirus	-	<i>Borrelia burgdorferi sensu stricto</i> strain B31	-	<i>Chlamydia trachomatis</i> Genovar F	-
SARS-CoV-2 strain 2019nCoV/USAWA1/2020	-	<i>Borrelia garinii</i> Type strain	-	<i>Chlamydia trachomatis</i> LGV	-
Enterovirus 68 and 71	-	<i>Borrelia japonica</i>	-	Cytomegalovirus AD-169	-
Enterovirus Echovirus 30	-	<i>Borrelia miyamotoi</i>	-	<i>Enterobacter aerogenes</i> Serotype Cloaca B	-
Enterovirus Coxsackievirus A24, A9 and B3	-	<i>Borrelia spielmanii</i>	-	<i>Enterobacter cloacae</i> Serotype Cloaca A	-
<i>Haemophilus influenzae</i> MinnA	-	Chikungunya virus strain F24	-	<i>Enterococcus faecalis</i>	-
Influenza A/New Caledonia/20/99(H1N1) virus	-	Chikungunya virus strain Martinique	-	<i>Enterococcus faecium</i> Serotype 11	-
Influenza A/California/7/2009(H1N1)pdm09 virus	-	Chikungunya virus strain S27 Petersfield	-	<i>Escherichia coli</i> 0.1285;O18:H7:K1	-
Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	Chikungunya virus WHO IS 11785/16 (from isolate R91064)	-	<i>Gardnerella vaginalis</i>	-
Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Coxiella burnetii</i> strain Nine Mile Q	-	<i>Haemophilus ducreyi</i> class 1	-
Influenza A/Victoria/210/2009 (H3N2) virus	-	Dengue virus serotype 1, strain Hawaii A	-	Hepatitis A	-
Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1)	-	Dengue virus serotype 2, strain New Guinea C	-	Herpes simplex virus 1(HSV-1) Strain <i>MacIntyre</i>	-
Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	Dengue virus serotype 3, strain H87	-	Herpes simplex virus 2(HSV-2) MS	-
Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	Dengue virus serotype 4, strain H241	-	Human DNA (male)	-
Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	Japanese Encephalitis virus strain Nakayama	-	Human papillomavirus 16 and 18	-
Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	<i>Leptospira</i>	-	<i>Klebsiella oxytoca</i>	-
Influenza A/Anhui/1/2013 (H7N9) virus	-	<i>Rickettsia conorii</i> strain Moroccan	-	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> serotype capsular 2	-
Influenza B/Brisbane/60/2008 virus	-	Rift Valley Fever Virus strain AR21229	-	<i>Listeria innocua</i> serotype 6a	-
Influenza B/Florida/04/06 virus	-	Rift Valley Fever Virus strain MP12	-	<i>Listeria ivanovii</i> subsp. <i>Ivanovii</i> serovar 5	-
Influenza B/Phuket/3073/2013 virus	-	St Louis Encephalitis virus	-	<i>Listeria monocytogenes</i> Serovar 4b	-

<i>Legionella bozemanii</i> serovar 1	-	<i>Plasmodium falciparum</i> 3D7	-	<i>Mycoplasma genitalium</i>	-
<i>Legionella dumoffii</i>	-	<i>Theileria annulata</i>	-	<i>Mycoplasma hominis</i>	-
<i>Legionella longbeachae</i>	-	Tick-Borne Encephalitis virus strain Neudorfl	-	<i>Neisseria gonorrhoeae</i> St 49226	-
<i>Legionella micdadei</i>	-	<i>Treponema phagedenis</i>	-	<i>Neisseria gonorrhoeae</i> Lvl Ng PorA	-
<i>Legionella pneumophila</i> serogroup 1	-	<i>Trypanosoma cruzi</i>	-	<i>Neisseria meningitidis</i> Serogroup A	-
Human metapneumovirus A and B	-	Usutu virus	-	<i>Proteus mirabilis</i>	-
<i>Moraxella catarrhalis</i>	-	West Nile virus strain NY99	-	<i>Pseudomonas aeruginosa</i>	-
<i>Mycoplasma pneumoniae</i>	-	West Nile virus strain Heja	-	<i>Serratia marcescens</i> subsp. <i>marcescens</i>	-
<i>Mycobacterium tuberculosis</i> not rifampin resistant	-	West Nile virus strain Ug37	-	<i>Stenotrophomonas maltophilia</i>	-
Human parainfluenza 1, 2, 3 and 4 viruses	-	Yellow Fever virus strain 17D	-	<i>Streptococcus agalactiae</i> Z019	-
<i>Pneumocytis jirovecii</i> Type A1 and g885652	-	Yellow Fever virus strain French Neurotropic	-	<i>Treponema pallidum</i>	-
Human rhinovirus type C	-	Zika virus Asian strain PF13/251013-18	-	<i>Trichomonas vaginalis</i>	-
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Zika virus African strain	-	<i>Ureaplasma parvum</i>	-
<i>Staphylococcus epidermidis</i>	-	Zika virus strain FB-GWUH-2016	-	<i>Ureaplasma urealyticum</i>	-
<i>Streptococcus pneumoniae</i> Z022	-	Zika Virus French Polynesian, strains 11468/16 and 11474/16	-	Varicella Zoster Virus OKA	-
<i>Streptococcus pyogenes</i>	-	<i>Acinetobacter baumannii</i>	-	Varicella Zoster Virus Ellen	-
<i>Streptococcus salivarius</i>	-				

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 <sup>(1)</sup>
AriaDx Real-Time PCR System	<b>Agilent Technologies</b>
<b>Applied Biosystems</b>	Mx3000P™ Real Time PCR System
7500 Fast Real-Time PCR System <sup>(1) (2)</sup>	Mx3005P™ Real Time PCR System
7500 Fast Dx Real-Time PCR System <sup>(1) (2)</sup>	<b>Applied Biosystems</b>
QuantStudio™ 12K Flex 96-well Fast	7300 Real-Time PCR System <sup>(1) (4)</sup>
QuantStudio™ 6 Flex 96-well Fast	7500 Real-Time PCR System <sup>(1)</sup>
QuantStudio™ 7 Flex 96-well Fast	7900 HT Real-Time PCR System <sup>(3)</sup>
QuantStudio™ 5 Real-Time PCR System	ABI PRISM 7000 <sup>(3)</sup>
QuantStudio™ 5 Fast Real-Time PCR System	ABI PRISM 7700 <sup>(3)</sup>
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
StepOne Plus™ Real-Time PCR System <sup>(3)</sup>	QuantStudio™ 6 Flex 96-well
StepOne™ Real-Time PCR System <sup>(3) (4)</sup>	QuantStudio™ 7 Flex 96-well
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 5 Fast Real-Time PCR System
<b>Azure Biosystems</b>	QuantStudio™ 5 Real-Time PCR System
Azure Cielo 3 <sup>(5)</sup>	QuantStudio™ 3 Real-Time PCR System
Azure Cielo 6	ViiA™ 7 Real-Time PCR System
<b>BIONEER</b>	<b>Analytik Jena</b>
Exicycler™ 96 Fast	qTOWER <sup>(6)</sup>
<b>Bio-Rad</b>	<b>BIONEER</b>
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
CFX96™ IVD Real-Time PCR Detection System	<b>BIOER</b>
Mini Opticon™ Real-Time PCR Detection System <sup>(5)</sup>	QuantGene 9600
CFX Opus 96	<b>Bio-Rad</b>
<b>Roche</b>	CFX96™ Deep Well Real-Time PCR System
LightCycler®480 Real-Time PCR System <sup>(1) (6)</sup>	CFX96™ Deep Well IDV Real-Time PCR System
LightCycler®96 Real-Time PCR System	iCycler iQ™ Real-Time PCR Detection System
Cobas z480 Analyzer <sup>(1) (6)</sup>	iCycler iQ™5 Real-Time PCR Detection System
	My iQ™ Real-Time PCR Detection System <sup>(5)</sup>
	My iQ™ 2 Real-Time PCR Detection System <sup>(5)</sup>
<b>Formatos especiales <sup>(7)</sup></b>	<b>DNA-Technology</b>
<b>Bio Molecular Systems</b>	DTlite Real-Time PCR System <sup>(8)</sup>
Mic Real Time PCR Cyclor	DTprime Real-time Detection Thermal Cyclor <sup>(8)</sup>
<b>Cepheid</b>	<b>Eppendorf</b>
SmartCycler®	Mastercycler™ ep <i>realplex</i>
<b>Qiagen</b>	<b>Qiagen</b>
Rotor-Gene® Q	QIAquant 96 <sup>(6)</sup>

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## Intended use

Vitassay qPCR MonkeyPox allows the qualitative detection of genomic DNA specific for Monkeypox virus by real-time PCR in lesion swabs of patients suspected of Monkeypox virus infection. This product is intended to aid in the diagnosis of Monkeypox virus infections, alongside the patient's clinical data, epidemiological risk factors and other laboratory tests outcomes.

## References

Vitassay qPCR MonkeyPox 4x8-well strip, low profile	7041080
Vitassay qPCR MonkeyPox 4x8-well strip, high profile	7042080

## Reagents provided

In references 7041080 and 7042080:

Reference	Reagent/Material	Colour	Amount
7041S080/ 7042S080	MonkeyPox strips low/high profile	-	4 x 8-well strip
7C080	MonkeyPox Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8-cap strip

## Material and equipment required, not provided

- Collection and transport system.
- Laboratory freezers (- 30°C to - 10°C and/or ≤ -70°C).
- DNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortex
- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

## Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated on the label.
- The resuspended positive control should be stored at -20°C. To avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

## Summary

Human monkeypox virus, *Monkeypox* o MPXV, is one of the four viruses in the *Orthopoxvirus* genus of the *Poxviridae* family with capacity of infecting humans. The other three are the virus causing smallpox (VARV, or *variola major* virus), the *variola minor* virus, and the cowpox virus. This zoonotic double-stranded virus was not recognized as a distinct infection in humans until the 1970s, when the virus was isolated from a 9-month-old patient with smallpox suspicion in the Democratic Republic of Congo.

Most of the clinical manifestations caused by human monkeypox virus infection are indistinguishable to the smallpox. Before the onset of symptoms, there is an incubation period that usually lasts from 7 to 14 days. After that, there is an initial febrile prodrome which is accompanied by generalized headache and fatigue. The symptomatology usually evolves to a maxillary, cervical, or inguinal lymphadenopathy in many patients, prior to a typical rash development after 1 to 3 days after the fever appears. The rash progress

to macules, papules, vesicles, pustules, and scabs before falling off. The symptoms usually last 2 to 4 weeks and the fatality rate has been reported to be between 1-10%, being higher in children.

Since the first description in humans in 1970, the monkey pox has been reported in 11 African countries and different outbreaks have been experienced. Specially in Nigeria, where since 2017 over 500 suspected and 200 confirmed cases have been registered (WHO).

This is a disease of global public health concern. In 2003, the first monkeypox outbreak outside of Africa was in the United States of America and was linked to contact with infected pet prairie dogs. Since 13 May 2022, cases of monkeypox have been reported to WHO from 12 Member States that are not endemic for monkeypox virus, across three WHO regions.

Primary animal-to-human infection is believed to occur when handling monkeypox-infected animals, via direct or indirect contact, although the exact mechanisms are not yet defined. Human-to-human transmission between humans is thought to be through aerosols as well as direct contact with body fluids, lesions on the skin or mucosal tissue.

There are no specific treatments for monkeypox, but symptomatic management and support, as well as drug administration for the secondary bacterial infections are advisable. Furthermore, there is not any vaccine specifically designed against this virus that has been yet licensed. However, it has been observed that the smallpox vaccine grants a cross-protection.

Diagnosis is carried out with specimens from skin lesions as swabs of vesicular lesions, exudate, or crusts. In case of viral cultures, samples must be obtained from the oropharyngeal or nasopharyngeal cavities with a swab, however molecular methods, such as PCR or sequencing are preferred, especially for the differentiation of strains.

## Test principle

Vitassay qPCR MonkeyPox is based on the real-time amplification of a conserved region of *G2R* and *F3L* genes of MonkeyPox virus. After DNA extraction, the MonkeyPox virus presence is detected by an increase in observed fluorescence during the reaction, after hydrolysis of the fluorescent probe.

The assay is based on 5' exonuclease activity using two primers and a fluorogenic hydrolysis probe to detect accumulation of the amplified target sequence during the PCR reaction. When the polymerase begins to spread the primers, the probe is hydrolyzed by its exonuclease 5'-3' activity causing the spatial separation of the fluorophore and the quencher. The increase in the resulting fluorescent signal is proportional to the amount of amplified product in the sample and is detected by means of real-time PCR equipment.

Vitassay qPCR MonkeyPox is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, a set of primers and probe for an endogenous internal control detection allows the extraction process control and the detection of a possible inhibition reaction. The assay uses a human housekeeping gene as an endogenous internal control (IC) (*HGB* gene, hemoglobin subunit beta gene present in human DNA) that is expected to be present in all nucleated human cells. The detection channels of the target sequences are described below:

Target	Detection Channel
<i>G2R</i> and <i>F3L</i> genes	FAM
Endogenous Internal Control (IC)	HEX, VIC, or JOE (depending on the equipment used)

## Precautions

- For professional *in vitro* diagnostic use (use by qualified and trained clinical laboratory personnel).
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or is broken or if the foil has been broken or damaged.
- Do not remove desiccant from reagent pouches once is open.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles, and mask.
- Do not eat, drink, or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment, and reagents to the area in which the previous step was performed. Use separate areas for the patient samples and controls preparation.
- Avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile disposable aerosol resistant pipette tips or RNase/DNase positive displacement pipette tips is recommended.

- An appearance of the reaction mixture in stabilized format, normally found at the tube bottom, different from the usual one (without conical shape, inhomogeneous, smaller/larger in size and/or color different from whitish) does not alter the test functionality.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious and they must be managed according to the national safety legislation and national health waste legislation. Take necessary precautions during the collection, transport, storage, treatment, and disposal of samples.
- Use personal protective equipment (PPE) and biological safety cabinet for handling of potentially infectious samples and reagents according to current recommendations.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

## Procedures

### Specimen collection, transport, and conservation

For specimen collection, conservation, and transport, user-validated conditions must be followed. The Vitassay qPCR MonkeyPox has been tested with lesion swabs. Other sample types must be validated by the user.

Overall, clinical samples should be collected appropriately in clean containers with or without transport media (depending on sample type), correctly labelled and promptly processed to ensure the test's quality. It is recommended to use fresh specimens for the test. Transport must always be carried out in accordance with local and national regulations for the transport of biological samples. The specimens should be transported at 2 to 8°C within an hour of collection. For long term transport (more than 1 hour), we recommend shipping at -20°C or lower. The samples can be stored at 2 to 8°C or frozen at -20°C or lower (at -80°C ideally) for long-term conservation. Repeated freeze-thaw cycles should be avoided to prevent sample and nucleic acids degradation.

The clinical specimens must be collected, transported, and stored according to appropriate laboratory guidelines. For more information, refer to the following guidelines:

- CDC (Specimen collection guidelines. Website: <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf>)
- IDSA (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94)

### DNA extraction

For nucleic acid isolation from clinical specimens, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. For the sample's pretreatment follow the instructions for use of the extraction kit used. The assay has been validated with the following extraction kits:

- MagDEA Dx SV Kit, using the magLEAD® 12gC instrument (Precision System Science Co.).
- MagCore® Viral Nucleic Acid Extraction Kit (Low PCR Inhibition) using the MagCore® System (RBC bioscience).

### Positive control preparation

Reconstitute the lyophilized MonkeyPox Positive Control (red tube) with 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After first use, dispense into aliquots to avoid multiple freeze-thaw cycles, and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend open and manipulate it in a separate laboratory area away from the other components and samples.

### Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run for each assay.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add them into each well.
- Pipette 5 µL of DNA sample, negative control (yellow tube) or reconstituted positive control (red tube) and add them into the corresponding wells.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

## Programme your thermocycler

Set your thermocycler following the conditions below:

Step	Temperature	Time	Cycles
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (\*) through the FAM (Monkeypox virus) and HEX, VIC, or JOE (Endogenous Internal Control) channels. Depending on the equipment used select the proper detection channel (Attached II).

## Analysis and interpretation of results

The results analysis is done by the software itself of the used real-time PCR system following manufacturer's instructions. To verify the correct operation of the amplification mix, and the extraction procedure, check the endogenous internal control (IC) signal emission.

Use the positive control amplification curve as a starting point during reaction validation (prior to interpretation of sample results), to ensure that the threshold falls within the exponential phase of the amplification curves and above any background noise signals. The threshold value may vary between different instruments due to different signal intensities. It is recommended to set the threshold values for each channel (target) independently by the end user.

For a valid diagnostic test run, the following control conditions must be met:

### Positive control

The positive control used in each run must show an amplification curve (Ct ≤40) in FAM and HEX, VIC, or JOE (depending on equipment used) channels, which validates the reaction.

The positive control includes the human *housekeeping HHB* gene target; therefore, amplification signals are observed in all target channels, including the Endogenous Internal Control.

### Negative control

The negative control included in each run must show signal' absence (Ct >40 or no signal) in FAM and HEX, VIC, or JOE (depending on equipment used) channels, which validates the reaction.

In Negative Control, it is possible that non-specific amplifications in late Ct values (Ct > 35) may appear due to human manipulation of the process, since human *HBB* gene is a human housekeeping gene that is expected to be present in all nucleated human cells.

The experiment seems to be failed if there is amplification signal in negative control or signal absence in the positive control for any channel. The assay should be repeated.

The clinical samples test results assessment should be performed once controls' results have been validated. The result interpretation is summarized in the following table:

MonkeyPox (FAM)	Endogenous Internal Control (HEX, VIC, or JOE)	Interpretation	
+	+/- <sup>1</sup>	Valid	MonkeyPox DNA detected
-	+ <sup>2</sup>	Valid	MonkeyPox DNA not Detected <sup>2</sup>
-	- <sup>2</sup>	Invalid	Test failed <sup>2</sup>

(+) **Positive:** Amplification signal (Ct ≤40)

(-) **Negative:** No amplification signal (Ct >40 or no signal)

<sup>1</sup> Sometimes, the Endogenous Internal Control (IC) detection is not necessary since a high copies number of the target can cause a preferential amplification of target-specific nucleic acids. The Endogenous Internal Control (IC) shows or not an amplification signal (Ct ≤40 or no signal).

<sup>2</sup> In the case of negative MonkeyPox virus target genes detection, endogenous IC must show an amplification signal with Ct ≤35. The Ct value could be very variable due to the Endogenous Internal Control is a human *housekeeping* gene that should be present in all human nucleated cells

in the original sample. If there is a signal' absence or Ct value > 35 of Endogenous Internal Control, the result is considered 'Invalid', and retesting is required. It is recommended to repeat the qPCR by diluting the DNA sample 1:10 and/or 1:100, or to re-extract and retest to check for possible failure in the extraction procedure and/or inhibition problems.

If the obtained result is confusing or doubtful, it is necessary to check that all the steps have been carried out correctly, to verify the correct performance of each qPCR steps and to review all the parameters, the sigmoid shape of the curve and the fluorescence intensity. It is also recommended to repeat the assay, preferably in duplicate, depending on the available material (repeat qPCR with the same isolated DNA sample, or re-extract and retest another aliquot of the same specimen or, obtain a new specimen and retest).

The test results must be evaluated by a health professional, together with the medical history, clinical symptoms and/or the results obtained in other diagnostic tests.

## Quality Control

To confirm the appropriate performance of the molecular diagnostic technique, an endogenous Internal Control (IC) is included in each reaction. Moreover, a positive and a negative control must be included in each assay to interpret the results correctly.

## Performance evaluation

### Clinical sensitivity and specificity

To evaluate the clinical sensitivity and specificity of Vitassay qPCR MonkeyPox, lesion swab samples were used. These samples were collected from patients with lesions compatible with Monkeypox virus infection.

These samples were tested using the RealStar Orthopoxvirus PCR kit 1.0 (Altona Diagnostics) and Vitassay qPCR MonkeyPox, their results were compared, and they were confirmed by sequencing. A total of 132 lesion swabs were evaluated. After comparison and discrepancy analysis, 40 samples were positive for Monkeypox virus in both kits and 92 samples were negative. All samples were considered true positives and no discrepancies were found.

Overall, the clinical sensitivity and specificity of Monkeypox virus were 1 (0.95-1) and 1 (0.98-1), respectively.

Sample	TP	TN	FP	FN	SE	SP	PPV	NPV	LR+	LR-	OA%
Monkey Pox	40	92	0	0	1 (0.95-1)	1 (0.98-1)	1 (0.90-1)	1 (0.95-1)	183.7 (11.57-2916)	0.012 (0.001-0.193)	100 (97.2-100)

TP: True Positive; TN: True Negative; FP: False Positive; FN: False Negative; SE: Sensitivity; SP: Specificity; PPV: Positive Predictive Value; NPV: Negative Predictive Value, OA: Overall agreement, LR: Likelihood Ratio.

### Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of a standard from the Monkeypox virus ranging from  $10^7$  to  $10^1$  copies/reaction. This assay has a detection limit of 8 DNA copies per reaction for Monkeypox virus in lesion swabs (positive rate 95%).

### Analytical specificity

The analytical specificity for the Monkeypox virus detection was tested within the panel of following microorganisms, where no cross-reactivity was observed with any of the species:

Pruebas de reactividad cruzada							
Human Adenovirus types 1-5, 8, 15, 31, 40 and 41	-	Respiratory syncytial virus (RSV) A and B	-	<i>Aspergillus fumigatus</i>	-		
Bocavirus	-	<i>Anaplasma marginale</i>	-	<i>Atopobium vaginae</i>	-		
<i>Bordetella bronchiseptica</i>	-	<i>Bartonella henselae</i> strain Houston-1	-	<i>Bacteroides fragilis</i>	-		
<i>Bordetella holmesii</i>	-	<i>Borrelia hermsii</i>	-	<i>Candida albicans</i>	-		
<i>Bordetella parapertussis</i>	-	<i>Borrelia lusitanae</i>	-	<i>Candida dubliniensis</i>	-		
<i>Bordetella pertussis</i>	-	<i>Borrelia valaisiana</i>	-	<i>Candida glabrata</i>	-		
<i>Chlamydia caviae</i>	-	<i>Borrelia azfeli</i> strain P-Ko/1984	-	<i>Candida krusei</i>	-		
<i>Chlamydia psittaci</i> genotypes A and C	-	<i>Borrelia bavariensis</i>	-	<i>Candida parapsilosis</i>	-		
<i>Chlamydia pneumoniae</i> CM-1	-	<i>Borrelia bisetti</i>	-	<i>Candida tropicalis</i>	-		

Human coronavirus 229E, OC43, NL63 and HKU1	-	<i>Borrelia burgdorferi</i> strain IRS	-	<i>Chlamydia trachomatis</i> SW	-
MERS-Coronavirus	-	<i>Borrelia burgdorferi sensu stricto</i> strain B31	-	<i>Chlamydia trachomatis</i> Genovar F	-
SARS-CoV-2 strain 2019nCoV/USAWA1/2020	-	<i>Borrelia garinii</i> Type strain	-	<i>Chlamydia trachomatis</i> LGV	-
Enterovirus 68 and 71	-	<i>Borrelia japonica</i>	-	Cytomegalovirus AD-169	-
Enterovirus Echovirus 30	-	<i>Borrelia miyamotoi</i>	-	<i>Enterobacter aerogenes</i> Serotype Cloaca B	-
Enterovirus Coxsackievirus A24, A9 and B3	-	<i>Borrelia spielmanii</i>	-	<i>Enterobacter cloacae</i> Serotype Cloaca A	-
<i>Haemophilus influenzae</i> MinnA	-	Chikungunya virus strain F24	-	<i>Enterococcus faecalis</i>	-
Influenza A/New Caledonia/20/99(H1N1) virus	-	Chikungunya virus strain Martinique	-	<i>Enterococcus faecium</i> Serotype 11	-
Influenza A/California/7/2009(H1N1)pdm09 virus	-	Chikungunya virus strain S27 Petersfield	-	<i>Escherichia coli</i> 0.1285;O18:H7:K1	-
Influenza A/Michigan/45/2015 (H1N1)pdm09 virus	-	Chikungunya virus WHO IS 11785/16 (from isolate R91064)	-	<i>Gardnerella vaginalis</i>	-
Influenza A/Singapore/GP1908/2015, IVR-180 (H1N1)pdm09 virus	-	<i>Coxiella burnetii</i> strain Nine Mile Q	-	<i>Haemophilus ducreyi</i> class 1	-
Influenza A/Victoria/210/2009 (H3N2) virus	-	Dengue virus serotype 1, strain Hawaii A	-	Hepatitis A	-
Influenza A/Thüringen/5/2017 (H3N2) virus (Clade 3C2a.1)	-	Dengue virus serotype 2, strain New Guinea C	-	Herpes simplex virus 1(HSV-1) Strain <i>MacIntyre</i>	-
Influenza A/Switzerland/9715293/2013 (H3N2) virus	-	Dengue virus serotype 3, strain H87	-	Herpes simplex virus 2(HSV-2) MS	-
Influenza A/Hong Kong/4801/2014, NYMC X-263B (H3N2) virus	-	Dengue virus serotype 4, strain H241	-	Human DNA (male)	-
Influenza A/South Australia/55/2014, IVR-175 (H3N2) virus	-	Japanese Encephalitis virus strain Nakayama	-	Human papillomavirus 16 and 18	-
Influenza A/DE-SH/Reiherente/AR8444/ 2016 (H5N8) virus	-	<i>Leptospira</i>	-	<i>Klebsiella oxytoca</i>	-
Influenza A/Anhui/1/2013 (H7N9) virus	-	<i>Rickettsia conorii</i> strain Moroccan	-	<i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> serotype capsular 2	-
Influenza B/Brisbane/60/2008 virus	-	Rift Valley Fever Virus strain AR21229	-	<i>Listeria innocua</i> serotype 6a	-
Influenza B/Florida/04/06 virus	-	Rift Valley Fever Virus strain MP12	-	<i>Listeria ivanovii</i> subsp. <i>Ivanovii</i> serovar 5	-
Influenza B/Phuket/3073/2013 virus	-	St Louis Encephalitis virus	-	<i>Listeria monocytogenes</i> Serovar 4b	-
<i>Legionella bozemanii</i> serovar 1	-	<i>Plasmodium falciparum</i> 3D7	-	<i>Mycoplasma genitalium</i>	-
<i>Legionella dumoffii</i>	-	<i>Theileria annulata</i>	-	<i>Mycoplasma hominis</i>	-
<i>Legionella longbeachae</i>	-	Tick-Borne Encephalitis virus strain Neudorfl	-	<i>Neisseria gonorrhoeae</i> St 49226	-
<i>Legionella micdadei</i>	-	<i>Treponema phagedenis</i>	-	<i>Neisseria gonorrhoeae</i> Lvl Ng PorA	-
<i>Legionella pneumophila</i> serogroup 1	-	<i>Trypanosoma cruzi</i>	-	<i>Neisseria meningitidis</i> Serogroup A	-
Human metapneumovirus A and B	-	Usutu virus	-	<i>Proteus mirabilis</i>	-
<i>Moraxella catarrhalis</i>	-	West Nile virus strain NY99	-	<i>Pseudomonas aeruginosa</i>	-
<i>Mycoplasma pneumoniae</i>	-	West Nile virus strain Heja	-	<i>Serratia marcescens</i> subsp. <i>marcescens</i>	-
<i>Mycobacterium tuberculosis</i> not rifampin resistant	-	West Nile virus strain Ug37	-	<i>Stenotrophomonas maltophilia</i>	-
Human parainfluenza 1, 2, 3 and 4 viruses	-	Yellow Fever virus strain 17D	-	<i>Streptococcus agalactiae</i> Z019	-
<i>Pneumocystis jirovecii</i> Type A1 and g885652	-	Yellow Fever virus strain French Neurotropic	-	<i>Treponema pallidum</i>	-

Human rhinovirus type C	-	Zika virus Asian strain PF13/251013-18	-	<i>Trichomonas vaginalis</i>	-
<i>Staphylococcus aureus</i> subsp. <i>aureus</i>	-	Zika virus African strain	-	<i>Ureaplasma parvum</i>	-
<i>Staphylococcus epidermidis</i>	-	Zika virus strain FB- GWUH-2016	-	<i>Ureaplasma urealyticum</i>	-
<i>Streptococcus pneumoniae</i> Z022	-	Zika Virus French Polynesian, strains 11468/16 and 11474/16	-	Varicella Zoster Virus OKA	-
<i>Streptococcus pyogenes</i>	-	<i>Acinetobacter baumannii</i>	-	Varicella Zoster Virus Ellen	-
<i>Streptococcus salivarius</i>	-				

## Analytical reactivity

The Vitassay qPCR MonkeyPox kit's reactivity was evaluated against Monkeypox virus, showing positive results.

## Compatibles real-time PCR equipment

Vitassay qPCR MonkeyPox has been validated on the following equipment:

- Cobas z480 Analyzer (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) <sup>I</sup>
- CFX96™ Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTPPrime Real Time Detection Thermal Cycler (DNA-Technology)
- DTlite Real-Time PCR System (DNA-Technology)
- CFX Opus 96 Touch™ Real-Time PCR Detection System (Bio-Rad)
- LightCycler 480 Instrument II (Roche)
- NEOS-96 Real Time PCR System (Linear)
- QuantStudio™ 5 Real-Time PCR System (Applied Biosystems)
- QuantGene 9600 (BioEr Technology)
- Rotor-Gene® Q (Qiagen) <sup>II</sup>

<sup>I</sup>: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

<sup>II</sup>: For Rotor-Gene® Q thermocycler the product should be reconstituted following the procedure and transferred into specific tubes of the equipment.

## Limitations

- All obtained results must be interpreted by a specialist in conjunction with available clinical information and laboratory findings.
- This assay has been validated with DNA extracted from lesion swabs. The use of other samples has not been established.
- The correct test performance depends on the sample's quality; proper DNA from clinical specimens must be extracted.
- This test does not provide quantitative values nor indicate the number of organisms present. This is a qualitative test.
- Extremely low levels of target below the limit of detection might be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination either by MonkeyPox Positive Control during its reconstitution, by samples containing high concentrations of target DNA or contamination due to PCR products from previous reactions.
- Detection may be affected by several factors and their combinations which may lead to false negative results, including a) inadequate specimen sampling, shipping, storage, handling; b) procedural errors (including DNA extraction); c) DNA degradation during specimen shipping, storage, and/or preparation; d) pathogen load below the limit of detection for the assay; e) the presence of Real-Time amplification inhibitors or other types of interference (an interference study evaluating the effect of vaccines, antiviral therapeutics, antibiotics, chemotherapeutics or immunosuppressant drugs used to prevent the infection or used during the treatment of the infection was not performed); f) failure to follow the manufacturer's instructions and procedures; g) mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown Monkeypox virus clades.
- Viral DNA detection may not indicate the presence of viable and/or infectious virus or that Monkeypox virus is the causative agent for clinical symptoms.
- The negative results do not preclude Monkeypox virus infection and should not be the sole basis of a patient treatment/management decision. Optimum specimen types for Monkeypox virus identification and/or stage of infection most suitable for its collection have not been established yet. Consider the collection of multiple specimens from the same patient at different time points, which may increase the probability of detecting the virus.

- If the patient clinical data, laboratory tests and epidemiological studies suggest possible Monkeypox virus infection, and other diseases have been discarded, a false negative result might not be discarded, and additional tests should be performed.
- Fluorescence values may vary due to multiple factors such as PCR equipment, extraction system, sample type, and its previous treatment, among others.
- Some samples may fail to exhibit *HBB* amplification curves due to low human cell numbers in the original clinical sample. A negative endogenous IC signal does not preclude the Monkeypox virus presence in a clinical specimen.

## Attached I. Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following guidance table separated by tube type. Please refer to the table and check the equipment and its specifications before running the assay. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
<b>Agilent Technologies</b>	<b>Abbott</b>
AriaMx Real-Time PCR System	Abbott m2000 <sup>(1)</sup>
AriaDx Real-Time PCR System	<b>Agilent Technologies</b>
<b>Applied Biosystems</b>	Mx3000P™ Real Time PCR System
7500 Fast Real-Time PCR System <sup>(1) (2)</sup>	Mx3005P™ Real Time PCR System
7500 Fast Dx Real-Time PCR System <sup>(1) (2)</sup>	<b>Applied Biosystems</b>
QuantStudio™ 12K Flex 96-well Fast	7300 Real-Time PCR System <sup>(1) (4)</sup>
QuantStudio™ 6 Flex 96-well Fast	7500 Real-Time PCR System <sup>(1)</sup>
QuantStudio™ 7 Flex 96-well Fast	7900 HT Real-Time PCR System <sup>(3)</sup>
QuantStudio™ 5 Real-Time PCR System	ABI PRISM 7000 <sup>(3)</sup>
QuantStudio™ 5 Fast Real-Time PCR System	ABI PRISM 7700 <sup>(3)</sup>
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
StepOne Plus™ Real-Time PCR System <sup>(3)</sup>	QuantStudio™ 6 Flex 96-well
StepOne™ Real-Time PCR System <sup>(3) (4)</sup>	QuantStudio™ 7 Flex 96-well
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 5 Fast Real-Time PCR System
<b>Azure Biosystems</b>	QuantStudio™ 5 Real-Time PCR System
Azure Cielo 3 <sup>(5)</sup>	QuantStudio™ 3 Real-Time PCR System
Azure Cielo 6	ViiA™ 7 Real-Time PCR System
<b>BIONEER</b>	<b>Analytik Jena</b>
Exicycler™ 96 Fast	qTOWER <sup>(6)</sup>
<b>Bio-Rad</b>	<b>BIONEER</b>
CFX96™ Real-Time PCR Detection System	Exicycler™ 96
CFX96™ IVD Real-Time PCR Detection System	<b>BIOER</b>
Mini Opticon™ Real-Time PCR Detection System <sup>(5)</sup>	QuantGene 9600
CFX Opus 96	<b>Bio-Rad</b>
<b>Roche</b>	CFX96™ Deep Well Real-Time PCR System
LightCycler® 480 Real-Time PCR System <sup>(1) (6)</sup>	CFX96™ Deep Well IDV Real-Time PCR System
LightCycler® 96 Real-Time PCR System	iCycler iQ™ Real-Time PCR Detection System
Cobas z480 Analyzer <sup>(1) (6)</sup>	iCycler iQ™ 5 Real-Time PCR Detection System
	My iQ™ Real-Time PCR Detection System <sup>(5)</sup>
<b>Special Formats <sup>(7)</sup></b>	My iQ™ 2 Real-Time PCR Detection System <sup>(5)</sup>
<b>Bio Molecular Systems</b>	<b>DNA-Technology</b>
Mic Real Time PCR Cyclor	DTlite Real-Time PCR System <sup>(8)</sup>
<b>Cepheid</b>	DTprime Real-time Detection Thermal Cyclor <sup>(8)</sup>
SmartCycler®	<b>Eppendorf</b>
<b>Qiagen</b>	Mastercycler™ ep <i>realplex</i>
Rotor-Gene® Q	<b>Qiagen</b>
	QIAquant 96 <sup>(6)</sup>

(1) A special grid is needed to fit these real-time PCR kits.

(2) Select Ramp Speed "Standard" in New Experiment/Advanced Set-up/Experiment Properties.

(3) No Cy5 caption.

(4) No ROX caption.

(5) Only FAM and HEX caption.

(6) Specific compensation color is required.

(7) The product must be reconstituted following the appropriate procedure (see Test procedure) and transferred to the specific tubes for Mic, SmartCycler®, or Rotor-Gene® Q.

(8) See Attached III to configure exposure settings.

## Attached II. Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real-Time PCR thermocyclers are specified in the following table:

Thermocycler	Vitassay Channel	Detection Channel	Observations
<b>Bio-Rad CFX96™</b>	FAM	FAM	During the first few cycles of an analysis, some wells may exhibit abnormal RFU values and show a non-sigmoid rising line. To correct this effect, select the Apply Fluorescence Drift Correction option in the Settings menu for the baseline configuration.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Applied Biosystems ABI 7500</b>	FAM	FAM	Passive reference option ROX is none. During the first few cycles of an assay, some wells may exhibit abnormal RFU values and show a non-sigmoidal rising line. To correct for this effect, select the Start cycle and End cycle values so that the baseline ends before the detection of significant fluorescence.
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Roche Lightcycler®480II</b>	FAM	465/510	Color Compensation required.
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
<b>Roche Cobas z 480</b>	FAM	465/510	Color Compensation required.
	HEX	540/580	
	ROX	540/610	
	Cy5	610/670	
<b>Cepheid Smartcycler®</b>	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
<b>Abbott m2000rt</b>	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
<b>Agilent Technologies Mx3000P™ Mx3005P™</b>	FAM	FAM	Passive reference option ROX is none.
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Bio Molecular Systems Mic Real Time PCR Cycler</b>	FAM	Green	In the "Run Profile" menu, introduce the correct parameters for "Temperature Control" (Standard TAQ (v3)), Volume (20 ul) and the thermal profile. In the "Cycling" window, select the "Acquire on" option for all the channels by clicking on them. Use the default "Gain" values for each channel (Green = 3, Yellow = 10, Orange = 10, Red = 10).
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>Agilent Technologies AriaMx</b>	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
<b>Qiagen Rotor-Gene® Q</b>	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range must be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
<b>BIONEER Exicycler™ 96</b>	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

**Attached III. Optical measurement exposure setting**

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

<b>Thermocycler</b>	<b>Vitassay channel</b>	<b>Exposure values</b>
<b>DTlite Real-Time PCR System (DNA-Technology)</b>	FAM	500
	HEX	500
	ROX	500
	Cy5	500
<b>DTprime Real-time Detection Thermal Cycler (DNA-Technology)</b>	FAM	500*
	HEX	1000
	ROX	1000
	Cy5	1000

\*If result in FAM channel is not as expected, there are no amplifications or high background noise is observed, please lower the exposure values indicated above to 150.

---


## Bibliography/Bibliografía

- Berhanu A, Prigge JT, Silvera PM, Honeychurch KM, Hruby DE, Grosenbach DW. (2015). Treatment with the smallpox antiviral tecovirimat (ST-246) alone or in combination with ACAM2000 vaccination is effective as a postsymptomatic therapy for monkeypox virus infection. *Antimicrob Agents Chemother.* Jul;59(7):4296-300.
- Centers for Disease Control and Prevention (CDC). Monkeypox. Available at: <https://www.cdc.gov/poxvirus/monkeypox/> Accessed May 2022
- Edghill-Smith Y, Golding H, Manischewitz J, King LR, Scott D, Bray M, Nalca A, Hooper JW, Whitehouse CA, Schmitz JE, Reimann KA, Franchini G. (2005). Smallpox vaccine-induced antibodies are necessary and sufficient for protection against monkeypox virus. *Nat Med.* Jul;11(7):740-7.
- Ladnyj ID, Ziegler P, Kima E. (1972) A human infection caused by monkeypox virus in Basankusu Territory, Democratic Republic of the Congo. *Bull World Health Organ.* 46(5):593-7.
- Li Y, Zhao H, Wilkins K, Hughes C, Damon IK. (2010). Real-time PCR assays for the specific detection of monkeypox virus West African and Congo Basin strain DNA. *J Virol Methods.* Oct;169(1):223-7.
- McCollum AM, Damon IK. (2014). Human monkeypox. *Clin Infect Dis.* 2014 Jan;58(2):260-7. Erratum in: *Clin Infect Dis.* Jun;58(12):1792.
- Petersen E, Kantele A, Koopmans M, Asogun D, Yinka-Ogunleye A, Ihekweazu C, Zumla A. (2019). Human Monkeypox: Epidemiologic and Clinical Characteristics, Diagnosis, and Prevention. *Infect Dis Clin North Am.* Dec;33(4):1027-1043.
- World Health Organization (WHO). MonkeyPox. Available at: <https://www.who.int/news-room/fact-sheets/detail/monkeypox>. Accessed May 2022
- World Health Organization (WHO). "Multi-Country Monkeypox Outbreak in Non-Endemic Countries.". *Disease Outbreak News.*, 21 May 2022, Available at: <https://www.who.int/emergencies/disease-outbreak-news/item/2022-DON385> Accessed May 2022

## Trademarks

All trademarks that may appear in this package insert are the property of their respective owners.

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