## **MASTAFLUOR™ Toxoplasma**

MASTAFLUOR™ Toxoplasma Screen	<b>REF</b>	631182	UDI-DI: 4250729700095	10 x 10 Tests
MASTAFLUOR™ Toxoplasma IgG	<b>REF</b>	631184	UDI-DI: 4250729700118	10 x 10 Tests

### **Gebrauchsanweisung / Instructions for Use / Notice d'utilisation**

**Nur zur Verwendung durch Fachpersonal /  
For professional use only /  
Exclusivement pour un usage professionnel**

**CE** 2797

	Deutsch	Seiten	02–05
	English	Pages	06–09
	Français	Pages	10–13

## MASTAFLUOR™ Toxoplasma

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## MASTAFLUOR™ Toxoplasma

Immunofluorescence assay for the detection of antibodies to *Toxoplasma gondii* in human serum.

### Intended Purpose

Semi-quantitative (titratable) immunofluorescence assay for the detection of IgG antibodies (and additionally IgM antibodies using the Screen assay) directed against *Toxoplasma gondii* tachyzoites in human serum as an aid to diagnosis of immunity against toxoplasmosis (presence of specific antibodies).

The assay is suitable for manual use and is intended for professional in-vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

### Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

### Principle of the Test

Diluted specimens, the negative and positive controls are applied to the wells of microscope slides and incubated. All wells are coated with purified antigens of the very pathogen. If specific antibodies are present in the serum they will bind to the fixed pathogen antigens forming a stable antigen-antibody complex. Slides are then washed to remove any unbound material. Complexed antibodies are detected by the addition of a fluorescence labelled anti-human IgG, IgM, or IgA immunoglobulin conjugate. After a further washing step to remove any unbound conjugate, slides are viewed under a fluorescence microscope.

A green fluorescence is observed if pathogen specific antibodies are present in the sample material.

## Kit Contents

1.	<b>SLIDE</b>	<b>Slides</b>	10 x 10 Wells
		coated with <i>Toxoplasma gondii</i>	
2.	<b>CONTROL</b> +	<b>Positive Control</b>	0,5 mL
		IgG positive serum, human, ready-to-use, contains < 0,1 % NaN <sub>3</sub> (potentially infectious, see warnings)	
3.	<b>CONTROL</b> -	<b>Negative Control</b>	0,5 mL
		IgG and IgM negative serum, human, ready-to-use, contains < 0,1 % NaN <sub>3</sub> (potentially infectious, see warnings)	
4.	<b>CONJ</b> <b>G</b>	<b>FITC-conjugate IgG</b>	3 mL
		FITC-labelled anti-human IgG (γ-chain) ready-to-use, contains < 0,1 % NaN <sub>3</sub>	
		or	
	<b>CONJ</b> <b>POLY</b>	<b>FITC-conjugate Screen</b>	3 mL
		FITC-labelled anti-human IgG/IgM/IgA ready-to-use, contains < 0,1 % NaN <sub>3</sub>	
5.	<b>DIL</b> <b>5x</b>	<b>Diluent</b>	2 x 10 mL
		5 x concentrated, to be diluted 1:5 before use	
6.	<b>EVANS BLUE</b>	<b>Evans Blue</b>	3 mL
		ready-to-use, contains < 0,1 % NaN <sub>3</sub>	
7.	<b>MOUNTING MEDIUM</b>	<b>Mounting Medium</b>	3 mL
		ready-to-use, contains < 0,1 % NaN <sub>3</sub>	
8.	<b>BUFFER</b>	<b>PBS</b>	2 Sachets
		Dissolve 1 sachet in 1 L distilled water, pH 7,2 ± 0,2	

### Further Abbreviations

1. **RTU** ready to use

### Materials Required but not Provided

1. Sterile test tubes.
2. Micropipettes and tips.
3. Staining dish or Coplin jar.
4. Moist chamber.
5. Volumetric flask for PBS.
6. Distilled or deionised water.
7. Forceps.
8. Cover slips.
9. Wash bottle.
10. Fluorescence microscope with a filter combination suitable for FITC (e.g. 490 nm excitation filter and a 510nm barrier filter).

## Warnings

1. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
2. Do not cross-contaminate reagents or interchange caps on bottles. Do not use contaminated reagents. Do not use reagents with damaged bottles or packaging. Use a separate pipette or pipette tips for each sample and reagent.
3. Do not use haemolysed, lipaemic or icteric serum.
4. The control/calibrator materials of human origin provided have been tested for antibodies to HIV, HCV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the samples are free of infections or microbial contamination.
5. The kit reagents may contain bovine components. To protect against infections with prions or zoonotic pathogens, standard occupational safety regulations (GLP) must be observed.
6. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose of by flushing to drain with plenty of water.
7. The MASTAFLUOR™ Toxoplasma Screen detects both IgG and IgM antibodies directed against *T.gondii*. It is not possible to differentiate between IgG and IgM.
8. If a positive result is obtained from a pregnant woman it is recommended to carry out further tests to determine the identity of the antibodies.  
  
If a negative result is obtained from a pregnant woman it is advisable to repeat the test during the period of pregnancy to detect a possible seroconversion. In case of doubt, it is recommended to perform another test or to consult the physician.

## Precautions / Notes

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
3. Do not use reagents beyond the expiry date.
4. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
5. Use disposable plasticware wherever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
6. Always use a new pipette tip for each pipetting step.
7. Do not pipette with your mouth.
8. Do not allow wells to dry out during the assay procedures.
9. Do not expose slides to intense sunlight or similar adverse conditions while incubating.

10. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.

## Disposal

Samples and contaminated disposables should be disposed according to relevant disposal directives and regulations for infectious materials.

For the decontamination of surfaces use suitable surface disinfectants. Contaminated glassware should be autoclaved at 121 °C.

## Stability and Storage

MASTAFLUOR™ Toxoplasma can be used until the end of expiry date as indicated on kit label. All kits components and reagents should be stored at 2–8 °C and brought to room temperature before use.

After first opening reagents are stable for at least 3 months. A slide shall be used shortly after opening the pouch.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Samples (serum) may be stored according to general recommendation in literature. In general samples can be kept at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

## Sample Material

- a) Serum
- b) Samples can be stored at 2–8 °C for up to 3 days. Serum should be aliquoted immediately after sampling and kept at -20°C for longer storage.  
  
The samples should not be frozen and thawed repeatedly.  
  
After thawing samples should be briefly vortexed carefully before being used in the assay.  
  
Lipaemic, haemolytic, icteric or bacterially contaminated samples can cause false positive or false negative results (Ref. 5, 6, 7).

## Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or deionised water.  
  
It is recommended to use the diluent to prepare the sample dilutions. Dilute the 5x concentrated diluent 1:5 with distilled or deionised water (e.g. 10 mL of diluent + 40 mL of water).
3. **A) Qualitative:** Dilute all serum specimens with diluent or PBS buffer as indicated below:  
  
**IgG assays: 1:64**  
e.g. 10 µL serum + 630 µL diluent or PBS  
**Screen assays: 1:64**  
e.g. 10 µL serum + 630 µL diluent or PBS

## B) Semi-quantitative: For the IgG and screen assay:

For the IgG and screen approach: To determine the titer, the samples can be further diluted starting from the addition titer (1:64). To do this, prepare a serial dilution series with Diluent or PBS, e.g.: prepare 50 µl Diluent or PBS and titrate with 50 µl 1:64 diluted sample (1:2).

4. All **kit controls are ready-to-use** and **should not be diluted** like specimen.
5. Carefully remove the slides from the pouch. Do not touch wells.
6. Apply 20–25 µL of ready-to-use controls and diluted specimens to respective wells on the slides. Ensure that all the wells are covered. Do not scratch the surface of the wells.
7. Add slides to a moist chamber and incubate at room temperature

### IgG assays 30 min

### Screen assays 30 min

8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the center of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min at room temperature with a change of PBS after 5 min to increase washing strength.
10. Remove slides from the staining jar and drain off any excess buffer. Using a blotter, dry the area outside wells. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of the respective FITC-conjugate to each well. The conjugate shall cover the whole well. Do not allow the wells to dry.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 min.
13. After incubation wash the slides with PBS as in steps 9 and 10.  
Optionally: Add 3 to 5 drops of Evans Blue for counterstaining to the last wash step. Counterstaining may improve differentiation between negative and positive reactions.  
After the Evans Blue counterstain dip the slide in fresh PBS to remove excess Evans Blue.
14. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.  
Note: Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip. Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering.
15. Examine reactions under a fluorescence microscope at a total magnification of 400x or at 800x.

## Interpretation of Results

**IgG and Screen positive reaction:** Positive sera show a homogenous green fluorescence of the cytoplasmic membrane of *Toxoplasma gondii* tachyzoites. A reaction can be rated from very strongly fluorescent (4+), strongly fluorescent (3+), fluorescent (2+) to slightly fluorescent (1+).

Cytoplasm may show green staining or even fluorescence together with the cytoplasmic membrane. Cytoplasmic reactions alone must be read as negative.

**Negative reaction:** The *Toxoplasma* membrane does not show any fluorescence. Fluorescence of cytoplasm only is considered as negative.

### Controls

Controls have to fit the lot specific reactions stated in the QC certificate. The figures in the table below are for orientation only and do not reflect the data of the current lot.

Specimen	Expected Results
Positive control	3 + bis 4 +
Negative control	negative

For a correct interpretation, the results should be compared with positive and negative controls.

### Interpretation of Specimen Results

**Qualitative:** A patient sample is considered positive if it shows a positive reaction (1+ - 4+) in the 1:64 dilution.

**Semi-quantitative:** A patient sample is considered positive if the 1:64 dilution or a higher dilution shows fluorescence of the toxoplasma membranes (1+ - 4+).

**IgG antibody:** ≥ 1:64

**Screen / polyvalent:** ≥ 1:64

The final titer is the highest dilution with a minimum reaction of 1+.

For result interpretation recommendations of (inter-) national reference centers should be followed (e.g. [www.rki.de](http://www.rki.de)).

### Limitations / Interferences

1. Autoantibodies (ANA type) may lead to cytoplasmic fluorescence of *Toxoplasma gondii*.
2. The MASTFLUOR™ *Toxoplasma* test has been demonstrated to be highly sensitive and specific. Like for other lab based assays final diagnosis shall not be made on this assay alone, always interpret the data in connection with the patient's medical history and the clinical picture.

### Performance characteristics

#### Analytical Sensitivity and Specificity, Cut-off

A panel of 140/153 clinically positive and negatively evaluated left over samples was selected to determine the performance data (cut-off ≥ 1:64). The TOXOREAGENT

assay (MAST Diagnostica GmbH, Reinfeld, Germany), Toxoreagent assay (EIKEN Chemicals Ltd, Tokyo, Japan) were used as comparable products. The VIR-ELISA ANTI-TOXO-IgG/IgM (Viro-Immun Labor-Diagnostics GmbH, Oberursel, Germany) were used to control the sample.

The MASTFLUOR™ Toxoplasma obtained the following results in this panel:

	Formula	IgG	Screen
Sensitivity (Sens)	$\frac{TP}{(TP + FN)}$	90,57 %	95,61 %
Specificity (Spec)	$\frac{TN}{(TN + FP)}$	94,12 %	94,87 %
Positive predictive value (PPV)	$\frac{TP}{(TP + FP)}$	97,96 %	98,20 %
Negative predictive value (NPV)	$\frac{TN}{(TN + FN)}$	76,19 %	88,10 %
Efficiency (EFF)	$\frac{(TP + TN)}{total}$	91,43 %	95,42 %

(Abbreviations: TP: True positive, TN: True negative, FN: False negative, FP: False positive)

No differences in intra-assay and inter-assay variance was found.

## Cross-Reactivity

47 sera containing antinuclear antibodies were tested for cross-reactivity. No false positive or negative results were obtained but the samples showed a high background (see Limitations).

In one study, 23 sera from pregnant women were used to determine possible cross-reactions. No false positive or negative results were detected. Nevertheless, a warning was included (see warnings).

## Traceability

The NIBSC sample 01/600 or TOXM can be used in the qualitative evaluation to monitor the positive control of the Toxoplasma IgG / Screen Kit. The NIBSC sample TOXM can be used in the semi-quantitative evaluation to monitor the positive control.

## Clinical Performance

The clinical performance was demonstrated by proficiency testing.

For the MASTFLUOR™ Toxoplasma IgG/Screen assay, proficiency tests (INSTAND e.V. Ringversuche, Section 452, Toxoplasmosis) are performed 1-2 times per year. Data on the interlaboratory comparisons will be provided on request.

## Availability of the summary of safety and performance (Art. 29, IVDR)

The report will be available in the EUDAMED database (as soon as the module is available).

<https://ec.europa.eu/tools/eudamed/#/screen/home>.

The SSP report will be provided on request.

## Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

## References and Change History

You can find the references and the Change History at the end of the instruction for use.

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Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Kit Contents	A warning has been added and a reference to the section 'Warnings' has been included.
Warnings and precautions	Subdivision into warnings, precautionary measures and disposal
Warnings	Addition of warnings against bovine components, Toxoplasma Screen antibodies and sera of pregnant women.
Stability and Storage	Specification of the use of the slide
Sample Material	Addition of references 5-7.
Test procedure	Addition of semi-quantitative dilution serious
Interpretation of Results	Addition of the fluorescent intensities. Addition of semi-quantitative interpretation (titer).
Performance characteristics	Addition of a detailed description of analytical sensitivity/specificity, cut-off, cross-reactivity (ANA/ENA samples and sera of pregnant women), traceability, clinical performance.
References	Addition of references 5-7.

**Mast Diagnostica GmbH,**

Feldstraße 20,  
23858 Reinfeld,  
Deutschland

Tel: +49 (0)4533 2007 0

Fax: +49 (0)4533 2007 68

email: mast@mast-diagnostica.de

Web: www.mast-group.com

**Mast Group Ltd.**

Mast House, Derby Road  
Bootle, Merseyside, L20 1EA  
United Kingdom

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**Icons are used according to DIN EN ISO 15223-1**

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MD V.4.0

2024-11-11