



**Mast Group Ltd.**  
Mast House, Derby Road,  
Bootle, Merseyside, L20 1EA  
United Kingdom  
Tel: + 44 (0) 151 472 1444  
Fax: + 44 (0) 151 944 1332  
email: sales@mast-group.com  
Web: www.mast-group.com



**Mast Diagnostica GmbH**  
Feldstrasse 20  
DE-23858 Reinfeld  
Germany  
Tel: + 49 (0) 4533 2007 0  
Fax: + 49 (0) 4533 2007 68  
email: mast@mast-diagnostica.de  
Web: www.mast-group.com

**Mast Diagnostic**  
12 rue Jean-Jacques Mention  
CS91106, 80011 Amiens, CEDEX 1  
France  
Tél: + 33 (0) 3 22 80 80 67  
Fax: + 33 (0) 3 22 80 99 22  
email: info@mast-diagnostic.fr  
Web: www.mast-group.com



**Mast  
Group**

## MASTDISCS® *Combi* AmpC and ESβL Detection Set

### D68C

#### Intended use

For the detection of AmpC and/or extended spectrum beta-lactamase (ESβL) enzyme production in Enterobacterales.

FOR IN VITRO DIAGNOSTIC USE ONLY

#### Contents and Formulation\*

4 cartridges per pack, each cartridge containing approximately 50 discs:

<b>Cartridge A</b>	Cefpodoxime 10 µg
<b>Cartridge B</b>	Cefpodoxime 10 µg + ESβL inhibitor
<b>Cartridge C</b>	Cefpodoxime 10 µg + AmpC inhibitor
<b>Cartridge D</b>	Cefpodoxime 10 µg + ESβL inhibitor + AmpC inhibitor

#### Storage and shelf life

Store at 2 to 8°C in the containers provided until the expiry date shown on the pack label. Allow to equilibrate to room temperature before opening.

#### Precautions

For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet.

#### Materials required but not provided

Standard microbiological supplies and equipment such as loops, MAST® culture media, Mueller-Hinton agar, swabs, forceps, callipers etc., as well as an incubator capable of maintaining 35 ± 1 °C.

#### Procedure

- Using a pure, fresh culture of the test organism, prepare a suspension equivalent in density to a 0.5 McFarland standard in physiological saline.
- Using a sterile swab, spread the suspension uniformly across the surface of a single Mueller Hinton Agar plate in accordance with the European Committee on Antimicrobial Susceptibility Testing (EUCAST) procedure.
- Using a MAST® DISCMASTER Dispenser, or alternatively a sterile needle or forceps, place one of each disc onto the plate of inoculated medium, ensuring sufficient space between the discs to allow formation of clearly defined zones of inhibition.
- Incubate at 35 ± 1 °C for 18 ± 2 hours.
- Measure and record the diameter of any zones of inhibition, to the nearest whole millimetre. Discs showing no zone of inhibition should be recorded as 6 mm.

#### Interpretation of results

To interpret results based on observed zones of inhibition, use the D68C calculator. The calculator is available for download and can be accessed via [www.mast-group.com](http://www.mast-group.com), in the registered members section. Alternatively, results can be interpreted manually by comparing inhibition zone diameters in the sequence described below:-

**Step 1** - Compare the zone of inhibition of the cefpodoxime disc (**A**) to the zones of inhibition of each of the cefpodoxime plus inhibitor discs (**B, C, and D**). If all zones are within 2mm of each other, record the organism as demonstrating neither ESBL nor AmpC activity.

**Step 2** - Subtract **A** from **B**, and **C** from **D**. If each of **B – A** and **D – C** ≥ 5mm **AND** when compared, the differences in zone diameter between discs **B** and **D** and discs **A** and **C** are 4mm or less then the organism is demonstrating ESBL activity alone.

**Step 3** – Subtract **B** from **D** and **A** from **C**. If each of **D – B** and **C – A** ≥ 5mm **AND** when compared the differences in zone diameter between discs **A** and **B** and discs **C** and **D** are 4mm or less then the organism is demonstrating AmpC activity alone

**Step 4** – Subtract **C** from **D**. If **D – C** ≥ 5mm **AND** when compared the differences in zone diameter between discs **A** and **B** is 4mm or less then the organism is demonstrating ESBL and AmpC combined activity.

#### Quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction. Zones of inhibition obtained against a negative control organism *E. coli*†ATCC® 25922 should be equal or show no greater difference in diameter than ±2mm. Any greater difference implies malfunction or deterioration. Do not use the product if the reactions with the control organisms are incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain.

Test Organism	Result
<i>Escherichia coli</i> NCTC 13351	ESβL Positive
<i>Escherichia coli</i> NCTC 13352	ESβL Positive
<i>Escherichia coli</i> NCTC 13353	ESβL Positive
<i>Enterobacter cloacae</i> NCTC 13406	AmpC Positive
<i>Escherichia coli</i> ATCC® 25922	ESβL and AmpC Negative

#### Limitations

D68C is not suitable for use with *Pseudomonas* spp. or *Acinetobacter* spp. To avoid potentially erroneous results do not mix cartridges from different batches of D68C and ensure all discs in the set are tested on the same plate. Organisms producing a fully resistant profile i.e. no zone of inhibition on all discs could indicate demonstration of an MBL or KPC carbapenemase production, which could also be masking concurrent ESβL or AmpC expression. Users are obliged to always use the latest version of the D68C calculator.

#### References

Bibliography available on request.