



Mast Group Ltd. Mast House, Derby Road, Bootle, Merseyside, L20 1EA United Kingdom

Tel: + 44 (0) 151 472 1444 Fax: + 44 (0) 151 944 1332 email: sales@mast-group.com Web: www.mast-group.com



Mast Diagnostica GmbH Feldstrasse 20 DE-23858 Reinfeld

Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 email: mast@mast-diagnostica.de Web: www.mast-group.com

Germany



12 rue Jean-Jacques Mention CS91106, 80011 Amiens, CEDEX 1 France

Tél: + 33 (0) 3 22 80 80 67 Fax: + 33 (0) 3 22 80 99 22 email: info@mast-diagnostic.fr Web: www.mast-group.com



# MASTDISCS® Combi Cefotaxime ESβL ID Disc Set

## **D62C**

#### Intended use

For the detection of extended spectrum beta-lactamases (ES\( \mathbb{L} \mathbb{S} \)) in Enterobacterales.

# FOR IN VITRO DIAGNOSTIC USE ONLY

#### Contents and Formulation\*

3 x paired sets of cartridges per pack, each cartridge containing approximately 50 discs:

CTX30 Cefotaxime 30 µg discs (x3)

CTXCV Cefotaxime 30 µg + clavulanic acid 10 µg

discs (x3)

### Storage and shelf life

Store at 2 to 8°C in the containers provided until the expiry date shown on the pack label. Allow to equilibrate to room temperature before opening.

#### **Precautions**

For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet.

#### Materials required but not provided

Standard microbiological supplies and equipment such as loops, MAST® culture media, Mueller-Hinton agar, swabs, forceps, callipers etc., as well as an incubator capable of maintaining 35±2 °C.

# Procedure

- Using a pure, fresh culture of the test organism, prepare a suspension equivalent in density to a 0.5 McFarland standard.
- Using a sterile swab, spread the suspension uniformly across the surface of a single Mueller Hinton Agar plate in accordance with the Clinical and Laboratory Standards Institute (CLSI) procedure.
- 3. Using a MAST® DISCMASTER Dispenser, or alternatively a sterile needle or forceps, place one of each type of disc onto the plate of inoculated medium, ensuring sufficient space between the discs to allow formation of clearly defined zones of inhibition.
- 4. Incubate at 35±2 °C for 17±1 hours.
- Measure and record the diameter of any zones of inhibition, to the nearest whole millimetre. Discs showing no zone of inhibition should be recorded as 6 mm.

# Interpretation of results

Compare the zone of inhibition for the cefotaxime disc to that of the cefotaxime plus clavulanic acid combination disc. An increase in zone diameter of  $\geq 5$  mm in the presence of clavulanic acid indicates the presence of ES $\beta$ L in the test organism.

# **Quality control**

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction. Zones of inhibition obtained using the combination disc plus clavulanic acid and corresponding cefotaxime only disc against ES $\beta$ L-negative control organism *E. coli* ATCC $^{\odot}$  25922, should be equal or show no greater difference in diameter than  $\pm 2$  mm. Any greater difference implies malfunction or deterioration. Do not use the product if the reactions with the control organisms are incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain:

Test Organism	Result
Escherichia coli	Positive
NCTC 13351	
Escherichia coli	Positive
NCTC 13353	
Escherichia coli	Negative
ATCC® 25922	

#### Limitations

D62C is not suitable for testing *Pseudomonas* spp. or *Acinetobacter* spp. D62C should always be used in conjunction with **MAST**DISCS® *Combi* Ceftazidime ES $\beta$ L ID disc set (D64C); a positive result using one or both tests indicates the presence of an ES $\beta$ L in the test organism. To avoid potentially erroneous results do not mix cartridges from different batches of D62C and ensure both discs in the set are tested on the same plate.

# References

Bibliography available on request.