

Mas Mast Boot

Mast Group Ltd. Mast House, Derby Road, Bootle, Merseyside, L20 1EA United Kingdom

Tel: + 44 (0) 151 472 1444 Fax: + 44 (0) 151 944 1332 email: sales@mast-group.com Web: www.mast-group.com



Mast Diagnostica GmbH Feldstrasse 20 DE-23858 Reinfeld

Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 email: mast@mast-diagnostica.de Web: www.mast-group.com

Germany

Mast Diagnostic

12 rue Jean-Jacques Mention CS91106, 80011 Amiens, CEDEX 1 France

Tél: + 33 (0) 3 22 80 80 67 Fax: + 33 (0) 3 22 80 99 22 email: info@mast-diagnostic.fr Web: www.mast-group.com



CAMP (Preston Blood Free) MAST® SELECTAVIAL

SV18 Series

Intended Use

For the selective culture of *Campylobacter* spp.

FOR IN VITRO DIAGNOSTIC USE ONLY

Contents

10 vials of supplement.

Formulation

Material:	Concentration in medium:
Cefoperazone	32mg/L

Storage and shelf life

Store unopened at 2 to 8°C until expiry date shown on pack label. Once reconstituted use immediately.

Precautions

For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet.

Materials required but not provided

Standard microbiological supplies and equipment such as loops, MAST® culture media, swabs, applicator sticks, incinerators and incubators, etc., as well as serological and biochemical reagents, and additives such as blood.

Procedure

- 1. Sterilise appropriate volume of MAST® Preston Campylobacter Agar Base (DM251D) cool to 50 to 55°C and hold at this temperature.
- Reconstitute the contents of one vial using the diluent specified on the pack label. The best method is to aseptically add the diluent using a sterile needle and syringe. Draw the diluent into the syringe and after removing the plastic cap, inject through the rubber stopper of the vial. The lyophilised supplement will rapidly dissolve and may be withdrawn into the syringe.
- 3. Add the antibiotic supplement to the volume of medium specified on the pack label and discard the needle into an approved container.
- Mix gently but thoroughly to evenly distribute the selective agents. Pour culture plates (15 to 20 mL per plate) and allow to set.
- 5. Prepared culture plates may be used immediately or stored in plastic bags at 2 to 8°C for up to one week before use.
- 6. Directly inoculate the surface of a dried plate from the specimen to achieve isolated colonies.

7. Plates should be incubated in an atmosphere of reduced oxygen of 5% oxygen, 10% carbon dioxide and 85% nitrogen, and examined after 48 hours at 37°C.

Interpretation of results

The use of CAMP (Preston Blood Free) MAST® SELECTAVIAL suppresses the growth of normal flora, thus allowing *Campylobacter* spp. to be readily identified. Colonies of *C. jejuni* will appear grey, moist and spreading. Strains of *C. coli* as creamy grey raised, moist and often discrete colonies. Further identification tests should be carried out to confirm identification e.g. MAST® ID CAMP IDENTIFICATION SYSTEM (CAMP ID).

Quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction. Do not use the product if the reactions with the control organisms are incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain.

Test Organisms	Result
Proteus mirabilis	No growth
ATCC® 43071	
Escherichia coli	No growth
ATCC® 25922	
Enterococcus faecalis	No growth
ATCC® 29212	
Campylobacter jejuni	Growth
ATCC® 29428	

References

Bibliography available on request.