



Mast Group Ltd.
Mast House, Derby Road,
Bootle, Merseyside, L20 1EA

United Kingdom Tel: + 44 (0) 151 472 1444 Fax: + 44 (0) 151 944 1332

Fax: + 44 (0) 151 944 1332 email: sales@mast-group.com Web: www.mast-group.com



Mast Diagnostica GmbH
Feldstrasse 20
DE-23858 Reinfeld

Germany

Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 email: mast@mast-diagnostica.de Web: www.mast-group.com **Mast Diagnostic**

12 rue Jean-Jacques Mention CS91106, 80011 Amiens, CEDEX 1 France

Tél: + 33 (0) 3 22 80 80 67 Fax: + 33 (0) 3 22 80 99 22 email: info@mast-diagnostic.fr Web: www.mast-group.com



Neomycin MAST® SELECTAVIAL

SV8 Series

Intended Use

For the selective isolation of Clostridia and other anaerobes.

FOR IN VITRO DIAGNOSTIC USE ONLY

Contents

10 vials of MAST® SELECTAVIAL.

Formulation

Material:	Concentration in medium:
Neomycin	75mg/L

Storage and shelf life

Store unopened at 2 to 8°C until expiry date shown on pack label. Once reconstituted use immediately.

Precautions

For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet.

Materials required but not provided

Standard microbiological supplies and equipment such as loops, MAST® culture media, swabs, applicator sticks, incinerators and incubators, etc., as well as serological and biochemical reagents, and additives such as blood.

Procedure

- Sterilise appropriate volume of MAST® Blood Agar Base Special (DM101D), Columbia Agar (DM115D) or Brucella Medium (DM107D), cool to 50 to 55°C and hold in a water bath at this temperature.
- 2. Reconstitute the contents of one vial using the diluent specified on the pack label. The best method is to aseptically add the diluent using a sterile needle and syringe. Draw the diluent into the syringe and after removing the plastic cap of the vial, inject through the rubber-stopper of the vial. The lyophilised supplement will rapidly dissolve and may be withdrawn into the syringe.
- Add the antibiotic solution to appropriate amount of media and discard the needle into an approved container.
- 4. Mix gently but thoroughly to evenly distribute the selective agents.
- Supplement the medium with 5 to 7% defibrinated horse blood. Other growth factors such as haemin and menadione may also be added as required. Mix well and pour culture plates of normal thickness (15 to 20 mL per plate). Allow to set.

- 6. Prepared culture plates may be used immediately or stored in plastic bags at 2 to 8°C for up to one week before use.
- Inoculated plates should be incubated at 35 to 37°C in an anaerobic jar or cabinet. Examine the plates after 48 hours incubation but continue incubation for up to 5 days.

Interpretation of results

Neomycin blood agar will allow the growth of clostridia, most *Bacteroides fragilis* and some anaerobic cocci whilst suppressing the growth of most Gram negative bacteria.

Quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction. Do not use the product if the reactions with the control organisms are incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain.

Test Organisms	Result
Staphylococcus aureus	No growth
ATCC® 25923	
Proteus mirabilis	No growth
ATCC® 43071	
Pseudomonas aeruginosa	Growth
ATCC® 27853	
Enterococcus faecalis	No growth
ATCC® 29212	
Bacteroides fragilis	Growth
ATCC® 25285	
Clostridium perfringens	Growth
ATCC® 13124	
Clostridium sporogenes	Growth
ATCC® 19404	

References

Bibliography available on request.