



Mast Group Ltd.
Mast House, Derby Road,
Bootle, Merseyside, L20 1EA
United Kingdom
Tel: + 44 (0) 151 472 1444
Fax: + 44 (0) 151 944 1332
email: sales@mast-group.com
Web: www.mast-group.com



Mast Diagnostica GmbH
Feldstrasse 20
DE-23858 Reinfeld
Germany
Tel: + 49 (0) 4533 2007 0
Fax: + 49 (0) 4533 2007 68
email: mast@mast-diagnostica.de
Web: www.mast-group.com

Mast Diagnostic
12 rue Jean-Jacques Mention
CS91106, 80011 Amiens, CEDEX 1
France
Tél: + 33 (0) 3 22 80 80 67
Fax: + 33 (0) 3 22 80 99 22
email: info@mast-diagnostic.fr
Web: www.mast-group.com



**Mast
Group**

MAST® ID PYR Strips

ET07

Intended use

A rapid strip test for the detection of pyrrolidonyl amino peptidase activity in streptococci and enterococci.

FOR IN VITRO DIAGNOSTIC USE ONLY

Contents

25 strips (ET07)

Formulation*

Filter paper strips 5.7cm by 0.6cm, which are printed to identify the test, positive control and negative control areas. The strips are impregnated with L-pyrrolidonyl-β-naphthylamide.

Storage and shelf life

Store at 2 to 8°C in the containers provided until the expiry date shown on the pack label. Allow to equilibrate to room temperature before opening.

Precautions

For *in vitro* diagnostic use only. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet.

Materials required but not provided

Standard microbiological supplies and equipment such as loops, MAST® culture media, swabs, applicator sticks, incinerators and incubators, etc., as well as serological and biochemical reagents and additives such as 4-(Dimethylamino) cinnamaldehyde (DMACA), CAS No. 6203-18-5.

Procedure

- Using a sterile needle or forceps, place a PYR Strip onto a suitable surface e.g. a clean microscope slide or empty Petri dish.
- Aseptically add a small drop of sterile de-ionised or distilled water to each area of the strip. The paper should be moist but not saturated.
- Using a pure, fresh culture of the test organism, remove several colonies by using either a wooden applicator stick or a loop and rub onto the test area of the strip.
- Incubate at 35 to 37°C for 5 minutes.
- Apply approximately 10µl DMACA reagent (as prepared below) to each organism applied.
- Dissolve 1g of 4-(Dimethylamino) cinnamaldehyde (DMACA), CAS No. 6203-18-5, in 100ml of 10% v/v concentrated hydrochloric acid. The solution can be stored in the dark for up to 2 months at room temperature. DMACA reagent is classified as Irritant
- Observe any colour change occurring within 30 seconds.

Interpretation of results

Positive - Fuschia pink colour development.

Negative - No colour change

Quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate a positive reaction and at least one organism to demonstrate a negative reaction, these should be applied to the appropriate areas of the strip. Do not use the product if the reactions with the control organisms are incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain.

Test Organisms	Result
<i>Streptococcus pyogenes</i> ATCC® 19615	Positive
<i>Enterococcus faecalis</i> ATCC® 29212	Positive
<i>Streptococcus agalactiae</i> ATCC® 13813	Negative

Limitations

It is recommended that biochemical and/or serological tests are performed on colonies from pure culture to confirm identification e.g. MAST® STREP - Enzyme Extraction (RST201).

Organisms other than streptococci or enterococci e.g. certain *Klebsiella* spp. may appear PYR positive. Only catalase negative, Gram positive cocci should be tested.

A small proportion of group A streptococci are unable to hydrolyse PYR.

Certain non enterococcal group D streptococci have been reported to be PYR positive.

References

Bibliography available on request.