CHROMagar™ MH Orientation

Instructions For Use
Available in several languages

NT-EXT-035

Version 3.0

ENGLISH

English Version

FRANCAIS


Version Française

ESPAÑOL

Version Español

DEUTSCH

Deutsch Version

MEDIUM PURPOSE

Chromogenic medium for the differentiation and direct antibiotic susceptibility testing in the same plate.

COMPOSITION

The product is composed of a powder base (B) and a supplement (S).

| Product | = | Base (B) | Supplement (S) |
|-----------------|---|---|--------------------|
| Total g/L | | 39.4 g/L | 4.0 ml/L |
| Composition g/L | | Agar 17.0 Peptone and growth factors 21.0 Chromogenic mix 1.4 | Growth factors mix |
| Aspect | | Powder Form | Liquid Form |

| | |
|---------|----------|
| STORAGE | 15-30 °C |
|---------|----------|

| | |
|----------------|-------------|
| FINAL MEDIA pH | 7.1 +/- 0.2 |
|----------------|-------------|

PREPARATION (Calculation for 1L)

| | |
|---|--|
| Step 1 Preparation of the mix | <ul style="list-style-type: none"> Disperse slowly 39.4 g of powder base in 1L of purified water. Add 4ml of supplement (S) into slurry. Stir until agar is well thickened. Heat and bring to boil (100 °C) while swirling or stirring regularly. <p style="color: green;">Advice 1: For the 100°C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).</p> <ul style="list-style-type: none"> AUTOCLAVE at 121 °C during 15 min. |
| Step 2 Pouring | <ul style="list-style-type: none"> Cool in a water bath to 45-50 °C, swirling or stirring gently. Pour into sterile Petri dishes. Let it solidify and dry. |
| Storage | <ul style="list-style-type: none"> Store in the dark before use. Prepared media plates can be kept for one day at room temperature. Plates can be stored for up to 2 months under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration. |

INOCULATION

Spread a bacterial suspension according to classical procedure. Alternatively perform a direct spreading of the sample on the plate allowing a confluent growth.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Incubate in aerobic conditions at 37 °C for 18-24 hours. (In some cases a first reading can be done as soon as 8 hours after incubation).

Typical Samples

e.g. respiratory samples,
urine

Possible enrichment step
Spreading technique


PRE-IDENTIFICATION INTERPRETATION

| Microorganism | Typical colony appearance |
|--|--|
| Gram (-) | |
| <i>E.coli</i> | → dark pink to reddish |
| <i>Klebsiella, Enterobacter, Citrobacter, Serratia</i> | → metallic blue (+/- reddish halo) |
| <i>Proteus, Morganella, Providencia</i> | → brown halo |
| <i>Proteus vulgaris</i> | → blue with brown halo |
| <i>Pseudomonas</i> | → translucent (+/- natural pigmentation cream to green) |
| <i>Acinetobacter</i> | → cream |
| <i>Stenotrophomonas</i> | → colourless |

| Microorganism | Typical colony appearance |
|------------------------|---------------------------|
| Gram (+) | |
| <i>Enterococcus</i> | → turquoise blue |
| <i>S.aureus</i> | → golden, opaque, small |
| <i>S.saprophyticus</i> | → pink, opaque, small |
| <i>StrepB</i> | → light blue |

CHROMagar™ MH Orientation

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE & LIMITATIONS

- Do not supplement with blood as it will affect the chromogenic coloration.
- Final identification may require additional testing such as biochemical or immunological test.

QUALITY CONTROL

Perform Quality Control according to the use of the medium and the local QC regulations and norms.

Good preparation of the medium can be tested, isolating the ATCC strains below:

| Microorganism | Typical colony coloration |
|------------------------------------|---------------------------|
| <i>E.coli</i> ATCC® 25922 | → reddish |
| <i>K.pneumoniae</i> ATCC® 13883 | → metallic blue |
| <i>E.faecalis</i> ATCC® 29212 | → turquoise blue |
| <i>S.aureus</i> ATCC® 12600 | → golden yellow |
| <i>S.saprophyticus</i> ATCC® 15305 | → pink |

WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For *in vitro* diagnostic use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121°C for at least 20 minutes.

REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity

Pack Size

| Pack Size | Ordering References |
|-------------------------------------|---------------------|
| 5000 ml <i>250 Tests of 20ml</i> | MH482 |
| 25 L <i>1250 Tests of 20ml</i> | MH483-25 |

| Base (B) | Supplement (S) |
|-------------------------------|----------------------------------|
| MH482(B) Weight: 197 gr | MH482(S) Volume: 20ml |
| MH483-25(B) Weight: 985 gr | 5 x MH482(S) Volume: 5 x 20ml |

Need some Technical Documents?

- Available for download on www.CHROMagar.com
- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach
ATCC® is a registered trademark of the American Type Culture Collection
NT-EXT-035 V3.0 / 01-Apr-16

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