# CHROMagar™ ECC

## Instructions For Use

Available in several languages

**NT-EXT-016**

Version **8.0**

Click below for:


**ENGLISH**

**FRANCAIS**

**ESPAÑOL**

**DEUTSCH**

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## MEDIUM PURPOSE

Chromogenic medium for the detection and enumeration of  $\beta$ -glucuronidase positive *E. coli* and coliforms in food, water and environmental samples.

Coliforms, *Enterobacteriaceae* able to ferment lactose (lactose positive *Enterobacteriaceae*), are bacteria present in human and warm blooded animals intestinal flora, in the soil and water. Coliforms are proof of organic, environmental or faecal contamination. Faecal contamination, due to coliforms coming from animal waste, consists mainly of *Escherichia coli* and thermotolerant *Klebsiella*. Strict regulations exist for *E. coli*/Coliform presence in water and food samples.

This can be explained by the importance of these germs in determining water and food safety.

## COMPOSITION

The product is composed of a single powder medium.

Product	=	Pack
Total g/L		32.8 g/L
Composition g/L		Agar 15.0 Peptone and yeast extract 8.0 NaCl 5.0 Chromogenic mix 4.8
Aspect		Powder Form
STORAGE		15/30 °C
FINAL MEDIA pH		7.2 +/- 0.2

## PREPARATION (Calculation for 1 L)

### Step 1

Preparation

- Disperse slowly 32.8 g of powder in 1 L of purified water.
- Stir until agar is thickened.
- Heat and bring to boil (100° C) while swirling or stirring regularly.

**Advice 1:** For the 100 °C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

**Advice 2:** in case of product samples containing a high load of *Pseudomonas* and/or *Aeromonas*, cefsulodin can be added at 7.5 mg/L.

### Step 2

Pour plates

- Cool in a water bath to 45-50 °C.
- Swirl or stir gently to homogenize.
- Pour medium into Petri dishes.
- Let it solidify and dry.

*If using pouring technique procedure, please refer to Inoculation part.*

### Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to 2 months under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

*Some tiny white crystals can appear after storage at 2/8 °C but **do not interfere** with the performances of the media and will vanish as soon as the plates are warmed at room temperature.*

## INOCULATION

Related samples (e.g. Processed food, raw materials, water, milk & environment) can be processed by direct streaking on the plate, as well as prior appropriate enrichment step.

### --> IF USING SURFACE TECHNIQUE PROCEDURE:

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak the sample or place the inoculated membranes on plate surface.
- Incubate in aerobic conditions at 37 °C during 24 h.

### --> IF USING POURING TECHNIQUE PROCEDURE:

- Prepare 90 mm Ø sterile Petri dishes and add 1 mL of inoculum in each.
- Then pour 10 mL of melted medium. Mix and let it solidify.
- Incubate in aerobic conditions at 37 °C during 24 h.

### Advice 3: Incubation possibilities:

If research is focused on faecal coliform bacteria

Incubate at 44 °C, 24 h


If research is targeted to maximise total coliform detection

Incubate at 30 °C, 24 h

## INTERPRETATION

Microorganism	Typical colony appearance
<i>E. coli</i>	→ blue
Other (faecal*) coliform bacteria	→ mauve
Other gram (-) bacteria	→ colourless

### Typical colony appearance

The pictures shown are not contractual.

## PERFORMANCE & LIMITATIONS

- Sensitivity for *E. coli* is 97 % (Ogden et al. 1991).
- Rare β-glucuronidase negative *E. coli* strains are false negative on this medium (typically O157 *E. coli*). *If research is focused on rare pathogenic strains such as O157 E. coli : please refer to CHROMagar™ O157 product.*
- If your research is focused on total coliform, few *Hafnia* are false negative and have a colourless appearance.

## QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
<i>E. coli</i> ATCC® 25922	→ blue
<i>C. freundii</i> ATCC® 8090	→ mauve
<i>E. cloacae</i> ATCC® 35030	→ mauve
<i>K. pneumoniae</i> ATCC® 13883	→ mauve
<i>S. aureus</i> ATCC® 25923	→ inhibited
<i>E. faecalis</i> ATCC® 29212	→ inhibited

## WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- Collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

## DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

## REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

## IFU/LABEL INDEX

- Catalogue reference
- Consult instructions for use
- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity
- Protect from light
- Manufacturer

### Need some Technical Documents?

Available for download on [www.CHROMagar.com](http://www.CHROMagar.com)

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

### Σ Pack Size

5000 mL

250 Tests of 20 mL

=

### Ordering References

EF322

Weight: 164 g

25 L

1250 Tests of 20 mL

=

EF323-25

Weight: 820 g

10 Kg

15243 Tests of 20 mL

=

EF323-10Kg

Weight: 10 Kg

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach  
ATCC® is a registered trademark of the American Type Culture Collection  
NT-EXT-016 V8.0 / EN 28-May-19

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