CHROMagar™ **Candida**

Instructions For Use
Available in several languages

NT-EXT-001

Version **9.0**


Click below for:

ENGLISH

FRANCAIS

ESPAÑOL

DEUTSCH

REFERENCES

Pack Size	Ordering References	Weight
5000 mL = 250 Tests of 20 mL	CA222	Weight: 238.5 g
25 L = 1250 Tests of 20 mL	CA223-25	Weight: 1192.5 g
10 kg = 10 450 Tests of 20 mL	CA223-10kg	Weight: 10 kg

MEDIUM PURPOSE

Chromogenic medium for the detection of *Candida* spp.

More commonly, *Candida* species are involved in oropharyngeal, digestive, urogenital infections and wet skin wounds. "Early diagnosis is essential for effective management of the patients." WHO Guidelines on Standard Operating Procedures for Laboratory Diagnosis of HIV-Opportunistic Infections.

Although *C. albicans* remains the major species involved, other types such as *C. tropicalis*, *C. krusei* or *C. glabrata* have increased proportionally as new antifungal agents have worked very effectively against *C. albicans*.

COMPOSITION

The product is composed of a powder base.

Product	Pack
Total g/L	47.7 g/L
Composition g/L	Agar 15.0 Peptone 10.2 Chloramphenicol 0.5 Chromogenic mix 22.0
Aspect	Powder Form
STORAGE	15-30 °C
FINAL MEDIA pH	6.1 +/- 0.2

Need some
Technical Documents?

Available
for download on
www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

PREPARATION (Calculation for 1 L)

Step 1

Preparation of the mix

- Disperse slowly 47.7 g of powder base in 1 L of purified water.
- Stir until agar is well thickened.
- Heat and bring to boiling (100 °C) while swirling or stirring regularly. DO NOT HEAT TO MORE THAN 100 °C. DO NOT AUTOCLAVE AT 121 °C.

Warning 1: If using an autoclave, do so without pressure.

Advice 1: For the 100 °C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

Step 2

Pouring

- Cool in a water bath to 45-50 °C, swirling or stirring gently.
- Pour medium into sterile Petri dishes.
- Let it solidify and dry.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to two months under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

CHROMagar™ Candida

SPECIMEN COLLECTION AND HANDLING

CHROMagar™ Candida can be used with the following specimens: Blood, skin, sputum, throat, urine, ears, stools and vaginal specimens.

Use of transport devices approved for collection of such specimens is recommended.

MATERIAL REQUIRED BUT NOT PROVIDED

Standard microbiological laboratory material for culture media preparation, control, streaking, incubation and waste disposal.

INOCULATION



Related samples can be processed by direct streaking on the plate.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate in aerobic conditions at 35-37 °C for 20-48 hours.

INTERPRETATION

Microorganism	Typical colony appearance
<i>C. albicans</i>	→ green
<i>C. tropicalis</i>	→ metallic blue
<i>C. krusei</i>	→ pink, fuzzy
<i>C. kefyr</i> , <i>C. glabrata</i>	→ mauve-brown
Other species	→ white to mauve

Typical colony appearance

C. tropicalis

The pictures shown are not contractual.

PERFORMANCE

In the following study, 382 pure strains and 339 samples (vaginal, oral and anorectal swabs) were tested and read after 48 h incubation at 37 °C.

	CHROMagar™ Candida
Sensibility	> 99 % *
Specificity	> 98 % *
NPV	100 % *
PPV	95 % *

* Data obtained from the study «CHROMagar™ Candida, a new differential isolation medium for presumptive identification of clinically important candida species» Odds F. et al. 1994, Journal of Clinical Microbiology, 32 : 1923-1929

LIMITATIONS AND COMPLEMENTARY TESTS

- Definite identification requires additional testing.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms.

Good preparation of the medium can be tested, isolating the following ATCC strains:

Microorganism	Typical colony appearance
<i>C. krusei</i> ATCC® 14243	→ pink and fuzzy
<i>C. albicans</i> ATCC® 60193	→ green
<i>C. tropicalis</i> ATCC® 1369	→ metallic blue
<i>C. glabrata</i> ATCC® 2001	→ mauve
<i>E. coli</i> ATCC® 25922	→ inhibited

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- This laboratory product should be used only by trained personnel (healthcare professional, etc). Wear appropriate protective clothing, gloves and eye/face protection and handle appropriately with procedures and good laboratory practices.
- Use of the medium may be difficult for people who have problems recognising colours.
- For a good microbial detection, collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.
- Culture media should not be used as manufacturing material or components.
- Do not ingest or inhale the product.
- Do not use the product after the expiry date.
- Do not use the product if it show any evidence of contamination or any sign of deterioration.
- Do not use the product if the packaging is damaged.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- Reading and interpretation should be performed using isolated colonies.
- Interpretation of the test results should be made taking into consideration colonial and microscopic morphology and if necessary, the results of any other tests performed.
- Laboratory, chemical or biohazardous wastes must be handled and discarded in accordance with all local and national regulations.
- For hazard and precaution recommendations related to some chemical components in this medium, please refer to the pictogram(s) mentioned on the labels. The Safety Data Sheet (SDS) is available on www.chromagar.com

CHROMagar™ Candida

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

LITERATURE REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX

-  Catalogue reference
-  Consult instructions for use
-  Quantity of powder sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

REVISION HISTORY

This is version V9.0 of this document

Changing version is related to the new 3 pages format of the IFU.

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