FecalSwab[®]

Instructions for Use

CE 0123

Copan FecalSwab™ Collection, Transport and Preservation System of Enteric Bacteria Instructions for use

INTENDED USE

The Copan FecalSwab™ Collection, Transport and Preservation System is intended for the collection of rectal swabs and fecal specimens and to preserve the viability of enteric pathogenic bacteria during transport from the collection site to the testing laboratory. In the laboratory, FecalSwab™ specimens are processed using standard clinical laboratory operating procedures for culture.

SUMMARY AND PRINCIPLES

Enteric infections can be caused by different types of bacteria. With such a wide array of pathogens and the need for cost containment, physician input and practice guidelines can help the laboratory determine which tests are appropriate for detecting the etiological agent of diarrhea. Microbiology laboratories should review the local epidemiology of bacterial enterocolitis and implement routine stool culture methods that will allow recovery and detection of all the major pathogens causing most of the cases in their geographic area. All microbiology laboratories should routinely test for the presence of *Salmonella* spp., *Shigella* spp., and *Campylobacter* spp. on all stool cultures⁽²⁾. One of the routine procedures in the diagnosis of enteric infections involves the collection and safe transportation of rectal swab samples or stool samples. This can be accomplished using the Copan FecalSwab™ Collection, Transport and Preservation System. Copan FecalSwab™ incorporates a modified Cary-Blair medium which is a non-nutritive transport and preservation medium containing chloride salts, sodium salts, phosphate buffer, L-Cysteine, agar and water. The medium is designed to maintain the viability of enteric pathogenic bacteria during transit to the testing laboratory.

Copan FecalSwab™ Collection, Transport and Preservation System is supplied in a sterile collection kit format comprising a package containing a tube filled with 2 ml of FecalSwab™ transport and preservation medium and a regular size nylon flocked specimen collection swab or a flocked specimen collection swab with plastic ring marker on the shaft which helps to visualize the maximum depth for the rectal sampling. The flocked swab applicator can be used to collect the clinical rectal specimen or as a transferring tool for stool specimens.

Once the sample is collected, it should be placed immediately into the FecalSwab™ transport tube where it comes into contact with the transport medium. Swab specimens for bacterial investigations collected using the FecalSwab™ should be transported directly to the laboratory, preferably within 2 hours from the collection^(1-5,16,17) to maintain optimum organism viability. If immediate delivery or processing is delayed, then specimens should be refrigerated at 2-8°C and processed within 72 hours or stored at room temperature (20-25°C) and processed within 48 hours. In case of *C. difficile* culture investigation, specimens should be refrigerated at 2-8°C and processed within 48 hours or stored at room temperature (20-25°C) and processed within 24 hours. Independent scientific studies on swab transport systems have shown that for certain bacteria viability is superior at refrigerated temperatures compared with room temperature^(7,9-10).

REAGENTS

FecalSwab™ Transport and Preservation Medium

Chloride salts
Sodium salts
Phosphate buffer
L-Cysteine
Agar
Distilled water

WARNINGS AND PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
3. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
4. All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimens, containers and media after their use. Observe other CDC Biosafety Level 2 recommendations^(8,12-15).
5. It must be assumed that all specimens contain infectious micro-organisms and materials used to process them should be considered potentially infectious. Therefore they must be handled with appropriate precautions and dispose according to laboratory regulations for infectious waste.
6. Directions should be read and followed carefully.
7. Tube label shows a maximum filling line. If sample collected exceeds maximum fill line discard the swab and the tube. A second specimen should be collected using a new Copan FecalSwab™ kit.
8. Do not re-sterilize unused products.
9. Do not re-pack.
10. Not suitable to collect and transport microorganisms other than enteric pathogenic bacteria.
11. Not suitable for any other application than intended use.
12. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
13. Do not use if the swab is visibly damaged (i.e., if the swab tip or swab shaft is broken).
14. Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
15. Do not ingest the medium.
16. Do not use the FecalSwab™ medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
17. The use by the patient of antacids, barium, bismuth, anti-diarrheal medication, antibiotics, histamine, nonsteroidal anti-inflammatory drug or oily laxatives prior the specimen collection might affected the results of the analysis.
18. Stool specimen should not contain urine or water.

STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 5-25°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual collection unit and the specimen transport tube label.

PRODUCT DETERIORATION

Copan FecalSwab™ should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed, (4) the package is open, or (5) there are other signs of deterioration.

SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Rectal swab specimens and stool specimens collected for microbiological investigations which comprise the isolation of enteric pathogenic bacteria should be collected and handled following published manuals and guidelines⁽¹⁻⁵⁾.

To maintain optimum organism viability, transport specimens collected using FecalSwab™ directly to the laboratory, preferably within 2 hours of collection^(1-5,16,17). If immediate delivery or processing is delayed, then specimens should be refrigerated at 2-8°C and processed within 72 hours or stored at room temperature (20-25°C) and processed within 48 hours. In case of *C. difficile* culture investigation, specimens should be refrigerated at 2-8°C and processed within 48 hours or stored at room temperature (20-25°C) and processed within 24 hours.

Specific requirements for the shipment and handling of specimens should be in full compliance with local regulations⁽¹²⁻¹⁵⁾. Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

MATERIALS SUPPLIED

Fifty (50) FecalSwab™ collection kits are contained in a “Vi-Pak” shelf pack and 10 x 50 units are contained in a box. Each collection kit consists of a package containing a plastic screw-cap tube with conical or round shaped bottom filled with 2 ml of FecalSwab™ transport and preservation medium and a specimen collection swab that has a tip flocked with soft nylon fiber

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolating and culturing enteric pathogenic bacteria. These materials include culture media plates or tubes and incubation systems. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques for enteric pathogenic bacteria from clinical swab samples ⁽¹⁻⁵⁾.

DIRECTIONS FOR USE

The Copan FecalSwab™ Collection, Transport and Preservation System is available in product configurations indicated in the table below.

Table 1

Catalog No.	Copan FecalSwab™ Product Descriptions	Pack Size	Suitable for Automation	Capture Cap Feature
470CE	Sterile single use sample collection pack containing: - Green screw-cap tube with internal conical shape filled with 2ml of modified Cary-Blair medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per ‘Vi-pack’ 10x50 units per box	NO	YES
470CE.A	Sterile single use sample collection pack containing: - Green screw-cap tube with internal conical shape filled with 2ml of modified Cary-Blair medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per ‘Vi-pack’ 10x50 units per box	YES	YES
4E020S.A	Sterile single use sample collection pack containing: - Green screw-cap tube with round shape filled with 2 ml of Modified Cary-Blair medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per ‘Vi-pack’ 10x50 units per box	YES	YES
4E048S	Single use sample collection pack containing: - Polypropylene green screw-cap tube with internal conical shape filled with 2ml of FecalSwab™ Medium. - One applicator swab with flocked nylon fiber tip and a plastic ring marker on the shaft	50 units per ‘Vi-pack’ 10x50 units per box	NO	YES
4E049S	Single use sample collection pack containing: - Polypropylene orange screw-cap tube with internal conical shape filled with 2ml of FecalSwab™ Medium. - One applicator swab with flocked nylon fiber tip and a plastic ring marker on the shaft	50 units per ‘Vi-pack’ 10x50 units per box	NO	YES


Specimen Collection

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. The patient should be cautioned against the use of antacids, barium, bismuth, anti-diarrheal medication, antibiotics, histamine, nonsteroidal anti-inflammatory drug or oily laxatives prior to collection of the specimen ⁽²²⁾.

For Rectal Swab Collection:

1. Peel open the kit package and remove the tube of medium and the flocked swab applicator (see Figure 1a).
2. Use the flocked swab to collect the clinical specimen. The operator must touch the swab applicator only above the marked breakpoint line (the area from the line to the end of the swab shaft), as illustrated in Figure 2, which is the opposite end to the nylon fiber tip. At all times when handling the swab applicator, the operator must not touch the area below the marked breakpoint line as this will lead to contamination of the applicator shaft and the subsequent culture thus invalidating the test results.
3. Insert the flocked swab through the rectal sphincter 2.5 to 3.5 cm and gently rotate ⁽¹⁷⁾. **NOTE: the rectal flocked swab has a plastic ring on the shaft which is a marker for the maximum depth for the rectal sampling. Do not insert the rectal swab beyond this marker.**
4. Withdraw and examine to make sure there is fecal material visible on the tip ⁽¹⁸⁾.
5. After collection, transfer the swab into the tube with the preservation medium and visually check that the maximum filling line (“MAX. FILL”) indicated on the label is not exceeded. **NOTE: If sample collected exceeds maximum fill line discard the swab and the tube. A second specimen should be collected using a new Copan FecalSwab™ kit.**
6. Holding the swab shaft between thumb and finger, mash and mix the stool specimen against the side of the tube to evenly disperse and suspend the specimen in the preservation medium making sure not to smear stool in the outside of the tube.
7. Break the swab off into the tube as follows:
 - Hold the tube in one hand pointing the opening of the tube away from your face.
 - With the other hand grasp the swab shaft at the very end with the thumb and first finger.
 - Lean the part of the shaft with the breaking point against the rim of the tube.
 - Bend the swab shaft at a 180 degrees angle to break it off at the breakpoint. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
 - Discard the broken upper part of the swab shaft.
8. Place the screw cap on the tube and tighten.
9. Shake the vial until the sample appears homogeneous.

Fig 1a. Operation of the FecalSwab™ Collection, Transport and Preservation System For Rectal Swab Collection

For Stool Specimen Collection:

1. Have the patient obtain stool specimen. Stool specimen should not contain urine or water. Patient should pass the stool into a clean, dry pan or a special container mounted on the toilet for this purpose. **NOTE: Toilet paper is not suitable to collect stool because it may be impregnated with barium salts, which are inhibitory to some fecal pathogens** ⁽¹⁸⁾.
2. Peel open the kit package and remove the tube of medium and the flocked swab applicator (see Figure 1.b).
3. The operator must touch the swab applicator only above the marked breakpoint line (the area from the line to the end of the swab shaft), as illustrated in Figure 2, which is the opposite end to the nylon fiber tip. At all times when handling the swab applicator, the operator must not touch the area below the marked breakpoint line (the area from the line to the tip of the nylon flocked swab) as this will lead to contamination of the applicator shaft and the subsequent culture thus invalidating the test results.
4. Collect a small amount of stool by inserting all the tip of the flocked swab into stool sample and rotate it ⁽¹⁸⁾. Bloody, slimy or watery area of stools should be selected and sampled ^(16, 19, 20, 21).
5. After collection examine the swab to make sure there is fecal material visible on the tip ⁽¹⁶⁾. In case it is not, insert again the flocked swab into stool sample and rotate taking care all the area of the swab tip is in contact with the sample. **NOTE: the swab should not be used as a paddle or spoon but as a probe. DO NOT try to collect and transfer an excessive amount of fecal sample into the transport medium tube. The swab tip only needs to be coated with sample material.**
6. After collection, transfer the swab into the tube with the preservation medium and visually check that the maximum filling line ("MAX. FILL") indicated on the label is not exceeded. **NOTE: If sample collected exceeds maximum fill line discard the swab and the tube a second specimen should be collected using a new Copan FecalSwab™ kit.**
7. Holding the swab shaft between thumb and finger, mash and mix the stool specimen against the side of the tube to evenly disperse and suspend the specimen in the preservation medium making sure not to smear stool in the outside of the tube.
8. Break the swab off into the tube as follows:
 - Hold the tube in one hand pointing the opening of the tube away from your face.
 - With the other hand grasp the swab shaft at the very end with the thumb and first finger.
 - Lean the part of the shaft with the breaking point against the rim of the tube.
 - Bend the swab shaft at a 180 degrees angle to break it off at the breakpoint. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
 - Discard the broken upper part of the swab shaft.
9. Place the screw cap on the tube and tighten.
10. Shake the vial until the sample appears homogeneous.
11. Write patient's name and demographics on the tube label and send sample to laboratory

Fig 1.b Operation of the FecalSwab™ Collection, Transport and Preservation System For Stool Specimen Collection



Fig 2. Collection swab showing breakpoint indication line and area for holding the applicator


REGULAR FLOCKED SWAB

Do not touch the applicator in the area below the breakpoint indication line

During specimen collection hold applicator above the breakpoint indication line, in this area

FLOCKED SWAB WITH PLASTIC RING

The operator must only handle the part of the swab applicator shaft above the breakpoint indication line as shown in Fig 2. After the swab sample is taken from the patient or from the stool specimen, the maximum filling line checked, the swab applicator shaft is broken off at the marked breakpoint indication line into the FecalSwab™ tube of transport medium. The operator then discards the handle part of the swab. The tube's screw cap is then replaced and secured tightly. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap (see Fig 3). This molded funnel shape captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip.

Fig 3. Capture of broken swab applicator stick by FecalSwab™ tube cap

In the testing laboratory when the FecalSwab™ cap is unscrewed and removed, the swab applicator stick is securely attached to the cap. This feature allows the operator to conveniently remove the swab and perform various microbiology analyses using the tube cap as a handle to hold and manipulate the swab.

Processing FecalSwab™ Specimens in the Laboratory-Bacteriology

FecalSwab™ samples should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines (1-5).

Culture investigations of swab specimens for the presence of bacteria routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for inoculation of FecalSwab™ samples onto solid agar in Petri dishes is as follows.

Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (8, 12-15).

Vortex mix the FecalSwab™ tube containing the swab sample for 5 seconds to evenly disperse and suspend the patient specimen in the medium

1. Unscrew the FecalSwab™ cap and remove the swab applicator.
2. Roll the tip of the FecalSwab™ applicator onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
3. If it is necessary to inoculate additional culture media plates, return the FecalSwab™ applicator to the transport medium tube for two seconds and recharge the swab applicator tip with transport medium/patient sample suspension then repeat Step No. 2.

The procedure described above utilizes the FecalSwab™ applicator like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum (see Fig 4-1).

Alternatively, the operator can vortex mix the FecalSwab™ tube with the swab inside for 5 seconds and then transfer 100µl volumes of the suspension onto each culture plate using a volumetric pipettor and sterile pipet tips (see Fig 4-2). Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 5).


Note that following approved internal laboratory procedures, different volumes of suspension can be transferred onto each culture plate using inoculating loops, volumetric pipettor and sterile pipet tips.

Fig 4. Procedures for inoculation of FecalSwab™ specimens onto solid agar in Petri dishes

1. Using swab to inoculate specimen
2. Using pipettor and sterile pipet tips to inoculate 100µl of specimen

Fig 5. Procedure for streaking FecalSwab™ specimens on agar Petri dishes for primary isolation (15)

Seed a primary inoculum of FecalSwab™ specimen onto the surface of an appropriate agar culture plate in the first quadrant.

Use a sterile bacteriology loop to streak the primary inoculum across the surface of the second, third and fourth quadrants of the agar culture plate.

Processing FecalSwab™ specimens with automatic systems

Copan FecalSwab™ code 470CE.A and 4E020S.A is compatible and can be processed with automated systems for planting and streaking on traditional culture media plates such as Copan WASP® and WASPLab®. Please refer to automation manufacturer's instructions for FecalSwab™ processing.

LIMITATIONS

1. For recovery of *C. difficile* FecalSwab™ specimens should be refrigerated at 2-8°C and processed within 48 hours or stored at room temperature (20-25°C) and processed within 24 hours.
2. In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations (8, 12-15) when handling or analyzing patient samples.
3. Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection (1-5, 7, 9, 10).
4. FecalSwab™ is intended for use as a collection and transport medium for enteric pathogenic bacteria. FecalSwab™ cannot be used as enrichment, selective or differential medium.
5. Performance testing with Copan FecalSwab™ was conducted using laboratory strains spiked onto a swab following the test protocols based on those described in CLSI M40-A2(6).
6. The rectal collection and the stool specimen collection shall be performed using Copan flocked swab included into FecalSwab™ pouch.

RESULTS

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

PERFORMANCE CHARACTERISTICS

The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in CLSI M40-A2⁽⁶⁾. FecalSwab™ system has an intended use limited to enteric pathogenic bacteria, therefore its field application is restricted to fecal/stool specimen. For this reason the bacterial recovery studies were conducted with the following microorganisms, described and defined in CLSI M40-A2, *Quality Control of Microbiological Transport Systems: Approved Standard* and included only the following enteric pathogenic bacteria strains from the paragraph 9.3.1 of the CLSI M40-A2 document, in particular:

<i>Escherichia coli</i>	ATCC® 25922
<i>Salmonella typhimurium</i>	ATCC® 14028
<i>Shigella flexneri</i>	ATCC® 12022

In addition, Copan included testing of additional clinically relevant enteric pathogenic bacteria. The specific bacterial strains used in these studies are here listed:

<i>Campylobacter jejuni</i>	ATCC® 33291
<i>Yersinia enterocolitica</i>	ATCC® 9610
<i>Vibrio parahaemolyticus</i>	ATCC® 17802
<i>Enterococcus faecalis vancomycin resistant (VRE)</i>	ATCC® 51299
<i>Escherichia coli</i> O157:H7	ATCC® 700728
<i>Clostridium difficile</i>	ATCC® 9689

In addition, Copan included an evaluation of overgrowth as per CLSI M40-A2 document with the following strain:

<i>Pseudomonas aeruginosa</i>	ATCC® BAA-427
-------------------------------	---------------

All cultures of bacteria were ATCC® (American Type Culture Collection) and were obtained commercially.

The selection of these organisms also reflects those enteric pathogenic bacteria of potential interest from rectal and stool specimens that can be present in specimens collected and analyzed in a typical clinical microbiology laboratory.

Bacterial viability studies were performed on inoculated Copan FecalSwab™ hold at two different ranges of temperature, 2-8°C and 20-25°C, corresponding to cold temperature and controlled room temperature, respectively. FecalSwab™ collection swabs were inoculated with 100µl of specific concentrations of organism suspension. Swabs were then placed in transport medium tubes and were held for 0 hrs, 6hrs, 24 hrs, 48 hrs and 72 hrs (72 hours is applicable only for the 2-8°C cold temperature range). At the appropriate time intervals, each transport tube was processed swabbing onto appropriate agar plate.

Additional viability tests were performed on *E. coli* O157:H7 ATCC® 700728, *Vibrio parahaemolyticus* (ATCC® 17802) and *Salmonella typhimurium* (ATCC® 14028) with fecal matrix in order to test the device under the conditions of its intended use.

TEST RESULTS

Viability performance is measured for each test organism at the 48 hrs time point when tubes were held at controlled room temperature (20-25°C), or at the 72 hrs time point when held at cold temperature (2-8°C), and compared with the acceptance criteria.

Copan FecalSwab™ System was able to maintain acceptable recovery of all organisms evaluated at both cold temperature (2-8°C) and controlled room temperature (20-25°C). In case of *C. difficile* culture investigation, specimens should be refrigerated at 2-8°C and processed within 48 hours or stored at room temperature (20-25°C) and processed within 24 hours. Acceptable recovery is defined as a plate count that remains within 2 log₁₀ of the initial microorganism concentration for each microorganism tested.

SUMMARY OF RESULTS FOR BACTERIAL RECOVERY IN PBS*

Organism*	COLD TEMPERATURE (2-8°C)			CONTROLLED ROOM TEMPERATURE (20-25°C)		
	Average CFU recovered Time 0 hrs	Average CFU recovered Time 72hrs	T=72 hrs Log reduction (-) or Log increase (+)	Average CFU recovered Time 0 hrs	Average CFU recovered Time 48 hrs	T=48 hrs Log reduction (-) or Log increase (+)
<i>Escherichia coli</i> ATCC® 25922	1.07E+02	7.29E+01	-0.17	1.07E+02	5.72E+03	1.73
<i>Escherichia coli</i> O157:H7 ATCC® 700728	8.99E+01	1.04E+02	0.06	8.99E+01	4.12E+03	1.66
<i>Salmonella typhimurium</i> ATCC® 14028	1.38E+02	9.70E+01	-0.15	1.38E+02	9.72E+03	1.85
<i>Shigella flexneri</i> ATCC® 12022	1.27E+02	1.16E+02	-0.04	1.27E+02	9.67E+03	1.88
<i>Vibrio parahaemolyticus</i> ATCC® 17802	2.00E+02	1.54E+02	-0.11	2.00E+02	1.58E+04	1.90
<i>Enterococcus faecalis</i> <i>vancomycin resistant (VRE)</i> ATCC® 51299	1.68E+02	1.16E+02	-0.16	1.68E+02	2.26E+03	1.13
<i>Yersinia enterocolitica</i> ATCC® 9610	1.17E+02	1.04E+02	-0.05	1.17E+02	9.78E+03	1.92
<i>Campylobacter jejuni</i> ATCC® 33291	2.14E+02	3.42E+01	-0.80	2.14E+02	4.28E+00	-1.70

* the organism was diluted in PBS and the pure suspension tested in studies based on CLSI M40-A2

Organism*	COLD TEMPERATURE (2-8°C)			CONTROLLED ROOM TEMPERATURE (20-25°C)		
	Average CFU recovered Time 0 hrs	Average CFU recovered Time 48 hrs	T=48 hrs Log reduction (-) or Log increase (+)	Average CFU recovered Time 0 hrs	Average CFU recovered Time 24 hrs	T=24 hrs Log reduction (-) or Log increase (+)
<i>Clostridium difficile</i> ATCC® 9689	4.42E+01	6.30E-01	-1.85	4.42E+01	5.30E-01	-1.92

* the organism was diluted in PBS and the pure suspension tested in studies based on CLSI M40-A2

SUMMARY OF RESULTS FOR BACTERIAL RECOVERY IN FECAL MATRIX

Organism*	COLD TEMPERATURE (2-8°C)			CONTROLLED ROOM TEMPERATURE (20-25°C)		
	Average CFU recovered Time 0 hrs	Average CFU recovered Time 72 hrs	T=72 hrs Log reduction (-) or Log increase (+)	Average CFU recovered Time 0 hrs	Average CFU recovered Time 48 hrs	T=48 hrs Log reduction (-) or Log increase (+)
<i>Escherichia coli</i> ATCC® 25922	1.23E+02	1.60E+02	0.11	1.23E+02	7.54E+03	1.79
<i>Salmonella typhimurium</i> ATCC® 14028	9.46E+01	1.43E+02	0.18	9.46E+01	6.38E+03	1.83
<i>Vibrio parahaemolyticus</i> ATCC® 17802	1.11E+02	1.36E+02	0.09	1.11E+02	7.36E+03	1.82

* the organism was diluted in fecal matrix and the suspension tested in studies based on CLSI M40-A2

SUMMARY OF RESULTS FOR BACTERIAL OVERGROWTH STUDY AT 4-8°C

Organism*	COLD TEMPERATURE (4-8°C)		
	Average CFU recovered Time 0 hrs	Average CFU recovered Time 48hrs	T=48hrs Log reduction (-) or Log increase (+)
<i>Pseudomonas aeruginosa</i> ATCC® BAA-427	3.79E+01	3.22E+01	-0.07

* the organism was diluted in PBS and the pure suspension tested in studies based on CLSI M40-A2

Please refer to symbol table at the end of the instructions for use

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BIBLIOGRAPHY

1. James Versalovic, Karen C. Carroll, Guido Funke, James H. Jorgensen, Marie Louise Kandry, David W. Warnock, eds. *Manual of Clinical Microbiology*. 10th edition. Washington, DC: ASM; 2011.
2. Lynne S. Garcia. (Editor in Chief). *Clinical Microbiology Procedures Handbook*. 3rd edition. American Society for Microbiology, Washington, DC, 2007.
3. Isenberg HD. (Editor in Chief). *Essential Procedures for Clinical Microbiology*. American Society for Microbiology, Washington DC, 1998.
4. C.H. Collins, Patricia M. Lyne, J.M. Grange. *Microbiological Methods*. Sixth Edition. Butterworths, 1989.
5. Miller JM. *A Guide to Specimen Management in Clinical Microbiology*. Second Edition. American Society for Microbiology. Washington, DC. 1999.
6. Quality Control of Microbiology Transport Systems; Approved Standard-Second Edition. M40-A2 Vol 34 No. 9. June 2014. Clinical and Laboratory Standards Institute (formerly NCCLS).
7. Kevin Gough, MD, Michelle Alfa, PhD, and Godfrey Harding, MD. Evaluation of routine enteric pathogens in hospitalized patients: A Canadian perspective. *The Canadian Journal of Infectious Diseases* Pulsus Group.
8. Baron EJ, Miller JM, Weinstein MP, Richter SS, Gilligan PH, Thomson RB, Jr., et al. *A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2013*
9. Public Health England. UK Standard for Microbiology Investigations: Investigation of Faecal Specimens for Enteric Pathogens. Standard Unit, Microbiology Services, PHE. *Bacteriology B30*, Issue no: 8.1, 24 April 2014.
10. M. Wasfy, B. Oyoyo, A. Elgindy, A. Churilla. Comparison of Preservation Media for Storage of Stool Samples. *Journal of Clinical Microbiology*, Aug. 1995.
11. W. Nugent, M.J. Russell, S. Beck and G.P. Leonard. Nassau University Medical Center (NY). Evaluation of the physical and diagnostic nature of swabs. 23rd Annual Clinical Virology Symposium, Clearwater Beach (FL), 2007.
12. Fleming D. *Biological Safety: Principles and Practices*. January 2000. ASM, Washington DC.
13. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/jrtext.htm>.
14. Richardson JH. *Biosafety in Microbiological and Biomedical Laboratories*. December 1994. Diane Publishing Company.
15. Hansen DJ. *Healthcare, Laboratories and Biosafety*. Vol 2., 1992. CRC Press.
16. Peter H. Gilligan, J. Michael Janda, Mohamed A. Karmali, J. Michael Miller, *Laboratory Diagnosis of Bacterial Diarrhea CUMITECH 12°*, ASM April 1992
17. www.cdc.gov/foodsafety/outbreaks/investigating-outbreaks/specimen-collection.html#1
18. JOB AIDS: How To Collect a Fecal Specimen And Transfer To Transport Medium Centers for Disease Control and Prevention
19. Bailey and Scott, *Diagnostic Microbiology* Chapter 11 01st edition of Bailey and Scott's *Diagnostic Microbiology* (1962)
20. Patzi-Vargas S¹, Zaidi M, Bernal-Reynaga R, León-Cen M, Michel A, Estrada-Garcia TJ. Persistent bloody diarrhoea without fever associated with diffusely adherent *Escherichia coli* in a young child. *Med Microbiol*. 2013 Dec;62(Pt 12):1907-10.
21. Marchou B. Traveller's diarrhea: epidemiology, clinical practice guideline for the prevention and treatment *Presse Med*. 2013 Jan;42(1):76-81.
22. Diagnostic Procedure for Stool Collection, Centers for Disease Control and Prevention Richard K. Root, *Clinical Infectious Disease* 1999 p. 18

INDEX OF SYMBOLS

Symbol / Simbolo / Símbolo / Symbol / Symbole / Símbolos / Symbol / Szimbólum / Symbol / Symbol / Símbol / Символ	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Betyder / Jelentés / Betydning / Znaczenie / Značenje / Betydelse / Anlami / Значення
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Producent / Gyártó / Produsent / Producent / Proizvođač / Tillverkare / Üretici / Виробник
	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro / Dispositivo de diagnóstico in vitro / Diagnostisk medicinsk anordning in vitro / In vitro diagnosztikai eszköz / In vitro diagnostisk medisinsk utstyr / Urządzenie do diagnostyki in vitro / In vitro dijagnostički uređaj / Medicinsk anordning för in vitro-diagnostik / In vitro tani cihazı / Медицинский прибор для диагностики in vitro
CE 0123	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Identifikationsnummer for det bemyndigede organ / A bejelentett szervezet azonosító száma / Ateranvänd inte / Tekrar kulanman / Numer identyfikacyjny jednostki notyfikowanej / Identifikacioni broj prijavljenog tela / Det anmälda organets identifieringskod / Onaylı kurumun tanımlama numarası / Ідентифікаційний номер уповноваженого органу
	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiação ionizante / Sterilisationsmåde: Bestråling / Ionizáló sugárzással sterilizált / Sterilisert ved bruk av ioniserende stråling / Sterylizacja za pomocą promieniowania jonizującego / Sterilisanu upotrebom jonizujućeg zračenja / Steriliserad med strålning / İşinlana kullanılan yöntem / Стерилизований завдяки опроміненню
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Brug det ikke igen / Ne használja újra / Må ikke gjenbrukes / Nie używać ponownie / Nemojte ponovo koristiti / Ateranvänd inte / Tekrar kulanman / Не використовувати повторно
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Katalognummer / Katalog száma / Katalognr / Numer katalogowy / Kataloški broj / Katalognummer / Katalog numarası / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Límites de temperatura / Límites de temperatura / Temperaturgrænser / Hőmérséklet határok / Temperaturgrenser / Limity temperature / Opseg temperature / Temperaturgränser / Sıcaklık limitleri / Температурні обмеження
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Anvendes før / Felhasználható / Må brukes innen / Limity temperature / Uputrebiti do / Sista förbrukningsdag / Son kullanma tarihi / Використати до
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultar as instruções de utilização / Se brugsvejledningen / Olvassa el a használati utasítást / Se instruksjoner for bruk / Zapoznać się z instrukcją użycia / Pročitajte uputstvo za upotrebu / Läs bruksanvisningen / Kullanma talimatlarına bakınız / Користуватись інструкціями по використанню
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacar / Riv op for at åbne / Tépje le a nyitáshoz / Riv opp for at åpne / Oderwać / Ogulite / Dra för att öppna / Açmak için yırtın / Відкривати тут
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Seriennummer (parti) / Kötégkód (kötég) / Kod partii (Seria) / Lot-nummer (parti) / Serijski (lot) broj / Partinummer (lotto) / Lot (seri) kodu / Код серії (партії)
	Contains sufficient for <n> tests / Contenu suffisant pour <n> tests / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Inhold tilstrækkelig til <n> prøver / A tartalom elegendő az <n> tesztekhez / Innhold tilstrekkelig for <n> prøver / Zawiera materiał pozwalający na wykonanie <n> testów / Sadržaj dovoljan za testova / Innehåller tillräckligt för <n> tester / <n> test için yeterli içerik / Придатний для проведення <n> кількості випробувань
	Handle with care / Maneggiare con cura / Tratar con cuidado / Mit Vorsicht behandeln / Manipuler avec soin / Manusear com cuidado / Håndtag med omhu / Kezelje körültekintően / Behandles forsigtig / Ostrożnie / Rukovati pažljivo / Hanteras varsamt / Dikkatlice muamele edin / Обережно поводитися
	Do not use if package is damaged / Non utilizzare in caso di confezionamento danneggiato / No utilizar en caso de paquete dañado / Bei Beschädigung der Verpackung nicht verwenden / Ne pas utiliser si l'emballage est abîmé / Não utilizar se a embalagem estiver danificada / Må ikke anvendes, hvis emballagen er beskadiget / Sérült csomagolás esetén ne használja / Må ikke brukes hvis emballasjen er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Nemojte koristiti ako je ambalaža oštećena / Använd inte vid skadad förpackning / Ambalajı hasar görmüşse kullanmayın / Не використовувати, якщо упаковка пошкоджена

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia, Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com