IVD solutions through partnership



CHROMagar[™] Candida

For isolation and differentiation of major clinical-significant *Candida* species



CHROMagar[™] Candida

Plate Reading

Candida albicans \rightarrow areen

 Candida tropicalis → metallic blue

 Candida krusei → pink. fuzzv



Manufacturer:

CHROMagar 4 place du 18 juin 1940 75006 Paris - France e-mail: CHROMagar@CHROMagar.com www.CHROMagar.com

Distributed by: Mast Diagnostica GmbH Feldstraße 20 DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0 Fax: +49 (0)4533 2007 68 e-mail: mast@mast-diagnostica.de www.mast-group.com

For isolation and differentiation of major clinical-significant Candida species

Background

Perhaps as a result of progress in medical therapy with a related decrease in patient resistance and increase in the number of susceptible hosts, the incidence of bloodstream Candida infections rose through the 1990s but has stabilised in recent years. More commonly, Candida species are involved in superficial oropharyngeal and urogenital infections particularly for immuno-depressed people such as the elderly and AIDS victims. "Early diagnosis is essential for early effective management of the patients." (WHO Guidelines on Standard Operating Procedures for Laboratory Diagnosis of HIV-Opportunistic Infections)

Although C.albicans remains the major species involved, other types such as C.tropicalis, C.krusei or C.glabrata have increased proportionately as new antifungal agents have worked very effectively against C.albicans. This shows the importance of an accurate detection for a proper antifungal therapy choice.

Medium Performance

High contrast

Intense colony colouration helping to differentiate species.

Fast

3

4

5

20-48 h incubation.

Panoramic view

Panoramic view of the Candida population mix with the ability to detect major species but also the presence of minor populations.

Powerful tool

It displays high specificity and sensitivity for 3 of the major Candida species at a time.

Specificity and sensitivity:

Candida albicans → Exceed 99 %* Candida tropicalis \rightarrow Exceed 99 %* Candida krusei \rightarrow Exceed 99 %^{*}

*Specificity and sensitivity from scientific study: "CHROMagar Candida, a new differential isolation medium for the presumptive identification of clinically important Candida species." Odds F. et al 1994, Journal of Clinical Microbiology, 32: 1923-1929.

Less workload

With traditional media such as Sabouraud agar, identification of mixed cultures is difficult, whereas CHROMagar™ Candida, just by the colour of the colony, will instantly differentiate various Candida species and help locating their colonies in mixed cultures.

Medium Description

Powder Base	Total
Usual Samples	Blood, skin, sputum, throat, urine, ears, stools and vaginal specimens.
Procedure	Direct Streaking. Incubation 18-48 h at 35-37 °C. Aerobic conditions.
Scientific Publications on this product: available on www.CHROMagar.com Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com	

Ordering Information P

Product	Order Code	
CHROMagar™ Candida dry media, 5 liter	15CA222	
CHROMagar™ Candida ready to use plates, 20 pcs.	201400	

CHROMagar is a trademark of Dr A. Rambac

Version 6.0 / Nov-19 (Mast)

LF-EXT-001