
IVD solutions through partnership



CHROMagar[™] MRSA

For isolation and differentiation of Methicillin Resistant *Staphylococcus aureus* (MRSA) including low level MRSA





● CHROMagar[™] MRSA

For isolation and differentiation of Methicillin Resistant *Staphylococcus aureus* (MRSA) including low level MRSA

Background

Leading cause of nosocomial infections, especially in intensive care units, the MRSA sources are either endogenous (the patient) or through cross contamination (environmental or by person to person contact).

The major issue with this pathogen is its resistance to a large panel of antibiotics, among them beta-lactam antibiotics, limiting the therapeutic options for clinicians.

Early detection is essential for controlling the spread of MRSA, providing appropriate care, and avoiding complex and expensive treatments. Pre-admission screening for MRSA has proved to be an effective method for reducing the hospital burden of MRSA-colonised patients. The savings due to consistent decolonisation before elective admission outweigh the costs of screening. Today, in the US, the extra-expenses linked to difficult treatments of MRSA infections are estimated at \$2.4 billion for about 370,000 hospital stays. (Genetic Engineering and Biotechnology News, August 2009).

In the UK, the estimation of the additional cost of discharging every hospital patient who acquires MRSA is £9,000.

Medium Performance

Absolutely reliable

CHROMagar[™] MRSA, introduced in 2002, was the <u>first</u> chromogenic medium for MRSA detection. It lead to such significant reductions in both, the response time and laboratory workload, that it allowed an absolutely necessary wide-scale patient screening.

Efficient

The medium exhibits sensitivity and specificity values close to 100%. CHROMagar[™] MRSA allows an accurate detection of MRSA with a higher level of sensitivity than oxacillin containing media.

Fast & easy interpretation

Intense mauve \rightarrow colour in 18-24h.

Medium Description

Total 82.5 g/L Agar 15.0 Peptones and yeast extract. 40.0 Salts 25.0 Chromogenic mix 2.5 Storage at 15/30 °C - pH: 6.9 +/-0.2 5 Shelf Life > 18 months	
Powder form qsf 20 L20 mLStorage at 2/8 °CShelf Life > 18 months	
nasal, perineal, throat, rectal specimens	
Direct Streaking. Incubation 18-24 h at 35-37 °C. Aerobic conditions	

Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Plate Reading

- Methicillin Resistant
- Staphylococcus aureus (MRSA) \rightarrow rose to mauve
- → IOSE to IIIauve
- Methicillin Susceptible
 Staphylococcus aureus (MSSA)
 → inhibited
- Other bacteria
 → blue, colourless or inhibited





Manufacturer:

CHROMagar 4 place du 18 juin 1940 75006 Paris - France e-mail: CHROMagar@CHROMagar.com www.CHROMagar.com

Distributed by: Mast Diagnostica GmbH Feldstraße 20 DE-23858 Reinfeld Tel: +49 (0)4533 2007 0 Fax: +49 (0)4533 2007 68 e-mail: mast@mast-diagnostica.de www.mast-group.com

Ordering Information

	Product	Order Code
	CHROMagar™ MRSA dry media, 5 liter	15MR502
	CHROMagar™ MRSA ready to use plates, 20 pcs.	201402