**BC+™**

## Package insert and How to use guide

# BC+™

ENGLISH

## Copan BC+™ Information Leaflet & User Guide

See the glossary of symbols at the end of the information leaflet

### INTENDED USE

Copan BC+™ is a system for the transfer and transport of positive blood cultures for a subsequent use for culturing or microscopic investigation.

### SUMMARY AND PRINCIPLES

#### Blood culture

The culturing of blood is the recommended method for detecting the presence of organisms in the blood. The culturing of blood is essential for the microbiological diagnosis of bacteremia, fungemia, endocarditis and many other diseases that occur with a fever of unknown origin. It's also essential for the diagnosis of infections of prosthetic materials (joint or vascular) as well as sepsis associated with central venous catheters. Blood cultures may also detect bloodstream infections (BSI) associated with other diseases such as pneumonia, septic arthritis or osteomyelitis<sup>1,2</sup>. The recent increase in antibiotic resistance, especially with Gram negative antibiotics, is often the cause of empirical treatment failure in BSI. A fast and correct microbiological diagnosis is crucial and provides valuable information on correct patient management. A reduction in the turn around time in each stage of the process from transportation to reporting is a fundamental objective for optimal blood culture management<sup>1,2,3</sup>. If vials incubated in continuous monitoring instruments give positive results, a microscopic examination (Gram stain) is generally prepared along with seeding onto solid culture medium<sup>1,2,3,5</sup>. The investigation process of a positive blood culture is however not easy to manage and is difficult to automate. Copan BC+™ allows the blood culture to be transported and transferred in a more manageable and automated format. In order to ensure a proper diagnosis, the sample transferred to the BC+™ must be processed as soon as possible after arriving in the laboratory.

### PRODUCT DESCRIPTION

Copan BC+™ is composed of a plastic vacuum tube and a screw cap with membrane (Fig. 1).

The vacuum in the Copan BC+™ test tubes is calibrated to aspirate a volume of 2.5 ml of a liquid that has the density of water.


Fig.1 BC+™

Copan BC+™ is supplied in a re-sealable bag containing 50 test tubes. 12 re-sealable bags containing 50 test-tubes each are packaged in one cardboard box.

Copan BC+™ is available in the product configurations indicated in the table below:

Catalogue No.	Product Description	Packaging
2U026S	12 x 80 sterile vacuum test tube with round bottom and red screw cap with pierceable membrane	Re-sealable bag containing 50 test-tubes 12 x 50 test tubes/box

### REQUIRED MATERIALS THAT ARE NOT INCLUDED

Vials of blood culture to be tested. Suitable perforators. Copan has tested the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. Qualification tests of the BC+™ were performed by Copan using the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. The instructions provided here were developed for the use of the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. The use of any other perforator must be validated by the user before use

# BC+™

and may require different or additional operating procedures and safety precautions than the ones listed in this leaflet.

Appropriate material for isolation and culturing of microorganisms, such as plates with culture medium or test tubes and incubation systems.

Appropriate material for the preparation of a microscopic examination.

The user must refer to the laboratory manuals for the protocols relating to the microorganism culturing, microscopy and identification techniques of clinical samples.

## STORAGE OF THE PRODUCT

This product is ready to use and does not require additional preparation. It must be kept in its original package at temperatures between 5 and 25°C until the time of use. Do not overheat. Do not incubate or freeze prior to use. Improper storage could have an influence on product performance. Do not use after the expiry date printed on the external package and on the label of the test tube.

## LIMITATIONS

1. Copan BC+™ is designed for the transfer and transport of positive blood cultures for a subsequent use for microbiological culture on solid culture medium or microscopic investigation.
2. Copan BC+™ contains no preservatives, anticoagulants or other additives.
3. The sample volume of blood drawn can vary with sample density, altitude, ambient temperature, barometric pressure, residual life of the product and the transfer technique.
4. For use in automated WASP™ systems the amount of blood culture that is drawn should be within the window on the label of the BC+™ test-tube.
5. Manual filling of the test-tube is not recommended.
6. The blood culture may be a very critical sample from a biosafety point of view. Follow the BSL recommendations given in the internal laboratory procedures for the handling of samples (transfer into a Copan BC+™ test tube and subsequent seeding on a plate). According to the WHO/CDC a minimum safety level of BSL 2 must be ensured. If an infection with organisms belonging to risk group 3 or higher is suspected (such as *Mycobacterium tuberculosis*, *Brucella* species, *Francisella* species, *Y. pestis*, *B. mallei*, *B. pseudomallei*, *S. typhi* and *S. paratyphi*) and based on the origin of the sample or the patient's clinical condition, the handling operations may require biosafety conditions that are higher than level 2<sup>2,6,7</sup>.
7. The conditions, times and volume of the sample collected are important variables for obtaining reliable results. Follow the recommended guidelines for sample collection<sup>1,2,3</sup>.

## WARNINGS AND PRECAUTIONS

1. This product is intended for *in vitro* diagnostic use only.
2. Qualification tests of the BC+™ were performed by Copan using the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. The instructions provided here were developed for the use of the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. **The use of any other perforator must be validated by the user before use and may require different or additional operating procedures and safety precautions than the ones listed in this leaflet.**
3. Take the necessary precautions regarding the biological risk and commonly recognized aseptic techniques. This product must only be used by properly trained and qualified personnel<sup>4</sup>.
4. All samples and materials used for product analysis are to be considered potentially infected and must therefore be handled in such a way as to prevent infection of laboratory personnel. After use, sterilize all biohazardous waste including samples, containers and culture media. Follow the appropriate biosafety recommendations required by current regulations<sup>4</sup>.
5. Carefully read and follow the instructions.
6. This product is for single use only; reuse may cause cross-contamination and/or give inaccurate results.
7. The product must only be handled by trained personnel<sup>4</sup>.
8. Do not repackage or relabel.
9. Not validated for transferring and transporting other fluids besides blood cultures.
10. Not validated for applications other than those for which it is intended.
11. Do not unscrew the cap or pierce the membrane before actually collecting the sample. The capacity of the vacuum in the test tube and therefore its sterility and effectiveness in drawing the blood culture are compromised after opening the cap or piercing the membrane.
12. Unscrewing the cap of the test tube under vacuum can cause the membrane to detach from the cap and the membrane may adhere to the threads of the test tube.
13. Do not use the same perforator for more than one vial of blood culture. This may cause cross-contamination and misdiagnosis.
14. Do not transfer the blood culture sample into the test tube using a syringe with needle (since the test tube is under vacuum, injecting the sample creates a positive pressure that acts on the syringe plunger and on the sample and could lead to spillage and expose the operator to potentially infectious material).

15. Do not use if the device is visibly damaged.
16. Do not use if foreign bodies are present inside the test tube.
17. Before vortexing, make sure that the cap is securely closed.
18. The BC+™ blood culture transfer procedure requires the use of a perforator: take appropriate precautions.

# BC+™

19. Do not ingest any component of the device.
20. The plastic components of the BC+™ device are not guaranteed for cryoresistance. Extremely low temperatures and/or exposure to nitrogen vapours may affect the strength of the plastic and the device could break and/or may not function properly (e.g. sample leakage).

## PRODUCT DETERIORATION

Do not use Copan BC+™ if: 1. the product has visible signs of damage or contamination, 2. the expiry date has passed, 3. the test tube is open, or 4. if there are other signs of deterioration.

## COLLECTION, STORAGE AND TRANSPORT OF THE SAMPLES

The samples for positive blood culture must be transferred into the Copan BC+™ following the instructions and be handled following the guidelines and procedures outlined in the published manuals<sup>5</sup>. The positive blood culture samples transferred into Copan BC+™ test tubes must be tested as soon as possible follow the laboratory SOPs. The transport and handling of the samples within the laboratory, within the building or between buildings must be done in accordance with current standards.

## INSTRUCTIONS FOR USE

### Sample Transfer, Transport and Storage with Copan BC+™

Follow the recommended procedures to process blood samples<sup>1,2,3</sup> and the instructions provided by the manufacturer of the blood culture vials used. Once the blood culture is positive, it must be processed to confirm the presence of microorganisms and to identify the microorganisms in the culture<sup>1,2,3,5</sup>. To make the performing of the confirmation tests easier and to automate the process, the blood culture can be easily transferred into the BC+™ using a perforator (not supplied).

The blood culture may be a very critical sample from a biosafety point of view. Follow the BSL recommendations given in the internal laboratory procedures for the handling of samples (transfer into a Copan BC+™ test tube and subsequent seeding on a plate). According to the WHO/CDC a minimum level of BSL 2 must be ensured. If an infection with organisms belonging to risk group 3 or higher is suspected (such as *Mycobacterium tuberculosis*, *Brucella* species, *Francisella* species, *Y. pestis*, *B. mallei*, *B. pseudomallei*, *S. typhi* and *S. paratyphi*) and based on the origin of the sample or the patient's clinical condition, the handling operations may require biosafety conditions higher than 2<sup>2,6,7</sup>.


Fig. 2: Procedure for the transfer of positive blood culture from the vial to the Copan BC+™ test tube (a, b, c, d, e)

**NOTE:** the instructions provided here were developed for the use of the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. The use of any other perforator must be validated by the user before use and may require different or additional operating procedures and safety precautions than the ones listed in this leaflet.

1. Put it in a hood having an appropriate biosafety level.
2. Unpack the perforator. **NOTE:** qualification tests of the BC+™ were performed by Copan using the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. The instructions provided here were developed for the use of the VACUETTE® QUICKSHIELD complete PLUS 21G perforator manufactured by Greiner Bio One International GmbH. The use of any other perforator must be validated by the user before use and may require different or additional operating procedures and safety precautions than the ones listed in this leaflet.  
**To avoid accidental punctures, be careful not to touch the needles in the perforator.**
3. Wipe the cap of the blood culture vial and the BC+™ with alcohol and allow to dry<sup>5</sup>.
4. Invert the blood culture vial several times<sup>5</sup>.
5. Place the vial in a vertical position on the work surface of the biosafety hood as shown in Figure 2a. Remove the needle cap from the perforator and position it as shown in fig. 2a.  
Pierce the membrane of the blood culture vial cap with the needle and fully insert it into the vial.  
**NOTE: a positive pressure may have formed in the blood culture vial due to the formation of gas during growth. If this is the case, it is recommended that the vial be vented beforehand. Refer to the instructions for use of the vial.**
6. Hold the vial containing the blood culture at an angle in one hand so that the needle is immersed in the sample while being careful not to slide the perforator down (Fig. 2b). **NOTE: at this point the blood culture could begin to leak out of the perforator.**
7. With the other hand take the Copan BC+™ test tube and insert it in the perforator.  
Pierce the membrane and press the test tube against the perforator, pushing it all the way down so that the needle is fully inserted in the BC+™ (fig. 2b). Keep the BC+™ test tube in this position throughout the operation. Do not twist the perforator or the test tube to force the perforation of the cap (the needle must be completely immersed in the blood culture): the blood culture will be drawn into the test tube. **NOTE: do not open or perforate the test tube before starting the transfer procedure.**
8. Once the vacuum in the test tube can no longer aspirate the sample, put the vial and test tube in a horizontal position and slowly remove the BC+™ test tube from the adapter while securely holding the perforator with one hand (Fig. 2c).  
**NOTE: For the correct processing of the BC+™ test tube in WASP™, the volume level of the blood culture that is drawn must be within the window on the label of the test tube.**
9. Place the BC+™ test tube containing the blood culture in a special test tube rack (Fig. 2d).
10. Rotate the vial to its initial position by placing it on the work surface under the hood (with the cap facing upwards) (Figure 2e).
11. Lift the perforator in order to slowly remove the needle from the blood culture vial; cover the needle with the appropriate protective cap and dispose

# BC+™

of the perforator in an appropriate container for sharp objects. **NOTE: lift the perforator slowly in order to avoid or limit the splashing of blood culture.**

12. Label the BC+™ test tube containing the transferred blood culture.
13. Use the blood culture transferred into the BC+™ test tube to seed the plate or to prepare a slide for microbiological tests<sup>1,2,3</sup>. The seeding of the plate and slide preparation can be automated by processing the BC+™ test tube in WASP™.

**NOTE: In the laboratory, wear protective gloves and other suitable protective equipment when opening sterile products and handling clinical samples. Handle the perforator with care to avoid getting accidentally pricked. Dispose of the perforator in the appropriate container<sup>4</sup>.**

## Use in the laboratory

The blood culture samples transferred into the Copan BC+™ test tube can be tested:

- for microbiological culture using appropriate culture media and laboratory techniques;
- by microscopic investigation after slide smear.

Follow the recommendations given for biosafety level 2 according to the WHO/CDC or that required by the current regulations for the handling of samples. If an infection with organisms belonging to risk group 3 is suspected (such as *Mycobacterium tuberculosis*, *Brucella* species, *Francisella* species, *Y. pestis*, *B. mallei*, *B. pseudomallei*, *S. typhi* and *S. paratyphi*) the transfer into a Copan BC+™ test tube and subsequent plate seeding operations must be done at biosafety 3 containment conditions<sup>2,5,7</sup>.

Follow the laboratory's internal standard operating procedures (SOPs) and the guidelines for the treatment of blood culture samples.

## Using BC+™ with automated systems

Positive blood cultures properly transferred into Copan BC+™ can be processed with automated systems such as Copan WASP™.

For the correct processing of the BC+™ test tube in WASP™, the volume level of the blood culture that is drawn must be within the window on the label of the test tube.

The Copan BC+™ can be used on other automated seeding systems after validation by the end user.

Consult the manufacturer's instructions for the processing of Copan BC+™.

## **RESULTS**

Copan has validated that BC+™ does not interfere with the positive blood culture sample after contact of up to 8 hours at 20-24°C using Bact/ALERT®, Bioréieux. The results of blood culture confirmation tests depend greatly on the quality of the blood culture and that the laboratory analysis is done in a timely manner after the blood culture in the vial has been determined to be positive by a continuous monitoring instrument. False positive results for blood cultures in vials have been reported for automated incubation and seroconversion detection systems. These results depend on many instrument factors as well as factors relating to the sampling quality, such as excessive sample volume<sup>1</sup>, or the presence of a number high of leukocytes<sup>2</sup>. If the plate seeding and/or microscopic examination is negative, perform the checks suggested in the guidelines, for example: verification of the leukocytes, verification of the instrumental growth curves of the vial, use of alternative stains or culture media.

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**BC+™****BIBLIOGRAPHY**

1. Mandell G., Bennett J.E., Dolin R.. Principles and Practice of Infectious Diseases . Sixth ed. 2005. Elsevier Churchill Livingstone.
2. UK Standards for Microbiology Investigations. Investigation of Blood Cultures (for Organisms other than *Mycobacterium* species). B 87, 2014.
3. L.S. Garcia. Clinical Microbiology Procedures Handbook, third ed. ASM, Washington, DC. 2010
4. Laboratory biosafety manual. World Health Organization. Third ed, 2004.
5. Clinical and Laboratory Standards Institute. Principles and Procedures for Blood Cultures ; Approved guideline M47-A. 2007.
6. Mangalgi S, Sajjan A. Comparison of three blood culture techniques in the diagnosis of human Brucellosis. J Lab Physicians 2014;6:14-7
7. Fact sheet on typhoid and paratyphoid fever A Central Committee on Infectious Diseases Nov 2003

**TABELLA DEI SIMBOLI / TABLE OF SYMBOLS / TABLEAU DES PICTOGRAMMES / TABELLE DER VERWENDETEN SYMBOLE / TABLA DE SIMBOLOS / TABELA DE SIMBOLOS**

Simbolo / Symbol / Pictogramme / Symbol	Significato / Meaning / Signification / Bedeutung / Significado
	Fabbricante / Manufacturer / Fabricant / Hersteller / Fabricante
	Dispositivo diagnostico in vitro / In vitro diagnostic device / Dispositif de diagnostic in vitro / In-vitro-Diagnostikum / Dispositivo de diagnóstico in vitro
	Sterilizzato usando radiazioni ionizzanti / Sterilized using ionizing radiation / Stérilisé à l'aide de rayonnements ionisants / Sterilisation durch Bestrahlung / Esterilizado por radiaciones ionizantes / Esterilizado por radiação ionizante
	Non riutilizzare / Do not reuse / Ne pas réutiliser / Nicht zur Wiederverwendung / No reutilizar / Não reutilizar
	Limiti di temperatura / Temperature limits / Limites de température / Temperaturbegrenzung / Límites de temperatura / Limites de temperatura
	Utilizzare entro / Use before / À utiliser avant / Verwendbar bis / Usar antes de / Usar antes de
	Consultare le istruzioni per l'uso / Consult the instructions for use / Consulter les instructions d'utilisation / Gebrauchsanweisung beachten / Consultar las instrucciones de uso / Consultar as instruções de utilização
	Codice del lotto (partita) / Batch code (lot) / Code du lot (lot) / Chargenbezeichnung / Código del lote (Lote) / Código do lote (remessa)
	Contenuto sufficiente per <n> test / Contents sufficient for <n> tests / Contenu suffisant pour <n> tests / Inhalt ausreichend für <n> Prüfungen / Contenido suficiente para <n> pruebas / Conteúdo suficiente para <n> testes
	Non utilizzare in caso di confezionamento danneggiato / Do not use if the packaging is damaged / Ne pas utiliser si l'emballage est endommagé / Nicht verwenden, wenn Verpackung beschädigt / No utilizar si el envase está dañado / Não utilizar se a embalagem apresentar danos

Copan Italia SpA  
Via Perotti, 10  
Brescia, Italy

Copan Italia SpA  
Via Perotti 10  
25125 Brescia Italy  
Tel: +39 030 2687211  
Fax: +39 030 2687250  
E-mail: info@copangroup.com  
Website: www.copangroup.com

North American Distributor:  
Copan Diagnostics Inc.  
26055 Jefferson Avenue  
Murrieta, CA 92562 USA  
Tel: 951-696-6957  
Fax: 951-600-1832

E-mail: customerservice@copanusa.net  
Website: www.copanusa.com

Innovating Together™