

# MAST® Culture Media and Supplements

## Technical Information Sheet

Product Code DM 160



## Mannitol Salt Agar

A selective medium for the isolation of pathogenic staphylococci.

### 1. Description

Following the observation by Koch<sup>1</sup> that usually only staphylococci are able to grow on 7.5% salt agar, Chapman<sup>2</sup> described the use of Mannitol Salt Agar for the differential isolation of pathogenic staphylococci. It is worth noting that Gunn *et al*<sup>3</sup> improved the medium by the addition of 2% egg yolk which identified 94.7% of coagulase producing strains. Recent work has shown that some strains of staphylococci may be partially inhibited by the high level of sodium chloride in Chapman's medium.<sup>4,5</sup> To overcome this, the medium has been

modified by reducing the sodium chloride concentration from 7.5 to 3.0%. Inhibition of Gram negative bacilli has been maintained by the inclusion of Lithium chloride. Pyruvate and glycine are included as growth promoters to enhance the recovery of *S.aureus* from samples containing small numbers of the organism.<sup>6,7</sup> The medium is suitable for use with Oxacillin Selectatab™ (MS29) for the detection of Methicillin Resistant *Staphylococcus aureus* (MRSA) using the method described by Lally *et al*<sup>8</sup>.

### 2. Technical Formula\*

Formula	grams per litre
Peptone	8.0
Yeast extract	2.0
Lactalbumin	3.0
Sodium chloride	30.0
Phenol red	0.0225
Lithium chloride	7.0
Glycine	1.0
Sodium pyruvate	3.0
Agar B	12.0
<b>pH approx.7.4</b>	

### 3. Directions

1. Suspend by swirling 76g of powder in 1 litre or the contents of the sachet in the stated volume of distilled or deionised water.
2. Autoclave at 121°C (15p.s.i.) for 15 minutes.
3. Mix well before pouring.

#### 4. In Use

Inoculate plates by spreading the sample over the surface of the medium and incubate at 37°C for up to 48 hours. Pathogenic staphylococci grow well and ferment mannitol, indicated by yellow zones around the colonies, whereas non-pathogenic staphylococci produce smaller white colonies indicating non-fermentation of mannitol.

#### 5. References

1. Koch FE. *Zentbl Bakt ParasitKde (Abt.1)* 1942; **149**: 122-124.
2. Chapman GH. *J Bact.* 1945; **50**: 202-203.
3. Gunn BA, Dunkelberg WE, Creitz JR. *Am J Clin Path.* 1972; **57**: 236-240.
4. Van Enk RA, Thompson DK. *J Clin Microbiol.* 1992; **30**: 504-505.
5. Faiers M, George R, Jolly J, Wheat P. *Multipoint Methods in the Clinical Laboratory - A Handbook.* p41 PHLS, Lon 6.
6. Baird Parker AC. *J Appl Bact.* 1962; **25**: 12.
7. Giolloti C, Cantoni C. *J Appl Bact.* 1966; **29**: 395-398.
8. Lally RT, Ederer MN, Woolfrey BF. *J Clin Microbiol.* 1985; **22**: 501-504



**Mast Group Ltd.**  
Mast House, Derby Road, Bootle  
Liverpool, Merseyside, L20 1EA  
United Kingdom  
Tel: + 44 (0) 151 472 1444  
Fax: + 44 (0) 151 944 1332  
email: sales@mastgrp.com  
Web: www.mastgrp.com

**Mast Diagnostica GmbH**  
Feldstrasse 20  
DE-23858 Reinfeld  
Germany  
Tel: + 49 (0) 4533 2007 0  
Fax: + 49 (0) 4533 2007 68  
email: mast@mast-diagnostica.de  
Web: www.mastgrp.com

**Mast Diagnostic**  
12 rue Jean-Jacques Mention  
CS91106, 80011 Amiens, CEDEX 1  
France  
Tél: + 33 (0) 3 22 80 80 67  
Fax: + 33 (0) 3 22 80 99 22  
email: info@mast-diagnostic.fr  
Web: www.mastgrp.com

