# CHROMagar™ Enterobacteria

**Instructions For Use**  
Available in several languages

**NT-EXT-099**

Version 3.0


Click below for:

**ENGLISH**

**FRANCAIS**

**ESPAÑOL**

**DEUTSCH**

## MEDIUM PURPOSE

Chromogenic medium for detection and enumeration of *Enterobacteriaceae*.

«The *Enterobacteriaceae* and coliform bacteria within this family represent two of the most common groups of indicator organism used by the food industry. In some countries, depending on regulatory requirements, the food industry has moved towards testing for *Enterobacteriaceae*». ILSI Europe (International Life Sciences Institute)

## COMPOSITION

The product is composed of one powder Base.

Product	=	Pack
Total		42.3 g/L
Composition		Agar 11.0 Peptones and yeast extract 22.0 Mineral salts 4.9 Chromogenic and selective mix 0.7 Growth factors 3.7
Aspect		Powder Form
STORAGE		15-30 °C
FINAL MEDIA pH		7.4 +/- 0.2

## PREPARATION (Calculation for 1 L)

### Step 1

Preparation of the mix

- Disperse slowly 42.3 g of powder base in 1 L of purified water.
- Stir until agar is well thickened.
- Heat and bring to boil (100 °C) while swirling or stirring regularly.

**Advice 1:** For the 100 °C heating step, mixture may also be brought to a boil in a microwave oven: after initial boiling, remove from oven, stir gently, then return to oven for short repeated bursts of heating until complete fusion of the agar grains has taken place (large bubbles replacing foam).

DO NOT HEAT TO MORE THAN 100 °C. DO NOT AUTOCLAVE AT 121 °C.

**Warning:** If using an autoclave, do so without pressure.

### Step 2

Cooling

- Cool in a water bath to 45-50 °C.
- Swirl or stir gently to homogenize.

### Step 3

Pouring in plates or bottles

- Pour medium into Petri dishes or 200 mL bottles.
- Let solidify.
- If pouring in plates, let them dry.

### Storage

- Store in the dark before use.
- Prepared media plates and bottles can be kept for one day at room temperature.
- Plates can be stored for up to 1 month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.
- Bottles can be stored for up to 6 months under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

## INOCULATION

### IF USING POURING TECHNIQUE PROCEDURE

• Prepare 90 mm Ø sterile Petri dishes and add 1 mL of inoculum in each. Then pour 10 mL of melted medium. Mix and let solidify. Add an overlay of the same medium and let solidify. Invert the plate and incubate at 37 °C for 24 hours. If searching for psychotropic *Enterobacteriaceae*, incubate at 30 °C.

### IF USING SURFACE TECHNIQUE PROCEDURE:

• Pour into sterile Petri dishes and allow to solidify and dry. If the agar plate has been refrigerated, allow to warm to room temperature before inoculation. Streak the sample into plate and incubate in aerobic conditions at 37 °C for 24 h. If searching for psychotropic *Enterobacteriaceae*, incubate at 30 °C.


## Typical Samples

e.g. food products, animal feed, environment samples  
\*\*\*  
Pouring, streaking or spreading techniques

## INTERPRETATION

Microorganism	Typical colony appearance
<i>E. coli</i>	→ blue with or without blue halo
<i>Proteus</i>	→ red with swarming
Other <i>Enterobacteriaceae</i>	→ pink to red
Gram (+) bacteria	→ inhibited
Other Gram (-) bacteria	→ mostly inhibited

### Typical colony appearance

The pictures shown are not contractual.

## PERFORMANCE & LIMITATIONS

- Definite identification may require additional testing.
- Some strains of *Pseudomonas* may appear as pink colonies. However, they can be differentiated by an Oxidase test.

## QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the following ATCC strains:

Microorganism	Typical colony appearance
<i>E. coli</i> ATCC® 25922	→ blue
<i>K. pneumoniae</i> ATCC® 13883	→ pink
<i>P. aeruginosa</i> ATCC® 9027	→ inhibited
<i>E. faecalis</i> ATCC® 29212	→ inhibited

## WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.

## REFERENCES

Σ Pack Size

5000 mL

=

250 Tests  
of 20 mL

Ordering References

EB042

=

EB042

Weight: 211.5 g

- For Laboratory use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.
- Interpretation of the test results should be made taking into consideration colonial and microscopic morphology and if necessary, the results of any other tests performed.
- Laboratory, chemical or biohazardous wastes must be handled and discarded in accordance with all local and national regulations.
- For hazard and precaution recommendations related to some chemical components in this medium, please refer to the pictogram(s) mentioned on the labels. The Safety Data Sheet (SDS) is available on [www.chromagar.com](http://www.chromagar.com)

## DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

## LITERATURE REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <https://www.chromagar.com/en/product/chromagar-enterobacteria/>

## IFU/LABEL INDEX

- Catalogue reference
- Consult instructions for use
- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity
- Protect from light
- Manufacturer

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach  
ATCC® is a registered trademark of the American Type Culture Collection  
NT-EXT-099 V3.0 / EN 07-Nov-23

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