CHROMagar™ **StrepA**

Instructions For Use
Available in several languages

NT-EXT-103

Version **1.2**


Click below for:

ENGLISH


FRANCAIS

ESPAÑOL

DEUTSCH

REFERENCES

Pack Size	Ordering References	Base (B)	Supplement (S1)	Supplement (S2)
5000 mL = 	SP372	SP372(B) Weight: 271.5 g	SP372(S1) Volume: 10 mL	SP372(S2) Weight: 0.5 g
25 L = 	SP373-25	SP373-25(B) Weight: 1357.5 g	SP373-25(S1) Volume: 50 mL	SP373-25(S2) Weight: 2.5 g

MEDIUM PURPOSE

Chromogenic medium for the screening of Group A Streptococci in throat samples.

Group A *Streptococcus* (GAS) *S. pyogenes* cause many different infections (strep throat or pharyngitis, scarlet fever, impetigo) that can be transmitted in a variety of ways. Throat culture is essential to confirm bacterial pharyngitis prior to the application of antibiotic therapy, preventing complications and reducing the risk of transmission, particularly in children and elderly patient.

COMPOSITION

The product is composed of a powder base (B) and 2 supplements (S1 + S2).

Product	Base (B)	Supplement S1	Supplement S2
Total g/L	54.3 g/L	2 mL/L	0.1 g/L
Composition g/L	Agar 15.0 Peptones and yeast extract 27.0 Salts 9.0 Growth factors 3.0 Chromogenic and selective mix 2.2	Growth factors 2.0	Selective mix 0.1
Aspect	Powder Form	Liquid Form	Powder Form
STORAGE	15/30 °C	15/30 °C	2/8 °C
FINAL MEDIA pH	7.2 +/- 0.2		

Need some Technical Documents?

Available for download on www.CHROMagar.com

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

PREPARATION (Calculation for 1 L)

Step 1

Base + S1

- Disperse slowly 54.3 g of powder base in 1 L of purified water.
- Add 2 mL of supplement S1 into slurry.
- Stir until agar is well thickened.
- Autoclave at 121 °C during 15 min.
- Cool at 45/50 °C keeping on stirring.

Step 2

S2

- In a transparent vessel, add 100 mg of supplement S2 in 1 mL of purified water.
- Filter to sterilize at 0.45 µm.

Step 3

Base + S1 + S2

- Aseptically add 1 mL of S2 preparation into base + S1 slurry, cooled at 45/50 °C while mixing.
- Swirl or stir gently to homogenize.

Step 4

Pouring

- Pour into sterile Petri dishes.
 - Let it solidify and dry.
- Once dried, the appearance of the plate is **translucent grey**.

Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to 1 month under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.

SPECIMEN COLLECTION AND HANDLING

CHROMagar™ StrepA can be used with the following specimens: throat swabs.

Use of transport devices approved for collection of such specimens is recommended.

MATERIAL REQUIRED BUT NOT PROVIDED

Standard microbiological laboratory material for culture media preparation, control, streaking, incubation and waste disposal.

INOCULATION


Related samples can be processed by direct streaking on the plate.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate at 35-37 °C for 18-24 hours in a CO₂ atmosphere.

INTERPRETATION

Microorganism	Typical colony appearance
<i>Streptococcus pyogenes</i> (group A)	→ orange to red
Other <i>Streptococcus</i>	→ steel blue or colourless
Other Gram (+) bacteria	→ inhibited
Yeasts	→ inhibited
Gram (-) bacteria	→ inhibited

Typical colony appearance

The pictures shown are not contractual.

PERFORMANCE

In the following study, 159 throat samples were tested, being positive 120 after 20 h incubation at 35 °C in a CO₂ atmosphere.

	CHROMagar™ StrepA	Reference Method (blood agar)
Sensibility	96,7 % *	90,8 %
Specificity	100 % *	--

* Data obtained from the study «Validation and implementation of Colorex™ CHROMagar™ Strep A agar on WASP™/WASPLab™ for screening for Streptococcus pyogenes using the ESwab™» Mark Gaskin et Al. ASM 2019

LIMITATIONS AND COMPLEMENTARY TESTS

- Rare strains of Group A *Streptococcus* may require an additional 24 h of incubation for a satisfactory colony size.
- Final identification may require additional testing such as biochemical or immunological tests like PYR and agglutination tests.

QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms.

Good preparation of the medium can be tested, isolating the following ATCC strains:

Microorganism	Typical colony appearance
<i>S. pyogenes</i> ATCC® 19615	→ orange to red
<i>E. faecalis</i> ATCC® 29212	→ steel blue
<i>E. coli</i> ATCC® 25922	→ inhibited
<i>S. aureus</i> ATCC® 25923	→ inhibited

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- This laboratory product should be used only by trained personnel (healthcare professional, etc). Wear appropriate protective clothing, gloves and eye/face protection and handle appropriately with procedures and good laboratory practices.
- Use of the medium may be difficult for people who have problems recognising colours.
- For a good microbial detection, collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.
- Culture media should not be used as manufacturing material or components.
- Do not ingest or inhale the product.
- Do not use the product after the expiry date.
- Do not use the product if it show any evidence of contamination or any sign of deterioration.
- Do not use the product if the packaging is damaged.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- Reading and interpretation should be performed using isolated colonies.
- Some precipitates may be observed in the agar but these do not affect the performance of the product.
- Interpretation of the test results should be made taking into consideration colonial and microscopic morphology and if necessary, the results of any other tests performed.
- Laboratory, chemical or biohazardous wastes must be handled and discarded in accordance with all local and national regulations.
- For hazard and precaution recommendations related to some chemical components in this medium, please refer to the pictogram(s) mentioned on the labels. The Safety Data Sheet (SDS) is available on www.chromagar.com

CHROMagar™ StrepA

DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

LITERATURE REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

IFU/LABEL INDEX


-  Catalogue reference
-  Consult instructions for use
-  Quantity of powder sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

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