LBM[®] Thiol Broth

Package insert and How-to-use guide

CE IVD

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THIOL Broth - Product Package Leaflet and Use Instructions

INTENDED USE

THIOL broth is an enrichment medium which facilitates the growth of facultative and obligate aerobic and anaerobic microorganisms. COPAN Thiol Broth is used for the enrichment of anaerobic bacteria isolated from cultures of clinical samples. THIOL Broth can be used with automatic plate culture systems (Copan WASP™ or Copan WASPLab™).

SUMMARY

THIOL is an enrichment medium consisting of dextrose, peptone and L-cystine solution and yeast extract, which provide the essential growth factors for bacterial reproduction.

The essential ions are provided by sodium chloride. Sodium thioglycollate and L-cystine act as reducing agents, which contribute to maintaining a low redox potential, supporting anaerobiosis.

Thiol broth also contains growth supplements such as Hemin and Vitamin K1 which facilitate the reproduction of some particularly demanding anaerobic bacterial strains.

REAGENTS

Components of THIOL Broth:

| Names of components |
|-----------------------------|
| L-cystine |
| Sodium chloride |
| Glucose |
| Yeast extract |
| Pancreatic digest of casein |
| Sodium thioglycollate |
| Hemin |
| Vitamin K1 |
| Distilled water |

STORAGE

Thiol Broth is a ready-made product and no additional preparation operations are required. The batch packed in its original package can be kept at 5-25°C until it is used or until the expiry date indicated on the tube and on the box. Do not overheat. Do not incubate or freeze before use. Inappropriate storage may cause a loss of product effectiveness.

PRODUCT DETERIORATION

Do not use THIOL Broth if: (1) the product shows visible signs of damage or contamination; (2) there are indications of leakages from the tube; (3) it is past its expiry date; (4) there are other visible deterioration signs.

SUPPLIED MATERIALS

| Code | Product Description | Package Dimensions | Suitable for automation |
|--------|--|---|-------------------------|
| 4U002N | 12x80 mm PET test tubes with PE perforated screw cap. One test tube contains 4 ml of THIOL Broth. | 50 units per internal package. 50x6 units per box. | YES |

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for culture and isolation of bacteria; refer to laboratory reference manuals for the protocols recommended for the culture and identification techniques.

INSTRUCTIONS FOR USE

For specific guidelines on sampling procedures, consult the reference manuals published. The Copan eSwab system can be used as a system for sampling, transport and preservation of the clinical sample. Refer to Copan for information on the most appropriate code, and refer to the product package leaflet for use instructions.

USER INSTRUCTIONS

MANUAL OPERATIONS

- Unscrew the THIOL broth tube cap.
- Inoculate the sample into the open tube.
For liquid samples, use a 10 µl loop or micropipette to transfer the sample to the tube.

For samples in eSwab:

- Homogenize the tube for 10 sec with vortex at 2000/2500 RPM, unscrew the eSwab tube cap, immerse the swab in the Thiol Broth medium and

shake vigorously for 5/10 seconds to release the sample into the enrichment medium, dispose of the swab and re-close the Thiol Broth tube.

- b) Using a micropipette, transfer 30 µl of liquid Amies medium directly into the THIOIOL Broth. Fit the caps on the eSwab and Thiol Broth tubes.
3. Homogenize for 5-10 seconds at 2000/2500 rpm by a vortex.
4. Incubate the inoculated THIOIOL Broth tubes at 35°C ± 2°C.
5. Examine the THIOIOL Broth tube after 24-48 hours. The tube must be checked daily. Absence of turbidity in the tube does not indicate absence of the microbe growth.
6. Proceed with surface culture on appropriate medium in a plate. Surface culture on plate can be performed through surface culture of a sample portion taken by a micropipette or a loop, followed by distribution on plate. The minimum recommended quantity for surface culture in a plate after incubation of the Thiol Broth tube is 1-10 ul. Incubate the plates for 24-96 hours at 35°C ± 2°C in appropriate anaerobiosis conditions, or according to standard laboratory procedure.
7. At the end of incubation, proceed with plate reading.
8. You are advised to plate the inoculated Thiol Broth sample daily.
9. You are advised to check growth in the plate daily.

OPERATIONS WITH AUTOMATIC SYSTEM (WASP™)

Thiol Broth can be used on automatic culture systems like COPAN WASP™. For correct use, see the WASP™/WASPLab™ instruction manual.

LIMITATIONS

1. Some difficult samples with specific requirements may not grow in this medium.
2. The conditions and the volume of the sample collected for culture are significant variables for obtaining reliable culture results.
3. An appropriate sample collection from patients is an extremely critical factor for isolation and identification of the bacterial strains of interest.
4. Absence of turbidity in the Thiol Broth tube is not indicative of a negative result.
5. Check the tube daily until the end of incubation.
6. You are advised to plate both after incubation of 24 hours and after incubation of 48 hours in an appropriate growth medium.
7. Maintain the Thiol Broth tube in incubation for a maximum of 96 hours at 35°C ± 2°C.

NOTICES

1. For diagnostic use in vitro.
2. This product is designed for single use; reuse may lead to unreliable results.
3. Not appropriate for any application other than its intended use.
4. Follow the recommended guidelines for collection of the samples.
5. Work in a laminar flow cabinet after inoculation of the product, wear gloves and other protections in line with universally recognized precautions for handling of clinical samples.
6. The use of this product with any diagnostic test or with any diagnostic tool shall be assessed and if necessary validated by the user before use.
7. Do not use the product if it is visibly damaged.
8. Do not ingest.
9. Do not inhale.
10. Do not inject.
11. The manufacturer shall not be considered liable for any inappropriate or non-qualified use of the product.
12. All clinical samples are considered potentially infective and shall be handled with appropriate precautions.
13. The THIOIOL Broth shall only be used for in vitro diagnostics. Its use for therapeutic or prevention purposes is not allowed under any circumstances.
14. Take appropriate precautions against biological hazards and apply aseptic techniques. The product may only be used by properly trained and qualified operators.
15. The instructions must be read and observed carefully.
16. Do not keep the tube in anaerobiosis conditions during incubation.
17. Do not place in an autoclave.
18. Do not filter.

WASTE DISPOSAL

Unused reagents can be considered as non-hazardous waste and disposed of accordingly. For disposal, refer to the product MSDS.

Used reagents and any other contaminated waste material shall be disposed of according to the procedures established for infected or potentially infected products. The laboratory shall be responsible for managing waste materials by treatment and disposal, as established in all applicable regulations.

QUALITY CONTROL PROCEDURE

PERFORMANCE TEST

- > From a fresh culture, prepare a 0.5 McFarland suspension in PBS of every ATCC organism.
- > Prepare a serial 10 -4 dilution of 0.5 McFarland suspension.
- > Inoculate with 400 µl of the dilution in PBS directly in THIOIOL broth.
- > Put the caps back on the tubes.
- > Homogenize for 10 seconds at 2000/2500 RPM with vortex.
- > For the control at zero time, seed 100 µl of the bacterial suspension on nutrient medium (e.g. blood agar plate) and incubate in appropriate growth conditions.
- > Incubate THIOIOL broth tube at 35°C ± 2°C for 24-48 hours.
- > After incubation, seed 100 µl THIOIOL broth on appropriate nutrient medium (e.g. blood agar plate).
- > Incubate the plates at 35°C ± 2°C for 24-48 hours in appropriate growth conditions.

LIMITS OF ACCEPTABILITY: Time zero: CFU counting 30-300

CFU/100 µl

Time 24-48 h: growth

RESULTS

| STRAIN | ZERO DEAD TIME COUNTING: CFU/PLATE | CFU COUNTING AFTER 18-24 HOURS AT 35°C±2°C; CFU/PLATE |
|--|------------------------------------|---|
| <i>Peptostreptococcus anaerobius</i> ATCC 27337* | 216 | GROWTH |
| <i>Bacteroides Vulgatus</i> ATCC 8482* | 198 | GROWTH |
| <i>Clostridium Perfringens</i> ATCC 13124* | 134 | GROWTH |
| <i>Prevotella Melaninogenica</i> ATCC 25845 | 121 | GROWTH |
| <i>Propionibacterium acnes</i> ATCC 6919 | 166 | GROWTH |
| <i>Bacteroides fragilis</i> ATCC 25285 | 200 | GROWTH |
| <i>Staphylococcus aureus</i> ATCC 25923 | 210 | GROWTH |

The results obtained and presented are representative of the values obtained during general evaluation of the product.
The results obtained and presented are representative of tests performed at the COPAN laboratory with ATCC strains.
The results obtained and presented were NOT obtained using clinical samples.

* Bacterial performance tests conducted using strains established by the CLSI-M22-A3 reference standard, specific for quality control of enrichment mediums for the growth of anaerobic bacteria.

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
BIBLIOGRAPHY

1. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. *Manual of Clinical Microbiology*. 7th ed. ASM, Washington, D.C..
2. National Committee for Clinical Laboratory Standards (NCCLS). 2004. *Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3*
3. Miller, J. M. 1999. *A Guide to Specimen Management in Clinical Microbiology*, 2nd ed. ASM, Washington, DC.
4. Isenberg, H. D., 2004. *Clinical Microbiology Procedures Handbook*, 2nd ed. ASM, Washington, DC.
5. *Atlas, Handbook of Microbiological Media* 4th edition 2010 CRC Press, Boca Raton, F

Tabella dei Simboli / Index of Symbols / Index of Symbols / Index of Symbols / Table des Symboles / Tabela de símbolos / Tabel med symboler / Tabel med symboler / Tabell över symboler / Sembollerin Tablosu

| Simbolo / Symbol / Symbol / Symbole / Symbole / Símbolos / Symbol / Symbol / Symbol / Symbol / Sembol | Significato / Meaning / Signification / Signification / Signification / Significado / Betyder / Betydning / Betyder / Anlam |
|---|--|
|  | Fabbricante / Manufacturer / Fabricante / Hersteller / Fabricant / Fabricante / Producent / Produzenten / Tillverkare / Üretici |
|  | Marchio CE / CE marking / Marca CE / CE-Kennzeichnung / Marquage CE / CE-mærke / CE işareti |
|  | Dispositivo diagnóstico in vitro/ In vitro diagnostic device / Dispositif de diagnostic in vitro/ Diagnosegerät in vitro/ Dispositivo de diagnóstico in vitro/ Medicinsk anordning för in vitro-diagnostik/ In vitro diagnostisk medicinsk anordning/ In vitro tani cihazı |
|  | Non riutilizzare / Do not reuse / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a usar / Brug det ikke igen / Må ikke gjenbrukes / Använd inte den igen / Tekrar kullannamayın |
|  | Numero di catalogo / Catalogue number / Número de catálogo / Bestellnummer / Référence du catalogue / Referência do catálogo / Katalognummer / Katalognr / Katalognummer / Katalog numarası |
|  | Limiti di temperatura / Temperature limitation / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Temperaturgrænser / Temperaturgrænser / Temperaturgrænser / Sıcaklık limitleri |
|  | Utilizzare entro / Use by / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Anvendes før / Må brukes innen / Ska användas innan / Son kullanna tarihi |
|  | Consultare le istruzioni per l'uso / Consult Instructions for Use / Consulte las instrucciones de uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultar as instruções de utilização / Se vejledningen til brug / Se instruksjoner for bruk / Les bruksanvisningen / Kullanna talimatlarına bakınız |
|  | Strappare per aprire / Peel / Desprender / Abziehen / Décoller / Destacável / Riv opp for å åpne / Skall for åpen / Skal för öppen / Açmak için yırtın |
|  | Codice del lotto (partita) / Batch code (Lot) / 7 Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Seriennummer (parti) / Batch-nummer (parti) / Dra för att öppna / Lot (seri) kodu |
|  | Contenuto sufficiente per <n> test / Contains sufficient for <n> tests / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contémo suficiente para <n> testes / Indhold tilstrækkeligt til <n> prøver / Innhold tilstrekkelig for <n> test / Innehåller tillräckligt för <n> tester / <n> test için yeterli içerik |
|  | Non utilizzare in caso di confezionamento danneggiato/ Do not use if package is damaged/ Ne pas utiliser si l'emballage est abimé/ Bei Beschädigung der Verpackung nicht verwenden/ No utilizar en caso de paquete dañado/ Não utilizar em caso de embalagem danificada/ Ska inte användas vid skadade förpackningar/ Skal ikke brukes i tilfelle av skadet emballasje/ Må ikke brukes hvis pakningen er skadet/ Ambalajı hasar görmüşse kullannamayın |

| | | | |
|---|--|--|--|
|  <p>Copan Italia S.p.A. Via F. Perotti, 10 25125 - Brescia, Italy</p> | <p>Copan Italia S.p.A. Via F. Perotti, 10 25125 Brescia Italy Tel +39 030 2687211 Fax +39 030 2687250 Email: info@copangroup.com Website: www.copangroup.com</p> | <p>North American Distributor: Copan Diagnostics Inc. 26055 Jefferson Avenue Murrieta, CA 92562 USA Tel: 951-696-6957 Fax: 951-600-1832 E-mail: customerservice@copanusa.net Website: www.copanusa.com</p> |  <p>COPAN Innovating Together™</p> |
|---|--|--|--|