LBM™ TSB Kit

Product Insert & How to Use Guide

TSB KIT™ enrichment broth - Product Insert & How to Use Guide

INTENDED USE

The **COPAN TSB KIT™ enrichment broth** product is a ready to use enrichment medium that supports a wide variety of aerobic bacteria. The COPAN TSB KIT™ product can be used to inoculate clinical specimens. The broth can be used with a Flocked Swab to inoculate the primary sample and after appropriate incubation, plated on nutritive or selective agar.

SUMMARY AND PRINCIPLES

Copan TSB KIT™ is an enrichment medium that supports a wide variety of aerobic bacteria. The TSB is a general nutritive medium that can be used to support the growth of aerobic and facultative anaerobic bacteria. The pancreatic digest of casein or papaic digest of soybean meal flour are amino acids and other complex nitrogenous substances source that support the growth of microorganism; the Glucose (dextrose) is considered an energetic source and the sodium chlorine maintains the osmotic balance. The dibasic potassium phosphate is used as a buffer that controls the pH.

PRODUCT DESCRIPTION

The COPAN TSB KIT™ enrichment broth product is a ready to use tube filled with 2 ml of nutritive medium combined with 2 sterile regular Flocked Swabs. The tube is provided with screw cap. The COPAN TSB KIT™ enrichment broth tube can be used with automation streaking (WASP™ or WASPLAB™) or manual streaking to inoculate the original sample into the broth and after appropriate incubation it is possible to streak with 4 Quadrant protocol in automation or by manual methods as per internal defined and validated protocols.

The product is available in the following configuration:

Catalog No.	Product Descriptions	Pack Size
468CE02	TSB 12x80 mm polypropylene tube screw-cap with conical shape, filled with 2 ml of TSB medium plus 2 sterile regular Flocked Swabs in paper peel pouch	50 units per shelf pack 6x50 units per box

REAGENTS

TSB KIT™ enrichment broth product: components

Pancreatic digest of casein
Papaic digest of soybean meal
Sodium Chloride
Dibasic potassium phosphate
Glucose

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for the cultivation and isolation of bacteria. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques.

STORAGE

This product is ready to use and no further preparation is necessary. The unopened bulk pack or individual kit pack can be stored at 5 - 25°C until used or until the expiration date. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box.

LIMITATIONS

1. In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens.
2. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.
3. Performance testing with Copan TSB KIT™ product was conducted using laboratory ATCC strains spiked into the TSB tube and not using human specimens.
4. Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals. Specimens should be collected as soon as possible after the clinical onset of disease. Highest bacterial titres are present during the acute illness.

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use and in no way intended for a curative or prophylactic purposes.
2. Not suitable for any other application than intended use.
3.  Copan TSB is for single use only; reuse may cause a risk of infection and/or inaccurate results.
4. Do not re-pack.
5. The use of this product in association with any diagnostic assay or with any diagnostic instrumentation should be validated by the user before using.
6. Do not use if the product is visibly damaged.
7. Do not ingest the medium.
8. Directions for use must be followed carefully.
9. Product to be handled by trained personnel only.
10. It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes must be disposed of according to laboratory regulations for infectious waste.
11. Results obtained depend on adequate specimen collection and timely transport and processing in the laboratory.

INSTRUCTIONS FOR USE

1. Open the paper film pouch containing the tube of TSB.
2. Open the paper film pouch containing the regular flocked swabs.
3. If you need to take two samples from the same patient, take one of the swabs to collect the first specimen (i.e. nose, throat or perineum sample).
WARNING: Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shape may have breaking point and diameter change points that could result less resistant to pressure or bending. Do not use excessive force, pressure or bending when collecting swab sample in order to avoid accidental breakage. Unscrew the TSB tube.
4. Insert the swab into the TSB tube all the way to reach the bottom of the tube. Dip and gently stir the swab for 5 seconds.
5. Lift up the swab from the liquid medium and swirl the swab against the tube walls 5 times to allow release of the sample from the flocked fibre. Remove the swab from the tube and recap.

6. Discard the first used swab into Biohazard Waste.
7. Take the last swab to collect the second specimen (i.e. nose, throat or perineum sample) and then break the swab at the moulded breaking point to deposit the swab into the tube according to the following procedure (Fig.1):
 - Hold the tube in one hand pointing the opening of the tube away from your face
 - With the other hand grasp the swab shaft at the very end with the thumb and first finger
 - Lean the part of the shaft with the breaking point against the rim of the tube
 - Bend the swab shaft at a 180 degrees angle to break it off at the colored ink breakpoint mark. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
 - Discard the broken handle part of the swab shaft into an approved medical waste disposal container.


Fig.1: Swab breakage

8. Recap the TSB tube
9. Incubate the TSB tube at 35°C±2°C for minimum 18-24 hours or proceed as per internal defined and validated protocol.
10. After incubation proceed as per internal defined and validated protocol.

WASTE DISPOSAL

Unused reagents may be considered as non hazardous waste and disposed of accordingly. Dispose of used reagents as well as any other contaminated disposable materials following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to handle waste and effluents produced according to their nature and degree of hazardousness and to treat and dispose of them (or have them treated and disposed of) in accordance with any applicable regulations.

QUALITY CONTROL

Procedure for testing:

- Starting from a fresh culture prepare 0.5 Mc Farland suspension of each organism in PBS.
- From the 0.5 McF suspension prepare appropriate dilution to contain from 300 to 3000 CFU/100 ul
- Using a micropipette inoculate the TSB tube with 200 ul of each prepared bacteria suspension.Vortex inoculated TSB tube for 10 seconds to mix.
- Plate 100ul of inoculate TSB onto TSA plates and incubate at 35+/-2°C for 18-24 hours as zero time baseline colony count.
- Incubate the inoculate TSB tube at 35+/-2°C for 18-24 hours
- After 18-24 hours take-out from the incubator the TSB tube and vortex for 10 seconds.
- Plate 100 ul of the inoculated TSB onto TSA plate as 24 hours time
- Incubate the inoculated TSA at 35+/-2°C for 18-24 hours
- Read and record the turbidity of the TSB and the UFC.

ACCEPTABILITY LIMITS: GROWTH.

PERFORMANCE TEST RESULTS:

STRAIN*	ZERO TIME COUNT; CFU/PLATE	CFU COUNT AFTER 18-24 HOURS at 35°C±2°C; CFU/PLATE
<i>Escherichia coli</i> ATCC 25922	188	CONFLUENT GROWTH
<i>Staphylococcus aureus</i> ATCC 25923	201	CONFLUENT GROWTH
<i>C.frendii</i> ATCC 8090	233	CONFLUENT GROWTH

* Performance testing with Copan TSB was conducted using laboratory strains. Performance testing was not conducted using human specimens. The product TSB is in compliance with CLSI-M22-A3.

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BIBLIOGRAPHY

- Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. Manual of Clinical Microbiology. 7th ed. ASM, Washington, D.C..
- National Committee for Clinical Laboratory Standards (CLSI). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
- Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
- Clinical Microbiology Procedures Handbook - 3rd ed./editor in chief, third edition and 2007 update , Lynne S.Garcia . 2010 ASM Press Washington, DC
- Atlas,Fourth Ed. 2010 Handbook of Microbiological Media p. 1842. CRC Press, Boca Raton, FL

Tabella dei Simboli/ Index of Symbols/ Table des Symboles/ Symboltabelle/ Tabla de símbolos/ Tabela de símbolos/ Tabell över symboler/ Symboltabelle/ Tabel med symboler

Simbolo/Symbol/Symbole	Significato /Meaning/ Signification/ Bedeutung/ Significado/ Betyder/ Betydning
	Fabbricante/ Manufacturer/ Fabricant/ Hersteller/ Fabricante/ Tillverkare/ Fabrikant/ Producent
	Dispositivo diagnostico in vitro/ In vitro diagnostic device / Dispositif de diagnostic in vitro/ Diagnosegerät in vitro/ Dispositivo de diagnóstico in vitro/ Medicinsk anordning för in vitro-diagnostik/ In vitro diagnostisk medicinsk anordning/
	Non riutilizzare/ Do not reuse/Ne pas réutiliser/Nicht zur Wiederverwendung/No reutilizar/ Não volte a usar/ Må ikke brukes på nytt/ Må ikke genbruges/ Får inte återanvändas
	Numero di catalogo/ Catalogue number/Référence du catalogue/Bestellnummer/Número de catálogo/ Referência do catálogo/ Katalognummer/ Katalognr
	Limiti di temperatura/Temperature limitation/Limites de temperature/Temperaturbegrenzung/Limites de temperatura/ Limites de temperatura/ Temperaturgränser / Temperaturbegrensninger
	Utilizzare entro/Use by/Utiliser jusque/Verwendbar bis/Fecha de caducidad/ Prazo de validade/ Ska användas innan/ Må brukes innen/ Anvendes før
	Consultare le istruzioni per l'uso/Consult Instructions for Use/Consulter les instructions d'utilisation/Gebrauchsanweisung beachten/Consulte las instrucciones de uso/ Consultar as instruções de utilização/ Se instruksjoner for bruk/ Se vejledningen til brug/ Les bruksanvisningen
	Codice del lotto (partita)/ Batch code (Lot)/ Code de lot (Lot)/Chargencode (Chagenbezeichnung)/ Código de lote (Lote)/ Serienummer(parti)/ Lot nummer (parti)/ Batch-nummer (parti)
	Contenuto sufficiente per <n> test/ Contains sufficient for <n> tests/Contenu suffisant pour <n> tests/ Ausreichend für <n> Tests/ Contenido suficiente para <n> pruebas/ Contêmo suficiente para <n> testes/ Innehåller tillräckligt för<n> tester/ Inhold tilstrekkelig for <n> test/ Indhold tilstrækkeligt til <n> prøver

 <p>Copan Italia Spa Via Perotti, 10 25125 Brescia, Italy</p>	<p>Copan Italia Spa Via Perotti 10 25125 Brescia Italy Tel: +39 030 2687211 Fax: +39 030 2687250 E-mail: info@copangroup.com Website: www.copangroup.com</p>	 <p>COPAN Innovating Together™</p>
<p><u>North American Distributor:</u> Copan Diagnostics Inc. 26055 Jefferson Avenue Murrieta, CA 92562 USA Tel: 951-696-6957 Fax: 951-600-1832</p> <p>E-mail: customerservice@copanusa.net Website: www.copanusa.com</p>		