

IVD solutions through partnership




CHROMagar™ Candida Plus

For detection and differentiation of major clinical Candida species, including *C. auris*

A large circular inset image shows a petri dish with CHROMagar Candida Plus. The agar surface is covered with various colonies of Candida species. Some colonies are pinkish-purple, some are yellowish, and some are bright green. The colonies vary in size and shape, from small dots to larger, more elongated structures. The background of the image is a light, neutral color.

CHROMagar
The Chromogenic Media Pioneer

● CHROMagar™ Candida Plus



Back Front **Pate Reading**

- *C. auris*
→ Light blue with blue halo
Blue from the back side
- *C. albicans*
→ Green-blue
- *C. tropicalis*
→ Metallic blue with pink halo
- *C. glabrata*
→ Mauve
- *C. krusei*
→ Mauve and fuzzy



Manufacturer:
CHROMagar
4 place du 18 juin 1940 75006 Paris - France
e-mail: CHROMagar@CHROMagar.com
www.CHROMagar.com

For detection and differentiation of major clinical Candida species, including *C.auris*

Background

The Candida are yeasts involved in various infections called Candidiasis, which can affect damaged skin, respiratory tract, digestive and urogenital systems. These Candidiasis can be severe with significant morbidity for nosocomial infections or in immunocompromised patients. Although *C. albicans* is still the main species involved, the use of antifungal agents has given rise to other species such as *C. tropicalis*, *C. krusei* and *C. glabrata*.

In 2016, The World Health Organization added to this list *C. auris*, with a prevalence of over 90 % resistant to fluconazole. In addition, some strains are multidrug resistant to amphotericin B, voriconazole, and/or echinocandins.

It is recommended to carry out an early diagnosis of Candida in order to provide specific treatment as quickly as possible. Candida can be isolated by swabbing the skin, throat, rectum, or urogenital tract.

CHROMagar™ Candida Plus is the first chromogenic isolation medium to detect and differentiate *C. auris* in addition to other major clinical Candida species such as *C. albicans*, *C. tropicalis*, *C. glabrata* or *C. krusei*.

Medium Performance

1 High specificity:

Differentiation of the most common Candida species with very high specificity:

- C. albicans* ≈ 100 % *
- C. tropicalis* ≈ 100 % *
- C. krusei* ≈ 100 % *

2 Unique medium to differentiate *C. auris* from other candida species.

Owe to its high specificity, it can be used also as a screening tool in case of outbreaks.

For <i>C. auris</i>	Specificity ≈ 100 % *
	Sensitivity ≈ 100 % *

* Specificity and Sensitivity from scientific study: «Evaluation of a novel chromogenic medium for Candida spp. identification and comparison with CHROMagar™ Candida for the detection of *Candida auris* in surveillance samples» Juan V. Mulet et al., 2020.

3 Easy identification:

Identification by MALDI-TOF can be carried directly from a colony. No need of subculture.

Medium Description

Powder Base	Total	50.9 g/L
	Agar	15.0
	Peptones	11.0
	Chromogenic and selective mix	24.9
	Storage at 15/30°C - pH: 6.1 +/- 0.2	
	Shelf Life	12 months

Usual Samples	Skin, throat, armpits, urogenital tract and rectal swab.
Procedure	Direct Streaking. Incubation 36-48h at 30-37 °C Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Distributed by:
Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
e-mail: mast@mast-diagnostica.de
www.mast-group.com

Ordering Information

Product	Order Code
CHROMagar™ Candida Plus dry media, 5 liter	15CA242
CHROMagar™ Candida Plus ready to use plates, 20 pcs.	201406