**UTM<sup>®</sup>**

## **Instructions for Use**

**CE** **IVD**

English	3
Italiano	9
Español	15
Deutsch	21
Français	28
Português	34
БЪЛГАРСКИ	41
Česky	48
Dansk	55
ΕΛΛΗΝΙΚΑ	62
Eesti Keel	70
Hrvatski	77
Latviešu	83
Lietuvių K.	90
Magyar	97
Nederlands	104
Norsk	111
Polski	117
Română	125
Slovenčina	132
Slovenščina	139
Srpski	145
Suomi	152
Svenska	160
Türkçe	167
УКРАЇНСЬКА	173

## Copan Universal Transport Medium (UTM-RT®) System Instructions for use

### INTENDED USE

Copan Universal Transport Medium (UTM-RT®) System is intended for the collection and transport of clinical specimens containing viruses, chlamydiae, mycoplasma or ureaplasma from the collection site to the testing laboratory. UTM-RT® can be processed using standard clinical laboratory operating procedures for viral, chlamydial, mycoplasma and ureaplasma culture.

### SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of infections caused by viruses, chlamydiae, mycoplasmas or ureaplasmas involves collection and refrigerated transport of biological specimens. Using the UTM-RT® System, the collected specimen can be stored for up to 48 hours at 2-25°C.

The UTM-RT® consists of a Hanks' Balanced Salt Solution (HBSS) enriched with proteins and sugars with a neutral pH and pH indicator. The medium contains some antibiotics and antimycotics to inhibit overgrowth of bacteria and yeasts, maintain cellular integrity and encourage preservation of viruses and chlamydiae.

Copan UTM-RT® System medium is provided in labeled screw-cap tubes designed for transport of the clinical sample. Copan UTM-RT® System is also supplied as a sample collection kit that comprises a package which contains one screw-cap tube of UTM-RT® medium and a peel pouch incorporating one or two sterile specimen collection swabs. A range of UTM-RT® sample collection kits are available which incorporate different types of shaft swabs which facilitate the collection of specimens from different sites of the patient as described below in the Directions for Use section. Once a swab sample is collected it should be placed immediately into the transport tube where it comes into contact with transport medium. To maintain optimum microorganism viability, transport the specimens to the laboratory as soon as possible.

After collection, the specimen should be stored at 2-25°C and processed within 48 hours.

If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder.

### REAGENTS

The UTM-RT® formulation includes proteins for virus stabilization<sup>17</sup>, antibiotics and antimycotics to prevent overgrowth of bacterial and fungal flora and a buffer solution to maintain a neutral pH.

Components
Sucrose
HBSS solution
Bovine serum albumin
Buffered solution
Gelatin
Amino acids
Antibiotics
Phenol Red

pH 7,3 ± 0,2 a 2+25°C

### REQUIRED MATERIALS BUT NOT PROVIDED

Materials suitable for isolation, differentiation and culture of viruses, chlamydiae, mycoplasmas and ureaplasmas.

### STORAGE

The product must be stored in its original packaging at a temperature between 2 and 25°C until the time of use. Do not overheat or freeze prior to use.

### SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Specimens for virus, chlamydia, mycoplasma or ureaplasma investigation should be collected and handled following published manuals and guidelines<sup>2,3,4,7,9,10,11</sup>. To maintain optimum viability, transport the specimen to the laboratory as soon as possible. After the collection, the specimen should be stored at 2-25°C and processed within 48 hours. If there will be a long delay exceeding 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations<sup>8,11,12</sup>. Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

### MATERIALS SUPPLIED

Copan UTM-RT® System is available in the format of kit (tube with specimen collection swabs) or tube only, the system is contained in box and boxes are contained in an external carton (Pack size are provided in table one).


Copan UTM-RT® System includes a screw-cap tube with conical or round bottom containing 1.5ml, 1ml, 2ml, 2.5ml, 3ml or 10ml of light orange-red transport medium. For all the product codes, excluding 3U008N, 331C, 348C and 3E056N, the tube contains three glass beads that allow homogenizing the specimen by vortexing.

UTM-RT® System tubes of transport medium are supplied alone or in a kit format. Each kit unit consists of a package containing: a pre-labeled screw-cap tube and a peel pouch incorporating one or two sterile specimen collection swabs (e.g. applicator swabs polyester tipped, flexible size applicator swab with flocked nylon fiber). For more details on available configurations please refer to Table 1.

The collection swab applicators FLOQSwabs® and CLASSIQSwabs™ provided with UTM-RT® have a molded breakpoint in the shaft of the applicator. For the line of FLOQSwabs® applicator the molded breakpoint show a colored mark.

UTM-RT® tube capture caps have an internal moulded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a moulded docking receptacle in the cap (Fig.1). In the testing laboratory when the cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube.

Fig 1. Capture of broken swab applicator stick by UTM-RT® tube cap

**LIMITATIONS**

1. Because calcium alginate swabs are toxic for many enveloped viruses<sup>5</sup> and may interfere with immunofluorescence tests<sup>2</sup>, they should not be used for specimen collection.
2. Wooden shaft swabs may contain toxins and formaldehydes<sup>2,6</sup> and should not be used.
3. The product codes without beads (REF 3U008N, 331C, 348C and 3E056N) are not suitable for use with mucous or particularly viscous specimens.
4. UTM-RT® kits are intended to be used with the medium tubes and swabs provided in the kit. The use of tubes of medium or swabs from another source could affect the performance of the product.

**WARNINGS AND PRECAUTIONS**

1. Single-use device for professional in vitro diagnostic use.
2. UTM-RT® System is ready for use and requires no further preparation.
3. Do not use beyond the expiry date.
4. Do not use the UTM-RT® medium for premoistening or prewetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
5. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.<sup>7</sup>
6. Do not re-sterilize unused swabs.
7. Do not re-pack.
8. Not suitable to collect and transport microorganisms other than viruses, chlamydiae, mycoplasma and ureaplasma.
9. Do not ingest the medium.
10. Specimens for the search of viruses, chlamydiae, mycoplasmas and ureaplasmas must be collected and handled using personal protective equipment against biological risk according to published manuals and guidelines<sup>1,4,6,7,9,17</sup>.
11. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms.
12. Do not use UTM-RT® if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the color of the medium has changed from light orange-red, (4) the swab pouch is open, or (5) there are other signs of deterioration.
13. The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.
14. Do not bend or shape the swab before the collection of the specimen. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
15. Due to the design of the flexible minitip, the swab will coil when placed in the tube. Therefore, if necessary, remove the swab from the tube, use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.
16. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.

**INSTRUCTION FOR USE**

UTM-RT® System is ready for use and requires no further preparation. It is available in the various configurations listed in **Table 1**.

REF	PRODUCT DESCRIPTION		PACK SIZE	SAMPLING SITES*	PRODUCT DESIGNED FOR HAVING CAPTURE CAP FEATURE (FOR UTM KIT ONLY)
	TUBE	SWAB			
302C; 302C.LC.	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	Two regular size applicator swabs polyester tipped with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
305C; 305CMH;	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
328C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swabs polyester tipped with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
306C; CA302MH; 346C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO*

307C; CA303MH;	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO
357C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
321C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip size applicator swab with flocked nylon fiber tip with breaking point  One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Flexible Minitip: Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites  Regular: Nose, throat, vagina, rectum, faeces and wounds	NO*
323C,	2 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO*
338C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point  One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Minitip: eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites  Regular: Nose, throat, vagina, rectum, faeces and wounds	NO*
347C	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One mini flexible applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
355C; 355CW	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
356C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab polyester tipped with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	NO*
358C	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
359C; 359CE.A	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES

360C	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
361C	1 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
365C	2 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
366C; 3C004N	2 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Nose, throat, vagina, rectum, faeces and wounds	YES
3U006N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
3U008N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom – No Beads.	One flexible minitip applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
3C041N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	One Mid Turbinate applicator swab with flocked nylon fiber tip with breaking point	50 tubes per package 6 x 50 tubes per box	Midturbinate	NO*
3C042N	3 mL of UTM-RT® medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 tubes per package 6 x 50 tubes per box	Eye, ear, nasal, passages, nasopharynx, throat, urogenital tracts and pediatric sites	NO*
330C;	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box		NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service
331C; CA339MH	10 mL of UTM-RT® medium in 25x90 mm screw-cap tube with internal shaped conical bottom.	NA	15 tubes per package 6 x 15 tubes per box		NO
348C	2,5 mL of UTM-RT® medium in 16x100 mm screw-cap tube with internal shaped conical bottom – No Beads.	NA	50 tubes per package 6 x 50 tubes per box		NO

349C;	1,5 mL of UTM-RT <sup>®</sup> medium in 16x100 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box		NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service
3U044N	3 mL of UTM-RT <sup>®</sup> medium in 16x100 mm screw- blue cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box		NO
3E056N	3 mL of UTM-RT <sup>®</sup> medium in 16x100 mm screw-cap tube with internal shaped conical bottom – No Beads.	NA	50 tubes per package 6 x 50 tubes per box		NO
350C; 350CV	1 mL of UTM-RT <sup>®</sup> medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box		NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service
353C	3 mL of UTM-RT <sup>®</sup> medium in 12x80 mm screw-cap tube with internal shaped conical bottom.	NA	50 tubes per package 6 x 50 tubes per box		NOT APPLICABLE Capture cap feature depends on the type of swab. For information regarding the correct combined, please contact customer service

**Table 1: product description**

Not all the product codes (REF) are salable in all countries. Please contact Copan Customer care service for product codes availability for a specific country. NO\* = The swab could be accidentally and not firmly captured and might fall down.

*The sampling sites Table is just a suggested table. Performance testing with Copan UTM-RT<sup>®</sup> System was conducted using laboratory strains. Performance testing was not conducted using human specimens. Please refer to your internal procedures to choose the most appropriate device for the specific sampling site.*


**Specimen Collection**

Proper collection of the specimen from the patient is a crucial aspect for successful isolation and identification of infectious organisms. Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute illness.


**UTM-RT<sup>®</sup> in kit**

1. Open the UTM-RT<sup>®</sup> kit package and remove the medium test tube and the internal bag containing the sterile swab.
2. Take the sterile swab out of its bag and collect the clinical specimen; to prevent the risk of contamination, make sure that the swab tip comes into contact with the collection site only. **NOTE:** Do not bend the swab before the collection of the specimen. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
3. After collecting the specimen, insert the swab into the test tube until the breakpoint is level with the test tube opening.
4. Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
5. Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
6. Screw the cap back onto the test tube and hermetically seal it.

**Fig 1. Collection swab showing breakpoint indication line and area for holding the applicator**

Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens and care should be taken to avoid splashes and aerosols when breaking the swab stick into the tube of medium. During sample collection when handling the swab applicator, the operator must not touch the area below the colored breakpoint indication line; that is the area from the line to the tip of the nylon flocked swab, as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.

#### UTM-RT<sup>®</sup> in bulk

1. Aseptically remove cap from tube taking care not to spill the medium.
2. Aseptically place vesicle aspirates<sup>26</sup>, corneal or conjunctival scrapings<sup>28</sup>, small pieces of tissue or stool<sup>30</sup> samples into the tube with UTM-RT<sup>®</sup> medium.
3. Screw the cap onto the test tube and hermetically seal it.
4. Identify the tube containing the specimen.
5. Send to the laboratory for immediate analysis.

If processing is delayed (over 48 hours), the specimens must be frozen at -70°C or colder.

#### DISPOSAL

Waste must be disposed of in compliance with local legislation.

#### QUALITY CONTROL

The UTM-RT<sup>®</sup> lots are tested for microbial contamination, toxicity for the cellular lines used for the viral cultures and the ability to maintain the viability of viral, chlamydia and mycoplasma strains for up to 48 hours at 2-25°C in accordance with the methods described in CLSI M40-A2<sup>5</sup>.

#### RESULTS AND PERFORMANCE

The results obtained largely depend on proper and adequate specimen collection as well as the promptness with which the specimens are transported to the laboratory and analysed.

Viability studies were performed using Copan UTM-RT<sup>®</sup> with a panel of representative strains of the various families supported by the UTM-RT<sup>®</sup>. The swabs that accompany each transport system were directly inoculated in triplicate with 100 µl of organism suspension using laboratory strains. Performance testing was not conducted using human specimens. Subsequently, the swabs were inserted in the respective test tubes containing the transport medium and stored for 0 and 48 hours at 2-6°C and at controlled room temperature (20-25°C). At the time of processing, each swab was vortexed for 20 seconds and removed from its transport medium test tube, then, an aliquot of the suspension was inoculated into the cellular line (200 µl) or into the appropriate culture medium<sup>6,15</sup>. All the cultures were processed using the standard laboratory culture technique<sup>6,15</sup>. Organism viability was determined by fluorescent cell counting for viral and chlamydia strains and CFU counting for mycoplasma and ureoplasma strains. The acceptability limits for time zero and for 48 hours were defined in accordance with the regulations M40-A2<sup>5</sup>. Copan UTM-RT<sup>®</sup> System preserved the viability of all the organisms tested for 48 hours at both controlled room temperature and in the refrigerator in the above described test conditions. The organisms evaluated and the results obtained are given in the table below.

Organism	ATCC <sup>®</sup> number	% of reduction of fluorescing infected cells after 48 hours time 2-6°C	% of reduction of fluorescing infected cells after 48 hours time 20-25°C
Herpes Simplex Virus Type 1	ATCC <sup>®</sup> VR-539	≤50%	≤50%
Herpes Simplex Virus Type 2	ATCC <sup>®</sup> VR-734	≤50%	≤50%
Respiratory Syncytial Virus	ATCC <sup>®</sup> VR-1580	≤50%	≤50%
Coxsackie B1 Virus	ATCC <sup>®</sup> VR-28	≤50%	≤50%
Chlamydia trachomatis	ATCC <sup>®</sup> VR-880	≤50%	≤50%
Influenza A	ATCC <sup>®</sup> VR-1679	≤50%	≤50%
Cytomegalovirus	ATCC <sup>®</sup> VR-977	≤50%	≤50%
Varicella-zoster virus	ATCC <sup>®</sup> VR-1367	≤50%	≤50%
Chlamydia pneumoniae	ATCC <sup>®</sup> VR-1360	≤50%	≤50%
Adenovirus	DHI 20-4740010	45%	47%
Parainfluenza virus type 3	DHI 20-4770010	0%	41%
Echovirus type 30	ATCC <sup>®</sup> VR-322	13%	24%

Organism	ATCC® number	% of reduction of CFUs after 48 hours time 2-6°C	% of reduction of CFUs after 48 hours time 20-25°C
Mycoplasma hominis	ATCC® 23114	35%	25%
Ureaplasma urealyticum	ATCC® 27618	0%	6%
Mycoplasma pneumoniae	ATCC® 15531	50%	50%

#### TABLE OF SYMBOLS

See the table of symbols at the end of the instructions for use.

#### NOTES FOR THE PROFESSIONAL USER

In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

#### REVISION HISTORY

Last Revision No.*	Release date	Changes made
01	07-2022	Revision of IFU sections (first revision in IVDR)

\* Should you need earlier revisions, contact Copan Customer Service.

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## BIBLIOGRAPHY

1. Gary W. Procop and Elmer W. Koneman, 2016. Color Atlas and Textbook of Diagnostic Microbiology, Seventh edition. Wolters Kluwer Health.
2. James H. Jorgensen, Michael A. Pfaller, Karen C. Carroll, Guido Funke, Marie Louise Landry, Sandra S. Richter, David W. Warnock, 2015. Manual of Clinical Microbiology, 11<sup>th</sup> Edition. ASM, Washington, DC.
3. James Versalovic, Karen C. Carroll, Guido Funke, James H. Jorgensen, Marie Louise Landry, David W. Warnock, 2011. Cumitech 15A. Laboratory Diagnosis of Viral Infections. ASM, Washington, DC.
4. Patricia Tille. 2014. Bailey & Scott's Diagnostic Microbiology, 13<sup>th</sup> Edition. Laboratory Medicine.
5. Clinical and Laboratory Standards Institute (CLSI), 2014. M40-A2 Quality Control of Microbiological Transport Systems; Approved Standard-Second Edition.
6. Clinical and Laboratory Standards Institute (CLSI), 2006. M41-A Viral Culture; Approved guidelines.
7. Wardford, A., M. Chernesky, and E. M. Peterson, 1999. Cumitech 19A, Laboratory Diagnosis of Chlamydia trachomatis Infections. ASM, Washington DC.
8. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
9. J. Michael Miller, Shelley A. Miller, 2017. A Guide to Specimen Management in Clinical Microbiology, Third Edition. ASM, Washington DC.
10. Centers for Disease Control and Prevention (CDC), 2016. Guide for Shipping Infectious Substances.
11. Centers for Disease Control and Prevention (CDC), 2009. Biosafety in Microbiological and Biomedical Laboratories 5<sup>th</sup> Edition.
12. World Health Organization 2015. Guidance on regulations for the Transport of Infectious Substances 2015 – 2016.
13. Centers for Disease Control and Prevention (CDC), 2002. Screening Tests to Detect Chlamydia trachomatis and Neisseria gonorrhoeae Infections - 2002.
14. Centers for Disease Control and Prevention (CDC), 2014. Recommendations for the Laboratory-Based Detection of Chlamydia trachomatis and Neisseria gonorrhoeae - 2014.
15. European Collection of Authenticated Cell Cultures (ECACC). Fundamental Techniques in Cell Culture Laboratory Handbook, 3<sup>rd</sup> edition (Protocol 9 – Detection of Mycoplasma by Culture Isolation).
16. J.B. Mahony, M.A. Chernesky Effect of Swab Type and Storage Temperature on the Isolation of Chlamydia trachomatis from Clinical Specimens Journal of Clinical Microbiology, Nov. 1985, p. 865-867.
17. S. Specter, R. L. Hodinka, S. A. Young. Clinical Virology Manual, fifth edition, 2016.
18. Peter Daley, Santina Castriciano, Max Chernesky, Marek Smieja. Comparison of Flocked and Rayon Swabs for Collection of Respiratory Epithelial Cells from Uninfected Volunteers and Symptomatic Patients. Journal of Clinical Microbiology, June 2006.
19. K. Loens, L. Van Heirstraeten, S. Malhotra-Kumar, H. Goossens, and M. Ieven. Optimal Sampling Sites and Methods for Detection of Pathogens Possibly Causing Community-Acquired Lower Respiratory Tract Infections. Journal of Clinical Microbiology, Jan. 2009.
20. Richard Garceau MD, Danielle Leblanc RT, Louise Thibault MD, Gabriel Girouard MD, Manon Mallet PhD. Herpes simplex virus type 1 is the leading cause of genital herpes in New Brunswick. Can J Infect Dis Med Microbiol 2012.
21. Marek Smieja, Santina Castriciano, Susan Carruthers, Geoffrey So, Sylvia Chong, Kathy Luinstra, James B. Mahony, Astrid Petrich, Max Chernesky, Mario Savarese, and Daniele Triva. Development and Evaluation of a Flocked Nasal Midturbinate Swab for Self-Collection in Respiratory Virus Infection Diagnostic Testing. Journal of Clinical Microbiology, Sept. 2010.
22. Donghyok Kwon, Kyeongcheol Shin, Mihwa Kwon, Hee-Bok Oh, Chun Kang, Joo-Yeon Lee. Development and Evaluation of Rapid Influenza Diagnostic Test for the Pandemic (H1N1) 2009 Virus. JCM Accepts, published online ahead of print on 27 October 2010.
23. Bhupesh K. Prusty, Christine Siegl, Petra Hauck, Johannes Hain, Suvi J. Korhonen, Eija Hiltunen- Back, Mirja Puolakkainen, Thomas Rudel. Chlamydia trachomatis Infection Induces Replication of Latent HHV-6. Plos One, April 2013.
24. Ellen Vancutsem, Oriane Soetens, Maria Breugelmanns, Walter Foulon, Anne Naessens. Modified Real-Time PCR for Detecting, Differentiating, and Quantifying *Ureaplasma urealyticum* and *Ureaplasma parvum*. The Journal of Molecular Diagnostics, Vol. 13, No. 2, March 2011.
25. Turkiya Al-Siyabi, Khalifa Binkhamis, Melanie Wilcox, Sallene Wong, Kanti Pabbaraju, Raymond Tellier, Todd F. Hatchette, and Jason J. Le Blanc. A cost-effective real-time PCR for the detection of adenovirus from viral swabs. Al-Siyabi et al. Virology Journal 2013.
26. Eidgenössisches Departement für Verteidigung, Bevölkerungsschutz und Sport VBS; SAMPLING AND DISPATCH OF SAMPLES - VIROLOGY Ident-Nr/Vers. 10014156143/01 Aktenzeichen:323-FachinfoViro2016.
27. Matthew L. Faron et al (2016); Multicenter Evaluation of Meridian Bioscience HSV 1&2 Molecular Assay for Detection of Herpes Simplex Virus 1 and 2 from Clinical Cutaneous and Mucocutaneous Specimens; Journal of Clinical Microbiology Volume 54 Number 8
28. Al-Siyabi et al. (2013); A cost effective real-time PCR for the detection of adenovirus from viral swabs. Virology Journal 10:184.
29. Le Roy C, Pereyre S, Hénin N, Bébéar C. 2017. French prospective clinical evaluation of the Aptima Mycoplasma genitalium CE-IVD assay and macrolide resistance detection using three distinct assays. J Clin Microbiol 55:3194–3200. <https://doi.org/10.1128/JCM.00579-17>.
30. Minodier et al. (2017) Clinical and virological factors associated with gastrointestinal symptoms in patients with acute respiratory infection: a two-year prospective study in general practice medicine BMC Infectious Diseases 17:729 DOI 10.1186/s12879-017-2823-9

Index of Symbols / Tabella dei Simboli / Tabla de símbolos / Symboltabelle / Table des Symboles / Tabela de símbolos / Индекс на символите / Seznam symbolů / Symbolforklaring / Ευρετήριο Συμβόλων / Tingmärkide kirjeldused / Symbolihakemisto / Kazalo simbola / Szimbólumok jegyzéke / Simbolių rodyklė / Simbolu rādītājs / Index van symbolen / Symbolindeks / Indeks symboli / Index simboluri / Zoznam symbolov / Indeks simbolov / Indeks simbola / Förteckning över symboler / Semboller Dizini / Таблица символов

<p>Symbol / Simbolo / Símbolo / Symbol / Symbole / Símbolos / Символ / Symbol / Symbol / Σύμβολο / Tingmärk / Symboli / Símbol / Szimbólum/ Simbolis / Simbols / Symbool / Symbol / Symbol / Símbol / Símbol / Symbol / Símbol / Símbol / Símbol / Símbol / Символ</p>	<p>Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значение / Význam / Betydning / Σημασία / Tähendus / Merkitys / Značenje / Jelentés/ Reikšmė / Nozīme / Betekenis / Betydning / Znaczenie / Semnificație / Význam / Pomen / Značenje / Betydelse / Anlam / Значення</p>
	<p>Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante /Производитель / Výrobce / Producent / Κατασκευαστής / Tootja / Valmistaja / Proizvođač / Gyártó/ Gamintojas / Ražotājs / Fabrikant / Produsent / Producent / Producător / Výrobca / Proizvajalec / Proizvođač / Tilverkare / Üretici / Виробник</p>
	<p>In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / Инвитро диагностично изделие / Diagnostický prostředek pro in vitro diagnostiku / Udstyr til in vitro-diagnostik / Διαγνωστικό προϊόν in vitro / In vitro diagnostika seade / In vitro -diagnostikkalaitte / In vitro dijagnostički proizvod / In vitro diagnosztikai eszköz / „In vitro“ diagnostikos prietaisas / In vitro diagnostikas ierīce / Hulpmiddel voor in-vitrodiagnostiek / In vitro-diagnostisk enhet / Urządzenie do diagnostyki in vitro / Dispositiv de diagnostic in vitro / Diagnostická pomůcka in vitro / In vitro diagnostični pripomoček / In vitro dijagnostičko medicinsko sredstvo / Enhet för in vitro-diagnostik / In vitro tani cihaz / Виріб для діагностики in vitro</p>
	<p>Medical Device / Dispositivo Medico / Producto sanitario / Medizinprodukt / Dispositif médical / Dispositivo Médico / Медицинско изделие / Zdravotnický prostředek / Medicinsk udstyr / Ιατροτεχνολογικό προϊόν / Meditsiiniseade / Lääkinnällinen laite / Medicinski proizvod / Orvostechnikai eszköz / Medicinos prietaisas / Mediciniska ierīce / Medisch hulpmiddel / Medisinsk utstyr / Wyrób medyczny / Dispozitiv medical / Zdravotnícka pomůcka / Medicinski pripomoček / Medicinski aparat / Medicinteknisk enhet / Tibbi Cihaz / Медицинский прибор</p>
	<p>Unique Device Identifier / Identificativo unico del dispositivo / Identificador único de dispositivo / Einmalige Produktkennung / Identifiant dispositif univoque / Identificador de dispositivo único / Уникален идентификатор на изделието / Jedinečný identifikátor zařízení / Unik enheds-id / Ατομική/Επίσημη ταυτοποίηση ιατροτεχνολογικού προϊόντος / Seadme kordumatu identifitseerimistunnus / Yksilötöivä laitetunniste / Jedinstveni identifikator proizvoda / Egyedi eszközazonosító / Unikalus įrenginio identifikatorius / Unikālais ierīces identifikators / Unieke hulpmiddelidentificatie / Identifikator for Unique-enhet / Unikatosy numer urządzenia / Identificator unic dispozitiv / Unikátny identifikátor pomůcky / Edinstveni identifikator naprave / Jedinstveni identifikator uređaja / Unik enhetsidentifierare / Benzersiz Cihaz Tanımlayıcı / Унікальний ідентифікатор пристрою</p>
	<p>CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / CE маркировка / Označení CE / CE-mærkning / Σήμανση CE / CE-mārgis / CE-merkintä / Oznaka CE / CE-jelölés / CE ženklas / CE markėjums / CE-markering / CE-merkning / Oznakowanie CE / Marcaj CE / Označenie CE / Oznaka CE / CE oznaka / CE-märkning / CE işareti / Маркування CE</p>
	<p>Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Идентификационный номер на нотифицирующая орган / Identifikačný číslo notifikovaného subjektu / Identifikationsnummer af bemyndiget organ / Αριθμός ταυτοποίησης κοινοποιημένου οργανισμού / Teavitatud asutuse identifitseerimisnumber / Ilmoitetun laitoksen tunnusnumero / Identifikációsjski broj obaviještenog tijela / A bejelentett szervezet azonosító száma/ Notifikuotosios įstaigos identifikavimo numeris / Plinvarotās iestādes identifikācijas numurs / Identificatienummer van de aangeamelde instantie / Identifikasjonsnummeret til det tekniske kontrollorganet / Numer identyfikacyjny jednostki notyfikowanej / Numărul de identificare al organismului notificat / Identifikačné číslo notifikovanej osoby / Identifikacijska številka priglašeneга organa / Identifikacioni broj nadležnog tela / Identifikationsnummer för anmänt organ / Onaylı kuruluş tanımlama numarasi / Идентифікаційний номер уповноваженого органу</p>

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