LBM[®]
MycoTB[™]

Instructions for Use

Français _____ **3**

English _____ **9**

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











SAMPLE PROCESSING/ PREPARATION PROCEDURE/PROTOCOL

INTRODUCTION

Myco TB[™] is intended for the fluidization, decontamination, concentration and resuspension of clinical specimens for culture investigations on solid and liquid media for the detection of Mycobacteria.

NOTE: Use tubes of Myco TB[™] taken from the same finished product lot number.

> Myco TB[™]

Mycob TB™ – Instructions for use

INTENDED USE

Mycob TB™ is intended to be used for fluidification, decontamination, concentration and resuspension of clinical specimens for culture screening of solid and liquid culture media to screen for mycobacteria.

SUMMARY AND PRINCIPLES

The bacteria belonging to the *Mycobacterium spp* genus on the cellular wall level have a high lipid content that characterises the alcohol-acid-resistance phenomenon. Mycobacteria infections, in particular tuberculosis, are still responsible for high morbidity and mortality rates in developing countries but also in industrialised countries. Tuberculosis is caused by species belonging to the *Mycobacterium tuberculosis* complex and is still a clinical and public health emergency. Recently, the problem has grown worse because of the appearance of a once rare but now more frequent resistance, also involving a larger number of drugs. Nontuberculous mycobacteria are ubiquitous in the environment, to date their interhuman transmission has been ruled out, quite a few species have been described, and they may normally be the cause of infections (mycobacteriosis) in immunocompromised but also in normoergic subjects.

Mycobacterium kansasii is the species most commonly responsible for pulmonary disease in immunocompetent subjects. Because all the organs can be affected by both *Mycobacterium tuberculosis* complex and nontuberculous mycobacteria, microbiology laboratories should expect to receive pulmonary and a variety of extra-pulmonary specimens; in such samples false-positives are unlikely while false-negatives are very frequent⁽¹⁾. Both pulmonary and extra-pulmonary samples need to undergo to decontamination, fluidification and concentration steps before to be processed for diagnostic tests, such as culture on solid and liquid culture media and/or microscopic and molecular testing. One of the widely used decontamination and fluidification techniques is that with N-acetyl-L-cysteine–sodium hydroxide (NALC-NaOH)⁽²⁾.

PRODUCT DESCRIPTION


Fig. 1: Mycob TB™ constituent components

Mycob TB™ is a product consisting of 4 test tubes: (A) REACTION TUBE: empty tube; (B) Mycob TB™ SOLUTION TUBE: tube containing the fluidification and decontamination solution; (C) NEUTRALIZING TUBE: tube containing the neutralizing reagent (neutralization/washing); (D) RESUSPENSION TUBE: tube containing the pellet resuspension solution.

- A. REACTION TUBE: a graduated PP test tube with screw cap with a capacity of 50mL.
- B. Mycob TB™ SOLUTION TUBE: a PET test tube with screw cap filled with 5mL of specimen fluidification and decontamination solution. The PET material of which the tube is made performs the function of oxygen barrier and hence does NOT allow oxidation of the solution containing the active ingredient for fluidification.
- C. NEUTRALIZING TUBE: an HDPE test tube with screw cap filled with 45mL of buffered saline solution (PBS) suitable for neutralization of the fluidification and decontamination reagent.
- D. RESUSPENSION TUBE: a sterile PP test tube with screw cap filled with 2mL of sterile buffered saline solution (PBS) useful for pellet resuspension after centrifugation.

| Code | Product description | Content |
|--------|--|---|
| 0U012N | Mycob TB™ Solution Kit Configuration | The test tubes making up the kit are supplied in a sealed bag containing everything necessary for one analysis. 1 box contains everything necessary for 30 analyses. (6 boxes of 5 kits) |
| 0E025N | Mycob TB™ Solution Bulk Configuration | The test tubes are supplied in a single box containing: (A) REACTION TUBE: one bag of 30 pieces (B) Mycob TB™ SOLUTION TUBE: one box of 30 pieces (C) NEUTRALIZING TUBE: two boxes of 15 pieces (D) RESUSPENSION TUBE: one box of 30 pieces 1 box contains everything necessary for 30 analyses. |

REAGENTS

The Myco TB[™] formulation comprises: ⁽³⁾

Reagent components Myco TB[™] SOLUTION TUBE (B):

Sodium hydroxide solution
L-N-Acetyl-cysteine

Neutralizing and pellet resuspension reagent components NEUTRALIZING TUBE (C) and RESUSPENSION TUBE (D): ⁽³⁾

Sodium phosphate
Sodium chloride
Potassium phosphate

REQUIRED MATERIALS BUT NOT PROVIDED

Suitable materials for the culture and isolation of mycobacteria.

List of suggested materials to perform the analysis:

1. Single-use graduated pipettes
2. Timer
3. Vortexer or mixer with adjustable speed up to a maximum of 2500/3000 RPM.
4. Centrifuge with rack suitable to hold test tubes with a capacity of 50mL. It is advisable to use a refrigerated centrifuge ⁽³⁾.

STORAGE

Myco TB[™] is ready for use and requires no further preparation.

It can be stored at a temperature between 2 and 25°C until use or the expiry date.

Do not overheat.

Do not incubate or freeze prior to use.

Improper storage can cause a loss of efficacy.

Do not use after the expiry date printed on the product packaging.

LIMITATIONS

1. Not suitable for any other application different from its intended use.
2. Specimens with a volume of over 5mL or under 1mL cannot be processed.

WARNINGS AND PRECAUTIONS

1. Myco TB[™] has been tested with pulmonary clinical specimens (sputum, induced sputum, bronchial aspirate or bronchial lavage), extra-pulmonary (urine, cavity fluids, gastric aspirate, biopsy, lymph nodes, purulent exudate) and negative clinical matrices inoculated with a reference ATCC (urine, cephalorachidian fluid). Due to the particularity of the observed specimens and pathogens and the low incidence of tuberculosis in Europe⁽⁴⁾, it is however recommended to use the product after internal laboratory validation.
2. The methods of collection, transport and storage of the specimen represent significant variables in obtaining reliable results. Since microbiological diagnosis of tuberculosis typically requires non-routine specific methods and reagents, as well as more time and biosecurity requirements than common bacteriological diagnostics⁽¹⁾ it is appropriate to follow the recommended guidelines for the detection of mycobacteria⁽³⁾ in clinical specimens.
3. The product is intended for *in vitro* use only and may not be used for therapeutic or prophylactic purposes under any circumstances.
4. This product is intended for single use; if reused, the results may be unreliable.
5. Use only test tubes from the same finished product batch. Product performance is not guaranteed if individual Myco TB[™] components from different finished product batches are used.
6. For code 0E025N, the finished product batch number is indicated on the package of each component.
7. The Myco TB[™] SOLUTION (B) must be COLOURLESS. Any colour indicates that the Myco TB[™] SOLUTION may have deteriorated and cannot be used.
8. Myco TB[™] is supplied complete with all the components necessary for specimen processing (fluidification, decontamination, concentration and resuspension). Do not use Myco TB[™] with components different from those supplied. Any precipitates in the Myco TB[™] SOLUTION TUBE (B) do not compromise use of the reagent and do not interfere with use of the product.
9. Use a structure and/or instruments suitable for handling clinical specimens suspected of containing risk group 3 microorganisms (laboratory with level 3 biosafety/containment measures, class I, II or III laminar flow hoods where all clinical specimens and potentially infected containers are to be handled, centrifuges with baskets and anti-aerosol devices suitable for handling clinical specimens suspected of containing risk group 3 microorganisms)^(5,3).
10. Use of this product in non-culture techniques and/or in combination with diagnostic tests or diagnostic instruments must be validated in advance.
11. Do not ingest the reagents.
12. All clinical specimens treated with Myco TB[™] are to be considered potentially infectious and must be handled with appropriate precautions. ^(5,3) In the event of spillage, appropriately disinfect following the laboratory guidelines or internal procedures.
13. The product may only be used by adequately trained and qualified persons. ^(5,3)
14. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.
15. Carefully read and follow the instructions.

16.

Danger - Contains Sodium Hydroxide - Causes severe skin burns and eye damage.

Wear protective gloves/protective clothing/eye protection/face protection - IF SWALLOWED: rinse mouth. Do NOT induce vomiting - IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]. Immediately call a POISON CENTER/doctor - IF INHALED: Remove person to fresh air and keep comfortable for breathing - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor - Dispose of contents/container in accordance with local regulation.

PRODUCT DETERIORATION

Do not use one or more Myco TB™ components if:

- (1) the product shows visible signs of damage or contamination;
- (2) the product shows visible signs of leakage from one of its components;
- (3) the product has passed its expiry date which is printed clearly on the product package;
- (4) the product shows any other visible signs of deterioration.

INSTRUCTIONS FOR USE

Specimen collection

Adequate collection of specimens from the patient is an extremely critical factor for the successful isolation and identification of infectious organisms.^(2,3)

For the collection of original specimens follow the guidelines or internal laboratory procedures.

Laboratory use

All the operations are to be carried out under a laminar flow hood. ⁽²⁻³⁻⁵⁾

Take the test tubes out of each package and prepare the material listed below:

1. A test tube rack
2. The REACTION TUBE (A) (empty tube).
3. Prepare a NEUTRALIZING TUBE (C) in correspondence to each REACTION TUBE (A).
4. Prepare a Myco TB™ SOLUTION TUBE (B) in correspondence to (possibly in front of) the REACTION TUBE (A).
5. Prepare a RESUSPENSION TUBE (D) in correspondence to (possibly behind) the REACTION TUBE (A).

Assess the quantity and quality of the original specimen (minimum treatable volume 1mL, maximum treatable volume 5mL).

Transfer of the specimen to the REACTION TUBE (A) (Fig. 2):

1. Open the REACTION TUBE (A).
2. Using a sterile pipette transfer 5mL of the original specimen to the REACTION TUBE (A).
3. Discard the pipette.

If the available volume of the original specimen is under 5mL:

make sure that you transfer at least 1mL of specimen to the REACTION TUBE (A) using a sterile pipette and bring up to a volume of 5mL by transferring the neutralizing reagent from the NEUTRALIZING TUBE (C) using another sterile pipette.


Fig. 2 - Specimen transfer to the REACTION TUBE (A)

Decontamination and fluidification phase (Fig. 3): ^(2,3)

1. Open the Myco TB™ SOLUTION TUBE (B).
2. Pour the entire content (5mL) into the REACTION TUBE (A) already containing the specimen.
3. Close the Myco TB™ SOLUTION TUBE (B) and discard it.
4. Close the REACTION TUBE (A).
5. Vortex the REACTION TUBE (A) for 10-30 seconds at 2000-2500 RPM based on the consistency of the original specimen. It is advisable to place the agitator under a laminar flow hood.⁽³⁾
6. Leave the REACTION TUBE (A) to rest at room temperature under a laminar flow hood for at least 3 minutes and up to a maximum of 10 minutes.
7. When the 3-minute rest time has elapsed, mix the REACTION TUBE (A) well with a vortexer and check the state of specimen fluidification.
NOTE: If the specimen is not fluidified, extend the rest time of the REACTION TUBE (A) up to a maximum of 10 minutes.
8. Open the REACTION TUBE (A) and place it on a test tube rack.
9. Open the NEUTRALIZING TUBE (C).
10. Pour a part of the content of the NEUTRALIZING TUBE (C) into the REACTION TUBE (A) and bring up to a volume of 50 mL.
11. Close the NEUTRALIZING TUBE (C) and discard it.


Fig. 3 - Decontamination and fluidification

Concentration and resuspension phase (Fig. 4): ^(2,3)

1. Securely screw the cap onto the REACTION TUBE (A) and mix the specimen 5-10 times by inversion.
2. Transfer the REACTION TUBE (A) to a centrifuge and centrifuge at 3300g for 5 minutes.
3. After centrifugation, leave to rest for 5 minutes under a laminar flow hood.
4. Open the REACTION TUBE (A).
5. Eliminate the supernatant leaving about 1-2mL of pellet on the bottom of the tube. A sterile pipette can be used or transfer the supernatant into an appropriate waste container.
6. Place the REACTION TUBE (A) in a test tube rack and open it.
7. Open the RESUSPENSION TUBE (D) containing 2mL of reagent.
8. Pour the entire reagent volume (2mL) into the REACTION TUBE (A).
9. Close the RESUSPENSION TUBE (D) and discard it.
10. Close the REACTION TUBE (A) and mix the thus treated specimen with a vortexer at 2000-2500 RPM for 5-10 seconds. It is advisable to place the agitator under a laminar flow hood. ⁽³⁾
11. The specimen is now ready to be used for the culture on solid or liquid media for the search for *Mycobacterium spp.*


Fig. 4 - Concentration and resuspension

WASTE DISPOSAL

Consult the product safety data sheet for disposal.

Used reagents and any other contaminated waste material must be disposed of following the procedures for infected or potentially infected products. The laboratory is responsible for handling the waste and waste fluids produced according to their nature and their degree of hazardousness, treating and disposing of them (or having them treated and disposed of) as established in any applicable regulation.

QUALITY CONTROL**Performance characteristics: Decontamination phase test**

All the operations are to be carried out under a laminar flow hood.

1. Prepare a 0.5 McFarland suspension.
2. Prepare a serial dilution such that the reagent under test is inoculated with a number of bacterial cells between 10^8 and 10^7 CFU/mL.
3. Open a Myco TB[™] package.
4. Place the REACTION TUBE (A) in a dedicated test tube rack.
5. Open the REACTION TUBE (A) under a laminar flow hood.
6. Using a micropipette, transfer 500µl of inoculum to the REACTION TUBE (A).
7. Bring up to a volume of 5mL by adding sterile PBS following the graduation on the REACTION TUBE (A).
8. Close the REACTION TUBE (A).
9. Vortex at 2000-2500 RPM for 10 seconds.
10. Open the REACTION TUBE (A), take 100µl and plate on a dedicated culture medium identifying the plate as time **T=0**.
11. Take the Myco TB[™] SOLUTION TUBE (B) and empty out the content into the REACTION TUBE (A).
12. Close the REACTION TUBE (A) and mix well with a vortexer at 2000-2500 RPM for 30 seconds.
13. Close and incubate the REACTION TUBE (A) at room temperature for 3 minutes.
14. Open the REACTION tube (A) and open the NEUTRALIZING TUBE (C).
15. Pour the content of the NEUTRALIZING TUBE (C) into the REACTION TUBE (A) filling it to a volume of 50mL reaching the maximum graduation notch on the REACTION TUBE (A).
16. Discard the NEUTRALIZING TUBE (C).
17. Close the REACTION TUBE (A) and insert it into the centrifuge.
18. Centrifuge the REACTION TUBE (A) at 3300g for 5 minutes.
19. Very carefully retrieve the REACTION TUBE (A) from the centrifuge and place it under a laminar flow hood.
20. Once under the hood, prepare a sterile container to collect the supernatant.
21. Very carefully transfer the supernatant from the REACTION TUBE (A) (approximately 48mL, leaving no more than 2mL on the bottom of the tube) into the container.
22. Close the REACTION TUBE (A) and the waste container and dispose of it according to the current biological waste protocol of the laboratory.
23. Open the RESUSPENSION TUBE (D) containing 2mL of PBS.

24. Pour the entire content of the RESUSPENSION TUBE (D) (2mL PBS) into the REACTION TUBE (A) containing the "pellet".
25. Close the REACTION TUBE (A) and vortex it at 2000-2500 RPM for 10 seconds.
26. Open the REACTION TUBE (A).
27. Plate 100µl on a dedicated culture medium and identify the plate. Make 5 plates of 100µl.
28. Incubate the plates according to the requirements of each individual microorganism under test.

The procedure from step 1 to 28 is performed on every individual strain teste

RESULTS

| STRAIN TESTER | CFU/PLATE at time zero | CFU/PLATE after 3 minutes of contact with the decontamination reagent |
|---|------------------------|---|
| <i>P. aeruginosa</i> ATCC 27853 | CONFLUENT GROWTH | ≤ 10 CFU/plate |
| <i>H. influenzae</i> ATCC 10211 | CONFLUENT GROWTH | ≤ 10 CFU/plate |
| <i>C. Albicans</i> ATCC 10231 | CONFLUENT GROWTH | ≤ 10 CFU/plate |
| <i>S. aureus</i> ATCC 43300 (MRSA strain) | CONFLUENT GROWTH | ≤ 10 CFU/plate |
| <i>K. pneumoniae</i> ATCC 700603 | CONFLUENT GROWTH | ≤ 10 CFU/plate |

TABLE OF SYMBOLS

See the table of symbols at the bottom of the instructions for use.

NOTES FOR THE PROFESSIONAL USER

In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

REVISION HISTORY

| Last Revision No.* | Release date | Changes made |
|--------------------|--------------|---|
| 02 | 03-2023 | Revision of IFU sections (first revision in IVDR) |

* Should you need earlier revisions, contact Copan Customer Service.

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Email: info@copangroup.com
 Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
 26055 Jefferson Avenue
 Murrieta, CA 92562, USA
 Tel: 951-696-6957
 Fax: 951-600-1832

E-mail: customerservice@copanusa.net
 Website: www.copanusa.com

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Copan Italia S.p.A.