# CHROMagar™ **B.cepacia**

**Instructions For Use**  
Available in several languages

**NT-EXT-124**

Version 2.0


Click below for:

**ENGLISH**

**FRANCAIS**

**ESPAÑOL**

**DEUTSCH**

# CHROMagar™ B.cepacia

Instructions For Use  
NT-EXT-124 V2.0 / 07-Nov-23

Chromogenic medium for detection and enumeration of *Burkholderia cepacia* complex.

## REFERENCES

Σ Pack Size

5000 mL

250 Tests  
of 20 mL

Ordering References

BK992

Base (B)

BK992  
Weight: 164,5 g

## INTENDED USE

CHROMagar™ B.cepacia is a selective chromogenic culture medium that is intended for use in the qualitative direct detection, differentiation and presumptive identification of *Burkholderia cepacia* complex (BCC) bacteria colonization in cystic fibrosis patients.

The test is performed by direct, or diluted, plating of sputum, expectoration, throat swabs, tracheal aspirates or bronchoalveolar lavages.

CHROMagar™ B.cepacia results can be interpreted after 36-72 h of aerobic incubation at 35-37 °C. A lack of growth or the absence of colonies on CHROMagar™ B.cepacia does not preclude the presence of BCC bacteria. CHROMagar™ B.cepacia is not intended to diagnose infection nor to guide nor to monitor treatment for infections. Further identification, susceptibility testing and epidemiological typing is needed on suspect colonies. The test has to be used in parallel of usual cultures antibiotic performed in cystic fibrosis patients.

CHROMagar™ B.cepacia can also be used in the qualitative and/or quantitative determination of BCC bacteria in the analyses of non-sterile, water-based pharmaceutical products, such as cream or sirup, in accordance with USP <60>. In this case, CHROMagar™ B.cepacia results can be interpreted after 36-72 h of aerobic incubation at 30-35 °C.

## COMPOSITION

The product is composed of a powder base.

| Product         | = | Base (B)  |
|-----------------|---|---|
| Total g/L       | = | 32.9 g/L  |
| Composition g/L | = | Agar 15.0<br>Peptones 11.0<br>Growth factors 6.6<br>Chromogenic and selective mix 0.3 |
| Aspect          | = | Powder Form   |
| STORAGE         | = | 15-30 °C  |
| FINAL MEDIUM pH | = | 6.2 +/- 0.2   |

Need some  
Technical Documents?

Available  
for download on  
[www.CHROMagar.com](http://www.CHROMagar.com)

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

## PREPARATION (Calculation for 1 L)

### Step 1 Preparation

- Suspend the powder base BK992 in the proportion of 32.9 g into 1 L of purified water.
- Stir until agar is well thickened.
- Heat and bring to boil (100 °C) while swirling or stirring regularly.  
**DO NOT HEAT TO MORE THAN 100 °C. DO NOT AUTOCLAVE AT 121 °C.**
- Cool in a water bath to 45/50 °C, keep stirring.

### Step 2 Pouring

- Pour into sterile Petri dishes.
- Let it solidify and dry.

### Storage

- Store in the dark before use.
- Prepared media plates can be kept for one day at room temperature.
- Plates can be stored for up to one month under refrigeration (2/8 °C) if properly prepared and protected from light and dehydration.

# CHROMagar™ B.cepacia

## SPECIMEN COLLECTION AND HANDLING

CHROMagar™ B.cepacia can be used with the following specimens:

**Clinical field:** Respiratory tract samples (bronchoalveolar lavages, sputum, nasopharyngeal aspirations or oropharyngeal swabs).

**Pharma industry:** Non-sterile products and purified water.

Use of transport devices approved for collection of such specimens is recommended.

## MATERIAL REQUIRED BUT NOT PROVIDED

Standard microbiological laboratory material for culture media preparation, control, streaking, incubation and waste disposal.

## INOCULATION

Related samples can be processed by direct streaking on the plate.

• If the agar plate has been refrigerated, allow to warm to room temperature before inoculation:

**Clinical field:** Streak sample. Incubate at 35 °C to 37 °C during 36 to 72 h in aerobic atmosphere.


**Pharma industry:** Streak sample or place the inoculated membrane (if using filtration\*) onto plate. Incubate at 30 °C to 35 °C during 36 to 72 h in aerobic atmosphere.

\* Use polycarbonate membranes for optimal performance.

## INTERPRETATION

| Microorganism       | Typical colony appearance     |
|---------------------|-------------------------------|
| <i>Burkholderia</i> | → green-blue +/- blue halo    |
| Yeasts and mold     | → Limited growth or inhibited |
| Gram (+) bacteria   | → Inhibited                   |
| Gram (-) bacteria   | → Limited growth or inhibited |

### Typical colony appearance

The pictures shown are not contractual.

## PERFORMANCE

|             | Analytical data *       | Clinical data**         |
|-------------|-------------------------|-------------------------|
|             | CHROMagar™<br>B.cepacia | CHROMagar™<br>B.cepacia |
| Sensitivity | ≈ 100 %                 | ≈ 100 %                 |
| Specificity | ≈ 100 %                 | 95 %                    |

\*Data obtained after 24-48 h incubation at 35 °C in aerobic conditions in the study «Evaluation of a new selective chromogenic agar for detection of *Burkholderia cepacia* complex». Katja Seme et al., poster. ASM Microbe 2022.

\*\* Data obtained after 24-48 h incubation at 35 °C in aerobic conditions in the study «Evaluation du milieu CHROMagar™ B.CEPACIA pour la détection de *B. cepacia* chez les patients atteints de mucoviscidose.» Julia MASOTTI et al., poster. RICA1 2021

## LIMITATIONS AND COMPLEMENTARY TESTS

- Most of the non-BCC strains are strongly inhibited or colorless, however a few non-fermenting strains can form green-blue colonies.
- Some of these strains detected were reported as of diagnostics interest.
- Final identification may require additional testing, including antibiotic susceptibility tests.
- Rare slow growing strains may require an additional 24 h incubation for a satisfactory result.
- The medium may brown slightly it does not affect the results
- Colorless inorganic filaments may appear on plate after 15-day storage at 2/8 °C (but they don't affect the reading)
- Not recommended when diagnosis requires a high negative predictive value (NPV).

## QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms.

Good preparation of the medium can be tested, isolating the following ATCC strains:

| Microorganism                     | Typical colony appearance |
|-----------------------------------|---------------------------|
| <i>B. cepacia</i> ATCC® 25416     | → blue +/- blue halo      |
| <i>S. epidermidis</i> ATCC® 12228 | → inhibited               |
| <i>K. pneumoniae</i> ATCC® 13883  | → inhibited               |
| <i>P. aeruginosa</i> ATCC® 27853  | → inhibited               |

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- This laboratory product should be used only by trained personnel (healthcare professional, etc). Wear appropriate protective clothing, gloves and eye/face protection and handle appropriately with procedures and good laboratory practices.
- Use of the medium may be difficult for people who have problems recognising colours.
- For a good microbial detection, collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.
- Culture media should not be used as manufacturing material or components.
- Do not ingest or inhale the product.
- Do not use the product after the expiry date.
- Do not use the product if it show any evidence of contamination or any sign of deterioration.
- Do not use the product if the packaging is damaged.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Inappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- Reading and interpretation should be performed using isolated colonies.
- Interpretation of the test results should be made taking into consideration the origin of the sample, the macroscopic and microscopic aspects of the colonies and, if necessary, the results of any other tests performed.

# CHROMagar™ B.cepacia

- Laboratory, chemical or biohazardous wastes must be handled and discarded in accordance with all local and national regulations.
- For hazard and precaution recommendations related to some chemical components in this medium, please refer to the pictogram(s) mentioned on the labels. The Safety Data Sheet (SDS) is available on [www.chromagar.com](http://www.chromagar.com)

## DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121 °C for at least 20 minutes.

## LITERATURE REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <https://www.chromagar.com/en/product/chromagar-b-cepacia-2/>

## IFU/LABEL INDEX

-  REF Catalogue reference
-  Consult instructions for use
-  Quantity of powder sufficient for X liters of media
-  Expiry date
-  Required storage temperature
-  Store away from humidity
-  Protect from light
-  Manufacturer

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