UTM[®]

Instructions for Use

CE **IVD**

English	3
Italiano	7
Español	12
Deutsch	17
Français	22
Português	27
БЪЛГАРСКИ	32
Česky	37
Dansk	42
ΕΛΛΗΝΙΚΑ	47
Eesti Keel	52
Hrvatski	57
Latviešu	61
Lietuvių K.	66
Magyar	71
Nederlands	76
Norsk	81
Polski	85
Română	90
Slovenčina	95
Slovenščina	100
Suomi	105
Svenska	109
Türkçe	114

Copan Universal Transport Medium (UTM-RT®) System Instructions for use

INTENDED USE

Copan Universal Transport Medium (UTM-RT®) System is intended for the collection and transport of clinical specimen containing viruses, chlamydiae, mycoplasma or ureaplasma from the collection site to the testing laboratory. UTM-RT® can be processed using standard clinical laboratory operating procedures for viral, chlamydial, mycoplasma and ureaplasma culture.

UTM-RT® is intended for the stabilization and transportation of an unprocessed upper respiratory clinical specimen suspected of containing respiratory viruses' nucleic acids. UTM-RT® is intended for use with compatible molecular assays.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of infections caused by viruses, chlamydiae, mycoplasmas or ureaplasmas involves collection and refrigerated transport of biological specimens.

The UTM-RT® consists of a Hanks' Balanced Salt Solution (HBSS) enriched with proteins and sugars with a neutral pH and pH indicator. The medium contains some antibiotics and antimycotics to inhibit overgrowth of bacteria and yeasts, maintain cellular integrity and encourage preservation of viruses and chlamydiae.

Copan UTM-RT® System medium is provided in labeled screw-cap tubes designed for transport of the clinical sample. Copan UTM-RT® System is also supplied as a sample collection kit that comprises a package which contains one screw-cap tube of UTM-RT® medium and a peel pouch incorporating sterile specimen collection swabs. A range of UTM-RT® sample collection kits are available which incorporate different types of shaft swabs and facilitate the collection of specimens from different sites of the patient as described below in the Product Description section. Once a swab sample is collected, it should be placed immediately into the transport tube where it comes into contact with transport medium. To maintain optimum microorganism viability, transport the specimens to the laboratory as soon as possible.

If specimen transported in UTM-RT® will be processed for viral, chlamydial, mycoplasma and ureaplasma culture, the specimen should be stored at 2–25°C and processed within 48 hours to maintain microorganism viability. If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder (this storage condition applies only in case UTM-RT® will be processed for culture as reported above).

Specimens transported in UTM-RT® for respiratory viral nucleic acids investigations should be processed within 96 hours when stored at 2–25°C. When the sample is processed both for microorganism culture and nucleic acids investigations, the specimen should be stored at 2–25°C and processed within 48 hours.

For use of the Copan UTM-RT® System medium with a molecular assay, consult the manufacturer's instructions for use for compatibility and specific target analyte processing directions.

REAGENTS

The UTM-RT® formulation includes proteins for virus stabilization¹⁴, antibiotics and antimycotics to prevent overgrowth of bacterial and fungal flora and a buffer solution to maintain a neutral pH along with a pH indicator.

Components
Sucrose
HBSS solution
Bovine serum albumin
Buffered solution
Gelatin
Amino acids
Antibiotics
Phenol Red

pH 7.3 ± 0.2 at 20–25°C

PRODUCT DESCRIPTION

UTM-RT® System is ready for use and requires no further preparation. Copan UTM-RT® System includes a screw-cap tube with round bottom containing 3mL of light orange-red transport medium.

UTM-RT® System tubes of transport medium are supplied alone (in bulk) or in a kit format. Each kit unit consists of a package containing: a pre-labeled screw cap tube and a peel pouch incorporating sterile specimen collection swabs. For more details on available configurations please refer to Table 1.

The collection swab applicators FLOQSwabs® provided with UTM-RT® have a molded breakpoint in the shaft of the applicator.

Table 1: Product description

REF	PRODUCT DESCRIPTION		PACK SIZE	PRODUCT DESIGNED FOR HAVING CAPTURE CAP FEATURE (FOR UTM KIT ONLY)
	TUBE	SWAB		
3U090N	3 mL of UTM-RT® medium in 16x100 mm screw-cap tube with round bottom.	NA	50 tubes per package 6 x 50 tubes per box	NO


3U091N01	3 mL of UTM-RT [®] medium in 16x100 mm screw-cap tube with round bottom.	One regular size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	NO
3U092N01	3 mL of UTM-RT [®] medium in 16x100 mm screw-cap tube with round bottom.	One minitip size applicator swab with flocked nylon fiber tip with breaking point	50 kits per package 6 x 50 kits per box	NO

Not all the product codes (REF) are available for sale in all Countries. Please contact Copan Customer care service for product codes availability for a specific Country.

Performance testing with Copan UTM-RT[®] System was conducted using laboratory strains. Performance testing was not conducted using human clinical specimens. Please refer to your internal procedures to choose the most appropriate device for the specific sampling site.

TECHNICAL NOTES:

a) Each of the swab(s) included in the kits listed in the Table 1 has a molded breaking point.

The molded breaking point is present and detectable as a **narrowing** of the stick.

Take time to visually locate the position of breaking point onto the stick when using the swab for collecting samples and for subsequent transfer in the transport tube as applicable.

b) A slight yellowing of the tip is a well-known phenomenon. This could be due to many factors: the type of raw material, the product sterilization treatment, and/or the product natural aging. Therefore, product yellowing is not necessarily indicative of product deterioration.

PRODUCT STORAGE

The product must be stored in its original packaging at a temperature between 2 and 25°C until the time of use prior to the indicated date of expiry (not to exceed 18 months from the date of manufacturing). Do not overheat or freeze prior to use.

REQUIRED MATERIALS BUT NOT PROVIDED

Materials suitable for isolation, differentiation and culture of viruses, chlamydiae, mycoplasmas and ureaplasmas and test systems/kits and reagents for molecular testing of viruses.

SPECIMEN COLLECTION, TRANSPORT AND PRESERVATION

Specimens for virus, chlamydia, mycoplasma or ureaplasma investigation should be collected and handled following published manuals and guidelines^{2,3,4,7,9,10,11}. To maintain optimum viability, transport the specimen to the laboratory as soon as possible.

If specimen transported in UTM will be processed for viral, chlamydial, mycoplasma and ureaplasma culture, the specimen should be stored at 2–25°C and processed within 48 hours to maintain microorganism viability. If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder (this storage condition applies only in case UTM-RT[®] will be processed for culture as reported above).

Specimens transported in UTM-RT[®] for respiratory viral nucleic acids investigations should be processed within 96 hours when stored at 2–25°C. When the sample is processed both for microorganism culture and nucleic acids investigations, the specimen should be stored at 2–25°C and processed within 48 hours.

Consult the manufacturer's instructions for use of the molecular assay platforms for specific processing directions.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations^{8,11,12}. Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

LIMITATIONS

1. Because calcium alginate swabs are toxic for many enveloped viruses⁵ and may interfere with immunofluorescence tests², they should not be used for specimen collection.
2. Wooden shaft swabs may contain toxins and formaldehydes^{2,6} and should not be used.
3. The product codes without beads are not suitable for use with mucous or particularly viscous specimens when homogenization of the specimen by vortexing in presence of the beads is required.
4. UTM-RT[®] kits are intended to be used with the medium tubes and swabs provided in the kit. The use of tubes of medium or swabs from any other source could affect the performance of the product.
5. Do not freeze specimens intended to be processed for nucleic acids investigations. Stability of nucleic acids when frozen in UTM-RT[®] has not been evaluated and validated.
6. UTM-RT[®] must be used by healthcare professional only.

WARNINGS AND PRECAUTIONS

1. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
2. For in vitro diagnostic use only.
3. UTM-RT[®] System is ready for use and requires no further preparation.
4. Do not use beyond the expiry date.
5. Do not use the UTM-RT[®] medium for premoistening or prewetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
6. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.⁷
7. Do not re-sterilize unused swabs.
8. Do not re-pack.
9. Not suitable to collect and transport microorganisms other than viruses, chlamydiae, mycoplasma and ureaplasma.
10. Do not ingest the medium.
11. Specimens for the search of viruses, chlamydiae, mycoplasmas and ureaplasmas must be collected and handled using personal protective equipment against biological risk according to published manuals and guidelines^{1,4,6,7,9,14}.
12. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms.
13. Do not use UTM-RT[®] if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the color of the medium has changed from light orange-red, (4) the swab pouch is open, or (5) there are other signs of deterioration.
14. The use of this product in combination with diagnostic kits or instruments must be validated by the user prior to use.
15. Do not bend or shape the swab before the collection of the specimen. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
16. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format and can be identified by the e-IFU indicator on the packaging label.


INSTRUCTIONS FOR USE

Proper collection of the specimen from the patient is a crucial aspect for successful isolation and identification of infectious organisms. Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute illness.


UTM-RT[®] in kit

1. Open the UTM-RT[®] kit package and remove the medium tube and the sterile swab in peel pouch.
2. Take the sterile swab out of its peel pouch and collect the clinical specimen; to prevent the risk of contamination, make sure that the swab tip comes into contact with the collection site only. **NOTE:** Do not bend the swab before the collection of the specimen. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft.
3. After collecting the specimen, insert the swab into the medium tube until the breakpoint is at the same level as the medium tube opening.
4. Bend the swab shaft at a 180 degrees angle to break it off at the breaking point. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
5. Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
6. Screw the cap back onto the test tube and hermetically seal it.
7. Identify the tube containing the specimen.
8. Send to the laboratory for immediate analysis.

Fig 1. Collection swab showing breakpoint and area for holding the applicator

Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens and care should be taken to avoid splashes and aerosols when breaking the swab stick into the tube of medium. During sample collection when handling the swab applicator, the operator must not touch the area below the breakpoint; that is the area from the breakpoint to the tip of the nylon flocked swab, as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.

UTM-RT[®] in bulk

1. Remove cap from tube, taking care not to spill or to contaminate the medium.
2. Place sample into the tube with UTM-RT[®] medium. If the sample is taken with a swab, insert and snap/cut off the swab(s) into the tube, pay attention to not splash the medium.
3. Screw the cap back onto the test tube and hermetically seal it.
4. Identify the tube containing the specimen.
5. Send to the laboratory for immediate analysis.

LABORATORY

If specimen transported in UTM-RT[®] will be processed for viral, chlamydial, mycoplasma and ureaplasma culture the specimen should be stored at 2–25°C and processed within 48 hours to maintain microorganism viability. If delivery and processing exceed 48 hours, specimens should be transported in dry ice and once in laboratory frozen at -70°C or colder (this storage condition applies only in case UTM-RT[®] will be processed for culture as reported above).

Specimens transported in UTM-RT[®] for respiratory viral nucleic acids investigations should be processed within 96 hours when stored at 2–25°C.

When the sample is processed both for microorganism culture and nucleic acids investigations, the specimen should be stored at 2–25°C and processed within 48 hours.

Consult the manufacturer's instructions for use of the molecular assay platforms for specific processing directions.

DISPOSAL

Waste must be disposed of in compliance with local legislation.

QUALITY CONTROL

The UTM-RT[®] medium is verified by USP membrane filtration test; cytotoxicity test for no toxic effects on cell culture; the ability to maintain the viability of viral, chlamydia, mycoplasma strains for up to 48 hours at 2–25°C in accordance with the methods described in CLSI M40-A2⁵; optimal pH range at 20–25°C between 7.10 and 7.50; the absence of amplifiable respiratory viruses' nucleic acids.

PERFORMANCE CHARACTERISTICS

The results obtained largely depend on proper and adequate specimen collection as well as the promptness with which the specimens are transported to the laboratory and analyzed.

Viability studies were performed using Copan UTM-RT[®] with a panel of representative strains of the various families supported by the UTM-RT[®]. The swabs that accompany each transport system were directly inoculated in triplicate with 100 µL of organism suspension using laboratory strains. Performance testing was not conducted using human clinical specimens. Subsequently, the swabs were inserted in the respective test tubes containing the transport medium and stored for 0 and 48 hours at 2–6°C and at controlled room temperature (20–25°C). At the time of processing, each swab was vortexed for 20 seconds and removed from its transport medium test tube, then, an aliquot of the suspension was inoculated onto the appropriate cell line (200 µL) or into the appropriate culture medium^{6,13}. All the cultures were processed using the standard laboratory culture technique^{6,13}. Organism viability was determined by fluorescent cell counting for viral and chlamydia strains and CFU counting for mycoplasma and ureaplasma strains. The acceptability limits for time zero and for 48 hours were defined in accordance with the regulations M40-A2⁵. Copan UTM-RT[®] System preserved the viability of all the organisms tested for 48 hours both at controlled room temperature and in the cold in the above-described test conditions. The organisms evaluated and the results obtained are given in the tables 2 and 3 below.

Table 2: Virus and chlamydia viability studies at room temperature and in the cold

Organism	ATCC [®] number	% of reduction of fluorescing infected cells after 48 hours time 2–6°C	% of reduction of fluorescing infected cells after 48 hours time 20–25°C
Herpes Simplex Virus Type 1	ATCC [®] VR-539	≤50%	≤50%
Herpes Simplex Virus Type 2	ATCC [®] VR-734	≤50%	≤50%
Respiratory Syncytial Virus	ATCC [®] VR-1580	≤50%	≤50%
Coxsackie B1 Virus	ATCC [®] VR-28	≤50%	≤50%
Chlamydia trachomatis	ATCC [®] VR-880	≤50%	≤50%
Influenza A	ATCC [®] VR-1679	≤50%	≤50%
Cytomegalovirus	ATCC [®] VR-977	≤50%	≤50%
Varicella-zoster virus	ATCC [®] VR-1367	≤50%	≤50%
Adenovirus	DHI 20-4740010	≤50%	≤50%
Parainfluenza virus type 3	DHI 20-4770010	≤50%	≤50%
Echovirus type 30	ATCC [®] VR-322	≤50%	≤50%
Chlamydia pneumoniae	ATCC [®] VR-1360	≤50%	≤50%

Table 3: Bacterial viability studies at room temperature and in the cold

Organism	ATCC® number	% of reduction of CFUs after 48 hours time 2–6°C	% of reduction of CFUs after 48 hours time 20–25°C
<i>Mycoplasma hominis</i>	ATCC® 23114	35%	25%
<i>Ureaplasma urealyticum</i>	ATCC® 27618	0%	6%
<i>Mycoplasma pneumoniae</i>	ATCC® 15531	50%	50%

Copan UTM-RT® was tested to verify the stabilization and transportation of respiratory viruses' nucleic acids using compatible molecular assay. One molecular assay was selected for testing: Cepheid Xpert® Xpress CoV-2/Flu/RSV plus assay, run on GeneXpert® Dx.

For each virus tested, a single viral strain stock form ATCC was diluted into a synthetic nasopharyngeal matrix and 100 µL of viral suspension were added to 3 mL UTM-RT®.

RSV, Flu A and Flu B target nucleic acids were amplified acceptably with stable ΔCt values relative to baseline after storage in UTM-RT® for up to 96 hours at both room temperature (22–28°C) and refrigerated temperature (2–8°C), as shown in the table below.

Virus	Beads	Lot age of UTM-RT® at testing	ΔCt at 96 hrs. ($T_{96} - T_0$) PASS criterion: $\Delta Ct < 3$			
			2–8°C	Result	22–28°C	Result
Flu A1	UTM-RT tube With beads or without beads	Freshly produced to up to 18 months	0–0.9	PASS	0.3–1.4	PASS
Flu A2			0–0.8	PASS	0.3–1.1	PASS
Flu B			-1–0.4	PASS	-0.8–1	PASS
RSV			-0.4–1	PASS	-0.1–1.1	PASS

CHANGE HISTORY

Last revision N.*	Date of issue	Introduced changes
00	07-2024	First issue

* Should you need earlier revisions, contact Copan Customer Service.

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BIBLIOGRAPHY

1. Gary W. Procop and Elmer W. Koneman, 2016. Color Atlas and Textbook of Diagnostic Microbiology, Seventh edition. Wolters Kluwer Health.
2. James H. Jorgensen, Michael A. Pfaller, Karen C. Carroll, Guido Funke, Marie Louise Landry, Sandra S. Richter, David W. Warnock, 2015. Manual of Clinical Microbiology, 11th Edition. ASM, Washington, DC.
3. James Versalovic, Karen C. Carroll, Guido Funke, James H. Jorgensen, Marie Louise Landry, David W. Warnock, 2011. Cumitech 15A. Laboratory Diagnosis of Viral Infections. ASM, Washington, DC.
4. Patricia Tille. 2014. Bailey & Scott's Diagnostic Microbiology, 13th Edition. Laboratory Medicine.
5. Clinical and Laboratory Standards Institute (CLSI), 2014. M40-A2 Quality Control of Microbiological Transport Systems; Approved Standard-Second Edition.
6. Clinical and Laboratory Standards Institute (CLSI), 2006. M41-A Viral Culture; Approved guidelines.
7. Wardford, A., M. Chernesky, and E. M. Peterson, 1999. Cumitech 19A, Laboratory Diagnosis of Chlamydia trachomatis Infections. ASM, Washington DC.
8. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
9. J. Michael Miller, Shelley A. Miller, 2017. A Guide to Specimen Management in Clinical Microbiology, Third Edition. ASM, Washington DC.
10. Centers for Disease Control and Prevention (CDC), 2016. Guide for Shipping Infectious Substances.
11. Centers for Disease Control and Prevention (CDC), 2009. Biosafety in Microbiological and Biomedical Laboratories 5th Edition.
12. World Health Organization 2015. Guidance on regulations for the Transport of Infectious Substances 2015 – 2016.
13. European Collection of Authenticated Cell Cultures (ECACC). Fundamental Techniques in Cell Culture Laboratory Handbook, 3rd edition (Protocol 9 – Detection of Mycoplasma by Culture Isolation).
14. S. Specter, R. L. Hodinka, S. A. Young. Clinical Virology Manual, fifth edition, 2016.

INDEX OF SYMBOLS / INDICE DEI SIMBOLI / ÍNDICE DE SÍMBOLOS / SYMBOLVERZEICHNIS / INDEX DES SYMBOLES / ÍNDICE DE SÍMBOLOS / ИНДЕКС НА СИМВОЛИТЕ / REJSTŘÍK SYMBOLŮ / INDEKS ÖVER SYMBOLER / ΕΡΕΤΗΡΙΟ ΣΥΜΒΟΛΩΝ / TINGMÄRKIDE KIRJELDUSED / KAZALO SIMBOLA / IZMANTOTIE SIMBOLI / SIMBOLIŲ RODYKLĖ / SZIMBÓLUMOK JEGYZÉKE / INDEX VAN SYMBOLEN / SYMBOLINDEKS / INDEKS SYMBOLI / INDEX SIMBOLURI / ZOZNAM SYMBOLOV / KAZALO SIMBOLOV / SYMBOLIHAKEMISTO / LISTA ÖVER SYMBOLER / SEMBOLLER DIZINI

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Email: info@copangroup.com
 Website: www.copangroup.com

North American Distributor:
 Copan Diagnostics Inc.
 26055 Jefferson Avenue
 Murrieta, CA 92562, USA
 Tel: 951-696-6957
 Fax: 951-600-1832

E-mail: customerservice@copanusa.net
 Website: www.copanusa.com

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Copan Italia S.p.A.