MASTAZYME™ DIPHTHERIA

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11. Literatur

1. Borcic, B et al. Immunity to diphtheria in the Croatian population in 1994. Lijec. Vjesn., 118: 227 (1996).
2. CDC. Diphtheria Epidemic – New Independent States of the Former Soviet Union. 1990-1994; MMWR 44: 177 (1995).
3. Department of Health and Human Services, Food and Drug Administration. Biological Products; Bacterial Vaccines and Toxoids; Implementation of Efficacy Review; Proposed Rule. Federal Register Vol 50 (240): 51002 (1985).
4. Frank, JW et al. Diphtheria overimmunization in children. Can. Med. Assoc. J. 141: 1241 (1989).
5. Gupta, RK et al. Diphtheria antitoxin levels in blood and plasma donors. J. Infect. Dis. 173: 1493 (1996).
6. Kjeldsen, K et al. Immunity against diphtheria and tetanus in the age group 30-70 years. Scand. J. Infect. Dis. 20: 177 (1988).
7. Lau, RC. Detection of diphtheria toxin antibodies in human sera by ELISA. J. Hyg. 96: 415 (1986).
8. Mueller, JH et al. Production of diphtheria toxin of high potency (100 Lf) on a reproducible medium. J. Immunol. 40: 21 (1941).
9. Melville-Smith, M et al. Estimation of diphtheria antitoxin in human sera by ELISA. J. Med. Microbiol. 25: 279 (1988).
10. Nosov, SD. Infectious Diseases of Childhood, 2. Auflage; Mir Publishers, Moskau, 61 (1984).
11. Recommendations of the Advisory Committee of Immunization Practices (ACIP). Diphtheria, Tetanus, and Pertussis: Recommendations for vaccine use and other preventive measures; MMWR 40: 10 (1991).
12. Rieger, J et al. Diphtheria immunity in the German population; Gesundheitswesen 56: 667 (1994).

Contents	Page
1. Intended Use	11
2. Introduction	11
3. Principle of the Test	11
4. Kit Contents	12
5. Materials Required but not Provided	13
6. Warnings and Precautions	13
7. Storage and Stability	13
8. Specimen Collection and Handling	14
9. Assay Procedure	14
10. Results and Interpretation	15
11. Assay Performance	16
12. References	16

1. Intended Use

MASTAZYME™ DIPHTHERIA has been designed for the detection and the quantification of specific IgG antibodies against the toxin of *Corynebacterium diphtheriae* in serum and plasma. Further applications in other body fluids are possible and can be provided on request.

This assay is intended for *in-vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

2. Introduction

Diphtheria is a bacterial infectious disease which appears predominantly during the childhood. The disease leads particularly to an inflammation of the pharynx, larynx and nasal mucosa. Additionally, bacterial toxins cause via long-distance effect damages of the heart, circulation and CNS. Only the toxigenic strains are pathogenic.

The etiologic agent is the *Corynebacterium diphtheriae*. These gram-positive bacteria prefer a microaerophil to anaerobe environment. Its pathogenicity is based on the secretion of an exotoxin that is circulating in the blood and effecting the heart muscle, kidneys and CNS. The diphtheria toxin will be produced by lysogenic strains.

Depending on the stage of disease the three types 'slight, middle and serious' can be distinguished. The natural source of infection is the sick individual whereas a carrier not absolutely shows symptoms. The infection is spread both through the aerial-droplet route and rarely by milk or smear infection. The appearance of Diphtheria shows a seasonal prevalence with the greatest incidence in winter.

Especially non-vaccinated children will be infected. The incubation time is depending on the number of invasive germs.

The place of infection is the mucosa of the respiratory tract where an acute local infection is developing. The secreted toxin leads to a superficial inflammation of the mucosa associated with the formation of a brown film (pseudo-membrane) upon it consisting of bacteria, necrotic epithelial cells, fibrin, red and white cells. From this local inflammation the toxin reaches other organs by using the blood and lymphatic circulation. Here it may cause severe damages. The grade of disease depends on the immunostate of the child. Usually, a limited Diphtheria arises whereas in case of an immunosuppression, a severe Diphtheria is observed. As a result of this course of disease the patient may die.

In most cases children will be vaccinated (e.g. DTP = Diphtheria-Tetanus-Pertussis) after the third month of life. The state of immunity can be monitored by determining the antitoxin IgG.

3. Principle of the Test

The principle of the test reaction can be described in four stages.

3.1 Serum incubation

Specific antibodies bind to the antigens on the solid phase to form a stable immune complex. After a 60 minutes incubation at room temperature the wells are washed with prediluted washing buffer to remove all non-reactive serum components.

3.2 Conjugate incubation

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich. After a 30 minutes incubation at room temperature the excess conjugate is removed by washing all wells with washing buffer.

3.3 Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and subsequently the colour development is stopped after 20 minutes incubation at room temperature by adding 0.5 M H₂SO₄ to the wells. The change in pH also causes the chromogen to change colour from blue to yellow.

3.4 Reading and interpretation

The intensity of the colour is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

4. Kit Contents

The kits contains sufficient reagents for 12 x 8 = 96 determinations. The strips and solutions have to be stored at 2–8 °C. The expiry date is mentioned on the labels.

12	Microtiter strips	single strips each with 8 break-apart wells coated with the purified toxin of <i>Corynebacterium diphtheriae</i>
1 x	Frame holder	
5 x 2 mL	Calibrators 1–5	human serum containing antibodies against diphtheria toxin (concentrations listed below) diluted in PBS and stabilised with 0.01 % methylisothiazolone and 0.01 % bromonitrodioxane as preservatives, ready to use

		IgG
Calibrator 1	Concentration (IU/mL)	0
Calibrator 2		0.01
Calibrator 3		0.1
Calibrator 4		0.5
Calibrator 5		1.0

IgG is calibrated against WHO reference preparation 00/456. The concentration is stated in IU/mL

1 x 60mL	Serum diluent	PBS/BSA buffer solution, contains < 0.1 % sodium azide as preservative, ready to use
1 x 12 mL	Enzyme conjugate solution	HRP-labelled goat anti-human-IgG, ready to use
1 x 12 mL	TMB substrate	3,3',5,5' Tetramethylbenzidine, ready to use
1 x 12 mL	Stopping solution	0.5 M sulfuric acid, ready to use
1 x 60mL	Washing buffer 10 x concentrated	PBS/Tween buffer solution 10 x concentrated to be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 minutes to avoid any crystals
2 x	Plate sealers	to cover microtiter strips during incubation
1 x	Plastic bag	re-sealable for dry storage of non-used strips

5. Materials Required but not Provided

- 5 µL-, 100 µL- and 500-µL micro- and multichannel pipets
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm)
- Microtiter Plate Washer (in case of manual washing: wash bottle)
- Reagent tubes for the serum dilution
- Measuring cylinder
- Distilled water or water of higher quality

6. Warning and Precautions

- For *in-vitro* diagnostic use only! Do not ingest or swallow! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera and plasma or buffers based upon have been tested to HBsAg, HIV and HCV respectively with generally accepted methods and were found negative. Nevertheless, precautions like the use of latex gloves have to be taken.
- Serum and reagent spills have to be wiped off with a disinfecting solution (e.g. sodium hypochlorite 5 %) and have to be disposed of properly.
- All reagents have to be brought to room temperature (18–24 °C) before performing the test.
- Before pipetting all reagents should be mixed thoroughly by gentle tilting or swinging. Vigorous shaking with formation of foam should be avoided.
- It is important to pipet with constant intervals so that all the wells of the microtiter plate have the same conditions.
- When removing reagents out of the bottles care has to be taken that the stoppers are not contaminated. Further a possible mix-up has to be avoided. The content of the bottles is usually sensitive to oxidation so that they should be opened only for a short time.
- In order to avoid a carry-over or a cross-contamination separate disposable pipet tips have to be used.
- No reagents from different kit lots should be used and they should not be mixed with one another.
- All reagents have to be used within shelf life.
- In accordance with a Good Laboratory Practice (GLP) or following ISO 9001 all laboratory devices employed should be regularly checked regarding the accuracy and precision. This refers e.g. to microliter pipets and washing or reading (ELISA Reader) instrumentation.
- The contact of certain reagents especially the stopping solution and the substrate with skin, eye and mucosa has to be avoided because possible irritations and acid burns could arise and there exists a danger of intoxication.

7. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted washing buffer is stable for up to 4 weeks when stored at 2–8 °C.

The opened kit should be used within three months.

8. Specimen Collection and Handling

Both serum or plasma (EDTA, heparin) can be used for the determination. Serum is separated from the blood which is aseptically drawn by venipuncture after clotting and centrifugation. The serum or plasma samples can be stored at 2–8 °C for up to 3 days. They should be kept at -20°C for a longer storage. The samples should not be frozen and thawed repeatedly. Lipemic, hemolytic or bacterially contaminated samples can cause false positive or false negative results.

Patient sera must be prediluted 1:101 in serum diluent (e.g. 5 µL serum + 500 µL serum diluent) prior to testing.

Samples containing concentrations higher than the highest calibrator have to be diluted further with serum diluent.

In case of interference with rheumatic factors, serum preabsorption with RF absorbent (MASTSORB™ Order Code: 651003) is recommended. **Do not absorb the calibrators.**

9. Assay Procedure

9.1. Preparation of Reagents

Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

Washing buffer: Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 60 mL buffer concentrate + 540 mL distilled water). Mix thoroughly.

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.
- A standard curve should be established with each assay.
- Put the unused microtiter strips back in the plastic bag and store them dry at 2–8 °C.

9.2. Assay Steps

Prepare a sufficient amount of microtiter wells for calibrators, controls and samples.

Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance.

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use calibrators into the appropriate wells.
2. Cover plate with the enclosed plate sealing foil and incubate at room temperature for 60 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Cover plate with plate sealing foil and incubate for 30 minutes at room temperature.
6. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.

7. Dispense 100 µL of TMB substrate into each well.
8. Cover plate with the plate sealing foil and incubate for 20 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After thorough mixing and wiping the bottom of the plate, read the optical density at 450 nm and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended.

The developed colour is stable for at least 60 minutes. Read optical densities during this time.

10. Results and Interpretation

Example:

	OD 450 nm	corrected OD	Mean OD Value
Blank	0.016		
Calibrator 1	0.037/0.045	0.021/0.029	0.025
Calibrator 2	0.072/0.063	0.056/0.047	0.052
Calibrator 3	0.376/0.386	0.360/0.370	0.365
Calibrator 4	1.480/1.461	1.464/1.445	1.455
Calibrator 5	2.149/2.126	2.133/2.110	2.122

The table above should be considered as an example which was achieved under arbitrary temperature and environmental conditions. **These data do NOT describe reference values which have to be found in other laboratories in the same way!**

10.1 Quantitative Calculation

The ready to use calibrators of the MASTAZYME™ DIPHTERIA are defined and values expressed in International Units (IU/mL). This gives access to an exact and reproducible quantification and in consequence patient antibody titer monitoring is possible. Concentration values for calibrators are printed on the labels of the vials.

A standard curve is plotted by entering the mean absorbance value of the calibrators on the Y-axis and the corresponding concentration on the X-axis using graph paper. The concentration of the patient samples can then be read directly from the graph.

The calculation of the result can be performed using a computer and a suitable software program.

Interpretation:

< 0.01 IU/mL	Basic vaccination recommended
0.01–0.1 IU/mL	Booster vaccination recommended
> 0.1 IU/mL	Good Immunity

11. Assay Performance

The assay characteristics of the MASTAZYME™ DIPHTHERIA have been established and assessed according the European IVD directive. Detailed validation data can be provided on special request.

12. References

1. Borcic, B et al. Immunity to diphtheria in the Croatian population in 1994. *Lijec. Vjesn.*, 118: 227 (1996).
2. CDC. Diphtheria Epidemic – New Independent States of the Former Soviet Union. 1990-1994; *MMWR* 44: 177 (1995).
3. Department of Health and Human Services, Food and Drug Administration. Biological Products; Bacterial Vaccines and Toxoids; Implementation of Efficacy Review; Proposed Rule. *Federal Register* Vol 50 (240): 51002 (1985).
4. Frank, JW et al. Diphtheria overimmunization in children. *Can. Med. Assoc. J.* 141: 1241 (1989).
5. Gupta, RK et al. Diphtheria antitoxin levels in blood and plasma donors. *J. Infect. Dis.* 173: 1493 (1996).
6. Kjeldsen, K et al. Immunity against diphtheria and tetanus in the age group 30-70 years. *Scand. J. Infect. Dis.* 20: 177 (1988).
7. Lau, RC. Detection of diphtheria toxin antibodies in human sera by ELISA. *J. Hyg.* 96: 415 (1986).
8. Mueller, JH et al. Production of diphtheria toxin of high potency (100 Lf) on a reproducible medium. *J. Immunol.* 40: 21 (1941).
9. Melville-Smith, M et al. Estimation of diphtheria antitoxin in human sera by ELISA. *J. Med. Microbiol.* 25: 279 (1988).
10. Nosov, SD. *Infectious Diseases of Childhood*, 2. Auflage; Mir Publishers, Moskau, 61 (1984).
11. Recommendations of the Advisory Committee of Immunization Practices (ACIP). Diphtheria, Tetanus, and Pertussis: Recommendations for vaccine use and other preventive measures; *MMWR* 40: 10 (1991).
12. Rieger, J et al. Diphtheria immunity in the German population; *Gesundheitswesen* 56: 667 (1994).

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Laboratoriumspräparate GmbH
MAST GROUP Ltd.