MASTAZYME™ FSME IgG

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Contents	Page
1. Intended use	10
2. Principle of the Test	10
3. Kit Contents	10
4. Materials Required but not Provided	11
5. Warnings and Precautions	11
6. Storage and Stability	12
7. Specimen Collection and Handling	12
8. Assay Procedure	12
9. Results and Interpretation	13
10. Assay Performance	14
11. References	14

1. Intended Use

MASTAZYME™ FSME IgG (TBE Virus = Tick-Borne Disease Virus) Antibody ELISA has been designed for the detection and the quantification of specific IgG antibodies respectively against TBE Virus in serum and plasma. Further applications in other body fluids are possible and can be provided on request.

This assay is intended for *in vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

2. Principle of the Test

The ELISA test can be described in four stages:

1.1 Serum incubation

Specific antibodies bind to antigen on the solid phase to form a stable immune complex. After incubation for 60 minutes at room temperature the wells are washed with a prediluted wash buffer to remove all non-reactive serum components.

1.2 Conjugate incubation

The anti-human-IgG horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG antibodies on the solid phase antigen to form a stable sandwich. After incubating for 30 minutes at room temperature the excess conjugate is removed by washing all wells with wash buffer.

1.3 Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and consequent colour development is stopped after a 20 minutes incubation at room temperature by adding 0.5 N H₂SO₄ to the wells. The shift in pH also causes the chromogen to change colour from blue to yellow.

1.4 Reading and interpretation

The colour intensity is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

The results can be read from a calibration curve or with an electronic graphing package using a 4 parameter sigmoidal curve.

3. Kit Contents

The kits contains sufficient reagents for 12 x 8 = 96 determinations. Kit components should be stored at 2–8 °C. The expiry date is shown on the individual labels.

12	Microtiter strips	Single strips each with 8 break-apart wells coated with the antigen of TBE Virus respectively.
1 x	Frame holder	
4 x 1,5 mL	Calibrators 1–4	Stabilized human plasma containing antibodies against the above antigens and pre-diluted in buffer. Ready to use (RTU).

	Anti-FSME IgG concentration [VU/mL]
Calibrator 1	1
Calibrator 2 (Cut-off)	100
Calibrator 3	400
Calibrator 4	1600
Positive Control	800
Negative Control	50

1 x 1,5 mL	Negative Control	Stabilized human plasma, RTU, can be integrated into the calibration curve.
1 x 1,5 mL	Positive Control	Stabilized human plasma, RTU, can be integrated into the calibration curve.
2 x 60 mL	Sample diluent	PBS solution with preservative, ready to use.
1 x 12 mL	Enzyme conjugate solution	HRP-labelled goat anti-human-IgG, ready to use.
1 x 12 mL	TMB substrate	3,3',5,5' Tetramethylbenzidine, ready to use.
1 x 12 mL	Stopping solution	0.5 N sulfuric acid, ready to use.
2 x 50 mL	Wash buffer 10 x concentrated	PBS/Tween buffer solution 10x concentrated. To be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 min to dissolve any crystals

4. Materials Required but not Provided

- 5 µL-, 100 µL- and 500 µL micropipettes and a multichannel pipette (optional).
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm)
- Microtiter plate washer (if washing manually: wash bottle)
- Reagent tubes for serum dilution
- Measuring cylinder
- Distilled water or water of similar or better quality

5. Warning and Precautions

- For *in vitro* diagnostic use only! Do not ingest! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera, plasma and buffers containing human biological material were found to be negative for Hepatitis B, C and HIV. Nevertheless, precautions such as the use of latex gloves must be taken.
- Serum and reagent spills should be cleaned with a disinfecting solution (e.g. sodium hypochlorite 5 %) and the refuse disposed of appropriately.
- All reagents should be brought to room temperature (18 to 24 °C) before starting the test procedure.
- Before pipetting, all reagents should be mixed thoroughly by gentle agitation. Vigorous shaking leading to formation of foam should be avoided.

- It is important to pipette with constant time intervals, so that all the wells of the microtiter plate have the same reaction conditions.
- When pipetting reagents out of the bottles care has to be taken that the stoppers are not contaminated. Furthermore, disposable pipette tips are strongly recommended to avoid cross-contamination. The contents of the bottles are usually sensitive to oxidation so should be opened only for a short time.
- No reagents from different kit batches should be used and reagents should not be mixed with one another.
- All reagents should be used within the listed shelf life.
- In accordance with Good Laboratory Practice (GLP) or following ISO 9001 all laboratory equipment should be checked regularly regarding accuracy and precision. This refers e.g. to microliter pipets and washing or measurement (ELISA Reader) instrumentation.
- Certain reagents (especially the stop solution and substrate) are irritants- avoid contact with skin, eyes and mucosal membranes. In case of accident, rinse with water and seek medical attention. Clean all equipment after use to avoid secondary contact incidents.

6. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted wash buffer is stable for up to 4 weeks when stored at 2–8 °C. However ensure that this is at room temperature before testing.

Kits should be used within three months of opening.

7. Specimen Collection and Handling

Both serum and plasma (EDTA, heparin) can be used for testing. Samples can be stored at 2–8 °C for up to 3 days, although for longer storage periods they should be aliquoted and kept at -20 °C. Repeated freezing and thawing is contra-indicated. Lipaemic, haemolytic, icteric and bacterially contaminated samples can cause false positive or false negative results.

Patient sera should be prediluted 1:101 in sample diluent (e.g. 5 µL serum + 500 µL sample diluent) prior to testing.

8. Assay Procedure

8.1. Preparation of Reagents

Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.

Washing buffer: Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 50 mL buffer concentrate + 450 mL distilled water). Mix thoroughly.

- Strictly follow the instructions for reliable test performance. The user takes sole responsibility for any modifications to the test procedure.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.

- A standard calibration curve should be established with each assay.
- Replace any unused microtiter strips in the resealable aluminium bag provided and store under dry conditions at 2–8 °C.

8.2. Assay Steps

Prepare a sufficient amount of microtiter wells for the calibrators, controls and samples.

Note: Other incubation conditions are possible. However in case of modifications to the recommended test procedure (e.g. an incubation temperature of 37 °C instead of RT) assay performance should be validated by the user.

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use calibrators and controls into the appropriate wells.
2. Incubate the test plate at room temperature for 60 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted wash buffer. Afterwards remove residual wash buffer by gentle tapping of the microtiter plate onto a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Incubate the test plate for 30 minutes at room temperature.
6. Discard the contents of the microwells and wash 3 times with 300 µL of diluted wash buffer. Afterwards remove residual wash buffer by gentle tapping of the microtiter plate onto a paper towel.
7. Pipette 100 µL of TMB substrate into each well.
8. Incubate for 20 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After thorough mixing and wiping the bottom of the plate, read the optical density at 450 nm and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended. Concentrations can be plotted with an electronic graphing package or by hand against the calibration curve.

The developed colour is stable for at least 30 minutes. Read optical densities during this time.

9. Results and Interpretation

Example for IgG calibrator curve

	VU/mL*	OD 450 nm
Cal. 1	1	0.009
Cal. 2 (Cut-off)	100	0.835
Cal. 3	400	1.821
Cal. 4	1600	2.454
Positive Control	800	2.154
Negative Control	50	0.483

* VU = Vienna Units

The table above should be considered as an example which was achieved under arbitrary temperature and environmental conditions. **These data do NOT describe reference values which have to be found in other laboratories in the same way!**

Qualitative

The calculated OD values for patient sera as mentioned above are compared with the value for the cut-off calibrator. If the value of the sample is higher, then it should be read as positive.

A value below the cut-off calibrator should be read as negative. It seems reasonable to define a range of $\pm 10\%$ around the OD value of the cut-off as a grey zone. It is recommended to repeat results laying within the grey zone using the same serum or a new sample of the same patient, taken after 2–4 weeks. Both samples should be measured in parallel in the same run.

Quantitative

The ready to use calibrators of MASTAZYME FSME IgG are defined and values are expressed in Vienna Units (VU/mL). This gives access to an exact and reproducible quantification and in consequence patient antibody titer monitoring is possible. Concentration values for calibrators are printed on the labels of the vials.

The calculation of the result can be performed using a computer and a suitable software program.

10. Assay Performance

The assay characteristics of the MASTAZYME™ FSME IgG ELISA have been established and assessed according the European IVD directive. Detailed validation data can be provided on special request.

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**Mast Diagnostica GmbH
**Mast Group Ltd.