# MASTAZYME™ MASERN / MEASLES

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## 11. Literatur

1. Altintas DU, Evliyaoglu N, Kilinc B; Sen'an DI, Gunesser S. The modification in measles vaccination age as a consequence of the earlier decline of transplacentally transferred anti-measles antibodies in Turkish infants. *Eur J Epidemiol*, 12(6): 647-8 (1996).
2. Bayas JM, Vilella A, Vidal J, Nebot X, Carbo JM, Navarro G, Prat A, Asenjo MA, Salleras L.: Susceptibility to measles, rubella and parotitis in young adults. *Med Clin (Barc)*, 106(15): 561-4 (1996).
3. Bouche F, Ammerlaan W, Berthet F, Houard S, Schneider F, Muller CP. Immunosorbent assay based on recombinant hemagglutinin protein produced in a high-efficiency mammalian expression system for surveillance of measles immunity. *J Clin Microbiol*, 36(3): 721-6 (1998).
4. Chiu HH, Lee CY, Chih TW, Lee PI, Chang LY, Lin YJ, Hsu CM, Huang LM. Seroepidemiological study of measles after the 1992 nationwide MMR revaccination program in Taiwan. *J Med Virol*, 51(1): 32-5 (1997).
5. Dagan R, Slater PE, Duvdevani P, Golubev N, Mendelson E. Decay of maternally derived measles antibody in a highly vaccinated population in southern Israel. *Pediatr Infect Dis J*. 14(11): 965-9 (1995).
6. De Souza VA, Pannuti CS, Sumita LM, de Andrade Junior HF: Enzyme-linked immunosorbent assay-IgG antibody avidity test for single sample serologic evaluation of measles vaccines. *J Med Virol*, 52(3): 275-9 (1997).
7. Duvdevani P, Varsano N, Slepon R, Lerman Y, Shohat T, Mendelson E: Determination of immunity to measles virus in young adults: comparative evaluation of a commercial enzyme immunoassay and the hemagglutination inhibition techniques. *Clin Diagn Virol*, 7(1): 1-6 (1996).
8. Johnson CE, Kumar ML, Whitwell JK, Staehle BO, Rome LP, Dinakar C, Hurni W, Nalin DR: Antibody persistence after primary measles-mumps-rubella vaccine and response to a second dose given at four to six vs. eleven to thirteen years. *Pediatr Infect Dis*, 15(8): 687-92 (1996).
9. King SM, Saunders EF, Petric M, Gold R: Response to measles, mumps and rubella vaccine in paediatric bone marrow transplant recipients. *Bone Marrow Transplant*, 17(4): 633-6 (1996).
10. Matter L, Germann D, Bally F, Schopfer K: Age-stratified seroprevalence of measles, mumps and rubella (MMR) virus infections in Switzerland after the introduction of MMR mass vaccination. *Eur J Epidemiol*, 13(1): 61-6 (1997).
11. Mendelson E, Duvdevani P, Varsano N, Lerman Y, Slepon R, Dagan R, Cohen D, Danon Y, Shohat T: Measles immunity and response to revaccination of a young adult population in Israel. *J Med Virol*, 50(3): 249-53 (1996).
12. Metintas S, Etiz S, Akgun Y, Kalyoncu C, Sariboyaci MA, Isikli B: A serological survey of measles vaccine in a rural region of Eskisehir in Turkey. *Public Healt*, 111(6): 373-6 (1997).
13. Narita M, Yamada S, Matsuzono Y, Itakura O, Togashi T, Kikuta H: Measles virus-specific immunoglobulin G subclass response in serum and cerebrospinal fluid. *Clin Diagn Virol*, 8(3): 233-9 (1997).
14. Nates S, Rey G, Giordano M, Medeot S, Depetris A, Boshell J, de Wolff CD: Immunoglobulin M antibody response to measles virus following natural virus infection, primary vaccination and reexposure to the virus. *Viral Immunol*, 10(3): 165-73 (1997).
15. Pison Garces FJ, Galbe Sanchez-Ventura J, Arcauz Eguren P, Aguirre y Daban C, Mengual Gil J, Larrad Mur L: Immunity to measles, mumps and rubella in children vaccinated with triple viral vaccine. *Aten Primaria*, 15(4): 235-7 (1995).
16. Vardas E, Andre M, Kreis S, Zachariades N, Schoub B: Measles virus transmission from infected to immune individuals--implications for measles control and elimination. *S Afr Med J*, 87(12): 1709 (1997).

<b>Contents</b>	<b>Page</b>
1. Intended Use	11
2. Introduction	11
3. Principle of the Test	11
4. Kit Contents	12
5. Materials Required but not Provided	12
6. Warnings and Precautions	13
7. Storage and Stability	13
8. Specimen Collection and Handling	13
9. Assay Procedure	14
10. Results and Interpretation	15
11. Assay Performance	15
12. References	16

## 1. Intended Use

The MASTAZYME™ MEASLES (measles virus) has been designed for the detection and the quantification of specific IgG / IgM / IgA antibodies respectively against measles virus in serum and plasma. Further applications in other body fluids are possible and can be provided on request.

This assay is intended for *in-vitro* diagnostic use only.

All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

## 2. Introduction

Measles is a highly contagious viral disease characterised by a clinically distinct prodrome of fever, coryza, conjunctivitis, cough and a pathognomic exanthem (Koplik's spots). The disease is the result of infection with the Measles Virus, genus Morbillivirus of the family Paramyxoviridae.

Ten to twelve days after infection, the most prominent and characteristic prodromal symptoms appear: Coryza, a persistent barking cough, keratoconjunctivitis often with photophobia and fever. Generally lymphadenopathy and splenomegaly are also frequent. During this period, Koplik's spots appear on the bucal mucosa that rapidly spread involving the entire mucous membrane. These spots are usually gone by the time the skin rash reaches its peak. The rash of Measles appears after a 3 to 5

days prodrome some 14 days after exposure. The rash quickly becomes maculopapular and spreads rapidly over the face, neck, trunk and extremities during the next three days. At its height, the eruption has generally deepened to a redish purple and may be associated with edema of the skin. Complications are: Otitis media, pneumonia and encephalitis. Measles have a more severe expression in younger or undernourished children with a higher incidence of hemorrhage Measles with 5 to 10 % of lethal cases.

In people who have been vaccinated with inactive virus (before 1968) the infection can have severe manifestations as: Pneumonia, peripheral edema, pleural effusion and atypical rash. Measles are one of the most contagious infectious diseases. The virus spreads through droplets emanating from the respiratory tract of infected persons or by direct contact. The incidence of Measles has declined since the introduction of vaccination programs

## 3. Principle of the Test

The principle of the test reaction can be described in four stages.

### 3.1 Serum incubation

Specific antibodies bind to the antigens on the solid phase to form a stable immune complex. After a 60 minutes incubation at room temperature the wells are washed with prediluted wash buffer to remove all non-reactive serum components.

### 3.2 Conjugate incubation

The anti-human-IgG /-IgM /-IgA horseradish peroxidase conjugate is added to all wells. The conjugate binds to IgG / IgM / IgA antibodies on the solid phase antigen to form a stable sandwich. After a 30 minutes incubation at room temperature the excess conjugate is removed by washing all wells with washing buffer.

### 3.3 Substrate reaction and stopping

The TMB substrate is dispensed into each well and the peroxidase enzyme/substrate reaction forms a stable blue chromogen. The reaction and subsequently the colour development is stopped after 20 minutes incubation at room temperature by adding 0.5 N H<sub>2</sub>SO<sub>4</sub> to the wells. The change in pH also causes the chromogen to change colour from blue to yellow.

### 3.4 Reading and interpretation

The intensity of the colour is read in a microtiter plate reader at 450 nm (recommended reference wavelength for bichromatic measurement: 600–690 nm). The intensity of the colour (OD) is directly proportional to the concentration of the specific antibody in the patient sample.

#### 4. Kit Contents

The kits contains sufficient reagents for 12 x 8 = 96 determinations. The strips and solutions have to be stored at 2–8 °C. The expiry date is mentioned on the labels.

12	Microtiter strips	single strips each with 8 break-apart wells coated with antigen of measles virus
1 x	Frame holder	
4 x 2 mL	Calibrators 1–4	human serum containing antibodies against measles virus (concentrations listed below) diluted in PBS and stabilised with 0.01 % methylisothiazolone and 0.01 % bromonitrodioxane as preservatives, ready to use

		IgG	IgM	IgA
Cal. 1 (negative)	Concentration (U/mL)	1	1	1
Cal. 2 (Cut-off)		10	10	10
Cal. 3 (weak positive)		40	35	30
Cal. 4 (positive)		250	150	150

1 x 60mL	Serum diluent	PBS/BSA buffer solution, contains < 0.1 % sodium azide as preservative, ready to use
1 x 12 mL	Enzyme conjugate solution	HRP-labelled goat anti-human-IgG /-IgM /-IgA, ready to use
1 x 12 mL	TMB substrate	3,3',5,5' Tetramethylbenzidine, ready to use
1 x 12 mL	Stopping solution	0.5 N sulfuric acid, ready to use
1 x 60mL	Washing buffer 10 x concentrated	PBS/Tween buffer solution 10 x concentrated to be diluted 1:10 prior to use; the concentrate should be warmed up to 37 °C for 15 minutes to avoid any crystals
2 x	Plate sealers	to cover microtiter strips during incubation
1 x	Plastic bag	re-sealable for dry storage of non-used strips

#### 5. Materials Required but not Provided

- 5 µL-, 100 µL- and 500 µL micro- and multichannel pipets
- Microtiter plate reader with a 450 nm filter (reference filter 600–690 nm)
- Microtiter Plate Washer (in case of manual washing: wash bottle)
- Reagent tubes for the serum dilution
- Measuring cylinder
- Distilled water or water of higher quality

## 6. Warning and Precautions

- For *in-vitro* diagnostic use only! Do not ingest or swallow! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
- All sera and plasma or buffers based upon have been tested to HBsAg, HIV and HCV respectively with generally accepted methods and were found negative. Nevertheless, precautions like the use of latex gloves have to be taken.
- Serum and reagent spills have to be wiped off with a disinfecting solution (e.g. sodium hypochlorite 5 %) and have to be disposed of properly.
- All reagents have to be brought to room temperature (18–24 °C) before performing the test.
- Before pipetting all reagents should be mixed thoroughly by gentle tilting or swinging. Vigorous shaking with formation of foam should be avoided.
- It is important to pipet with constant intervals so that all the wells of the microtiter plate have the same conditions.
- When removing reagents out of the bottles care has to be taken that the stoppers are not contaminated. Further a possible mix-up has to be avoided. The content of the bottles is usually sensitive to oxidation so that they should be opened only for a short time.
- In order to avoid a carry-over or a cross-contamination separate disposable pipet tips have to be used.
- No reagents from different kit lots should be used and they should not be mixed with one another.
- All reagents have to be used within shelf life.
- In accordance with a Good Laboratory Practice (GLP) or following ISO 9001 all laboratory devices employed should be regularly checked regarding the accuracy and precision. This refers e.g. to microliter pipets and washing or reading (ELISA Reader) instrumentation.
- The contact of certain reagents especially the stopping solution and the substrate with skin, eye and mucosa has to be avoided because possible irritations and acid burns could arise and there exists a danger of intoxication.

## 7. Storage and Stability

Store all reagents at 2–8 °C.

The expiry date of each reagent is printed on the individual labels. Do not use any reagents after the expiry date has been exceeded.

The diluted washing buffer is stable for up to 4 weeks when stored at 2–8 °C.

The opened kit should be used within three months.

## 8. Specimen Collection and Handling

Both serum or plasma (EDTA, heparin) can be used for the determination. Serum is separated from the blood which is aseptically drawn by venipuncture after clotting and centrifugation. The serum or plasma samples can be stored at 2–8 °C for up to 3 days. They should be kept at -20°C for a longer storage. The samples should not be frozen and thawed repeatedly. Lipemic, hemolytic or bacterially contaminated samples can cause false positive or false negative results.

Patient sera must be prediluted 1:101 in serum diluent (e.g. 5 µL serum + 500 µL serum diluent) prior to testing.

Samples containing concentrations higher than the highest calibrator have to be diluted further with serum diluent.

In case of interference with rheumatic factors, serum preabsorption with RF absorbent (MASTSORB™ Order Code: 651003) is recommended. **Do not absorb the calibrators.**

## 9. Assay Procedure

### 9.1. Preparation of Reagents

**Allow all kit components and specimens to reach room temperature (RT, 18–24 °C) prior to use and mix well.**

**Washing buffer:** Dissolve any crystals which may be in the bottle by warming to 37 °C and then mix well.

Dilute the concentrated washing buffer 1:10 with distilled water (e.g. 60 mL buffer concentrate + 540 mL distilled water). Mix thoroughly.

- Strictly follow the instructions for reliable test performance. Any changes or modifications are within the responsibility of the user.
- All reagents and samples must be brought to room temperature before use, but should not be left at this temperature for longer than necessary.
- A standard curve should be established with each assay.
- Put the unused microtiter strips back in the plastic bag and store them dry at 2–8 °C.

### 9.2. Assay Steps

Prepare a sufficient amount of microtiter wells for calibrators, controls and samples.

**Note: Other incubation conditions might be possible. In case of modifications of the recommended test procedure (e.g. incubation temperature 37 °C instead of RT) the user has to validate assay performance.**

1. Pipette 100 µL each of the diluted (1:101) samples and the ready to use calibrators into the appropriate wells.
2. Cover plate with the enclosed plate sealing foil and incubate at room temperature for 60 minutes.
3. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
4. Pipette 100 µL of enzyme conjugate solution into each well.
5. Cover plate with plate sealing foil and incubate for 30 minutes at room temperature.
6. Discard the contents of the microwells and wash 3 times with 300 µL of diluted washing buffer. Afterwards remove residues of the washing solution by gentle tapping of the microtiter plate on a paper towel.
7. Dispense 100 µL of TMB substrate into each well.
8. Cover plate with the plate sealing foil and incubate for 20 minutes in the dark (e.g. drawer) at room temperature.
9. Add 100 µL of stopping solution to each well.
10. After thorough mixing and wiping the bottom of the plate, read the optical density at 450 nm and calculate the results. Blank against air. A bichromatic measurement using a reference wavelength of 600–690 nm is recommended.

**The developed colour is stable for at least 60 minutes. Read optical densities during this time.**

## 10. Results and Interpretation

### Example

	OD 450 nm	corrected OD	Mean OD Value
Blank	0.010		
Calibrator 1 (negative)	0.039 / 0.041	0.029 / 0.031	0.030
Calibrator 2 (cut-off)	0.744 / 0.698	0.734 / 0.688	0.711
Calibrator 3 (weak positive)	1.318 / 1.360	1.298 / 1.350	1.324
Calibrator 4 (positive)	2.035 / 2.123	2.025 / 2.113	2.069

The table above should be considered as an example which was achieved under arbitrary temperature and environmental conditions. **These data do NOT describe reference values which have to be found in other laboratories in the same way!**

### 10.1. Qualitative Calculation

The calculated OD values for patient sera as mentioned above are compared with the value for the cut-off calibrator. If the value of the sample is higher, then it should be read as positive.

A value below the cut-off calibrator should be read as negative. It seems reasonable to define a range of  $\pm 20\%$  around the value of the cut-off as a grey zone. It is recommended to repeat results laying within the grey zone using the same serum or a new sample of the same patient, taken after 2–4 weeks. Both samples should be measured in parallel in the same run.

The positive calibrator must show an absorption value at least double the value received by the cut-off calibrator.

### 10.2 Quantitative Calculation

The ready to use calibrators of the MEASLES antibody kit are defined and values expressed are in arbitrary units (U/mL). This gives access to an exact and reproducible quantification and in consequence patient antibody titer monitoring is possible. Concentration values for calibrators are printed on the labels of the vials.

A standard curve is plotted by entering the mean absorbance value of the calibrators on the Y-axis and the corresponding concentration on the X-axis using graph paper. The concentration of the patient samples can then be read directly from the graph.

The calculation of the result can be performed using a computer and a suitable software program.

## 11. Assay Performance

The assay characteristics of the MASTAZYME™ MEASLES IgG / IgM / IgA ELISAs have been established and assessed according the European IVD directive. Detailed validation data can be provided on special request.

## 12. References

1. Altintas DU, Evliyaoglu N, Kilinc B; Sen'an DI, Guneser S. The modification in measles vaccination age as a consequence of the earlier decline of transplacentally transferred anti-measles antibodies in Turkish infants. *Eur J Epidemiol*, 12(6): 647-8 (1996).
2. Bayas JM, Vilella A, Vidal J, Nebot X, Carbo JM, Navarro G, Prat A, Asenjo MA, Salleras L.: Susceptibility to measles, rubella and parotitis in young adults. *Med Clin (Barc)*, 106(15): 561-4 (1996).
3. Bouche F, Ammerlaan W, Berthet F, Houard S, Schneider F, Muller CP. Immunosorbent assay based on recombinant hemagglutinin protein produced in a high-efficiency mammalian expression system for surveillance of measles immunity. *J Clin Microbiol*, 36(3): 721-6 (1998).
4. Chiu HH, Lee CY, Chih TW, Lee PI, Chang LY, Lin YJ, Hsu CM, Huang LM. Seroepidemiological study of measles after the 1992 nationwide MMR revaccination program in Taiwan. *J Med Virol*, 51(1): 32-5 (1997).
5. Dagan R, Slater PE, Duvdevani P, Golubev N, Mendelson E. Decay of maternally derived measles antibody in a highly vaccinated population in southern Israel. *Pediatr Infect Dis J*. 14(11): 965-9 (1995).
6. De Souza VA, Pannuti CS, Sumita LM, de Andrade Junior HF: Enzyme-linked immunosorbent assay-IgG antibody avidity test for single sample serologic evaluation of measles vaccines. *J Med Virol*, 52(3): 275-9 (1997).
7. Duvdevani P, Varsano N, Slepon R, Lerman Y, Shohat T, Mendelson E: Determination of immunity to measles virus in young adults: comparative evaluation of a commercial enzyme immunoassay and the hemagglutination inhibition techniques. *Clin Diagn Virol*, 7(1): 1-6 (1996).
8. Johnson CE, Kumar ML, Whitwell JK, Staehle BO, Rome LP, Dinakar C, Hurni W, Nalin DR: Antibody persistence after primary measles-mumps-rubella vaccine and response to a second dose given at four to six vs. eleven to thirteen years. *Pediatr Infect Dis*, 15(8): 687-92 (1996).
9. King SM, Saunders EF, Petric M, Gold R: Response to measles, mumps and rubella vaccine in paediatric bone marrow transplant recipients. *Bone Marrow Transplant*, 17(4): 633-6 (1996).
10. Matter L, Germann D, Bally F, Schopfer K: Age-stratified seroprevalence of measles, mumps and rubella (MMR) virus infections in Switzerland after the introduction of MMR mass vaccination. *Eur J Epidemiol*, 13(1): 61-6 (1997).
11. Mendelson E, Duvdevani P, Varsano N, Lerman Y, Slepon R, Dagan R, Cohen D, Danon Y, Shohat T: Measles immunity and response to revaccination of a young adult population in Israel. *J Med Virol*, 50(3): 249-53 (1996).
12. Metintas S, Etiz S, Akgun Y, Kalyoncu C, Sariboyaci MA, Isikli B: A serological survey of measles vaccine in a rural region of Eskisehir in Turkey. *Public Healt*, 111(6): 373-6 (1997).
13. Narita M, Yamada S, Matsuzono Y, Itakura O, Togashi T, Kikuta H: Measles virus-specific immunoglobulin G subclass response in serum and cerebrospinal fluid. *Clin Diagn Virol*, 8(3): 233-9 (1997).
14. Nates S, Rey G, Giordano M, Medeot S, Depetris A, Boshell J, de Wolff CD: Immunoglobulin M antibody response to measles virus following natural virus infection, primary vaccination and reexposure to the virus. *Viral Immunol*, 10(3): 165-73 (1997).
15. Pison Garcés FJ, Galbe Sanchez-Ventura J, Arcauz Eguren P, Aguirre y Daban C, Mengual Gil J, Larrad Mur L: Immunity to measles, mumps and rubella in children vaccinated with triple viral vaccine. *Aten Primaria*, 15(4): 235-7 (1995).
16. Vardas E, Andre M, Kreis S, Zachariades N, Schoub B: Measles virus transmission from infected to immune individuals-implications for measles control and elimination. *S Afr Med J*, 87(12): 1709 (1997).

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## 12. Bibliographie

1. Altintas DU, Evliyaoglu N, Kilinc B, Sen'an DI, Guneser S. The modification in measles vaccination age as a consequence of the earlier decline of transplacentally transferred anti-measles antibodies in Turkish infants. *Eur J Epidemiol*, 12(6): 647-8 (1996).
2. Bayas JM, Vilella A, Vidal J, Nebot X, Carbo JM, Navarro G, Prat A, Asenjo MA, Salleras L.: Susceptibility to measles, rubella and parotitis in young adults. *Med Clin (Barc)*, 106(15): 561-4 (1996).
3. Bouche F, Ammerlaan W, Berthet F, Houard S, Schneider F, Muller CP. Immunosorbent assay based on recombinant hemagglutinin protein produced in a high-efficiency mammalian expression system for surveillance of measles immunity. *J Clin Microbiol*, 36(3): 721-6 (1998).
4. Chiu HH, Lee CY, Chih TW, Lee PI, Chang LY, Lin YJ, Hsu CM, Huang LM. Seroepidemiological study of measles after the 1992 nationwide MMR revaccination program in Taiwan. *J Med Virol*, 51(1): 32-5 (1997).
5. Dagan R, Slater PE, Duvdevani P, Golubev N, Mendelson E. Decay of maternally derived measles antibody in a highly vaccinated population in southern Israel. *Pediatr Infect Dis J*. 14(11): 965-9 (1995).
6. De Souza VA, Pannuti CS, Sumita LM, de Andrade Junior HF: Enzyme-linked immunosorbent assay-IgG antibody avidity test for single sample serologic evaluation of measles vaccines. *J Med Virol*, 52(3): 275-9 (1997).
7. Duvdevani P, Varsano N, Slepon R, Lerman Y, Shohat T, Mendelson E: Determination of immunity to measles virus in young adults: comparative evaluation of a commercial enzyme immunoassay and the hemagglutination inhibition techniques. *Clin Diagn Virol*, 7(1): 1-6 (1996).
8. Johnson CE, Kumar ML, Whitwell JK, Staehle BO, Rome LP, Dinakar C, Hurni W, Nalin DR: Antibody persistence after primary measles-mumps-rubella vaccine and response to a second dose given at four to six vs. eleven to thirteen years. *Pediatr Infect Dis*, 15(8): 687-92 (1996).
9. King SM, Saunders EF, Petric M, Gold R: Response to measles, mumps and rubella vaccine in paediatric bone marrow transplant recipients. *Bone Marrow Transplant*, 17(4): 633-6 (1996).
10. Matter L, Germann D, Bally F, Schopfer K: Age-stratified seroprevalence of measles, mumps and rubella (MMR) virus infections in Switzerland after the introduction of MMR mass vaccination. *Eur J Epidemiol*, 13(1): 61-6 (1997).
11. Mendelson E, Duvdevani P, Varsano N, Lerman Y, Slepon R, Dagan R, Cohen D, Danon Y, Shohat T: Measles immunity and response to revaccination of a young adult population in Israel. *J Med Virol*, 50(3): 249-53 (1996).
12. Metintas S, Etiz S, Akgun Y, Kalyoncu C, Sariboyaci MA, Isikli B: A serological survey of measles vaccine in a rural region of Eskisehir in Turkey. *Public Healt*, 111(6): 373-6 (1997).
13. Narita M, Yamada S, Matsuzono Y, Itakura O, Togashi T, Kikuta H: Measles virus-specific immunoglobulin G subclass response in serum and cerebrospinal fluid. *Clin Diagn Virol*, 8(3): 233-9 (1997).
14. Nates S, Rey G, Giordano M, Medeot S, Depetris A, Boshell J, de Wolff CD: Immunoglobulin M antibody response to measles virus following natural virus infection, primary vaccination and reexposure to the virus. *Viral Immunol*, 10(3): 165-73 (1997).
15. Pison Garces FJ, Galbe Sanchez-Ventura J, Arcauz Eguren P, Aguirre y Daban C, Mengual Gil J, Larrad Mur L: Immunity to measles, mumps and rubella in children vaccinated with triple viral vaccine. *Aten Primaria*, 15(4): 235-7 (1995).
16. Vardas E, Andre M, Kreis S, Zachariades N, Schoub B: Measles virus transmission from infected to immune individuals-implications for measles control and elimination. *S Afr Med J*, 87(12): 1709 (1997).

**Notizen / Notes / Note:**

**Notizen / Notes / Note:**

**Notizen / Notes / Note:**

**Hersteller / Manufactured by:**  
**Fabriqué par:**

**MAST DIAGNOSTICA**  
**Laboratoriumspräparate GmbH**  
**Feldstraße 20, D-23858 Reinfeld**  
**Deutschland**  
**Tel.: +49 4533 2007-0**  
**Fax: +49 4533 2007-68**  
**www.mastgrp.com**  
**mast@mast-diagnostica.de**

**Distributed by:**

**MAST GROUP Ltd.**  
**MAST House, Derby Road, Bootle**  
**UK-Mersey Side L20 1EA**  
**Great Britain**  
**Phone: +44 151 933 7277**  
**Fax: +44 151 944 1332**  
**www.mastgrp.com**  
**sales@mastgrp.com**

**Distribué par:**

**MAST DIAGNOSTIC**  
**115, Rue Jules Barni**  
**80000 Amiens**  
**France**  
**Tél.: + 33 3 22 80 80 67**  
**Fax: + 33 3 22 80 99 22**  
**www.mastgrp.com**  
**service-commercial@mast-diagnostic.fr**