MASTAZYME™ Chlamydia

MASTAZYME™ Chlamydia

REF 695010

96 Tests

UDI-DI: 4250729700170

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MASTAZYME™ Chlamydia

Intended Purpose

MASTAZYME™ CHLAMYDIA is a qualitative, sensitive enzyme immunoassay for detection of Chlamydia antigen in endocervical and urethral specimens as an aid to diagnosis of chlamydiosis.

The monoclonal antibody used in the assay does react with LPS antigen from *Chlamydia trachomatis*, *Chlamydothrix psittaci* and *Chlamydothrix pneumoniae*.

The assay is suitable for manual or automated use on open EIA processor systems and is intended for professional in-vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

Der MASTAZYME™ Chlamydia is based on the principle of a one step enzyme immunoassay. With this test system, *Chlamydia trachomatis* antigen can be detected from urogenital swabs, which are sent in the MASTAZYME™ transport medium.

The sample material is boiled in the transport medium for 10 minutes. After cooling, the sample material is pipetted into the wells of the microtiter plate together with the conjugate and the monoclonal antibody. Chlamydia-specific lipopolysaccharides (LPS) are recognized during the incubation step (60 min 37 ° C) by monoclonal antibodies (mouse), which form a binding site for the conjugate antibody. Nonspecific and unbound reaction components are removed in a washing step. In the presence of specific immune complexes, the TMB substrate reaction produces a blue colour in the wells, which changes to a yellow colour after stopping.

The reaction batches can now be evaluated with an ELISA plate reader at 450 nm (if available with a reference filter: 600–690 nm). The absorbance is proportional to the antigen concentration.

Content

- STRIPS Microtiter strips**
12 single strips in a holder, each with 8 break-apart wells coated; pretreated.
- AB Monoclonal antibody**
1 x 5.5 mL mouse monoclonal antibody specific for Chlamydia lipopolysaccharide, diluted in blue coloured solution, containing Proclin; ready to use.
- CONJ Conjugate**
1 x 3 mL sheep anti-mouse IgG antibody conjugated with horseradish peroxidase, red coloured solution containing Proclin; ready to use
- CONTROL+ Positive control**
4 x 1 mL positiv control, contains *C. trachomatis* elementary bodies, ready to use.
- CONTROL- Negative control**
4 x 1 mL negative control, ready to use.
- WASH CONC Washing buffer, 30 x concentrated**
1 x 50 mL, 30 x concentrated, containing Proclin.
- SUBS TMB substrate**
1 x 22 mL 3,3', 5,5' Tetramethylbenzidine, ready to use.
- STOP Stop solution**
1 x 12 mL 2 N H₂SO₄ (sulphuric acid), ready to use.
- Instructions for use.**

Further Abbreviation

- RTU ready to use**

Materials Required but not Provided

- Vortex mixer
- Heating block or water bath for heating the samples to 100 °C
- Precision pipettes with 25 µL, 50 µL and 200 µL volume or multichannel pipettes for pipetting the corresponding volumes
- Disposable pipette tips
- Aqua dest. or HPLC grade water
- Clean vessels (measuring cylinders, beakers) with a volume of 500 mL, 1000 mL or 2000 mL
- Centrifuge (2500 x g)
- 37 °C incubator
- Automatic washer or suitable manual washing systems
- ELISA microplate reader with a 450 nm filter (optional reference filter: 620 nm)
- MASTAZYME™ Chlamydia Transport medium (Ref. No. 695020, UDI-DI: 4250729700187)
- MASTAZYME™ Chlamydia Specimen Collection Swabs (Ref. No. 605030, UDI-DI: 05060488970986); female endocervical and male urethral Dacron® tipped swabs only should be used. Wooden shafted

swabs, alginate swabs and agar or charcoal containing swabs should not be used. Dacron® is a registered trade-mark of DuPont Inc.).

Warning and Precautions

1. For in vitro diagnostic use only! Do not ingest! Laboratory safety precautions should be followed. Do not eat, drink or smoke in the laboratory.
2. Read instructions carefully before conducting the assay. Do not modify the instruction procedure. Only laboratory personnel should be permitted to perform the procedures described in this package insert.
3. The MASTAZYME™ Chlamydia is used to examine urogenital specimens. Regulations for the treatment of potentially infectious material must be followed.
4. Always use clean pipette tips for pipetting conjugate, samples and antibodies. Use micro pipettes according to the manufacturer's instructions.
5. Avoid cross-contamination of the reagents.
6. Handle all samples carefully and avoid cross-contamination of the wells.
7. Avoid pipetting at side walls of the wells.
8. Do not allow the wells to dry during the test.
9. Check opened conjugate and antibody vials for possible microbial contamination.
10. Do not continue to use the kit after the expiry date.
11. Do not combine kit reagents from different batches.
12. Do not use reagents that show colour changes or cloudiness.
13. Use the microtiter plate only once.
14. Do not pipette with your mouth.
15. The positive control was found to be non-infectious when grown on appropriate cell cultures; nevertheless, it should be treated as potentially infectious, since infectivity cannot be excluded with the known test methods.
16. Only use laboratory water of the highest quality (e.g. Aqua dest., Aqua bidest., Aqua iniectionis, HPLC grade).
17. All samples or controls must be equilibrated to room temperature before heat inactivation.
18. The chromogenic TMB substrate is flammable and can cause irritation on contact with the skin. Keep away from open flame.
19. The stop solution, 2N H₂SO₄, is corrosive. Avoid contact with skin, mucous membranes and eyes.
(R: 36/38; S: 26-27-45)
20. It must be ensured that all safety regulations relevant to the handling of infectious material are observed when handling patient samples. All controls and samples must be disposed of after use in accordance with the safety regulations for potentially infectious material.

Storage and Stability

1. The MASTAZYME™ Chlamydia and all unused components can be used until the expiry date displayed on the label, if stored at 2–8 °C. All samples or controls must be equilibrated to room temperature before heat inactivation.
2. The reagents are stable for 3 months after opening.
3. Clinical swab specimens can be stored in the transport medium at 2–8 °C for up to 8 days. For longer storage up to 6 months, the samples should be frozen at ≤ -70 °C. However, repeated defrosting and freezing cycles should be avoided.
4. Heat-treated patient samples can be stored at 2-8 °C for 24 hours, frozen at -20 °C for up to 4 weeks.
5. After dilution to working concentration, the washing buffer is stable for up to 30 days at 15–30 °C in tightly capped containers.

For longer storage, samples should be frozen at -20 °C. However, repeated freeze / thaw cycles should be avoided.

Preparation for testing

1. All samples or controls must be equilibrated to room temperature.
2. Fill 50 mL of the wash buffer concentrate with distilled water to 1500 mL (= 1:30 dilution) and mix well. It is recommended to only use the volume of wash solution required

Test procedure

Preliminary note: The following running order is recommended when pipetting the reagents and samples:

1. Monoclonal antibody (blue coloured solution)
2. Conjugate (red coloured solution)
3. Sample (s) or controls

To ensure optimal reaction conditions, no drops of reagents or samples should get stuck on the edge of the wells. The reaction mixture must be mixed carefully before the respective incubation steps.

1. Thoroughly mix the tubes with the sample material and the positive and negative control solutions by shaking (15 sec) before the test.
2. Heat the sample and control tubes for 10 min at 100 °C in a heating block or similar. Open the lids of the tubes slightly (pressure compensation!).
3. After completion of the heat step, take the samples and allow them to cool to room temperature.
4. In the meantime, assemble the microtiter plate with the required number of strips/wells. Consider two wells for the controls (1 negative, 1 positive)!
5. Pipette 50 µL of blue stained antibody solution into the wells for the preparations.
6. Pipette 25 µL of the red stained conjugate.

7. Mix the cooled samples and controls by shaking (15 sec) again.
8. Pipette 200 µL sample or control solution into the respective wells. Record the position of each sample to avoid mix-ups.
9. Seal the wells with the supplied cover foil and mix well.
10. Incubate the reaction mixture at 37 °C for 1 hour (+/- 3 min).
11. After incubation, wash the reaction mixtures 5 x with 350 µL Washing buffer. Remove any remaining washing solution from the wells by gentle tapping. (When using an automatic washing device, the device must be equilibrated with Washing buffer). Prolonged exposure of the Washing buffer to the preparation may lead to inconclusive results!
12. Pipette 200 µL ready-to-use Substrate solution into each well.
13. Cover the reaction wells and incubate for 20 min in the dark at room temperature. Chlamydia-positive reactions show a blue colour.
14. Pipette 50 µL Stop solution into each well in the same order and time interval as for substrate addition.
After addition of Stop solution, Chlamydia-positive reactions turn yellow.
15. Interpret the reaction mixtures within 30 min with an ELISA plate reader at 450 nm. Calibrate the device against air.
16. Record the results of the evaluation.

Criteria for interpretation

Criteria for interpretation: After calibration of the instrument, the mean extinction value of the negative control should be less than 0.2 OD values. The extinction value of the positive control should be higher than 0.7 OD values. If negative and positive controls do not fit into these criteria the test should not be interpreted and has to be repeated. For batch specific OD values see "Certificate of Kit Performance" supplied in the kit.

1. Calculate the average absorbance value obtained for negative controls.
2. Calculate the cut-off value as follows:
Cut-off = mean absorbance of negative controls + 0.1 absorbance units.
3. Calculate the upper limit for the borderline area as follows: Upper limit = cut-off value + 0.05 absorbance units.

Interpretation of Results

1. A test specimen with an absorbance value greater than the upper limit of the borderline area is considered positive.
2. A test specimen with an absorbance value equal to or less than the cut-off level is considered negative.
3. A test specimen with an absorbance between the cut-off level and the upper limit is considered borderline and should be re-tested.

Interpretation of the results during retesting

The determination should be repeated for each sample whose absorbance value lies within the "grey zone". However, it is not necessary to boil the samples again.

2 wells for the negative control and 1 well for the positive control should also be included in the repeat test. If possible, the controls should be taken from the same vials as in the original test.

The cut-off value should be calculated from the repeated negative and positive controls. Absorbances of samples higher than the cut-off value are considered positive. Samples having the same or lower absorbance value than the cut-off are considered negative.

Limitations of the test

1. MASTAZYME™ Chlamydia has not been validated with rectally collected samples.
2. All test results should be interpreted in conjunction with other clinical data and medical judgement.
3. A negative result does not necessarily indicate the absence of chlamydial infection. Numbers of organisms below the detection limit of the kit or improper sampling may give false negative results.
4. These immunoassay detects the LPS antigen from *Chlamydia sp.* In addition to *Chlamydia trachomatis*, *Chlamydophila psittaci* and *Chlamydophila pneumoniae* can lead to positive results.
5. MASTAZYME™ Chlamydia has not been assessed for its use in the determination of the patients response to therapy.
6. Diagnosis of chlamydial infection in forensic cases. MASTAZYME™ CHLAMYDIA should not be used as the only test for diagnosis.

Performance Characteristics

Urogenital specimens

MASTAZYME™ CHLAMYDIA has been tested on 320 samples of men and 1209 samples of women at three independent clinical trial centres in the USA. Specimens were randomly selected from high prevalence populations attending STD clinics with disease incidences. Results were compared with swab culture using McCoy cells, with discrepant results further examined by direct immunofluorescence testing on the centrifuged deposit of the specimen collected for MASTAZYME™ Chlamydia evaluation. Overall accuracy or efficiency of MASTAZYME™ Chlamydia against culture was shown to be 87% in men and 86% in women.

There was little difference in performance between men and women. Combined performance figures were as follows:

Sensitivity	= 89.4 %
Specificity	= 99.0 %

Cross-reactivity

Use of reagents other than those included in the kit or a change in test procedure may result in erroneous results.

The following organisms were found to be non-reactive in the MASTAZYME™ Chlamydia test when present at approximately 10⁹ CFU/test:

Acinetobacter calcoaceticus, *Acinetobacter lwoffii*, *Bacteroides bluus*, *Bacteroides capillus*, *Bacteroides fragilis* (3 isolates), *Bacteroides melanogenicus*, *Bacteroides necrophorus*, *Bacteroides nucleatum*, *Bacteroides sphaericus*, *Candida albicans* (7 isolates), *Corynebacterium diphtheriae* (2 isolates), *Clostridium perfringens*, *Clostridium spp.*, *Cytomegalovirus*, *Enterobacter agglomerans*, *Escherichia coli* (4 isolates), *Gardnerella vaginalis*, *Herpes simplex Viruses types I & II*, *Klebsiella oxytoca*, *Micrococcus sp.*, *Moraxella sp.*, *Neisseria gonorrhoea* (3 isolates), *Neisseria meningitidis*, *Proteus mirabilis*, *Proteus vulgaris*, *Pseudomonas aeruginosa* (2 isolates), *Staphylococcus aureus* (6 isolates, one Protein A - producing Strain), *Staphylococcus epidermidis*, *Staphylococcus saprophyticus*, *Treptococcus sp.* (Group A, B and G), *Ureaplasma urealyticum*.

Collection of Specimens


Specimens from the female cervix and urethra and from the male urethra should contain as many epithelial cells as possible.

If samples are to be additionally tested for gonorrhoea, this should be done with a special swab before taking the chlamydia test.

Female endocervical specimens


1. Clean the exocervical area with sterile gauze or swab before sampling, to remove excess mucus.
2. Insert the appropriated swab into the endocervical canal and rotate vigorously for 5–10 seconds.
3. Move the swab to the portio region and rotate vigorously, then withdraw the swab without touching the vaginal walls.
4. For urethral swabs (approximately 30 % of infections), insert the swab 2-3 cm into the urethra while turning. Rotate the swab vigorously for about 10 seconds to ensure contact with the cylindrical epithelium.
5. Dip the tip of the swab into the transport medium, shake vigorously and break off the plastic shaft at the notch.

6. Close the tube with swab and transport medium and send it to the appropriate laboratory.

Male urethral specimens

1. Specimens should not be collected if the male has urinated within the previous hours.
2. Insert the appropriate male swab 2–3 cm into the urethra rotating the swab for 10 seconds and ensuring that all surfaces of the urethra are contacted.
3. Withdraw the swab and immerse into MASTAZYME™ Chlamydia Transport Medium. Agitate vigorously, break off the swab shaft.
4. Close the container with the swab and transport medium and send it to the appropriate laboratory.

Availability of the summary of safety and performance (Art. 29, IVDR)

The report will be available in the EUDAMED database (as soon as the module is available).

<https://ec.europa.eu/tools/eudamed/#/screen/home>.

The SSP report will be provided on request.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

References and Change History

The references and the change history can be found at the end of the instructions for use.

Referenzen / References:

1. Grayston JT, Kuo CL et al. Chlamydia pneumoniae sp. nov for Chlamydia sp. strain TWAR. Int J Syst Bacteriol 1989; 39: 88-90
2. Schachter J. Chlamydia (Psittacosis-Lymphogranuloma Venerum-Trachoma Group) ch 85 in Manual of Clinical Microbiology 4th edition Ed Lenette A.S.M. 1985
3. Taylor Robinson D, Thomas BJ. The role of Chlamydia trachomatis in genital tract and associated diseases J Clin Pathol 1980; 33: 205-233
4. Brunham RC, Maclean IW et al. Chlamydia trachomatis; it's role in tubal infertility. J Infect Dis 1985; 152: 1275-1282
5. Weström L. Incidence, prevalence and trends in acute pelvic inflammatory disease and its consequence in industrial countries. Am Obstet & Gynecol 1980; 138: 880-892
6. Goh B. Chlamydia trachomatis genital infection. The Practitioner 1988; 232: 813-818
7. Alexander ER, Harrison HR. Role of Chlamydia trachomatis in perinatal infection. Rev Infect Dis 1983; 5: 713-719
8. Wentworth BB, Judson FN. Laboratory Methods for Diagnosis of STD. Am Public Health Ass. 1984.
9. Chlamydia infection, Annual Epidemiological Report for 2017
https://ecdc.europa.eu/sites/portal/files/documents/AER_for_2017-chlamydia-infection.pdf
10. Facts about Chlamydia
<https://ecdc.europa.eu/en/chlamydia/facts>

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Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Intended Purpose (English Version)	Change the wording "intended use" to "intended purpose".
Contents	Self-adhesive plate sealers are no longer included in the assay.
Availability of the summary of safety and performance	A new section Availability of the summary of safety and performance (Art. 29) was added.
Change History	Adding a reference to the Change History at the end of the IFU.

**Mast Diagnostica GmbH,**

Feldstraße 20,
23858 Reinfeld,
Deutschland

Tel: +49 (0)4533 2007 0

Fax: +49 (0)4533 2007 68

email: mast@mast-diagnostica.de

Web: www.mast-group.com

Mast Group Ltd.

Mast House, Derby Road
Bootle, Merseyside, L20 1EA
United Kingdom

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Icons are used according to DIN EN ISO 15223-1

MD V.2.0

2022-11-16