# CHROMagar™ MRSA

## Instructions For Use

Available in several languages

**CHROMagar™ MRSA**  
**Product codes MR500 & MR502**

**NT-EXT-012**

Version 8

**Page 2**

**ENGLISH**

English Version

**FRANCAIS**

Version Française

**ESPAÑOL**

Version Español

**DEUTSCH**

Deutsch Version

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**CHROMagar™ MRSA**  
**Product code MR533-10Kg**

**NT-EXT-079**

Version 2.0

**Page 13**

**ENGLISH**

English Version

**FRANCAIS**

Version Française

**ESPAÑOL**

Version Español

**DEUTSCH**

Deutsch Version

# CHROMagar™ MRSA

Product codes: MR500 & MR502

Instructions For Use

ENGLISH

Instructions For Use

## MEDIUM PURPOSE

Chromogenic medium for the isolation and differentiation of Methicillin Resistant *Staphylococcus aureus* (MRSA). The major issue with this pathogen is its resistance to a large panel of antibiotics, among them beta-lactam antibiotics, limiting the therapeutic options for clinicians.

## COMPOSITION

The product is composed of a powder base and 1 supplement.

Product	=	Base	+	Supplement
Total g/L		82,5 g/L		
Composition g/L		Agar 15.0 Peptones and yeast extract 40.0 Salts 25.0 Chromogenic mix 2.5		
Aspect		Powder Form		Freeze dried vial
STORAGE		15-30°C		2-8°C
FINAL MEDIA pH		6.9 +/- 0.2		

## PREPARATION (Calculation for 1L)

### Step 1

Preparation of the base  
CHROMagar  
MRSA base (B)

- Disperse slowly 82,5g of powder base in 1L of purified water.
  - Stir until agar is well thickened.
  - Autoclave at 110°C during 5 min.
- DO NOT AUTOCLAVE AT 121°C. DO NOT HEAT LONGER THAN 5 MIN.
- Cool in a water bath to 45-50°C.
  - Swirl or stir gently to homogenize.

### Step 2

Preparation of the  
Supplement (S)

- For reconstitution, aseptically rehydrate CHROMagar MRSA supplement ref SU620 with 20.0ml of sterile water.
- Mix slowly.
- Add 1mL of CHROMagar MRSA supplement to the prepared CHROMagar MRSA medium (Step1).

### Step 3

Pouring

- Pour into sterile Petri dishes.
- Let it solidify and dry.

### Storage

- Store in the dark before use.
  - Prepared media plates can be kept for one day at room temperature.
  - Plates can be stored for up to one month under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.
- Advice 2: If not fully used, rehydrated CHROMagar MRSA supplement can be stored one week at 2-8°C or up to 2 months at -20°C.

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
## PERFORMANCE & LIMITATIONS

- Some studies showed sensitivity for MRSA of 100% and specificity of 100% (de Gialully et al. 2004 & Taguchi et al. 2004).
- Definite identification as MRSA requires, in addition, a final identification as *S.aureus*.

## QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms.  
Good preparation of the medium can be tested, isolating the ATCC strains below:

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### • CHROMagar MRSA MR500 & MR502

Σ Pack Size	Ordering References	Base (B)	Supplement (S)
1000 ml 50 Tests of 20ml	MR500	MR500 Weight: 82.5gr	+ SU620
5000 ml 250 Tests of 20ml	MR502	MR502 Weight: 412,5gr	+ SU620

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach - ATCC® is a registered trademark of the American Type Culture Collection  
NT-EXT-012 V8 / 08-Jul-15

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Web link: <http://www.chromagar.com/publication.php>

## IFU/LABEL INDEX

- Quantity of powder sufficient for X liters of media
- Expiry date
- Required storage temperature
- Store away from humidity

### Need some Technical Documents?

- Available for download on [www.CHROMagar.com](http://www.CHROMagar.com)
- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

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**CHROMagar**  
The Chromogenic Media Pioneer

CHROMagar 4 place du 18 juin 1940  
75006 Paris - France  
Eメール: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
電話番号 +33 (0)1.45.48.05.05. ウェブサイト: [www.CHROMagar.com](http://www.CHROMagar.com)

IVD

CE

# CHROMagar™ MRSA

Product code MR533-10kg

Instructions For Use

ENGLISH

Instructions For Use

## MEDIUM PURPOSE

Chromogenic medium for the isolation and differentiation of Methicillin Resistant *Staphylococcus aureus* (MRSA).

The major issue with this pathogen is its resistance to a large panel of antibiotics, among them beta-lactam antibiotics, limiting the therapeutic options for clinicians.

## COMPOSITION

The product is composed of a powder base (A) and 3 supplement (B+C+SU).

Product	=	Base (A)	+	Liquid B	+	Mix C	+	Supplement SU
Total g/L		80.4 g/L		2 g/L		0.25 g/L		2.6 mg/L
Composition g/L		Agar 15.0 Peptones and yeast extract 40.0 Salts 25.0 Chromogenic mix 0.4		Growth factors 2.0		Selective and chromogenic mix		Selective mix
Aspect		Powder Form		Liquid Form		Powder Form		Powder Form
STORAGE		15-30°C		15-30°C		2-8°C		2-8°C
FINAL MEDIA pH		6.9 +/- 0.2						

## PREPARATION (Calculation for 1L)

<b>Step 1</b> Preparation of CHROMagar MRSA base (A) + Liquid B	<ul style="list-style-type: none"> <li>Disperse slowly 80.4g of CHROMagar MRSA <b>Base A</b> powder in 1L of purified water.</li> <li>Add 2g of CHROMagar MRSA <b>Liquid B</b> to it, into slurry.</li> <li>Stir until agar is well thickened.</li> <li>Autoclave (Base A + Liquid B) at 121°C during 15 min.</li> <li>Cool to 45/50°C keeping stirring.</li> </ul>	<b>Final Media</b> <b>HELPING CALCULATION</b> 5 L   Base 402g / Liquid B : 10g 10 L   Base 804g / Liquid B : 20g 25L   Base 2010g / Liquid B : 50g
<b>Step 2</b> Preparation of CHROMagar MRSA Mix C	<ul style="list-style-type: none"> <li>In a transparent vessel, add 250 mg of Mix C powder to 10 ml of purified water.</li> <li>Place under agitation with a magnetic stirring until Mix C is solubilized. (for 10ml --&gt; 30 min of stirring)</li> <li>Filter sterilise at 0.45 µm. <b>DO NOT USE 0.22 µm filter.</b></li> <li>And aseptically add into slurry (Base A + Liquid B) cooled to 45/50°C while mixing.</li> </ul>	<b>Final Media</b> <b>HELPING CALCULATION</b> 5 L   Mix C : 1.25g in 50 ml 10 L   Mix C : 2.5g in 100 ml 25L   Mix C : 6.25g in 250 ml
<b>Step 3</b> Preparation of CHROMagar MRSA Supplement	<ul style="list-style-type: none"> <li>For reconstitution make a stock solution of CHROMagar MRSA supplement ref SU625-BA at 2mg/ml (add 20mg of SU625-BA powder to 10ml of purified water). Homogenize and filter sterilise at 0,45 µm.</li> <li>Add 1.3 ml into melted (Base A + Liquid B + Mix C) cooled to 45/50°C.</li> </ul>	<b>Final Media</b> <b>HELPING CALCULATION</b> 5 L   Add 6.5 ml into melted base 10 L   Add 13 ml into melted base 25L   Add 32.5 ml into melted base
<b>Step 4</b> Pouring	<ul style="list-style-type: none"> <li>Swirl or stir gently to homogenize. Pour into sterile Petri dishes.</li> <li>Let it solidify and dry.</li> </ul>	
<b>Storage</b>	<ul style="list-style-type: none"> <li>Store in the dark before use.</li> <li>Prepared media plates can be kept for one day at room temperature.</li> <li>Plates can be stored for up to 2 weeks under refrigeration (2/8°C) if properly prepared and protected from light and dehydration.</li> </ul>	

## INOCULATION

Related samples can be processed by direct streaking on the plate, as well as prior appropriate enrichment step.

- If the agar plate has been refrigerated, allow to warm to room temperature before inoculation.
- Streak sample onto plate.
- Incubate in aerobic conditions at 37°C for 18-24 hours.

## Typical Samples

e.g. nasal, perineal, throat, rectal specimens  
\*\*\*

Possible enrichment step  
Direct streaking  
or spreading technique

# CHROMagar™ MRSA

Product code MR533-10kg

Instructions For Use

ENGLISH

Instructions For Use

## INTERPRETATION

Microorganism	Typical colony appearance
MRSA	→ pink to mauve
MSSA	→ inhibited
Other bacteria	→ inhibited, colourless, blue

### Typical colony appearance

The pictures shown are not contractual.

## PERFORMANCE & LIMITATIONS

- Some studies showed sensitivity for MRSA of 100% and specificity of 100% (de Gialully *et al.* 2004 & Taguchi *et al.* 2004).
- Definite identification as MRSA requires, in addition, a final identification as *S.aureus*.

## QUALITY CONTROL

Please perform Quality Control according to the use of the medium and the local QC regulations and norms. Good preparation of the medium can be tested, isolating the ATCC strains below:

Microorganism	Typical colony appearance
MR <i>S.aureus</i> ATCC® 33592	→ mauve
MS <i>S.aureus</i> ATCC® 25923	→ inhibited
<i>P.aeruginosa</i> ATCC® 9027	→ inhibited
<i>E.faecalis</i> ATCC® 29212	→ inhibited
<i>E.coli</i> ATCC® 25922	→ inhibited
<i>C.albicans</i> ATCC® 10231	→ inhibited

## WARNINGS

- Do not use plates if they show any evidence of contamination or any sign of deterioration.
- Do not use the product beyond its expiry date or if product shows any evidence of contamination or any sign of deterioration.
- For *in vitro* diagnostic use. This laboratory product should be used only by trained personnel in compliance with good laboratory practices.
- Any change or modification in the procedure may affect the results.
- Any change or modification of the required storage temperature may affect the performance of the product.
- Unappropriate storage may affect the shelf life of the product.
- Recap the bottles/vials tightly after each preparation and keep them in a low humidity environment, protected from moisture and light.
- For a good microbial detection: collection and transport of specimen should be well handled and adapted to the particular specimen according to good laboratory practices.

## DISPOSAL OF WASTE

After use, all plates and any other contaminated materials must be sterilized or disposed of by appropriate internal procedures and in accordance with local legislations. Plates can be destroyed by autoclaving at 121°C for at least 20 minutes.

## REFERENCES

Please refer to our website page «Publications» for scientific publications about this particular product.

Web link: <http://www.chromagar.com/publication.php>

## IFU/LABEL INDEX

Quantity of powder sufficient for X liters of media

Expiry date

Required storage temperature

Store away from humidity

Need some  
Technical Documents?

Available  
for download on  
[www.CHROMagar.com](http://www.CHROMagar.com)

- Certificate of Analysis (CoA) --> One per Lot
- Material Safety Data Sheet (MSDS)

Pack Size

10KG

6200 Tests  
of 20ml

=

Ordering References

MR533-10KG

=

Base (A)

MR533-10Kg

Liquid B

MR513-125(B)

+

Mix C

MR513-125(C)

Supplement SU

SU625-BA

CHROMagar™ and Rambach™ are trademarks created by Dr A. Rambach - ATCC® is a registered trademark of the American Type Culture Collection  
NT-EXT-079 V2.0 / 14-Feb-18

**CHROMagar**  
The Chromogenic Media Pioneer

CHROMagar 4 place du 18 juin 1940  
75006 Paris - France  
Email: [CHROMagar@CHROMagar.com](mailto:CHROMagar@CHROMagar.com)  
Tel +33 (0)1.45.48.05.05. Website: [www.CHROMagar.com](http://www.CHROMagar.com)

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