TOXOREAGENT

TOXOREAGENT

REF

RST7001

60 Tests

Art.-Nr. 487001

UDI-DI: 4250729700002

**Gebrauchsanweisung /
Instructions for Use /
Notice Technique**

**Nur zur Verwendung durch Fachpersonal /
For professional use only /
Exclusivement à un usage professionnel**

CE 2797

Deutsch

Seiten

02-06

English

Pages

07-11

Français

Pages

12-16

TOXOREAGENT

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TOXOREAGENT

Agglutination assay for the detection of *Toxoplasma gondii* antibodies in serum

Intended Purpose

Semi-quantitative (titratable) agglutination assay for the detection of antibodies directed against *Toxoplasma gondii* in human serum as an aid to diagnosis of immunity against toxoplasmosis (presence of specific antibodies).

The assay is suitable for manual use and is intended for professional *in-vitro* diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

Note that the CE-mark affixed to this product applies only to the use of the product with human samples, not with animal samples.

Important note for use of these kit instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page of the IFU. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

The Mast Toxoreagent Kit is an agglutination test based on latex particles. These are coated and stabilized with ultrasonically treated *Toxoplasma gondii* antigen (RH strain). The latex reagent is available as a liquid reagent component in the kit ready for use. The sample diluent is used for the dilution of the samples or the positive control. The positive control shows the reactivity of the test. The final titre to be obtained should be taken from the batch specific QC document. The test is performed in microtiter plates with U-bottom, which are not part of the kit. The test result is an agglutination pattern which is interpreted visually or by corresponding reader-controlled evaluation programs. The presence of antibodies directed against *Toxoplasma gondii* antigens leads to agglutination of the latex particles. This positive reaction depends on the concentration of antibodies in the serum sample. With negative results, the particles sediment to a button at the bottom of the microtiter well. The strength of the antibody reaction is determined by serial titration. Starting from the serum dilution of 1:8, the sample is titrated over a further 7 wells. The final agglutination reaction of this sample yields the titre. The maximum titration step to be achieved in this titration is a dilution step of 1:2048.

Contents

2. **AG/BEADS** Latex reagent or latex particles
1 x 12 mL, ready to use. Contains polystyrene latex particles coated with inactivated *Toxoplasma gondii* antigen and proclin (0.1% v/v) as preservative.
3. **CONTROL+** Positive Control
1 x 0.5 mL, ready to use; contains human serum with anti-*Toxoplasma gondii* antibodies and proclin (0.1% v/v) as preservative. See titre concentration on QC document supplied in the kit. (potentially infectious, see Warning and Precautions)
4. **DIL** Diluent
1 x 50 mL. Contains proclin (0.1% v/v) as preservative.
5. **Instructions for use**

Further abbreviations

6. **RTU** ready to use

Materials required but not provided

Standard microbiological supplies and equipment such as micropipettes, small dilution tubes and rigid U-shaped bottom microwell plates.

Warnings and precautions

- For *in vitro* diagnostic use only.
- Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
- Observe approved biohazard precautions and aseptic techniques when dealing with specimens of human origin.
- The positive control contains diluted human serum. All sera, plasma and buffers—containing human biological material were found to be negative for Anti-HIV 1 & 2, HIV Ag, Anti-HCV and HBsAg. However, precautions such as the use of gloves must be taken. As the presence of these viruses cannot be completely ruled out, precautions such as the use of personal protective equipment should be undertaken.
- Sterilise all biohazard waste before disposal.
- Agglutination tests are sensitive to the effects of heat, direct sunlight and vibration. Keep away from such sources during test incubation periods.
- Some kit components contain bovine serum albumin. As an infection risk after ingestion, inhalation or embrocation of reagents cannot be excluded entirely, care when handling the reagents is advised.
- The assay is based on latex particles. Latex may cause allergic reactions with skin. Wear gloves, coats, protection glasses etc. to reduce the risk of contact.
- Latex particles must be resuspended thoroughly before use. In case of insufficiently resuspended Latex reagent, reactions may be hard to read due to latex button size.

- The presence of anti-latex antibodies in the patient sample can lead to false positive results. Patients with a known latex allergy should be tested with different methods and all results should be interpreted together with clinical observations.
- Use a separate disposable tip for each sample to prevent cross contamination.
- Reseal caps on all reagents immediately after use.
- Do not use damaged or contaminated kit components.
- Kit components are matched and should not be interchanged between batches.
- The assay contains no hazardous substances according to Regulation (EC) 1907/2006. For further safety information see Material Safety Data Sheet, available on request or from the Mast website.

Stability and storage

- Store unopened at 2-8 °C in an upright position until the expiry date shown on the pack label.
- Once opened, the latex reagent and diluent should be stored at 2-8 °C with the lids tightly closed and may be used until the expiry date given on the label.
- The positive control should be stored at 2-8 °C with the lid tightly closed and may be used until the expiry date indicated on the label.
- **Do not freeze any of the reagents.**
- Serum samples may be stored at 2-8°C for up to 72 hours prior to testing. If longer storage is required, store at -20°C. Thawed samples must be mixed well prior to testing.

Sample Material

- c) Serum
- d) The serum samples can be stored at 2-8 °C for up to 3 days. They should be aliquoted and kept at -20 °C for longer storage.

The samples should not be frozen and thawed repeatedly.

After thawing samples should be briefly vortexed carefully before being used in the assay.

Lipemic, icteric, hemolytic or bacterially contaminated samples can cause false results.

Test Procedure

Allow the **TOXOREAGENT** components to equilibrate to **room temperature** (15-25 °C) before use.

A. Sample preparation

1. Use fresh serum samples obtained after blood clot removal by centrifugation. Serum samples may be stored at 2-8 °C for up to 72 hours or frozen for long-term storage.
2. Dilute all samples (except pos. control, see below) 1:8 by adding 50 µL of serum sample to a small tube and then add 350 µL of diluent. Mix well before use. Heat

inactivation of specimens is not required prior to testing.

B. Positive Control

The positive control is supplied ready to use

C. Procedure

1. Gently shake the latex reagent to suspend the particles. Avoid frothing.

Note: Make sure that all latex particles are thoroughly resuspended before use (see warnings and precautions).

2. Place 25 µL of diluent into wells 1-8 of row 1 of a rigid U-shaped bottom microwell plate using a micropipette. Repeat this operation for as many serum samples that are to be tested. All wells of column 9 are to be used as a control.
3. Dispense 25 µL of the reconstituted positive control to well 1 of the first row then add 25 µL of 1:8 diluted serum samples to well 1 of subsequent rows according to a pre-planned template.
4. Mix the contents of all wells in column 1 by pipetting up and down a few times and transfer 25 µL to the corresponding wells of column 2.
5. Continue mixing and transferring 25 µL amounts as in step 4 until you reach column 8. Finally mix and discard 25 µL from all wells of column 8. This will make a series of doubling dilutions from 1:16 in column 1 to 1:2048 in column 8. See diagram below.

a	1	2	3	4	5	6	7	8
b	1:16	1:32	1:64	1:128	1:256	1:512	1:1024	1:2048
c	25	25	25	25	25	25	25	25
d	25							
e	→	25 →	25 →	25 →	25 →	25 →	25 →	25 → dis- card
f	Add 25 µL latex reagent to all wells.							
	Mix contents of wells.							
	Leave tray at room temperature for at least 12 hours.							
	Read agglutination patterns.							


Key to Table

- a – Well number.
- b – Serum dilution.
- c – µL diluent.
- d – µL of positive control or 1:8 diluted serum samples.
- e – Mix diluted serum + diluent and transfer.
- f – Add latex reagent to all wells.

- Add 25 µL of latex reagent to all wells. **Ensure that the latex reagent is fully resuspended before use.**
- Mix the contents of the wells by gently tapping all four sides of the plate or by using a plate shaker set on high speed for at least 30 seconds.
- Cover the plate to avoid evaporation and incubate at room temperature (15-25 °C) for at least 12 hours, **on a horizontal surface, in an undisturbed position that is free from vibration and away from heat and direct sunlight.**
- Read the wells for agglutination patterns.

Reading and interpretation of results

Read the agglutination patterns according to the following criteria:

Result	Pattern	Interpretation
Intense agglutination with outer edge irregular and showing occasional tearing.		3+
Agglutination spread throughout the well.		2+
Agglutination is clearly present but does not cover the complete surface of the well.		1+
Weak agglutination pattern		0,5+
Negative reaction. A small and distinct circular sedimentation at the centre.		-

The titre is counted as the highest dilution of the test serum showing agglutination of 1+ or above.

The assay results can be considered as valid if the positive control matches the criteria given in the QC document.

For human serum

- An antibody titre (serum dilution factor) of 1:32 or above indicates the presence of toxoplasma antibodies. It is suggested that for all positive samples the patient should be monitored every 2 weeks to assess for antibody levels. A rising titre level or maintenance of a high level indicates toxoplasmosis.
- An antibody titre of 1:16 should be reported as positive but indicates the presence of a borderline level of toxoplasma antibodies. The probability of toxoplasmosis is low. Re-examination of the patient may be required by taking further blood samples to assess for a possible rise in titre value. If no rise in antibody titre is observed in subsequent samples the patient may be considered to be negative to toxoplasmosis.
- Antibody titres of less than 1:16 may indicate the absence of clinically significant levels toxoplasma antibodies. In case of any doubt the test shall be repeated using another method.

- If a negative result is obtained from a pregnant woman it is advisable to repeat the test during the period of pregnancy to detect a possible seroconversion.

In case of doubt it is recommended to request another serum sample and to consult the physician.

For pigs and cats

- An antibody titre (serum dilution factor) of 1:64 or above indicates the presence of toxoplasma antibodies.
- It is suggested that for all positive samples further samples should be taken and monitored every 2 weeks to assess for antibody levels.
- An antibody titre of 1:32 indicates the presence of a borderline level of toxoplasma antibodies. The probability of toxoplasmosis is low. Re-examination of the animal by taking further blood samples may be required to assess for a possible rise in titre value. If no rise in antibody titre is observed in subsequent samples the animal may be considered to be negative to toxoplasmosis.
- An antibody titre of 1:16 or less indicates the absence of clinically significant levels toxoplasma antibodies.

Limitations of use

- TOXOREAGENT detects IgG and IgM antibodies. The assay cannot differentiate between the two immuno globulin classes.
- TOXOREAGENT is useful for the serological diagnosis of toxoplasmosis. Human serum samples with titres of 1:16 are counted as weak positive (borderline level of toxoplasma antibodies) and should be retested after approximately 2 weeks.
- When testing samples from immunocompromised, congenitally infected or renal dialysis patients serological test results should be considered along with clinical history and other aspects of patient management to be considered diagnostically significant. It is advisable to confirm the patient's anti-toxoplasma status with an alternative test.
- False positive reactions have been noted in patient samples positive for antinuclear antibodies, cytomegalovirus IgM antibodies, or HBe antigens. Results from any serological test should always be considered along with the patient's clinical symptoms and history (Lit. 8, 9).
- Each diagnostic assay may give false-positive or false-negative results. In case of any doubt it is advisable to confirm the patient's anti-toxoplasma status with an alternative test.
- Serum samples only should be used in this test.
- Samples showing visibly lipaemic, haematic, icteric or other contaminations should be avoided.
- Samples up to the following concentrations were tested not showing any strange results, but nevertheless they should not be used:

Haemoglobin: 560 mg/dL
 Bilirubin: 24 mg/dL
 Triglycerides: 1500 mg/dL

- Occasionally a prozone effect may be observed with strongly positive sera. Please note that to prevent the prozone effect all samples should be fully titrated out. Screening is not recommended.

Performance characteristics

Sensitivity and Specificity

In a study of 147 human sera samples which were tested semi-quantitatively, the sensitivity of TOXOREAGENT was found to be 93.8 % and the specificity was found to be 97.1 %. Positive results taken as $\geq 1:32$ for TOXOREAGENT.

The analytical performance was calculated as follows:

	Formula	Value
Sensitivity (Sens)	$\frac{TP}{(TP + FN)}$	93.81%
Specificity (Spec)	$\frac{TN}{(TN + FP)}$	97.06%
Positive predictive value (PPV)	$\frac{TP}{(TP + FP)}$	99.07%
Negative predictive value (NPV)	$\frac{TN}{(TN + FN)}$	82.50%
Efficiency (Eff)	$\frac{(TP + TN)}{total}$	94.56%
Positive likelihood ratio (LR+)	$\frac{(Sens)}{(1 - Spec)}$	31.98
Negative likelihood ratio (LR-)	$\frac{(1 - Sens)}{(Spec)}$	0.06

Assay comparison

In a study with 141-147 sera, sensitivity and specificity of TOXOREAGENT is in line with performance data of comparable assays:

	TOXO-REAGENT	Toxoreagent (EIKEN / Mast Group Ltd.)	MASTA-FLUOR™ Toxoplasma Screen (Mast Diagnostica GmbH)
Art.No.		RST 701	631181
Sens	93.81%	95.58%	99.09%
Spec	97.06%	64.71%	83.87%
PPV	99.07%	90.00%	95.61%
NPV	82.50%	81.48%	96.30%
Eff	94.56%	88.44%	95.74%
LR+	31.98	2.71	6.14
LR-	0.06	0.07	0.01

Positivity or negativity of samples used in the above mentioned comparison were confirmed by enzyme immunoassay technology (Lit. 10, 11).

Traceability - Quality control and reference material

It is recommended that quality control should be performed with the positive control provided to verify that the latex reagent is working correctly, and should be repeated at regular intervals. The antibody titre of the positive control is shown in the QC document supplied with the kit. If the positive control result falls outside of this range the test run should be considered invalid.

Also periodically check that latex reagents give correct results with known reference serum samples. Check for signs of deterioration. Do not use reagents if they are contaminated or cloudy.

Control materials used in the kit were tested against International Standards:

NIBSC 01/600:	1:2048
TOXM:	1:2048
NIBSC 09/B588:	1:1024

Assay Precision

When positive control titration (serial dilution 1:16 to 1:2048) was tested 10 times from different lots according to the test procedure, all results were found to be within ± 1 titer level.

High dose hook effect (= prozone effect)

Occasionally a prozone effect may be observed with strongly positive sera. To prevent the prozone effect, all samples should be titrated according to kit instruction. Screening is not recommended (see also Limitations of use).

Measuring range

The positive control and samples were titrated (1:16-1:32768) to evaluate agglutination pattern (negative to 3+). Illustrative pictures can be found in section Reading and interpretation of results.

Cut-off

The cut-off ($\geq 1:32$) was defined on a panel consisting of 147 positive and negative sera samples corresponding to the mentioned sensitivity and specificity.

Availability of the summary of safety and performance (Art. 29, IVDR)

The report will be available in the EUDAMED database (as soon as the module is available).

<https://ec.europa.eu/tools/eudamed/#/screen/home>.

The SSP report will be provided on request.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

References and Change History

You can find the references and the Change History at the end of the instruction for use.

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1. Rodrigues IMX, Castro AM, Gomes MBF, Amaral WN, Avelino MM; Congenital toxoplasmosis: evaluation of serological methods for the detection of anti-Toxoplasma gondii IgM and IgA antibodies, Mem. Inst. Oswaldo Cruz, Rio de Janeiro, 2009, Vol. 104 (3), 434-440
1. Kobahashi A, Hirai N, Suzuki Y, Nishikawa H, Watanabe N. Evaluation of a commercial Toxoplasma latex agglutination test. Jap J Parasitol. 1977; 26: 175-180.
2. Balfour AH, Fleck DG, Hughes PA, Sharp D. Comparative study of three tests (dye test, indirect haemagglutination, latex agglutination test) for the detection of antibodies to *Toxoplasma gondii* in human sera. J Clin Pathol. 1982; 35: 228-232.
3. Hollimann, RE. Investigation of HIV-positive patients for Toxoplasmosis using the latex agglutination test. Serodiagn Immunother Inf Dis 1990; 4: 249-253.
4. <http://www.cdc.gov/parasites/toxoplasmosis/index.html>
5. <http://ecdc.europa.eu/en/healthtopics/toxoplasmosis/pages/index.aspx>
6. <http://www.rki.de/DE/Content/InfAZ/T/Toxoplasmosis/Toxoplasmosis.html>
7. <http://www.bfr.bund.de/en/toxoplasms-54399.html> BgVV-Hefte^{SEP} Herausgegeben von M. Hartung, Bericht über die epidemiologische Situation der Zoonosen in Deutschland für 2001, 249 – 251
8. Balfour AH, Prestage ES, Noel I. False positive results in the Toxoreagent test for *Toxoplasma gondii* in immunocompetent patients, J. Clin. Pathol. 1988; 35 (1): 1135-1136.
9. Holliman RE, Johnson J, Duffy K, New L; Discrepant toxoplasma latex agglutination test results, J. Clin. Pathol. 1989; 42: 200-203.
10. Liu Q, Wang Z-D, Huang S-Y, Zhu X-Q; Diagnosis of toxoplasmosis and typing of *Toxoplasma gondii*, Parasites & Vectors 2015; 8 (292): 13071-015.
11. Ybañez D, Ybañez A P and Nishikawa Y; Review on the Current Trends of Toxoplasmosis Serodiagnosis in Humans, Front. Cell. Infect. Microbiol. 2020; 10 :204.

Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Cover page	Addition of the UDI-DI.
Principle of the Test	Description included.
Materials required but not provided	Included.
Warnings and precautions	Addition of a warning that anti-latex antibodies could lead to false positive results. A note inserted that the instructions for use must be read carefully.
Test procedure	Changing the order of the sections in the IFU. The information that the reagents must be brought to RT before use has been moved directly under the test procedure. Multiple mention that the positive control is ready to use, for better understanding.
Limitations of use	Literature reference have been added. Triglyceride concentration updated.
Performance Characteristics	Addition of the sections Assay comparison, Precision, Measuring range, High hose hook effect, Clinical performance. Update of the sections Analytical sensitivity and specificity / cut-off (number of samples, abbreviations, formula, values, cut-off level). Update of the section Interference / Cross-reactions (concentration of triglycerides). A new section Availability of the summary of safety and performance (Art. 29) was added. A new section Reporting serious incidents was added.
References	Addition of four references.

**Mast Diagnostica GmbH,**

Feldstraße 20,
23858 Reinfeld,
Deutschland

Tel: +49 (0)4533 2007 0

Fax: +49 (0)4533 2007 68

email: mast@mast-diagnostics.de

Web: www.mast-group.com

Mast Group Ltd.

Mast House, Derby Road
Bootle, Merseyside, L20 1EA
United Kingdom

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Icons are used according to DIN EN ISO 15223-1

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MD V.2.0

2022-11-16