## **MASTAFLUOR™ FTA-ABS IgG**

MASTAFLUOR™ FTA-ABS IgG

**REF** 630522

10 x 10 Tests

**UDI-DI** 4250729700064

**Gebrauchsanweisung /  
Instructions for Use /  
Notice d'utilisation**

**Nur zur Verwendung durch Fachpersonal /  
For professional use only /  
Exclusivement pour un usage professionnel**

**CE** 2797

Deutsch

Seiten

02–06

English

Pages

07–11

Français

Pages

12–16

**MASTAFLUOR™ FTA-ABS IgG**

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## MASTAFLUOR™ FTA-ABS IgG

**Immunofluorescence assay for the detection of IgG antibodies to *Treponema pallidum* in human serum and plasma.**

### Intended Use

Semi-quantitative (titratable) immunofluorescence assay for the detection and confirmation of IgG antibodies directed against *Treponema pallidum* in human serum and plasma as an aid to diagnosis of syphilis.

The assay is suitable for manual use and is intended for professional in-vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

### Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

### Principle of the Test

Diluted specimens, the negative and positive controls or calibrators are applied to the wells of microscope slides and incubated. All wells are coated with purified antigens of the very pathogen. If specific antibodies are present in the serum they will bind to the fixed pathogen antigens forming a stable antigen-antibody complex. Slides are then washed to remove any unbound material. Complexed antibodies are detected by the addition of a fluorescence labelled anti-human IgG immunoglobulin conjugate. After a further washing step to remove any unbound conjugate, slides are viewed under a fluorescence microscope.

A green fluorescence is observed if pathogen specific antibodies are present in the sample material.

### Kit Contents

- SLIDE Slides**  
Ten slides with 10 wells each coated with fixed *Treponema pallidum* cells (Nichols strain).
- CONTROL+ Positive control**  
1 x 1 mL human serum with *Treponema pallidum* specific antibodies (IgG). Ready to use. (potentially infectious, see Warning and Precautions)  
Contains < 0.1 % (w/v) sodium azide as preservative.
- CONTROLAC Absorption control (serum)**  
1 x 1 mL of non-specific control serum (human). Ready to use. (potentially infectious, see Warning and Precautions)  
Contains < 0.1 % (w/v) sodium azide as preservative.
- CONJG FITC Conjugate IgG**  
1 x 3 mL of anti-human immunoglobulin ( $\gamma$  chain, rabbit) conjugated to fluorescein isothiocyanate (FITC). Ready to use.  
Contains Proclin (0.1 % v/v) as preservative.
- BUFFER PBS**  
2 sachets with 10 g PBS powder. Dissolve 1 sachet (10 g) in 1 litre distilled or deionised water to make a solution of phosphate buffered saline (PBS) at pH 7.2  $\pm$  0.2.
- MOUNTING MEDIUM Mounting Medium**  
1 x 3 mL of a buffered glycerol mounting medium. Contains Proclin (0.1 % v/v) as preservative. Ready-to-use.
- SORBENT Sorbent/Ultrasonicate**  
1 x 4 mL culture extract of *Treponema reiteri*, ultrasonicate, liquid, ready to use.  
Contains <0.1 % (w/v) sodium azide as preservative.
- Instructions for use**

### Further Abbreviations

- RTU** ready to use

### Materials Required but not Provided

- Sterile test tubes.
- Micropipettes and tips.
- Staining dish or Coplin jar.
- Moist chamber.
- Volumetric flask.
- Distilled water or deionised.
- Forceps.
- Cover slips.
- Wash bottle.
- Fluorescence microscope with a filter combination suitable for FITC (e.g. 490 nm excitation filter and a 510nm barrier filter).

## Warning and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
3. Do not use reagents beyond the expiry date.
4. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
5. Do not mouth pipette.
6. Use disposable plasticware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. The control/calibrator materials of human origin provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the samples are free of infections or microbial contamination.
8. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
9. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
10. Do not allow wells to dry out during the assay procedures.
11. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
12. Samples and contaminated disposables should be disposed according to relevant disposal directives and regulations for infectious materials. For the decontamination of surfaces use suitable surface disinfectants. Contaminated glassware should be autoclaved at 121 °C.
13. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
14. Microbial contaminated samples should not be used.
15. Do not use haemolysed, lipaemic or icteric samples.
16. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose of by flushing to drain with plenty of water.
17. Use distilled or deionised water.

## Stability and Storage

MASTAFLUOR™ FTA-ABS IgG can be used until the end of expiry date as indicated on kit label. All kits components and reagents should be stored at 2–8 °C and brought to room temperature before use.

After first opening reagents are stable for 90 days. An opened slide shall be used the same day.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Samples (serum, plasma) may be stored according to general recommendation in literature. In general samples can be kept at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

## Sample Material

- a) Serum / Plasma material is collected by professional personnel according to current standards / best practice guidelines. According to current understanding an infection by contaminated blood cannot be ruled out entirely. Any sample material should therefore be treated as potentially infectious.
- b) Samples should be taken and can be stored at 2–8 °C for up to 3 days according to general recommendation in literature. Serum should be aliquoted immediately after sampling and kept at -20°C for longer storage.

The samples should not be frozen and thawed repeatedly.

After thawing samples should be briefly vortexed carefully before being used in the assay.

Lipaemic, haemolytic, icteric or bacterially contaminated samples can cause false positive or false negative results.

## Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or deionised water.
3. The Sorbent is ready to use.  
Dilute all serum specimens 1:5 in Sorbent.  
Example: 10 µL serum + 40 µL Sorbent
4. Dilute control sera as indicated below:
  - a. Positive control + Sorbent (1:5)
  - b. Positive control + PBS (1:5)
  - c. Absorption control + Sorbent (1:5)
  - d. Absorption control + PBS (1:5)Optionally other controls can also be diluted with sorbent and PBS. After absorption samples / controls can be further diluted with PBS to determine the antibody content.
5. Incubate all samples at 37 °C for 30 min.
6. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
7. Apply 20–25 µL of treated specimens and controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells. Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.
8. Cover the moist chamber and incubate at 37 °C for 30 min.

9. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the center of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
10. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min with a change of PBS after 5 min to increase washing strength.
11. Remove slides from the staining jar and drain off any excess buffer. Using a blotter, dry the area outside wells. Do not touch the well surface.
12. Immediately transfer the slide to a moist chamber and add 20–25 µL of the respective FITC-conjugate to each well. The conjugate shall cover the whole well. Do not allow the wells to dry.
13. Incubate slides at 37 °C in a covered moist chamber in the dark for 30 min.
14. After incubation wash the slides with PBS as in steps 10 and 11.
15. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.  
Note: Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip. Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering.
16. Examine reactions under a fluorescence microscope at a total magnification of 400x or at 800x.
17. After examination used slides are to be discarded according to national law.

## Interpretation of Results

Controls have to fit the lot specific reactions stated in the QC certificate. The figures in the table below are for orientation only and do not reflect the data of the current lot.

Specimen	Expected Results*
Positive control + Sorbent	1 + to 4 +
Positive control + PBS	1 + to 4 +
Absorption control + Sorbent	negative**
Absorption control + PBS	1 + to 4 +

\*see measuring range for examples

\*\* negative: Treponemes can be visible but do not fluoresce.

For a correct interpretation the results should be compared with positive and negative controls.

If further dilution on the samples or controls are requested on the basis of the screening titer, this can be achieved by serial dilution with PBS on screening titer.

## Interpretation of Specimen Results

**Positive Samples:** Treponema show a green fluorescence. The fluorescence intensity depends on the concentration of the antibodies.

The intensity can be scored for fluorescence on a scale from 1+ to 4+. For titration a 1+ intensity is counted as minimally reactive.

**Borderline Samples:** Displaying a faint but nevertheless perceptible fluorescence (+/-) and should be repeated.

**Negative Samples:** In most cases no treponemes are visible. Treponemes which are stained green but are not fluorescent must be read NEGATIVE.

In some cases an inhomogeneous pattern is observed in the 1:5 screening dilution. Fluorescence is obtained in some treponemes or in parts of them, other treponemes in the same well are just stained without any fluorescence. These wells should be read as negative.

## Limitations / Interferences

1. The FTA-ABS IgG test may occasionally give a false positive result under certain patient conditions or diseases e.g. pregnancy, leprosy and systemic lupus erythematosus (see references 5, 6).
2. Sera with a high anti-Borrelia antibody titer (> 1:640) may show a green fluorescence of Treponema.
3. The FTA-ABS IgG test has been demonstrated to be highly sensitive and specific, however, results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.
4. Pathogenic treponemes include *T. pallidum*, the cause of syphilis, *T. pertenue*, the cause of yaws and *T. carateum*, the cause of pinta. The only treponeme of importance in Europe is *T. pallidum*. In addition many commensal species of treponemes exist and it is important to differentiate these from *T. pallidum* before a diagnosis of primary syphilis is given.
5. Commensal treponemes can be cultured in artificial media whereas attempts to culture pathogenic treponemes *in vitro* failed.
6. The clinical diagnosis of syphilis is confirmed in the laboratory by either demonstrating the presence of *T. pallidum* in the exudates from the lesions, or demonstrating the presence of serum antibodies against the organism. Methods used to measure antibody response to treponemal infection can be divided into two major categories:
  - a) Tests to measure antibodies against non-specific treponemal antigens i.e. cardiolipin or lipoidal antigen tests. These were formerly called 'Reagin' tests.
  - b) Tests to measure antibodies against antigens specific for pathogenic treponemes i.e. *T. pallidum* particle agglutination assay (TPPA) and the fluorescent treponemal antibody absorbance test (FTA-ABS).

## Performance characteristics

### Sensitivity and Specificity

The assay was validated on a panel consisting of 366 sera and plasma samples. The panel was derived from positive and negative clinical samples. MASTABLOT™ TP IgG (Mast Diagnostica, REF 663G08) and SERODIA®-TP•PA (Fujirebio Inc., Tokyo, REF 226414) were used as reference assays.

The analytical performance was calculated as follows:

	Formula	Value
Sensitivity	$\frac{TP}{(TP + FN)}$	99.52 %
Specificity	$\frac{TN}{(TN + FP)}$	94.63 %
Positive predictive value	$\frac{TP}{(TP + FP)}$	96.30 %
Negative predictive value	$\frac{TN}{(TN + FN)}$	99.30 %
Efficiency	$\frac{(TP + TN)}{total}$	97.49 %
Positive likelihood ratio	$\frac{(Sens)}{(1 - Spec)}$	18.54
Negative likelihood ratio	$\frac{(1 - Sens)}{(Spec)}$	0.01

(Abbreviations: TP: True positive, TN: True negative, FN: False negative, FP: False positive)

### Trueness

Sensitivity and specificity of MASTAFLUOR™ FTA-ABS IgG is in line with performance data of the reference assays.

### Cross-Reactivity

20 *Borrelia* samples were tested for cross-reactivity. Sera with a high anti-*Borrelia* antibody titer sometimes gave false positive results (see Limitations).

### Traceability

NIBSC QC 2, Lot: 17/B713 (Plasma TP 377) was used as a reference sample.

### Assay Precision


When positive control titration (serial dilution 1:5 to 1:640) was tested 10 times from different lots according to the test procedure, all results were found to be within  $\pm 1$  level of fluorescence.

### High dose hook effect (= prozone effect)

The MASTAFLUOR™ FTA-ABS IgG assay is an indirect immunofluorescence detection assays using an FITC-coupled antibody (2-step assay). A washing step has been added in which all unspecific or unbound antibodies are removed. For this reason, the high-dose hook effect is prevented.

## Measuring range

When positive control titration (IgG: 1:5-1:640) and sample dilution (1:5) was tested according to the test procedure, following results for measuring range were documented:

### Cut-off

The cut-off (1:5 dilution, minimum 1+ fluorescence intensity) was defined on a panel consisting of 366 positive and negative sera and plasma samples corresponding to the mentioned sensitivity and specificity.

### Clinical Performance

For the MASTAFLUOR™ FTA-ABS IgG assay, interlaboratory comparisons (INSTAND e.V. Ringversuche, Section 311, *Treponema pallidum*) are performed 1-2 times per year. Data on the proficiency tests are available on request.

The Zeus Scientific FTA assay (REF FA7001) was identified as an equivalence product. The following values of clinical performance were reported for the equivalence product in a comparative study of 303 serum samples (Binnicker et al., 2011):

Clinical sensitivity: 100%

Clinical specificity: 98.6%

k value: 0.98

## **Availability of the summary of safety and performance (Art. 29, IVDR)**

The report will be available in the EUDAMED database (as soon as the module is available).

(<https://ec.europa.eu/tools/eudamed/#/screen/home>).

The SSP report will be provided on request.

## **Reporting serious incidents**

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

## **References and Change History**

You can find the references and the Change History at the end of the instruction for use.

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1. Schöfer H, Brockmeyer NH, Hagedorn HJ, Hamouda O, Handrick W, Krause W, Marcus U et al. (2006): Syphilis. Leitlinie der Deutschen STD Gesellschaft zur Diagnostik und Therapie der Syphilis. *J. Dtsch. Dermatol. Ges.* 4 (2): 160-177.
2. Rath PM, Marsch WC, Brade V, Fehrenbach FJ (1994): Serological distinction between syphilis and Lyme borreliosis. *Zentralbl. Bakteriol.* 280 (3).
3. Hans-Jochen Hagedorn, Syphilis; in Lothar Thomas, Labor und Diagnose 8. Auflage (2012), TH-Books; Kapitel 42.14: 2006 - 2017
4. Annual epidemiological report 2014 - Antimicrobial resistance and healthcare-associated infections; Chapter: Syphilis, p. 51-53; <http://ecdc.europa.eu/>
5. Smikle MF, James OB, Prabhakar P (1990): Biological false positive serological tests for syphilis in the Jamaican population. *Sexually Transmitted Infections*, 66, p. 76-78
6. Wright DJM (1973): The significance of the fluorescent treponemal antibody (FTA-ABS) test in collagen disorders and leprosy *Journal of Clinical Pathology*, 26, p. 968-972
7. Syphilis, in Surveillance Report, Sexual Transmitted Infections in Europe, 2013, p. 25-31; [www.ecdc.europa.eu](http://www.ecdc.europa.eu)
8. Binnicker MJ, Jespersen DJ, Rollins LO (2011): Treponema-Specific Tests for Serodiagnosis of Syphilis: Comparative. Evaluation of Seven Assays. *J. Clin. Microb.*, 49 (4), p. 1313–1317.

Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Cover page	Addition of the number of the Notified Body
Stability and Storage	Change in-use stability from 3 months to 90 days
Limitations / Interferences	Addition of References Smikle et al. (Ref. 5) and Wright et al. (Ref. 6)
Performance Characteristics – Prozonenphänomen (german version)	Addition of the naming High dose hook effect
Performance Characteristics – Clinical performance	Proficiency tests added
Availability of the summary of safety and performance (Art. 29)	Addition of the webpage and sending on request
Change History	Adding a reference to the Change History at the end of the IFU.
References	Addition of new references for cross reactivity (Ref. 5, 6)

**Mast Diagnostica GmbH,**

Feldstraße 20,  
23858 Reinfeld,  
Deutschland

Tel: +49 (0)4533 2007 0

Fax: +49 (0)4533 2007 68

email: mast@mast-diagnostica.de

Web: www.mast-group.com

**Mast Group Ltd.**

Mast House, Derby Road  
Bootle, Merseyside, L20 1EA  
United Kingdom

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**Icons are used according to DIN EN ISO 15223-1**

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