



HPE buffer

Specification sheet

Art.no	M1940
Product specification	High Performance ELISA buffer, five fold concentrated.
Code	HPE
Application	General: ELISA using Poly-HRP-conjugated streptavidin to exclude an increased background. Especially: ELISA for determination of cytokines in plasma or sera.
Packing	Each bottle contains 55 ml five fold concentrated buffer.
Preservative	Merthiolate (0.001%).
Storage	High Performance ELISA buffer should be stored below -18°C to -32°C. Do not use this product beyond the stated expiration date.
Directions for use	Dilute the High Performance ELISA buffer five fold in distilled water.
Stability	After diluting to working strength buffer, it can be stored for 1 week at 2-8°C or 3 months below -18°C to -32°C.