

Mast Group Ltd. Mast House, Derby Road, Bootle Liverpool, Merseyside, L20 1EA United Kingdom

Tel: + 44 (0) 151 472 1444 Fax: + 44 (0) 151 944 1332 email: sales@mast-group.com Web: www.mast-group.com

## Mast Diagnostica GmbH

Feldstrasse 20 DE-23858 Reinfeld Germany

Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 email: mast@mast-diagnostica.de Web: www.mast-group.com

## **Mast Diagnostic**

12 rue Jean-Jacques Mention CS91106, 80011 Amiens, CEDEX 1 France

Tél: + 33 (0) 3 22 80 80 67 Fax: + 33 (0) 3 22 80 99 22 email: info@mast-diagnostic.fr Web: www.mast-group.com



# Acid Hydrolysate of Casein

#### **RM30**

## Intended use

A defined source of amino acids and salts.

#### Contents

See pack label.

## Storage and shelf life

All dehydrated culture media containers should be kept tightly closed and stored in a dry place at 10 to 25°C until the expiry date shown on the pack label.

### **Precautions**

For in vitro diagnostic use only. Observe approved hazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet (available on request or via MAST® website).

Since the late 1980s MAST® has sourced all its animalderived culture media ingredients from non-bovine animals wherever possible. All animal-derived materials, including bovine materials, are from BSE-free regions of the world, are from animals certified as disease-free by qualified veterinarians and have been heat-treated in accordance with European regulations. Despite these precautions, MAST® Media Raw Materials must not be used in the manufacture of vaccines or food ingredients, or in the manufacture of any other high-risk products involving culture processes such as those destined for in-vivo or agricultural use.

## Materials required but not provided

Standard microbiological supplies and equipment such as loops, MAST® selective supplements, swabs, applicator sticks, incinerators and incubators, etc., as well as serological and biochemical reagents and additives such as blood.

# Description

MAST® Acid Hydrolysate of Casein is rich in amino acids and essential mineral salts, but low in vitamin content. This preparation is used in MAST® Mueller Hinton Agar and Broth (DM170D and DM171D). The product is formed by hydrolysis of casein with hydrochloric acid under pressure, followed by neutralisation with sodium hydroxide. The digest, when made up in a 2% solution, is ideally suited for the preparation of chemically defined media.

## **Procedure**

Acid Hydrolysate of Casein should be added at the appropriate concentration to culture media with additional ingredients e.g. agar, meat extracts and peptones as specified in the formulation being prepared.

## Quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate expected performance. Do not use the product if the result with the control organism is incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain.

Test Organisms	Result
Escherichia coli ATCC® 25922	Growth*
Staphylococcus aureus ATCC® 25923	Growth*

\*2% w/v sterile solution Acid Hydrolysate of Casein with 0.5% w/v sodium chloride.

#### References

Bibliography available on request.