

MAST[®] DISCMASTER 5

Antibiotic Disc Dispenser Superior and accurate performance

- Guaranteed in-use disc stability of 4 weeks
- Optimal moisture protection
- Further enhanced to maintain superior performance



MAST[®] DISCMASTER 5

The **NEW MAST[®] DISCMASTER 5** dispenser has been developed to optimise performance and provides a reliable and robust solution for the dispensing of AST, ID and combination discs within the **MASTDISCS[®]** range. Due to its novel and innovative design, the **MAST[®] DISCMASTER 5** dispenser guarantees ***an in-use shelf life of 4 weeks for all antimicrobials including the carbapenem group*** when stored using charged desiccant.

The **MAST[®] DISCMASTER 5** is available in a convenient 6-place format, with an easy to use single-handed operation, permitting smooth integration into routine laboratory workflows and high throughput laboratories.

Guaranteed in-use stability

All **MASTDISCS[®] AST** including the carbapenem group are stable for 4 weeks, ensuring accurate and reliable results and convenience for stock control processes.

Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

Optimal moisture protection

New interlocking seals on the canister, providing improved protection against moisture. Designed to be used with **MAST[®] Silica Gel capsules**, available from stock **SILICA63**.

Positive disc placement

MASTDISCS[®] discs positively placed onto the agar surface permitting precise, centred positioning on each dispensing cycle. Stainless steel base plate and ejector pins provide robust and long lasting performance.

Lock-off indicated empty cartridge

Automatic lock-off, alerting user to an empty cartridge.



Ordering Information

Order Code	Product	Pack Size
MDD65	MAST[®] DISCMASTER 5 Dispenser - 6 place	1
SHD5	Single Cartridge Hand Dispenser	5
SILICA63	Silica Gel Capsule for Discmaster (MDD63 and above)	4
CANISTER65	Canister for MDD65 MAST[®] DISCMASTER Dispenser	1



IVD solutions through partnership

Mast Group Ltd.
Mast House, Derby Road, Bootle, Merseyside L20 1EA
Tel: + 44 (0) 151 933 7277 e-mail: sales@mast-group.com

Mast Diagnostica GmbH
Feldstrasse 20, DE-23858 Reinfeld
Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 e-mail: mast@mast-diagnostica.de

Mast Diagnostic
12 Rue Jean-Jacques Mention, CS 91106, 80011 Amiens CEDEX 1
Tél. + 33 (0) 322 80 80 67 Fax + 33 (0) 322 80 99 22 e-mail: info@mast-diagnostic.fr

