



## The PeliKine-Tool set

The PeliKine-Tool set has been developed to complement Sanquin's successful line of PeliKine-compact ELISA kits. The latter have been designed to combine the ease of a ready-to-use kit with the economy of a do-it-yourself assay. Therefore only a few simple buffer solutions and ELISA reagents need to be prepared.

The PeliKine-Tool set contains the chemicals and solutions to be used in the PeliKine-compact kits: coating and washing buffers, as well as substrate and stop solutions. The Tool set comprises of high quality ready prepared solutions, it is therefore time saving, reduce errors and ensures efficient highly reliable ELISA results.

*The PeliKine-Tool set: to make the most out of your Sanquin PeliKine-compact ELISA kits!*

### CONTENTS OF THE SET

The PeliKine-Tool set contains material sufficient for one PeliKine-compact kit (288 tests), reagents provided are:

Quantity	component	Format	Volume
1 bottle	Coating buffer capsules*	for 3 x 100 ml buffer	3 pcs
1 bottle	PBS tablets	for 3 x 200 ml buffer	3 pcs
2 bottles	Washing buffer	20-fold concentrated	2 x 50 ml
1 bottle	TMB substrate solution	ready-to-use	40 ml
1 bottle	Stop solution (0.18 M H <sub>2</sub> SO <sub>4</sub> )	ready-to-use	40 ml

*\*please note that this coating buffer should not be used for the following PeliKine-compacts™ ELISA kits: IL-6r, sCD27, Granzyme A and Granzyme B. Check your appropriate kit insert for the preparation of these coating buffers.*

Distributed by:

**Mast Group Ltd.**

Mast House, Derby Road, Bootle, Merseyside, L20

1EA, UK

Tel: +44 (0)151 933 7277

Fax: +44 (0)151 944 1332

[www.mastgrp.com](http://www.mastgrp.com)

[sales@mastgrp.com](mailto:sales@mastgrp.com)