

CRYOBANK[®]

Safe and Secure cryopreservation
of bacteria and fungi

**Full
traceability**

- Bar coded for full traceability
- Up to 9 years viable storage
- 25 beads per tube for multiple use
- Unique cryoprotectant
- Colour coded caps and beads

CRYOBANK®

CRYOBANK® is based on a cryovial system comprising chemically treated beads covered with a special cryogenic preserving solution. **CRYOBANK®** offers a reliable, convenient and versatile system for storing and preserving fastidious bacteria over long periods. **CRYOBANK®** permits the easy establishment of culture collections for laboratory accreditation of research purposes. Suitable for use in all laboratories dealing with bacteria including clinical, water, brewing, veterinary, food, pharmaceutical, dairy and biotechnology.

Features and benefits:

- Bar-coded for full traceability suitable for use with all commonly available bar code readers
- Stackable boxes
- Simple to use procedure
- Proven storage of fastidious bacteria, non-fastidious yeasts and moulds over long periods suitable for use down to -80 °C
- Supplied with picker for easy selection of cryovials
- Colour coded beads and caps for easy recognition
- External threaded tube reduce risk of contamination



Ordering information

Order Code	Product	Content
291800	Cryobank® starter set €€	1 Cryobank mixed colours, 1 Cryoblock
291801	Cryobank® - red €€	80 vials
291802	Cryobank® - green €€	80 vials
291803	Cryobank® - yellow €€	80 vials
291804	Cryobank® - blue €€	80 vials
291809	Cryobank® - mixed colours €€	80 vials
29180DRY	Cryobank® - mixed colours dry €€	80 vials
291690	Cryobank® storage box (bar-coded)	1 box for 81 vials
291610	Cryoblock	1 block for 18 vials
291611	Substitute polystyrene box for Cryoblock	1 pc

V 1.0 SJW 1217 (MD2018-03-19)

United Kingdom
Mast Group Ltd.
Mast House, Derby Road, Bootle
Merseyside L20 1EA

Tel: +44 (0)151 933 7277
 Fax: +44 (0)151 944 1332
 e-mail: sales@mastgrp.com

Germany
Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0
 Fax: +49 (0)4533 2007 68
 e-mail: mast@mast-diagnostica.de

France
Mast Diagnostic
12 rue Jean-Jacques Mention
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67
 Fax: +33 (0)322 80 99 22
 e-mail: info@mast-diagnostic.fr