

# Automaten-Adaption

	Reagenz 1 [µL]	Reagenz 2 [µL]	Probe [µL]	Anzahl Best.
Bayer Advia 1650	50	50	6	800
Hitachi 7700 R&D	140	140	2,8	285
Hitachi 911/917	150	150	3	265
Horiba Pentra 400	150	150	4	265
Konelab	100	100	2	400
Manuell	250	250	5	160
Olympus AU400/ AU640/AU2700/AU5400	150	150	3	265
Olympus AU5800	100 / 75	100 / 75	2 / 1,5	400 / 530
Randox Daytona	150	150	3	265
Roche Cobas	100	100	2	400
Roche Mod. P	100	100	2	400
Synermed	250	250	5	160