



Safety Data Sheet

Rev.01 Jan.06

Page 1 (4)

1. Product identification

Product name:

VIRO-SWAB T/LIQ.PLAS.

Product number:

147C

Supplier:

*COPAN ITALIA SPA 25125 Brescia -Italy- Via Perotti,10
Tel. ++ 39 0302687211 / Fax ++ 39 030 2687206*

2. Composition/information on ingredients

Hazardous ingredients:

Percentage on medium:

Classification:

Calcium Chloride

0,013%

R36

Sodium thioglycollate

0,05%

R21\22

3. Hazards identification

Skin contact:

Irritation

Eye contact:

Irritation

4. First-aid measures

Skin contact:

Immediately wash with water and soap and rinse thoroughly.

Eye contact:

Rinse opened eye for 15 minutes under running water. Then consult a doctor.

Inhalation:

In case of unconsciousness place patient on side position for transportation.

Swallowing:

Immediately call a doctor.

5. Fire-fighting measures



Safety Data Sheet

Page 2 (4)

Rev.01 Jan.06

Extinguishing media:	CO ₂ , extinguishing powder or water spray. Fight larger fires with water or alcohol resistant foam
Protective equipment:	Wear self-contained respiratory protective device Wear fully protective suit
6. Accidental release measures	
Precautions:	Wear protective equipment. Keep unprotected persons away
7. Handling and storage	
Handling:	No special measures required
Storage:	Store away from oxidizing agents Store in cool, dry conditions in well sealed receptacles.
8. Exposure controls/personal protection	
Respiratory protection:	N/A
Hand protection:	Latex gloves
Eye protection:	N/A
Skin protection:	Protective work clothing (lab coat)
9. Physical and chemical properties	
Appearance/odour:	Characteristic
pH:	7.00 ± 0.50
Density, kg/l:	Not determined
Melting point, °C:	Not determined
Boiling point, °C:	Not determined
Solubility:	Insoluble
Flash point, °C:	Not applicable
10. Stability and reactivity	
Materials and conditions to avoid:	No decomposition if used according to specifications
Hazardous decomposition products:	No dangerous decomposition products known



Safety Data Sheet

Page 3 (4)

Rev.01 Jan.06

11. Toxicological information		
Skin contact:	Irritant to skin and mucous membranes	
Eye contact:	Irritating effect	
Toxicological data:	Irritant	
	The products shows the following dangers according to internally approved calculation methods for preparations	
12. Ecological information		
Ecology:	The ecological effects have not been thoroughly investigated, but currently none have been identified.	
General Notes:	Not known to be hazardous to water.	
13. Disposal considerations		
Reccomendation:	Must adhere to state and federal regulations	
	Must not be disposed of with solid waste	
14. Transport information		
ADR/RID:	IMO/IMDG:	ICAO/IATA-DGR:
No	No	No
15. Regulatory information		
Product related hazard informations:	The product has been classified and marked in accordance with regulations on hazardous materials	
16. Other information		
<p><i>To the best of our knowledge, the information contained herein is accurate. However, neither Copan Italia or any of its subsidiaries assumes any liabilities whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user.</i></p> <p><i>All materials may present unknown hazards and should be used with caution.</i></p> <p><i>Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist.</i></p> <p><i>· Department issuing MSDS: Safety & Environment Department</i></p> <p><i>· Contact: Technical Service Representative</i></p>		

Safety Data Sheet
Rev.01 Jan.06



Page 4 (4)

The information in this safety data sheet has been obtained from current and reliable sources. However, the data is provided without warranty, express or implied, regarding its correctness or accuracy. It is the users' responsibility to determine safe conditions for use of this product and to assume liability for loss, injury, damage, or expense resulting from improper use of this product.
EN-2141-SDS-0198