

IVD solutions through partnership



Monoklonale Blutgruppenserien
AB0-Testserien / Rh-Untergruppen

Klonliste

Klonliste

AB0-Testseren / Rh-Untergruppen

Produkt	Best.-Nr.	Klon	Immunglobulin-Klasse	
Anti-A	5360103	BIRMA-1	IgM	CE 0344
Anti-B	5370103	LB-2	IgM	CE 0344
Anti-AB	5380103	ES-4 ES-15	IgM IgM	CE 0344
Anti-H	539010	10934C11	IgM	CE
Anti-D	521020	RUM-1	IgM	CE 0344
Anti-D mix (IgG/IgM)	521011	MS-26 TH-28	IgG IgM	CE 0344
Anti-C	520005	MS-24	IgM	CE 0344
Anti-C	520007	MS-273	IgM	CE 0344
Anti-c	520006	MS-33	IgM	CE 0344
Anti-c	520008	MS-35	IgM	CE 0344
Anti-E	520505	MS-258 MS-80	IgM IgM	CE 0344
Anti-E	520506	MS-260	IgM	CE 0344
Anti-e	520305	MS-21 MS-63	IgM IgM	CE 0344
Anti-e	520306	MS-62 MS-69	IgM IgM	CE 0344
Anti-C ^w	520512	MS-110	IgM	CE 0483
Anti-C ^w	520512SA	MS-110	IgM	CE

Weitere Testseren

Produkt	Best.-Nr.	Klon	Immunglobulin-Klasse	
Anti-Kell	530405	MS-56	IgM	CE 0483
Anti-k (Cellano)	520502	Lk1	IgM	CE
Anti-Le ^a	521502	LEA2	IgM	CE
Anti-Le ^b	521322	P3F234MD4	IgA	CE
Anti-M	520322	LM-110/140	IgG	CE
Anti-N	520102	1422C7	IgM	CE
Anti-N	520122	B03	IgG	CE
Anti-S	522502	MS-94	IgM	CE
Anti-S	522502SA	MS-94	IgM	CE
Anti-Fy ^a	522402	P3TIM	IgG	CE 0344
Anti-Jk ^a	522202	MS-15	IgM	CE 0344
Anti-Jk ^b	522302	MS-8	IgM	CE 0344
Anti-P1	521012	P3NIL100	IgM	CE
Anti-C ₃ d	545715	BRIC8	IgG1	CE

MD V.2.0 / 2022-02-17

United Kingdom
Mast Group Ltd.
Mast House, Derby Road, Bootle
Merseyside L20 1EA

Tel: +44 (0)151 933 7277
Fax: +44 (0)151 944 1332
e-mail: sales@mastgrp.com

Germany
Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
e-mail: mast@mast-diagnostica.de

France
Mast Diagnostic
12 rue Jean-Jacques Mention
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67
Fax: +33 (0)322 80 99 22
e-mail: info@mast-diagnostic.fr