



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	2U026S
GENERAL DESCRIPTION	VAC.CONT.12X80 RED CAP 1,5ML
DETAILED DESCRIPTION	12 envelopes each containing 50 vacuum tubes 1.5 ml with red screw cap (with rubber pad)

PRODUCT COMPONENTS

CAP	RED CAP	DIMENSIONAL SPECIFICATIONS
		EXTERNAL DIAMETER: 17.00 mm \pm 0.20 mm CAP LENGTH: 16.25 mm \pm 0.20 mm
PAD	PAD	DIMENSIONAL SPECIFICATIONS
		DIAMETER: 14.15 mm \pm 0.15 mm THICKNESS: 4.6 mm \pm 0.2 mm

This bulletin (Rev.00-28/04/2016) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

SAMPLING TUBE	TUBE	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm \pm 0.20 mm INTERNAL DIAMETER: 10.70 mm \pm 0.20 mm EXTERNAL DIAMETER: 13.40 mm \pm 0.20 mm

PACKAGING SPECIFICATION

PACKAGING	BULK
BOX QTY	600
ENVELOPE QTY	50

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	12 MONTHS