

MAST[®] DISCMASTER

Antibiotic Disc Dispenser Range
Optimal and superior performance



- 6 and 8-place suite of dispensers

- Guaranteed in-use disc stability of 4 weeks

- Optimal moisture protection

- Smooth, robust and superior performance

NEW!
8-place suite
now
available!

MAST® DISCMASTER

The **MAST®** DISCMASTER dispensers have been developed to optimise performance and provides a reliable and robust solution for the dispensing of Antimicrobial Susceptibility Test (AST), Identification (ID), and Antimicrobial Resistance (AMR) cartridge discs within the **MASTDISCS®** range. Due to its novel and innovative design, the **MAST®** DISCMASTER dispenser guarantees an in-use shelf life of 4 weeks for all antimicrobials including the carbapenem group when stored using charged desiccant.

The **MAST®** DISCMASTER range is available in both 6-place and 8-place formats, with an easy to use single-handed operation, permitting smooth integration into routine laboratory workflows and high throughput laboratories.

Guaranteed in-use stability

All **MASTDISCS®** including the carbapenem group are stable for 4 weeks, ensuring accurate and reliable results and convenience for stock control processes.

Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

Optimal moisture protection

Interlocking seals on the canister, providing improved protection against moisture. Designed to be used with **MAST®** Silica Gel capsules, available from stock (SILICA63).

Positive disc placement

MASTDISCS® positively placed onto the agar surface permitting precise, centred positioning on each dispensing cycle. Stainless steel base plate and ejector pins provide robust and long lasting performance.

Lock-off indicated empty cartridge

Automatic lock-off, alerting user to an empty cartridge.



Ordering Information

Order code	Product	Pack Size
773065	MAST® DISCMASTER Dispenser - 6 place (incl. storage container)	1
773080	MAST® DISCMASTER Dispenser - 8 place (incl. storage container)	1
773001	Single applicator for discs in magazines	5
760000	Desiccant for MAST® DISCMASTER Dispenser	2 x 2

ZT170 v.4.0 CEA 10/23 (V.1.1/MD-2024-08-26)

United Kingdom
Mast Group Ltd.
Mast House, Derby Road, Bootle
Merseyside L20 1EA

Tel: +44 (0)151 933 7277
Fax: +44 (0)151 944 1332
e-mail: sales@mastgrp.com

Germany
Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
e-mail: mast@mast-diagnostica.de

France
Mast Diagnostic
12 rue Jean-Jacques Mention
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67
Fax: +33 (0)322 80 99 22
e-mail: info@mast-diagnostic.fr

www.mast-group.com