



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	1U054S01
GENERAL DESCRIPTION	PLASTIC APPL.POLYESTER TIPPED SINGLE WRAPPED BP100
DETAILED DESCRIPTION	1 outer box contain 10 inner boxes each composed of 100 swabs individually packaged in paper-film pouch.

PRODUCT COMPONENTS

		SPECIFICATIONS
APPLICATOR	PLASTIC APPLICATOR	STICK DIAMETER: 2.50 mm ± 0.10 mm FIBRE TIP LENGTH: 15.00 mm ± 2.00 mm FIBRE TIP DIAMETER: 5.00 mm ± 0.50 mm TOTAL LENGTH: 144.00 mm ± 2.00 mm BREAK POINT: 98.00 mm ± 0.50 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.00-24/03/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.