

# MAST<sup>®</sup> DISCMASTER

## Antibiotic Disc Dispenser Range Optimal and superior performance

- 6 and 8-place suite of dispensers
- Guaranteed in-use disc stability of 4 weeks
- Optimal moisture protection
- Smooth, robust and superior performance



New  
8-place  
dispenser

# MAST<sup>®</sup> DISCMASTER

The **MAST<sup>®</sup>** DISCMASTER dispensers have been developed to optimise performance and provides a reliable and robust solution for the dispensing of Antimicrobial Susceptibility Test (AST), Identification (ID), and Antimicrobial Resistance (AMR) cartridge discs within the **MASTDISCS<sup>®</sup>** range. Due to its novel and innovative design, the **MAST<sup>®</sup>** DISCMASTER dispenser guarantees ***an in-use shelf life of 4 weeks for all antimicrobials including the carbapenem group*** when stored using charged desiccant.

The **MAST<sup>®</sup>** DISCMASTER range is available in both 6-place and 8-place formats, with an easy to use single-handed operation, permitting smooth integration into routine laboratory workflows and high throughput laboratories.

## Guaranteed in-use stability

All **MASTDISCS<sup>®</sup>** including the carbapenem group are stable for 4 weeks, ensuring accurate and reliable results and convenience for stock control processes.

## Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

## Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

## Optimal moisture protection

Interlocking seals on the canister, providing improved protection against moisture. Designed to be used with **MAST<sup>®</sup>** Silica Gel capsules, available from stock (SILICA63).

## Positive disc placement

**MASTDISCS<sup>®</sup>** positively placed onto the agar surface permitting precise, centred positioning on each dispensing cycle. Stainless steel base plate and ejector pins provide robust and long lasting performance.

## Lock-off indicated empty cartridge

Automatic lock-off, alerting user to an empty cartridge.



## Ordering Information

Order Code	Product	Pack Size
MDD65	<b>MAST<sup>®</sup></b> DISCMASTER 6-place Dispenser	1
MDD8	<b>MAST<sup>®</sup></b> DISCMASTER 8-place Dispenser	1
SHD5	Single Cartridge Hand Dispenser	5
SILICA63	Silica Gel Capsule for Discmaster	4
CANISTER65	Canister for <b>MAST<sup>®</sup></b> DISCMASTER Dispenser	1



IVD solutions through partnership

Mast Group Ltd.  
Mast House, Derby Road, Bootle, Merseyside L20 1EA  
Tel: + 44 (0) 151 933 7277 e-mail: sales@mast-group.com

Mast Diagnostica GmbH  
Feldstrasse 20, DE-23858 Reinfeld  
Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 e-mail: mast@mast-diagnostica.de

Mast Diagnostic  
12 Rue Jean-Jacques Mention, CS 91106, 80011 Amiens CEDEX 1  
Tél. + 33 (0) 322 80 80 67 Fax + 33 (0) 322 80 99 22 e-mail: info@mast-diagnostic.fr

