

# MAST<sup>®</sup> DISCMASTER

**Antibiotic Disc Dispenser Range**  
Optimal and superior performance



**NEW!**

8-place suite  
available  
soon!

- 6 and 8-place suite of dispensers
- Guaranteed in-use disc stability of 4 weeks
- Optimal moisture protection
- Smooth, robust and superior performance

# MAST® DISCMASTER

The **MAST®** DISCMASTER dispensers have been developed to optimise performance and provides a reliable and robust solution for the dispensing of Antimicrobial Susceptibility Test (AST), Identification (ID), and Antimicrobial Resistance (AMR) cartridge discs within the **MASTDISCS®** range. Due to its novel and innovative design, the **MAST®** DISCMASTER dispenser guarantees an in-use shelf life of 4 weeks for all antimicrobials including the carbapenem group when stored using charged desiccant.

The **MAST®** DISCMASTER range is available in both 6-place and 8-place formats, with an easy to use single-handed operation, permitting smooth integration into routine laboratory workflows and high throughput laboratories.

## Guaranteed in-use stability

All **MASTDISCS®** including the carbapenem group are stable for 4 weeks, ensuring accurate and reliable results and convenience for stock control processes.

## Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

## Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

## Optimal moisture protection

Interlocking seals on the canister, providing improved protection against moisture. Designed to be used with **MAST®** Silica Gel capsules, available from stock (SILICA63).

## Positive disc placement

**MASTDISCS®** positively placed onto the agar surface permitting precise, centred positioning on each dispensing cycle. Stainless steel base plate and ejector pins provide robust and long lasting performance.

## Lock-off indicated empty cartridge

Automatic lock-off, alerting user to an empty cartridge.



## Ordering Information

Order code	Product	Pack Size
773065	<b>MAST®</b> DISCMASTER Dispenser - 6 place (incl. storage container)	1
773080	<b>MAST®</b> DISCMASTER Dispenser - 8 place (incl. storage container)	1
773001	Single applicator for discs in magazines	5
760000	Desiccant for <b>MAST®</b> DISCMASTER Dispenser	2 x 2

ZT170 v.4.0 CEA 10/23 (V.1.0/MD-2024-05-06)

**United Kingdom**  
**Mast Group Ltd.**  
Mast House, Derby Road, Bootle  
Merseyside L20 1EA

Tel: +44 (0)151 933 7277  
Fax: +44 (0)151 944 1332  
e-mail: sales@mastgrp.com

**Germany**  
**Mast Diagnostica GmbH**  
Feldstraße 20  
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0  
Fax: +49 (0)4533 2007 68  
e-mail: mast@mast-diagnostica.de

**France**  
**Mast Diagnostic**  
12 rue Jean-Jacques Mention  
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67  
Fax: +33 (0)322 80 99 22  
e-mail: info@mast-diagnostic.fr

[www.mast-group.com](http://www.mast-group.com)