



MAST[®] DISCMASTER 5

Antibiotic Disc Dispenser
Superior and accurate performance



- In-use disc stability of 4 weeks

- Optimal moisture protection

- further enhanced to maintain superior performance

MAST® DISCMASTER 5

The NEW **MAST®** DISCMASTER 5 dispenser has been developed to optimise performance and provides a reliable and robust solution for the dispensing of AST, ID and combination discs within the **MASTDISCS®** range. Due to its novel and innovative design, the **MAST®** DISCMASTER 5 dispenser guarantees ***an in-use shelf life of 4 weeks for all antimicrobials including the carbapenem group*** when stored using charged desiccant.

The **MAST®** DISCMASTER 5 is available in a convenient 6-place format, with an easy to use single-handed operation, permitting smooth integration into routine laboratory workflows and high throughput laboratories.

Guaranteed in-use stability

All **MASTDISCS®** AST including the carbapenem group are stable for 4 weeks, ensuring accurate and reliable results and convenience for stock control processes.

Lock-off indicated empty cartridge

Automatic lock-off, alerting user to an empty cartridge.

Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

Optimal moisture protection

New interlocking seals on the canister, providing improved protection against moisture. Designed to be used with **MAST®** Silica Gel capsules, available from stock SILICA63.

Positive disc placement

MASTDISCS® discs positively placed onto the agar surface permitting precise, centred positioning on each dispensing cycle. Stainless steel base plate and ejector pins provide robust and long lasting performance.



Ordering Information

Order code	Product	Pack Size
773065	MAST® DISCMASTER 5 Dispenser - 6 place (incl. storage container)	1
773001	Single applicator for discs in magazines	5
760000	Desiccant for MAST® DISCMASTER Dispenser	2 x 2

V. 2.0 CEA 11/18 (MD2019-02-04)

United Kingdom
Mast Group Ltd.
Mast House, Derby Road, Bootle
Merseyside L20 1EA

Tel: +44 (0)151 933 7277
Fax: +44 (0)151 944 1332
e-mail: sales@mastgrp.com

Germany
Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
e-mail: mast@mast-diagnostica.de

France
Mast Diagnostic
12 rue Jean-Jacques Mention
CS91106, 80011 Amiens, CEDEX 1

Tél: +33 (0)322 80 80 67
Fax: +33 (0)322 80 99 22
e-mail: info@mast-diagnostic.fr

www.mast-group.com